

Protocol for REAP nuclear/cytoplasmic fractionation

REAGENTS

NP40 (Calbiochem, San Diego, CA, USA)

EQUIPMENT

Table-top centrifuge (Eppendorf, Barkhausenweg, Hamburg, Germany)

P1000 micropipette (Gilson, Middleton, WI, USA)

Micropipet tip for p1000 micropipet (Axygen, Union City, CA, USA)

1.5ml microcentrifuge tube (Axygen, Union City, CA, USA)

Cell scraper (Corning, Middleton, WI, USA)

Sonicator with microprobes (Misonix, Farmingdale NY, USA)

PROCEDURE

[Chill all reagents on ice.

Cell collection

1. Remove all culture medium from cell culture dish
2. Wash cell surface twice with ice-cold PBS
3. Add ice-cold PBS into cell culture dish (1 ml per 10 cm diameter dish)
4. Scrape cells with cell scraper
5. Collect cells with PBS into 1.5 ml microcentrifuge tube
6. Pop-Spin the tube for 10 sec. with table top centrifuge
7. Discard supernatant

[Keep cell dishes and tubes on ice

Fractionation: 0.1% NP40-PBS treatment

8. Triturate (pipet up & down) cell pellet 5 times with ice-cold 0.1% NP40-PBS (900-1000 μ L for cells from 10 cm diameter dish) using p1000 micropipette cut about 3 mm off at the end to enlarge the opening.

[Remove an aliquot into a fresh tube as the whole cell sample and keep on ice

9. Pop-Spin the remainder in the tube for 10 sec. with table top centrifuge
10. Transfer the supernatant to new tube and keep it on ice

[This is the cytoplasmic fraction.

11. Resuspend the pellet with ice-cold 0.1% NP40-PBS (1 ml for cells from 10cm diameter dish, triturate once)
12. Pop-Spin the tube for 10 sec. with table top centrifuge
13. Discard supernatant

[This is the nuclear pellet which should be white compared to the yellowish tone of the washed whole cell pellet at step 7.

Sample preparation for SDS-PAGE

- Whole cell lysate (from step 8)
 14. Mix 3 vol. of samples and 1 vol. of 4x Laemmli sample buffer
 15. Sonicate using microprobes at level 2, twice for 5 sec each on ice
 16. Boil for 1-4 min.
- Cytoplasmic fraction (from step 10)
 14. Mix 3 vol. of samples and 1 vol. of 4x Laemmli sample buffer
 15. Boil for 1-4 min.
- Nuclear fraction (from step 13)
 14. Resuspend the pellet with 1x Laemmli sample buffer (~200 μ L for pellet from 10cm diameter dish)
 15. Sonicate using microprobes at level 2, twice for 5 sec. each on ice
 16. Boil for 1-4 min.

Sample preparation for immunoprecipitation

- Whole cell lysate (from step 8)
 14. Sonicate using microprobes at level 2, twice for 5 sec each on ice
 15. Centrifuge at 8000 x g for 30 sec. at 4 °C, transfer supernatant to a fresh tube
 16. The supernatant can be immunoprecipitated directly, with or without addition of protease inhibitors and/or an equal volume of 2x RIPA buffer
- Cytoplasmic fraction (from step 10)
 14. Centrifuge at 8000 x g for 30 sec. at 4 °C, transfer to a fresh tube
 15. The supernatant can be immunoprecipitated directly, with or without addition of protease inhibitors and/or an equal volume of 2x RIPA buffer
- Nuclear fraction (from step 13)
 14. Resuspend the pellet with ice-cold 0.1%NP40-PBS (~500 μ L for pellet from 10 cm diameter dish)
 15. Sonicate using microprobes at level 2, twice for 5 sec each on ice
 16. Centrifuge at 8000 x g for 30 sec. at 4 °C, transfer to a fresh tube
 17. The supernatant can be immunoprecipitated directly, with or without addition of protease inhibitors and/or an equal volume of 2x RIPA buffer

RECIPES

10x PBS (pH7.4) (total volume: 1L)

Ingredient	Volume (g)	Final concentration (mM)
NaCl	80	1370
KCl	2	27

Na ₂ HPO ₄ ·12H ₂ O	29	80
KH ₂ PO ₄	2	15

0.1% NP40-PBS (total volume: 100mL)

10mL of 10x PBS

1mL of 10% (Vol. / Vol.) NP40

make vol. up to 100ml

4x Laemmli sample buffer (total volume: 10mL)

Ingredient	Volume	Final concentration
0.25M Tris-HCl (pH6.8)	10 ml	0.25%
β Mercaptoethanol	2 ml	20%
SDS	0.8 g	8%
Glycerol	2 ml	20%
Bromophenol Blue	0.8 mg	0.008%

2x IP (immunoprecipitation) buffer (total volume: 1L)

Ingredient	Volume	Final concentration (mM)
0.5 M Tris-HCl (pH 7.4)	20 ml	10
1 M NaCl	150 ml	150
1 M KCl	10 ml	10
1 M EDTA	1 ml	1
10% (V/V) NP-40	10 ml	0.1%
1 mg/ml pepstatin	1 ml	1 mg/ml
1 mg/ml aprotinin	1 ml	1 mg/ml
1 mg/ml leupeptin	1 ml	1 mg/ml