

Determination of protein concentration using a micro- Kjeldahl procedure¹

Reagents

(1) Digestion mixture: **CAUTION!** These reagents are highly corrosive and toxic.

0.4 ml selenium oxychloride (Aldrich # 336157 <http://www.sigmaaldrich.com/>)

8 gm potassium sulfate (Sigma-Aldrich P0772 <http://www.sigmaaldrich.com/>)

50 ml concentrated sulfuric acid (Fisher A468-250 (<https://www1.fishersci.com/>))

and water to 250 ml total

(2) Standard: Solutions of ammonium chloride (Acros AC42328-0010

<https://www1.fishersci.com/>) containing between 10-80 micrograms of nitrogen in 2 ml glass distilled or deionized water.

(3) hydrogen peroxide (30%) (Sigma Aldrich H1009).

(4) Color Reagent:

Nessler's Reagent (Fisher SN16I-500 250 <https://www1.fishersci.com/>)

PROCEDURE

- 1 Digest protein samples containing 20-80 micrograms of nitrogen with 0.5 ml of digestion mixture in 15 ml heavy walled Pyrex graduated tubes (Aldrich CLS808015) at 100-200 °C in a sand bath (tray filled with sand on a heavy duty hot plate) until the solvent evaporates.
- 2 Raise the temperature to 300-315 °C for 1-5 hours until the solution becomes colorless (the material will first turn, brown, then green, then clear).
- 3 Cool to room temperature and add approximately 0.2 ml of hydrogen peroxide (30%) (Sigma Aldrich H1009).

- 4 Heat the samples for an additional hour at 300 °C.
- 5 Cool to room temperature and bulk to 10 ml with glass distilled water.
- 6 Take two 1-ml aliquots for analysis. Add 1 ml of water and 1 ml of alkaline Nessler's Reagent (Fisher).
- 7 Let the color intensity develop for 10 minutes and measure the intensity at 420 nm.
- 8 The concentration of the protein is determined by calculating the percentage nitrogen from its amino acid composition.

REFERENCE

1. Lang, C. A.(1958) Simple microdetermination of Kjeldahl nitrogen in biological material. *Anal. Chem.* 30, 1692-1694.