

MICRO-C FOR MAMMALIAN CELLS (2018-v0p99)

I. PREPARE CROSSLINKED CHROMATIN FROM CELL CULTURE

1. Culture cells in recommended condition.
2. Harvest cells by centrifugation for 5 min at 800xg at room temperature.
3. Resuspend cells in **base media (w/o FBS)** in a concentration of **1x10⁶ cells/mL** (max. 30mL in 50mL tube).
4. Add **16% Formaldehyde** to a final concentration of **1%** (e.g. 2mL for 30mL sample).
Incubate for 10 min at room temperature with mixing.
5. Add **2M Tris pH7.5** to a final concentration of **0.75M** to quench the reaction (e.g. 18.75mL for 30mL sample).
Incubate for 5 min at room temperature.
Centrifuge for 5 min at >800xg at room temperature. Discard supernatant.
6. Wash cells twice by 1X PBS in a concentration of 1x10⁶ cells/mL.
Centrifuge for 5 min at >800xg at room temperature. Discard supernatant.
7. Freshly prepare long crosslinkers solution as described in table:

Crosslinkers	MW	Spacer (Å)	Stock	Working
DSG	326.26	7.7	300mM in DMSO	3.0mM in PBS
EGS	456.36	16.1	300mM in DMSO	3.0mM in PBS

Ex: 50mg DSG + 510.84 DMSO = 300mM DSG → +50mL PBS = 3mM DSG

8. Resuspend cell pellet in long crosslinker solution in a concentration of 1x10⁶ cells/mL.
Incubate for 45 min at room temperature with mixing.
9. Add **2M Tris pH7.5** to a final concentration of **0.75M** to quench the reaction.
Incubate for 5 min at room temperature.
Centrifuge for 5 min at >800xg at room temperature. Discard supernatant.
10. Wash cells twice by 1X PBS in a concentration of 1x10⁶ cells/mL.
Centrifuge for 5 min at >800xg at room temperature. Discard supernatant.
11. Snap freeze cell pellets by liquid nitrogen.

II. DIGEST CROSSLINKED CHROMATIN BY MICROCOCCAL NUCLEASE

1. Freshly complete **MB#1**:

Total	10 mL	5 mL	2.5 mL	Final
MB#1	9700 µL	4850 µL	2425 µL	50mM NaCl, 10mM Tris, 5mM MgCl ₂ , 1mM CaCl ₂ , 0.2% NP-40, 1X PIC
10% NP-40	200 µL	100 µL	50 µL	
100X PIC	100 µL	50 µL	25 µL	

2. Resuspend cell pellet in complete MB#1 in a concentration of 1x10⁶ cells/100µL.
Incubate for 20 min on ice.
Centrifuge for 5 min at ~10000xg for at 4°C. Discard supernatant.
3. Wash nuclei pellet in complete MB#1 in a concentration of 1x10⁶ cells/100µL.
Centrifuge for 5 min at ~10000xg for at 4°C. Discard supernatant.
4. Resuspend nuclei pellet in complete MB#1 in a concentration of 1x10⁶ cells/100µL.
5. Add appropriate amount of **MNase** to digest chromatin to **90% Monomer/10% Dimer**.
Incubate for 10 min at 37°C with shaking at 850 rpm.
6. Add **500mM EGTA** in a final concentration at **4mM** to stop reaction.
Incubate for 10 min at 65°C.
Centrifuge for 5 min at >10,000xg at 4°C. Discard supernatant.
7. Wash/Rinse nuclei pellet in 1mL of cold **MB#2** twice.
Centrifuge for 5 min at >10,000xg at 4°C. Discard supernatant.

III. REPAIR FRAGMENT ENDS

A. END-CHEWING			
Total	100 μ L	50 μ L	Final
Chromatin from 1-5 x 10 ⁶ cells	pellet	pellet	50mM NaCl, 10mM Tris, 10mM MgCl ₂ , 100 μ g/mL BSA 2mM ATP 5mM DTT
10X NEBuffer 2.1	10	5	
100mM ATP	2	1	
100mM DTT	5	2.5	
H ₂ O	68	34	
10U/μL T4 PNK	5	2.5	~10U/1 μ g DNA (NEB End repair module: 5U for 1 μ g DNA)
→ Incubate for 15 min at 37°C.			
5U/μL Klenow Fragment	10	5	~10U/1 μ g DNA (NEB End repair module: xxU for 1 μ g DNA)
→ Incubate for 15 min at 37°C.			

B. END-LABELING			
Total	150 μL	75 μL	Final
Chromatin	100 μL	100 μL	66mM dNTP / each 33mM NaCl, 23mM Tris, 10mM MgCl ₂ , 100μg/mL BSA 1.67mM ATP 6.67mM DTT
1mM Biotin-dATP	10	5	
1mM Biotin-dCTP	10	5	
10mM dTTP + dGTP	1	0.5	
10X T4 DNA Ligase Buffer	5	2.5	
20mg/mL BSA (200X)	0.25	0.125	
H ₂ O	23.75	11.875	
→ Incubate for 45 min at 25°C with interval mixing.			
→ Add 500mM EDTA in a final concentration at 30mM → 65°C for 20 min.			
→ Centrifuge for 5 min at > 10,000xg at 4°C. Discard supernatant.			
→ Rinse once with 1mL of cold MB#3 .			
→ Centrifuge for 5 min at > 10,000xg at 4°C. Discard supernatant.			

IV. PROXIMITY LIGATION & PURGE UNLIGATED ENDS

A. LIGATION			
Total	500 μ L	250 μ L	Final
Chromatin	pellet	pellet	50mM Tris, 10mM MgCl ₂ , 1mM ATP, 10mM DTT 100 μ g/mL BSA
Water	422.5	211.25	
10X T4 DNA ligase buffer <i>w/ATP</i>	50	25	
20mg/mL BSA (200X)	2.5	1.25	
400U/μL T4 DNA ligase	25	12.5	XXU/1 μ g DNA (optimal: ?)
→ Incubate for > 2.5 hours at room temperature with slow rotation.			
→ Centrifuge for 5 min at >16000xg at 4°C. Discard supernatant.			

B. REMOVE BIOTIN-dNTP FROM UNLIGATED ENDS			
Total	200 μ L	100 μ L	Final
Chromatin	pellet	pellet	10mM Bis-Tris-Propane-HCl, 10mM MgCl ₂ , 1mM DTT
10X NEBuffer#1	20	10	
Water	170	85	
100U/μL Exonuclease III	10	5	XXU/1 μ g DNA (optimal: ?)
→ Incubate at 37°C for 15 min with interval mixing.			

C. REVERSE CROSSLINKING			
Total	250 μ L	125 μ L	Final

Chromatin	200 μ L	100 μ L	-
80X Proteinase K	25	12.5	2mg/mL Proteinase K = 2X
10% SDS solution	25	12.5	1X SDS
→ Incubate at 65°C for overnight.			

V. DINUCLEOSOMAL DNA PURIFICATION

1. Phenol:Chloroform:Isoamylalcohol (PCI) extraction:

→ Add 1X volume of PCI → Vortex for 20 sec → Spin for 10min at 19800xg at room temperature → Keep upper layer.

2. Purify DNA by ethanol precipitation:

→ Add 0.1X volume of sodium acetate and 2.5X volume of 100% ethanol → Incubate for >1 hr at -80°C → Spin for 15min at 19800xg at 4°C → Wash pellet by 75% ethanol → Spin for 5 min at 19800xg at 4°C → Air dry pellet for 10min at room temperature → Resuspend pellet in 50 μ L of TE buffer (+1X RNase A) → Incubate for >30min at 37°C → ZymoClean.

3. Size-selection for di-nucleosomal DNA:

→ Separate **Monomer** and **Dimer** by running on **3.5% TAE / 3% TBE NuSieve agarose gel** → Cut band at size **250–400 bp (Avoid cutting < 200 bp to eliminate monomer)** → ZymoGel purification and resuspend DNA in 18 μ L of elution buffer → Quantify DNA by Qubit.

VI LIBRARY PREPARATION

1. END-IT

Total	50 μ L	
Input DNA	34	
10X End-it buffer	5	
10X ATP	5	
10X dNTP	5	
End-it enzyme mix	1	
→ 45 min @ 25°C		
→ 10 min @ 70°C		

2. STREPTAVIDIN PURIFICATION

- Wash 5 μ L beads/sample by 1X TBW.
Note: 10 μ g/1 μ L beads can bind 200 ng dsDNA or 5 pmol ssDNA. Input DNA is usually ~2 – 5 μ g DNA (20 – 100 pmol). 5 μ L can bind up to 1 μ g dsDNA or 100 pmol ssDNA.
- Resuspend in 150 μ L of 2X BW.
- Mix with 150 μ L of DNA sample for 20min at room temperature.
- Wash by 1X TBW @ 55°C w/ interval mixing **x2**.
- Wash by 10mM Tris.
- Resuspend in 25 μ L of EB buffer.

3. END REPAIR / A-TAILING

Total	30 μ L	
Input DNA	25	
End Prep Reaction Buffer	3.5	
End Prep Enzyme Mix	1.5	
→ 30 min @ 20°C w/ interval mixing		
→ 30 min @ 65°C		

4. ADPATER LIGATION

Total	46 μ L	
Input DNA	30	
Adapter for Illumina	0.5	
Ligation Master Mix	15	
Ligation Enhancer	0.5	
→ 30 min @ 20°C w/ interval mixing		
USER enzyme	1.5	
→ 15 min @ 37°C w/ interval mixing		

5. BEAD WASH

- Wash by 1X TBW @ 55°C w/ interval mixing.
- Wash by 10mM Tris.
- Resuspend in 24 μ L of EB buffer.

6. MINIMUM PCR & PURIFICATION

Total	50 μ L	
Streptavidin beads w/ EB	24	
2X KAPA HiFi Hot Start Mix	25	
10 μ M PE1.0 primer	0.5	
10 μ M PE2.0 primer	0.5	
Denaturation	98°C	45sec
8-12 cycles	98°C	15sec
	60°C	30sec
	72°C	30sec
Extension	72°C	1min
	4°C	Hold
→ Check library size by agarose gel		
→ 0.9X Ampure XP beads purification		
→ Quantify library by Qubit		

VII. DEEP SEQUENCING BY ILLUMINA PE-50 or PE-100

VIII. MATERIALS and REAGENTS

- **DSG (disuccinimidyl glutarate)** (ThermoFisher #20593)
- **EGS (ethylene glycol bis(succinimidyl succinate))** (ThermoFisher #21565)
- **16% Formaldehyde** (Fisher scientific #NC1040701)
- **2M Tris pH=7.5** (Sigma Aldrich #)
- **MBuffer#1:** 50mM NaCl, 10mM Tris-HCl pH 7.5, 5mM MgCl₂, 1mM CaCl₂, , 0.2% NP-40, 1X PIC.
- **MBuffer#2:** 50mM NaCl, 10mM Tris-HCl pH 7.5, 10mM MgCl₂
- **MBuffer#3:** 50mM Tris-HCl pH7.5, 10mM MgCl₂
- **Micrococcal Nuclease** (Worthington Biochem #LS004798)
Resuspended from lyophilized powder at 20 U/μl in Tris pH 7.4. Aliquot into tubes upon first use and freeze at -80°C.
- **0.5M EGTA** (Fisher scientific #)
- **0.5M EDTA** (Fisher scientific #)
- **T4 DNA Polymerase** (New England Biolabs #M0203)
- **DNA Polymerase I, Large (Klenow) Fragment** (New England Biolabs #M0210)
- **T4 Polynucleotide Kinase** (New England Biolabs #M0201)
- **T4 DNA Ligase** (New England Biolabs #M0202)
- **Exonuclease III (*E. coli*)** (New England Biolabs #M0206)
- **Biotin-14-dATP** (Jena Bioscience #NU-835-BIO14)
- **Biotin-11-dCTP** (Jena Bioscience #NU-809-BIOX)
- **20X Proteinase K solution** (Sigma Aldrich # 3115879001)
TE with 20 mg/ml proteinase K and 50% glycerol
- **Elution buffer:** 10 mM Tris-HCl pH 7.6
- **TE buffer:** 10 mM Tris-HCl pH 8.0, 1 mM EDTA
- **End-It DNA End-Repair Kit** (Lucigen # ER81050)
- **Dynabeads® MyOne Streptavidin C1** (Life Technologies # 65001)
- **KAPA HiFi HotStart ReadyMix** (KAPA Biosystems # KK2601)
- **NEBNext Ultra II** (New England Biolabs #E7645)