

	wild-type (n=30)	Fak56 <sup>N30/K24</sup> (n=30)	P-value
Single bouton area (μm²)	2.61 ± 0.55	2.60 ± 0.27	0.940
Bouton perimeter (μm)	6.02 ± 0.53	6.04 ± 0.35	0.976
Active zone / Bouton	4.15 ± 0.50	3.6 ± 0.29	0.349
Total active zone length (μm) / Bouton	1.95 ± 0.31	2.05 ± 0.18	0.776
Total active zone length (μm) /Bouton perimeter (μm)	0.30 ± 0.03	0.34 ± 0.02	0.279
T bar number / Bouton	1.45 ± 0.31	1.27 ± 0.20	0.623
Bouton area (μm²) /Active zone	0.60 ± 0.10	0.79 ± 0.08	0.150
T bar number /Active zone	0.32 ± 0.07	0.37 ± 0.05	0.642
Vesicle number /Bouton area (μm²)	81.86 ± 9.06	84.46 ± 7.62	0.827

Additional file 3. (A-B) Electron micrographs of cross-sections through a type I bouton of muscle 6/7 in wild-type and  $Fak56^{N30/K24}$  larvae. Squares show active zones with a T bar, which are enlarged in lower-left corner. Subsynaptic reticula (SSR), active zones (arrow), and mitochondria (Mt) are indicated. Scale bar, 2µm. (C) Features of the synaptic ultrastructure were quantified, and no significant difference in these parameters was found between wild-type and  $Fak56^{N30/K24}$ . Dissected larval body walls (including the CNS and motor axons) were fixed at RT for 30 minutes, followed by 4°C overnight in modified Trump's fixative (0.1 M sodium cacodylate buffer, 1% glutaraldehyde, and 4% formaldehyde). The fixed specimens were rinsed three times in 0.1 M sodium cacodylate buffer for 10 minutes, post-fixed for 30 minutes with 2% osmium tetroxide in

0.1 M sodium cacodylate buffer, rinsed three times for 10 minutes in 0.1 M sodium cacodylate buffer, and finally rinsed five times in ddH<sub>2</sub>O for 10 minutes. The muscle 6/7 in the A3 segment was knife-dissected out, and specimens were then stained *en bloc* in 2% aqueous uranyl acetate for 20 minutes, dehydrated in a graded ethanol series, and subsequently set into the Spurr's embedding medium. Thin sections (90 nm) were stained with uranyl acetate and lead citrate, and images were viewed on a Tecnai G2 Spirit TWIN electron microscope (FEI Company) and captured on a Gatan CCD Camera (794.10.BP2 MultiScan<sup>TM</sup>). TEM data were quantified by MetaMorph V6.3r7 (Molecular Devices).