

# A Successful Protocol for the Use of Pulse Oximetry to Classify Arterial Oxygenation into Four Fuzzy Categories

R. Matthew Sailors, ME, Thomas D. East, Ph.D., C. Jane Wallace, RN MS, Alan H. Morris, M.D.  
Pulmonary Division, LDS Hospital, 8th Ave. and C St., Salt Lake City, UT 84143.

## ABSTRACT

*Pulse oximetry is widely used in critical care medicine to noninvasively estimate arterial hemoglobin oxygen saturation. Despite the obvious benefits of using pulse oximetry to detect life threatening desaturations, it is unknown how well pulse oximetry is able to predict the finer graduations of arterial oxygenation needed for clinical decision making.*

*A computerized protocol was developed for the use of pulse oximetry to classify arterial oxygenation into four fuzzy categories and tested in a prospective clinical trial which compared the oxygenation category assigned by the protocol to one assigned by a respiratory therapist. In 3,742 classifications from 15 patients over a seven month period, the protocol showed 96% agreement with the therapists in the direction of therapy and 75% agreement with the oxygenation classes assigned by the therapists.*

## INTRODUCTION

Pulse oximetry is widely used in critical care medicine to noninvasively measure arterial hemoglobin saturation ( $\text{SaO}_2$ ). Although clinical practice guidelines exist for the use of pulse oximetry, they leave several important questions relating to the clinical use of pulse oximetry unanswered. [1]

1. How accurately does the pulse oximeter saturation ( $\text{SpO}_2$ ) measure  $\text{SaO}_2$  ?
2. How do these values relate to the arterial partial pressure of oxygen ( $\text{PaO}_2$ ) ?
  - a. How do you use  $\text{SpO}_2$  to determine arterial oxygenation?
  - b. How should  $\text{SpO}_2$  information be used in clinical decision making?
3. How do the "weasel words" used in the guidelines translate into specific values for the various clinical measurements and the intervals between reevaluation of the  $\text{SpO}_2$ – $\text{SaO}_2$  relationship? [2]

This paper discusses the development of an intensive care unit (ICU) computerized protocol that successfully addresses the aforementioned questions and classifies arterial oxygenation into four fuzzy categories (Table 2) by estimating the patient's  $\text{PaO}_2$  from the current  $\text{SpO}_2$ , past arterial blood gas (ABG)

data— $\text{PaO}_2$ ,  $\text{SaO}_2$ ,  $\text{SpO}_2$ , and the presence or absence of barotrauma.

The  $\text{SpO}_2$  protocol was designed to assess oxygenation as part of a set of protocols for managing mechanical ventilation. [3,4,5]

## METHODS

### Knowledge Engineering

The protocol logic was developed using a consensus committee to compile a list of common assumptions in the clinical use of  $\text{SpO}_2$  monitoring:

1.  $\text{SpO}_2$  and  $\text{PaO}_2$  are correlated.
2.  $\text{SpO}_2$  and  $\text{SaO}_2$  are not the same, but the gradient between them is small ( $< 4\%$ ) and stable within any given patient; therefore, single point calibration between  $\text{PaO}_2$  and  $\text{SpO}_2$  is adequate.
3. A patient's oxyhemoglobin dissociation curve may be shifted; however, the shift is stable.

Several different  $\text{SpO}_2$  protocols based on these common assumptions were developed and tested. Retrospective reviews and tests of these initial protocols revealed several fallacies in these assumptions:

1.  $\text{SpO}_2$  and  $\text{PaO}_2$  are not necessarily correlated. In 25% of cases,  $\text{PaO}_2$  and  $\text{SpO}_2$  do not even move in the same direction.  $\text{SpO}_2$  changes also reflect changes in parameters other than  $\text{PaO}_2$ . Carboxyhemoglobin (COHb), methemoglobin (MetHb), and peripheral blood flow also dramatically affect  $\text{SpO}_2$ .
2.  $\text{SpO}_2$  and  $\text{SaO}_2$  often differ widely (by  $> 6\%$ ) and their difference often change dramatically from the time of one blood gas to the next.
3. Shifts of the patient's apparent oxyhemoglobin dissociation curve ( $\text{SpO}_2$  vs.  $\text{PaO}_2$ ) may be quite large ( $> 20$  mmHg) and also vary dramatically from one blood gas to the next.

These problems indicate single point calibrations between  $\text{PaO}_2$  and  $\text{SpO}_2$  are not adequate, and indeed  $\text{PaO}_2$  and  $\text{SpO}_2$  may need to be recalibrated several times daily.

### $\text{SpO}_2$ Protocol

The current protocol (Figure 1) uses a combination of quality assurance rules, mathematical models,




Figure 1. SpO<sub>2</sub> Protocol in Flow Diagram Form  
 © 1995 IHC Hospitals, Inc. Reproduced with Permission

and thresholds to classify arterial oxygenation into one of the four oxygenation classes.

1. The entered SpO<sub>2</sub> is checked for a life threateningly low saturation (SpO<sub>2</sub> < 75%), so the ventilator management protocols can quickly respond without the need for an ABG or other data.
2. A quality assurance check is made to compare the difference between the SpO<sub>2</sub> measured at the time of the blood gas (SpO<sub>2</sub><sub>abg</sub>) and SaO<sub>2</sub>. If the difference is greater than 6 percentage points, the protocol exits and asks the clinician to try to reduce or eliminate the difference by changing, cleaning, or moving the SpO<sub>2</sub> probe.
3. Oxyhemoglobin curve shift parameters are calculated. ΔPaO<sub>2</sub>est (Eq. 1) adjusts for changes in temperature and acid balance (pCO<sub>2</sub>, [HCO<sub>3</sub><sup>-</sup>], pH), and Δ<sub>2</sub> (Eq. 2) adjusts the curve for factors, such as COHb and MetHb, which the pulse oximeter misinterprets as oxyhemoglobin.

$$\Delta PaO_{2est} = \frac{1}{0.070312 - 0.00060345 \cdot SpO_{2abg}} - PaO_2 \quad (1)$$

$$\Delta_2 = SpO_{2abg} - SaO_2 \quad (2)$$

Quality assurance checks are also made on ΔPaO<sub>2</sub>est to assure that it is within a reasonable range (-25 to 20 mmHg) and it is no more than 15 mmHg from the last value. Violations of these constraints cause the calculated ΔPaO<sub>2</sub>est to be replaced by the last stored value of ΔPaO<sub>2</sub>est if its magnitude is less than that of the calculated value. While ΔPaO<sub>2</sub>est and Δ<sub>2</sub> are recalculated every time a protocol SpO<sub>2</sub> is entered, they only change when a new PaO<sub>2</sub>, SaO<sub>2</sub>, SpO<sub>2</sub> set are available from an ABG.

4. PaO<sub>2</sub> is estimated from a shifted (ΔPaO<sub>2</sub>est and Δ<sub>2</sub>) hyperbolic curve (Eq. 3), fit to the standard oxyhemoglobin dissociation curve at 37° C, pH<sub>a</sub> = 7.40 between SaO<sub>2</sub> of 75% and 93%.

$$PaO_{2est} = \frac{1}{0.070312 - 0.00060345 \cdot (SpO_2 - \Delta_2)} - \Delta PaO_2 \quad (3)$$

Figure 2 illustrates the shifting effect of ΔPaO<sub>2</sub>est and Δ<sub>2</sub> on the oxyhemoglobin dissociation curve. [6] ΔPaO<sub>2</sub>est and Δ<sub>2</sub> model the changes in both position and shape of the dissociation curve resulting from SpO<sub>2</sub> measurement problems and oxygen-hemoglobin affinity changes due to variations in acid-base balance, temperature, 2,3-DPG, and molecules other than oxygen binding to the hemoglobin, such as in carbon monoxide poisoning. [7]




Figure 2. The effects of ΔPaO<sub>2</sub>est and Δ<sub>2</sub> on the oxyhemoglobin dissociation curve. [6,7]

Note: Scale and shape of curve have been exaggerated to emphasize the shift parameters.

5. Additional protocol logic was added to provide quality assurance for SpO<sub>2</sub> measurements. By limiting the interval between recalibrations, we were able to improve both the data quality and accuracy. We chose the level of patient support and O<sub>2</sub> Class as the guides for setting the recalibration interval. O<sub>2</sub> Class is assigned by a separate module using PaO<sub>2</sub>est and the membership functions listed in Table 2. Level of support is classified as either high, medium, or low and is based upon the Positive End Expiratory Pressure (PEEP) level and the inspired oxygen fraction (FiO<sub>2</sub>). The three levels of support, with the defining PEEP and FiO<sub>2</sub> levels, and the recalibration intervals are defined in Table 1. The recalibration interval is reevaluated every time a protocol SpO<sub>2</sub> is entered.

Table 1. Support Levels and Recalibration Intervals

O <sub>2</sub> Class (Table 2)	Support Level	Thresholds	Recalibration Interval
S	any		48.0 hours
T, M, A	low	FiO <sub>2</sub> < 0.4 PEEP < 5	24.0 hours
T, M, A	medium	0.4 < FiO <sub>2</sub> < 0.6 5 < PEEP < 15	12.5 hours
T, M, A	high	0.6 < FiO <sub>2</sub> 15 < PEEP	6.5 hours

S = Satisfactory, A = Acceptable, M = Marginal, T = Threatening

#### Retrospective Evaluation of SpO<sub>2</sub> Protocol

Computerized records of 14,424 paired ABG and SpO<sub>2</sub> measurements from 233 ARDS patients were used to retrospectively test the protocol. The first ABG was used to calculate ΔPaO<sub>2</sub>est and Δ<sub>2</sub>. The SpO<sub>2</sub> paired with the second ABG was used to calcu-

late PaO<sub>2</sub>est. PaO<sub>2</sub>est was then assigned to one of four clinically important oxygenation classes (O<sub>2</sub> Classes) based on its value and the presence or absence of barotrauma. The measured PaO<sub>2</sub> (PaO<sub>2</sub>meas) from the second ABG was also assigned to one of the four classes and the results of these classifications compared in a contingency table. [6]

Table 2. Fuzzy Arterial Oxygenation Classification Membership Functions

O <sub>2</sub> Class	Barotrauma	Nonbarotrauma
Threatening	PaO <sub>2</sub> < 50	PaO <sub>2</sub> < 50
Marginal	50 ≤ PaO <sub>2</sub> < 55	50 ≤ PaO <sub>2</sub> < 60
Acceptable	55 ≤ PaO <sub>2</sub> < 60	60 ≤ PaO <sub>2</sub> < 68
Satisfactory	60 ≤ PaO <sub>2</sub>	68 ≤ PaO <sub>2</sub>

### Prospective Clinical Trial of SpO<sub>2</sub> Protocol

The SpO<sub>2</sub> protocol was tested in a prospective trial comparing the O<sub>2</sub> Class assigned by the protocols to those assigned by respiratory therapists (RTs). The RTs scoring the patients were the staff of the Shock Trauma Intermountain Respiratory ICU (STRICU) at LDS Hospital. The STRICU staff has been using the ARDS protocols and assigning oxygenation classes to patients on a regular basis for the last seven years. [3,4,5]

The ventilator management protocols use the O<sub>2</sub> Class from the SpO<sub>2</sub> protocol unless the RT-assigned O<sub>2</sub> Class would have resulted in a change in therapy in the opposite direction (protocol would increase therapy, but RT would have reduced therapy). In the event of these disagreements, an ABG was drawn to assign an O<sub>2</sub> Class and recalibrate the shift parameters.

The SpO<sub>2</sub> protocol was subjected to a three month evaluation period, during which we monitored the protocol performance and responded to all real and perceived problems with the protocol. The SpO<sub>2</sub> protocol has been in daily clinical use for an additional four months in the ICUs of LDS Hospital and McKay-Dee Hospital Center, in Ogden, Utah.

## RESULTS AND DISCUSSION

### Retrospective Analysis Results

Table 3 compares oxygenation classifications from PaO<sub>2</sub>meas and PaO<sub>2</sub>est using the barotrauma oxygen classification functions.

The data was further lumped into two groups to assess the direction of subsequent therapy:

Increase therapy (↑ Tx) :

#### O<sub>2</sub> Classes

- marginal (see Table 2 for
- threatening for definitions )

Decrease or maintain therapy (↓/= Tx) :

#### O<sub>2</sub> Classes

- acceptable (see Table 2 for
- satisfactory for definitions )

Comparison of the appropriateness of the oxygenation classifications from PaO<sub>2</sub>est and PaO<sub>2</sub>meas shows that PaO<sub>2</sub>est drives therapy in the proper direction 84% of the time (Table 4).

Table 3. Retrospective Analysis of PaO<sub>2</sub>est using Barotrauma O<sub>2</sub> Class Thresholds

		O <sub>2</sub> Class from PaO <sub>2</sub> est				
		T	M.	A	S	Σ
O <sub>2</sub> Class from PaO <sub>2</sub> meas	T	1,581	351	111	104	2,147
	M	583	957	554	229	2,393
	A	178	606	985	871	2,640
	S	105	304	957	5,880	7,246
Σ		2,447	2,218	2,607	7,154	14,426

S=Satisfactory, A=Acceptable, M=Marginal, T=Threatening

Table 4. Retrospective Analysis of Therapy Changes

		PaO <sub>2</sub> est	
		↑ Tx	↓/= Tx
PaO <sub>2</sub> meas	↑ Tx	24.07%	7.4%
	↓/= Tx	8.27%	60.26%

↑ Tx : Increase Therapy ↓/= Tx : Decrease or Maintain Therapy

### Prospective Clinical Trial Results

The data from the prospective study, shown in Table 4, agree well with the results of the retrospective study. In 15 patients at two different hospitals, the protocol and RT O<sub>2</sub> Classes show 76.7 % agreement and are within one category in 98.9 % of the instances. Comparison of the direction of therapy change indicated a 96.1% agreement (Table 5).

Table 4. Prospective Study of SpO<sub>2</sub> Protocol

		O <sub>2</sub> Class from SpO <sub>2</sub>				
		T	M.	A	S	Σ
O <sub>2</sub> Class from RT	T	5	8	2	0	15
	M	20	246	76	18	360
	A	3	30	124	300	457
	S	0	17	394	2,499	2,910
Σ		28	301	596	2,817	3,742

S=Satisfactory, A=Acceptable, M=Marginal, T=Threatening

Table 5. Prospective Study of Therapy Changes

		SpO <sub>2</sub>	
		↑ Tx	↓/= Tx
RT	↑ Tx	7.46%	2.57%
	↓/= Tx	1.34%	88.64%

↑ Tx : Increase Therapy ↓/= Tx : Decrease or Maintain Therapy

Of the 146 instances where the RT and the SpO<sub>2</sub> protocol disagreed on the direction therapy should be adjusted, 96 were potentially inappropriate therapy reductions. The resolution of these 96 events is summarized in Table 6.

Table 6. Resolution of SpO<sub>2</sub> Protocol – RT Disagreement

Event	N
ABG drawn	37
Results agreed with Protocol	18
Results agreed with RT	19
No follow-up ABG	27
Protocol rerun with different data	32
	96

The largest single cause of these disagreements was data quality:

- bad or missing values for ABG-associated SpO<sub>2</sub>
- poor agreement between SpO<sub>2</sub> and SaO<sub>2</sub>
- use of last stored ΔPaO<sub>2</sub>est because calculated value does not meet quality assurance criteria (SpO<sub>2</sub> Protocol, item 3).

The next most common cause of disagreements was the calculated PaO<sub>2</sub>est. PaO<sub>2</sub>est was often incorrect in patients whose acid balance or temperature had changed by a clinically important amount since the last ABG. Misunderstandings of protocol function and incorrect "guesstimates" of the PaO<sub>2</sub> by clinical personnel also resulted in a large number of disagreements in assigning oxygen classifications.

Most of the disagreements on the direction of therapy were due to PaO<sub>2</sub>est overlapping, by less than 2 mmHg, into an O<sub>2</sub> Class immediately adjacent to the clinician assigned O<sub>2</sub> Class. For example, PaO<sub>2</sub>est = 54 in a patient with the presence of barotrauma which yields an O<sub>2</sub> Class of Marginal while the clinician estimates the O<sub>2</sub> Class to be Acceptable. This results in both different classifications and disagreement on the direction in which therapy should change.

None of the 96 disagreements resulted in negative patient outcomes or life-threatening hypoxemia.

Of the 32 instances where the protocol was rerun, more than half were rerun because of data entry and data quality errors.

## CONCLUSIONS

This project has lead to several interesting conclusions many of which can be generalized to other medical measurements and to any protocols or practice guidelines which use them.

1. SpO<sub>2</sub> and PaO<sub>2</sub> are not as well correlated as most clinicians think.
2. Protocols must include a great deal of quality assurance logic to be reliable.
3. Human decision making with SpO<sub>2</sub> is widely varied and has little formal quality assurance.
4. It is very difficult to extract and effect rules and standards for SpO<sub>2</sub>, a measurement that most clinicians feel is essential to the management of ICU patients.

Despite these problems, we have successfully created a protocol that standardizes the classification of arterial oxygenation assessed by pulse oximetry. This protocol may have application within pulse oximetry devices and ventilators in any clinical setting where SpO<sub>2</sub> is routinely monitored. Home care, assisting living centers, and emergency medical transport are only a few of the settings in which this protocol and others like it could be important assets.

## References

1. AARC clinical practice guideline: Pulse oximetry. *Resp Care* 1991;36:1406-1408.
2. MacDonald, Overnage Guidelines you can follow and can trust. *JAMA* 1994;271:872-87.
3. Kuperman, GJ, Gardner, RM, Pryor, TA HELP: A Dynamic Hospital Information System. Springer Verlag, New York, 1991.
4. Morris AH, Wallace CJ, Clemmer TP, et al. Final report: Computerized protocol-controlled clinical trial of new therapy which includes ECCO<sub>2</sub>R for ARDS. *Am Rev Respir Dis* 1992;145:4(2):A184.
5. Thomsen GE, East TD, Morris AH, et al. Clinical performance of a rule-based decision support system for mechanical ventilation of ARDS patients. In: Safran C, ed. *Proceedings of the Seventeenth Annual Symposium on Computer Applications in Medical Care* McGraw-Hill, New York, 1994.
6. East TJ, Morris AH, Wallace CJ, et al. Can pulse oximetry be used to reliably predict arterial oxygenation? *Crit Care Med* 1995;23:A27.
7. Guyton, AC Textbook of Medical Physiology, Eighth Edition. W.B. Saunders Company, Philadelphia, 1991.

**Acknowledgments:** This work was supported by AHCPR grant #HS 06594, the Deseret Foundation, Intermountain Health Care, Inc. (IHC), and LDS Hospital.