

Rapid and Sensitive Detection of Shiga Toxin-Producing *Escherichia coli* from Nonenriched Stool Specimens by Real-Time PCR in Comparison to Enzyme Immunoassay and Culture[†]

Thomas E. Gryns,^{1,2*} Lynne M. Sloan,¹ Jon E. Rosenblatt,^{1,3} and Robin Patel^{1,3}

Division of Clinical Microbiology, Mayo Clinic, Rochester, Minnesota¹; Molecular Microbiology Laboratory, Mayo Medical Laboratories, Andover, Massachusetts²; and Division of Infectious Diseases, Mayo Clinic, Rochester, Minnesota³

Received 17 October 2008/Returned for modification 12 January 2009/Accepted 5 May 2009

Shiga toxin (Stx)-producing *Escherichia coli* (STEC) bacteria are a frequent cause of food-borne gastroenteritis, hemorrhagic colitis, and hemolytic uremic syndrome. Because antimicrobial agents are generally contraindicated in patients infected with STEC, a sensitive and specific diagnostic test with rapid turnaround is essential. Current culture methods may fail to detect non-O157 STEC. We evaluated a Stx gene real-time PCR assay using hybridization probes and the LightCycler instrument with 204 prospectively collected stool specimens, which were also tested for Stx by enzyme immunoassay (EIA) (ProSpecT STEC; Remel, Lenexa, KS) and by culturing on chromogenic agar (Chromagar O157; BD BBL, Sparks, MD). In addition, 85 archived stool specimens previously tested for Stx (by EIA) and/or *E. coli* O157:H7 (by culture) were tested by PCR. Sample preparation for PCR included mixing the stool in sterile water and extraction of nucleic acid using the MagNA Pure LC instrument (Roche Diagnostics). The PCR assay had 100% sensitivity and specificity compared to EIA and culture for specimens collected prospectively (4 of 204 specimens were positive) and compared to culture and/or EIA for archival specimens (42 of 85 specimens were positive). Both the EIA and PCR produced positive results from a specimen containing an O103 serotype STEC in the prospective specimens, and the PCR test detected three positive specimens that contained nonviable STEC in the archived specimens. The PCR assay demonstrated 100% sensitivity and specificity compared to EIA and/or culture and more rapid turnaround than either EIA or culture.

Shiga toxin (Stx)-producing *Escherichia coli* (STEC) is a frequent cause of food-borne outbreaks of diarrhea (15). Disease caused by STEC is characterized by abdominal pain and bloody diarrhea, and 5 to 15% of those individuals infected with serotype O157:H7 develop hemolytic uremic syndrome (HUS), a potentially life-threatening condition consisting of hemolytic anemia, thrombocytopenia, and kidney failure caused primarily by Stx (8). STEC may carry genes for one or both types of Stx, Stx1 and Stx2 (17).

Although STEC strains are a diverse group of pathogens, up to the present, the most common serotype in the United States has been O157:H7. A common association is that of *E. coli* O157:H7 contaminating ground beef (3, 7), but recent large outbreaks have involved a variety of other foods, including leafy greens (6, 29). The diversity of potentially contaminated food means that patients may acquire STEC infection from many foodstuffs, far beyond the stereotypical risk of undercooked ground beef. The common denominator of tainted food products seems to be direct or indirect contamination from bovine feces. To best detect infected patients and potential outbreaks, clinical laboratories must have tools to quickly and accurately detect STEC in stool specimens. Culture on sorbitol MacConkey agar is an inexpensive, effective, and widely used method based on lack of sorbitol fermentation by *E. coli* O157:H7. Several drawbacks limit the utility of culture,

including slow turnaround, false-negative results in antibiotic-treated patients, and false-negative results due to emerging serotypes of non-O157 STEC that ferment sorbitol (1, 14, 16, 29). Alternatively, a method that is increasingly utilized is detection of Stx antigen from stool, either directly or after broth enrichment. Our experience concurs with enzyme immunoassay (EIA) product insert data that optimal sensitivity and specificity are achieved only when a broth enrichment step is employed; this results in slow turnaround.

Here, we describe a real-time PCR assay that can detect STEC using nucleic acid extracts of stool specimens. We evaluated the performance of this assay using both archived stool specimens and prospectively collected specimens and compared the results to those of culture and Stx antigen detection.

(This study was presented in part at the 48th Annual Inter-science Conference on Antimicrobial Agents and Chemotherapy, Washington, DC, 25 to 28 October 2008.)

MATERIALS AND METHODS

PCR assay. The PCR assay detects both the *stx*₁ and *stx*₂ genes by using primers designed for each of the genes. The master mixture (15 µl) containing 1× Roche LC FastStart DNA Master HybProbe (*Taq* DNA polymerase, reaction buffer, deoxyribonucleoside triphosphate mixture with dUTP instead of dTTP, and 10 mM MgCl₂), 3 mM MgCl₂, and 1× LC PCR primer-probe set (Table 1) (Stx1a, Stx1b, Stx1f, Stx1r, Stx2a, Stx2b, Stx2f, and Stx2r, kit no. 315; TIB MolBiol LLC, Adelphia, NJ) was added to the LightCycler cuvette. Extracted DNA (5 µl) was added to the reaction mixture. The cycling program was as follows: template denaturing at 95°C for 10 min; amplification of the template for 45 cycles of 10 s at 95°C, 15 s at 55°C (single acquisition), and 15 s at 72°C; and detection of the amplified product by melting analysis for 0 s at 95°C, 20 s at 59°C, 20 s at 40°C, ramp of 0.2°C/s, 0 s at 85°C, ramp of 0.2°C (continuous

* Corresponding author. Mailing address: Mayo Medical Laboratories New England, 160 Dascomb Road, Andover, MA 01810. Phone: (978) 253-6187. Fax: (978) 658-0439. E-mail: gryns.thomas@mayo.edu.

[†] Published ahead of print on 13 May 2009.

TABLE 1. Primers and probes used in the PCR assay^a

Gene	Name	Sequence (5' to 3')
<i>stx</i> ₁ ^b	Stx1a primer	CAAGAGCGATGTTACGGT
	Stx1b primer	AATTCTTCCTACACGAACAGA
	Stx1f probe ^d	CTGGGGAAGGTTGAGTAGCG
	Stx1r probe ^e	CCTGCCTGACTATCATGGACA
<i>stx</i> ₂ ^c	Stx2a primer	GGGACCACATCGGTGT
	Stx2b primer	CGGGCACTGATATATGTGTAA
	Stx2f probe ^d	CTGTGGATATACGAGGGCTTGATGTC
	Stx2r probe ^e	ATCAGGCGCGTTTTGACCATCT

^a Primer-Probe set, kit no. 315 (Tib MolBiol).

^b The *stx*₁ target corresponds to positions 2996451 to 2996243 of GenBank accession no. NC_002655 (21).

^c The *stx*₂ target corresponds to positions 1352455 to 1352659 of GenBank accession no. NC_002655 (21).

^d Labeled with fluorescein on the 3' end.

^e Labeled with LC640 on the 5' end and a phosphate on the 3' end.

acquisition), and 30 s at 40°C. *stx*₁ was differentiated from *stx*₂ amplification by melting curve analysis (Fig. 1). Positive and negative controls were included with each run. The assay was run successfully on two models of the LightCycler, 1.2 and 2.0.

Positive-control plasmids. Two positive control plasmids were constructed using the pCR 2.1 TOPO TA cloning kit (Invitrogen Corp., Carlsbad, CA) according to the manufacturer's instructions. The sources of the *stx*₁ and *stx*₂ sequences were *E. coli* O157:H7 (ATCC 43890) and O157:H7 (ATCC 43889), respectively. The sequences were amplified using primers Stx1a and Stx1b to clone 208 bp of *stx*₁ and Stx2a and Stx2b to clone 204 bp of *stx*₂. The sizes of the cloned sequences were confirmed by agarose gel electrophoresis, and the plasmids were purified using the High Pure Plasmid Isolation Kit (Roche Applied Science, Indianapolis, IN). The control plasmids were diluted in Tris-EDTA buffer (pH 8.0) and stored at 4°C.

Stool processing for PCR. A swab was used to transfer a pea-size amount of stool into 1 ml of sterile water. If the specimen was liquid enough to pipette, a 100- to 200-μl sample was mixed into 1 ml of sterile water. The stool-in-water samples were mixed by the use of a vortex apparatus and allowed to settle for at

least 2 min. Then, 200 μl of the supernatant was transferred to a sample cartridge for DNA extraction using a MagNA Pure LC Total Nucleic Acid Isolation Kit on a MagNA Pure LC instrument (Roche Diagnostics).

Analytical sensitivity and specificity, and inhibition. Stools sent to the clinical microbiology laboratory for *Clostridium difficile* testing were prepared as described above and used for the sensitivity and inhibition studies. For sensitivity studies, the stool-in-water samples were spiked with dilutions of *E. coli* O157:H7 before extraction. For inhibition studies, the stool-in-water sample was spiked with the plasmid control (final concentration, 100 copies/μl), and the mixture was extracted on the MagNA Pure LC instrument. To determine specificity, primer and probe sequences were used as subjects for a BLAST search on the National Center for Biotechnology Information website (<http://www.ncbi.nlm.nih.gov>). In addition, a panel of nucleic acid extracts from organisms found in stool (Table 2) was evaluated with the PCR assay.

Verification of the PCR results. A human subject protocol for clinical specimens included in this study was approved by the institutional review board prior to any testing. The PCR assay was compared to other methods using two sets of samples. The first was a prospective study including 204 consecutive stool specimens that were submitted to the clinical microbiology laboratory for Stx testing by the ProSpecT EIA (ProSpecT STEC; Remel, Lenexa, KS) with broth enrichment. In addition to the EIA, the specimens were cultured on Chromagar O157 (BD BBL, Sparks, MD). Isolates on Chromagar that were mauve were considered presumptive positive and were then confirmed on sorbitol MacConkey agar (BD BBL, Franklin Lakes, NJ). The EIA and culture plates were used in accordance with the manufacturers' instructions. The second set of stool samples included 85 archived clinical samples that had been frozen at -70°C. The set contained normal and diarrheal stools that were previously tested by routine stool culture (using sorbitol MacConkey agar) and/or tested for Stx by EIA (ProSpecT and/or Premier EHEC EIA [Meridian BioScience Cincinnati, OH]).

A positive PCR result was considered concordant with the comparison methods if a positive result was determined by either culture or EIA. Negative PCR results were considered concordant if all of the comparator methods used were negative. Samples with discordant results were retested by the same PCR assay, and nucleic acid extracts were sent to the Minnesota Department of Health (MDH) for analysis by a PCR method that used different primers (20).




FIG. 1. Melting curve analysis for the PCR assay. Amplification of *stx*₁ generates a curve with a peak at 55.77°C ± 2°C, while *stx*₂ generates a peak at 66.91°C ± 2°C. The lower intensity of the *stx*₂ peak is consistent whether detected alone or with *stx*₁.

TABLE 2. Organisms tested with the *stx* PCR assay (*n* = 66)

Organism	Accession no. or source ^a
<i>A. haemolyticus</i>	ATCC 17906
<i>A. haemolyticus</i>	ATCC 17977
<i>A. haemolyticus</i>	ATCC 19194
<i>Aeromonas hydrophila</i>	CAP-D-1-82
<i>Arcanobacterium pyogenes</i>	ATCC 19411
<i>Bacteroides distasonis</i>	ATCC 8503
<i>Bacteroides fragilis</i>	ATCC 25285
<i>Bacteroides thetaiotaomicron</i>	ATCC 29741
<i>Bacteroides vulgatus</i>	ATCC 29327
<i>Bifidobacterium adolescentis</i>	ATCC 15703
<i>Bifidobacterium bifidum</i>	ATCC 29521
<i>Bilophila wadsorthia</i>	NYS-01-Ed
<i>Campylobacter coli</i>	ATCC 33559
<i>Campylobacter jejuni</i>	ATCC 33560
<i>C. freundii</i>	ATCC 8090
<i>C. freundii</i>	NYS-2-01
<i>C. freundii</i>	NYS-2003-3 no. 5
<i>C. difficile</i>	ATCC 9689
<i>Clostridium perfringens</i>	ATCC 13124
<i>Clostridium ramosum</i>	ATCC 25582
<i>Collinsella aerofaciens</i>	ATCC 25986
<i>Cryptosporidium</i> sp.	Isolate from cat
<i>Dientamoeba fragilis</i>	ATCC 30948
<i>Eggerthella lenta</i>	ATCC 25559
<i>Encephalitozoon cuniculi</i>	ATCC 50602
<i>Encephalitozoon hellum</i>	ATCC 50451
<i>Encephalitozoon intestinalis</i>	ATCC 50651
<i>Entamoeba histolytica</i>	ATCC 30459
<i>Entamoeba moshkovskii</i>	ATCC 30042
<i>E. cloacae</i>	ATCC 13047
<i>E. cloacae</i>	CAP ID-13-B
<i>E. cloacae</i>	CAP-D-C-no.15
<i>E. cloacae</i>	CAP-D-3-84
<i>Enterococcus faecalis</i>	ATCC 19433-U
<i>Enterococcus faecium</i>	ATCC 19434
<i>E. coli</i> O142:K86(B):H6	ATCC 23985
<i>E. coli</i> O70:K:H42	ATCC 23533
<i>E. coli</i>	ATCC 25922
<i>Escherichia fergusonii</i>	ATCC 35469
<i>Escherichia hermanii</i>	ATCC 33650
<i>Escherichia vulneris</i>	ATCC 33821
<i>Eubacterium rectale</i>	ATCC 33656
<i>Fusobacterium nucleatum</i>	ATCC 25559
<i>Giardia lamblia</i>	ATCC 30957
<i>Klebsiella pneumoniae</i>	ATCC 700603
<i>Lactobacillus delbrueckii</i> subsp. <i>lacti</i>	ATCC 12315
<i>Lactobacillus rhamnosus</i>	ATCC 7469
<i>Mycobacterium avium</i>	ATCC 700398
<i>Peptostreptococcus magnus</i>	ATCC 29328
<i>Plesiomonas shigelloides</i>	ATCC 14029
<i>Proteus mirabilis</i>	ATCC 35659
<i>Pseudomonas aeruginosa</i>	ATCC 27853
<i>Salmonella enterica</i>	ATCC 35987
<i>Salmonella</i> group B	CAP-D-1-69
<i>Shigella boydii</i>	ATCC 29903
<i>Shigella dysenteriae</i>	ATCC 25931
<i>S. dysenteriae</i>	CDC 82-002-72
<i>Shigella flexneri</i>	CDC AB4-A04
<i>S. flexneri</i> serotype 2a	ATCC 29903
<i>Shigella sonnei</i>	ATCC 25931
<i>S. sonnei</i>	CDC 82-002-72
<i>Staphylococcus aureus</i>	ATCC 25923
<i>Staphylococcus epidermidis</i>	ATCC 14990
<i>Streptococcus bovis</i>	CAP-D-16-83
<i>Streptococcus sanguinis</i>	ATCC 10556
<i>Yersinia enterocolitica</i>	ATCC 9610

^a CAP, College of American Pathologists; NYS, New York State.TABLE 3. Detection of *stx* by PCR in stool compared to a combined standard^a

Real-time PCR result	No. of prospective specimens		No. of archived specimens		Total no. of specimens	
	Positive ^b	Negative ^c	Positive ^d	Negative ^e	Positive	Negative
Positive	4	0	42	0	46	0
Negative	0	200	0	43	0	243

^a Discordant results were resolved by a separate PCR method at the MDH. Results that were initially discordant with the real-time PCR results have been recalculated in this table to conform with the second PCR results from the MDH.

^b One positive by EIA only; three positive by both EIA and culture.

^c Negative by EIA and culture.

^d Positive by EIA and/or culture (*n* = 39). Three specimens were negative by EIA and/or culture but positive by the MDH PCR assay.

^e Negative by EIA and/or culture (*n* = 33). Specimens with discordant results were all negative by the MDH PCR assay (*n* = 10).

RESULTS

The analytical sensitivity of the PCR assay was 100% at 2 copies/μl of extracted specimen or 10 copies/reaction. This is roughly equivalent to 10⁴ CFU/g stool specimen. The inhibition studies demonstrated that none of the 53 specimen extracts tested contained substances that would cause a false-negative result in a low-positive sample. Melting curve analysis established that *stx*₁ melts at 55.77°C ± 2°C, while *stx*₂ melts at 66.91°C ± 2°C. Fifteen clinical samples yielded peaks between the two ranges. Six of these were sequenced, and all were found to be *Stx* sequences specific for *stx*₁ or *stx*₂. Therefore, a melting peak between 53.77 and 68.91 indicated a positive result for *stx* but did not allow differentiation of *stx*₁ versus *stx*₂.

The primers and probes were subjected to BLAST searches to identify potential cross-reacting sequences. Aside from the expected target, there were three organisms with exact matches to the *Stx*2a primer, *Acinetobacter haemolyticus*, *Enterobacter cloacae*, and *Citrobacter freundii*. At least three isolates of each of these species (strains acquired from ATCC or proficiency specimens from New York State or the College of American Pathologists) were tested in the cross-reactivity studies and were negative. Also, the BLAST search revealed no matches to the *Stx*2b primer, suggesting that no amplification would result in the event that *Stx*2a cross-reacted with an unintended sequence. In addition, 56 other organisms found in stool were tested in the assay, and none of them were positive (Table 2).

The accuracy experiments were divided into two parts. One was a series of 204 prospectively collected clinical stool specimens submitted for *Stx* testing by an EIA method, which we also cultured (Table 3). The other was a set of 85 archived stool specimens, some of which were previously positive for *Stx* (by EIA) and/or *E. coli* O157:H7 (by culture) and some of which were negative by toxin detection and culture. All of the prospective specimen results were concordant between PCR and EIA, including one positive specimen that was an O103 serotype STEC not detected by culture. Thirteen archived samples produced discordant results: six were previously positive by EIA, three were previously negative by EIA and culture, and four had previously discordant results of culture and EIA. Each of the discordant results was reextracted and retested by PCR. Each of the repeated results was the same as the initial result. To resolve the discordant results, the specimen extracts were sent to the MDH. The MDH tested the samples by

traditional PCR, which used primers for *stx*₁ and *stx*₂ different from those in the real-time PCR. The results from the MDH experiments were all concordant with the real-time PCR results. Archived specimen data and combined data from prospective and archived specimens are displayed in Table 3. All initially discordant results were recalculated using the data from MDH and are displayed accordingly in Table 3. The sensitivity and specificity of the assay were both 100% compared to the combined standard.

DISCUSSION

We developed and validated a real-time PCR assay that detects STEC directly from stool with 100% sensitivity and specificity and same-day results, a significant improvement in turnaround versus culture or antigen assays. The PCR assay has the sensitivity to detect as few as 1,000 CFU in a small sample of stool. In addition, two levels of sequence selection impart the high specificity of the assay. The first level involves the primers binding to and amplifying the target sequence. The second level is binding of the hybridization probes to sequences internal to the primer binding sites. Although some of the melting temperatures of the positive samples were between the typical *stx*₁ and *stx*₂ ranges, 100% of the positive specimens were detected.

The PCR assay has the ability to differentiate *stx*₁ and *stx*₂; however, the significance of this determination has not been resolved. Although some studies have shown that *stx*₂ tends to be associated with more serious disease (9, 10, 18), *stx*₁ has also been associated with severe cases. Complicating matters further, some STEC isolates carry genes for both *stx*₁ and *stx*₂ or more than one type of *stx*₂ (10). Notably, either type can cause serious disease, and the treatment for any STEC infection is primarily supportive—antibiotics should be withheld or discontinued (29, 31).

Interestingly, subtypes within the *stx*₂ group carry different disease risks. The five main subtypes are *stx*₂, *stx*_{2c}, *stx*_{2d}, *stx*_{2e}, and *stx*_{2f} (22, 26, 30). The BLAST searches revealed that the primer sequences for this assay matched *stx*₂ and *stx*_{2c} (GenBank accession no. Y10775 and M59432, respectively), which are the two types most commonly associated with the postinfection complication HUS, with 100% identity; the other subtypes have rarely or never been associated with HUS. The primer sequences also matched *stx*_{2d} (GenBank accession no. L11078), found commonly in non-O157:H7 strains, with 100% identity; this subtype is associated with diarrhea but rarely with HUS (23). Subtypes *stx*_{2e} and *stx*_{2f} (GenBank accession no. X81416 and M29153, respectively) matched the primer sequences for the assay with less than 90% identity, but these two subtypes are present in seldom-encountered serogroups of *E. coli* that are not typically human pathogens (10, 11, 24).

The PCR assay has the ability to detect non-O157 STEC pathogens, as well as culture-negative stools. Although non-O157 STEC strains are thought to be infrequent in the United States, culture methods using sorbitol MacConkey agar alone may not detect them, so their prevalence may be underestimated (16). Some states have found that non-O157 serotypes cause a substantial portion of STEC cases (5, 12, 19), and many studies from Europe have revealed that O157 serotypes are often the minority of STEC strains detected in diarrheal stools

(2, 14). Of the 204 prospective specimens included in this study, one of four positive results was due to a non-O157 serogroup, a rate of 25%. The isolate was an O103 serotype that fermented sorbitol and would have been missed if culture alone had been used. STEC was also detected by the PCR assay in five archived specimens that had been negative by culture (three of which were also negative by EIA). One limitation of the PCR method is that there is no isolate to characterize. Culture isolates are vital to public health efforts to detect outbreaks and track the epidemiology of STEC organisms.

One problem with Stx antigen assays is that a false positive may set off an inappropriate and potentially expensive investigation (4, 12). An additional benefit of the real-time PCR method described here is that the two sets of oligonucleotide primers and probes allow excellent selection of the desired target and minimize false-positive results. This was exemplified by the 100% specificity determined in this study. Another crucial aspect of the real-time platform is the ability to amplify and detect nucleic acid in a closed system, which minimizes the potential for contamination and false positives.

On the other hand, a false-negative result can be equally troubling, as a patient may receive antibiotic treatment, which could increase the risk of HUS. Some previous reports of PCR assays for Stx genes that use stool have demonstrated significant inhibition, and it was found necessary to use either a diluted extract or an enriched specimen (13). The assay described here found a 0% inhibition rate from stool, suggesting that false-negative results will be rare. Our laboratory has previously demonstrated 0% inhibition for 92 stools when they were prepared in the same manner for a *C. difficile* real-time PCR assay (27). In addition, the fact that the same extraction methodology preparatory to PCR can be used on whole blood without significant inhibition suggests that the presence of blood in stool will not impart a risk of false negativity (28).

Our assay is, to the best of our knowledge, the first comparison to culture and EIA of a LightCycler real-time PCR assay using a nonenriched stool specimen for STEC detection. Other groups have designed assays to detect STEC, either real-time PCR from isolated colonies or an enrichment broth or traditional PCR directly from a stool specimen, neither of which significantly improves turnaround time compared to culture or EIA. Schuurman et al. (25) described two real-time assays using stool specimens that were validated with a panel of STEC isolates, a group of 19 previously positive stools, and prospective specimens (3 of 115 stools were positive and confirmed by culture). Our assay compared favorably to the limit of detection but, in contrast, detected no PCR inhibition. We were able to include more than twice as many stool specimens, and our gold standard was a composite of both EIA and culture. Our assay will provide same-day turnaround with excellent sensitivity and specificity, which will provide the fastest possible result to the clinician.

ACKNOWLEDGMENTS

We thank Brian Duresko of the Mayo Clinic for help in gathering samples. We also appreciate the kind collaboration of the MDH, in particular, John Besser, Billie Anne Juni, Bonnie Koziol, Sharon Pen-dergrass, and Charlotte Taylor.

REFERENCES

- Bettelheim, K. A. 1998. Reliability of CHROMagar O157 for the detection of enterohaemorrhagic *Escherichia coli* (EHEC) O157 but not EHEC belonging to other serogroups. *J. Appl. Microbiol.* **85**:425–428.
- Bonnet, R., B. Souweine, G. Gauthier, C. Rich, V. Livrelli, J. Sirot, B. Joly, and C. Forestier. 1998. Non-O157:H7 Stx2-producing *Escherichia coli* strains associated with sporadic cases of hemolytic-uremic syndrome in adults. *J. Clin. Microbiol.* **36**:1777–1780.
- CDC. 2002. Multistate outbreak of *Escherichia coli* O157:H7 infections associated with eating ground beef—United States, June–July 2002. *JAMA* **288**:690–691.
- CDC. 2006. Importance of culture confirmation of Shiga toxin-producing *Escherichia coli* infection as illustrated by outbreaks of gastroenteritis—New York and North Carolina, 2005. *MMWR Morb. Mortal. Wkly. Rep.* **55**:1042–1045.
- CDC. 2007. Laboratory-confirmed non-O157 Shiga toxin-producing *Escherichia coli*—Connecticut, 2000–2005. *MMWR Morb. Mortal. Wkly. Rep.* **56**:29–31.
- CDC. 2006. Ongoing multistate outbreak of *Escherichia coli* serotype O157:H7 infections associated with consumption of fresh spinach—United States, September 2006. *MMWR Morb. Mortal. Wkly. Rep.* **55**:1045–1046.
- CDC. 1993. Update: multistate outbreak of *Escherichia coli* O157:H7 infections from hamburgers—western United States, 1992–1993. *MMWR Morb. Mortal. Wkly. Rep.* **42**:258–263.
- Clardy, T. G. 2004. The role of Shiga-toxin-producing *Escherichia coli* in hemorrhagic colitis and hemolytic uremic syndrome. *Semin. Pediatr. Infect. Dis.* **15**:260–265.
- Cornu, G., W. Proesmans, A. Dediste, F. Jacobs, J. Van De Walle, A. Mertens, J. Ramet, and S. Lauwers. 1999. Hemolytic uremic syndrome in Belgium: incidence and association with verocytotoxin-producing *Escherichia coli* infection. *Clin. Microbiol. Infect.* **5**:16–22.
- Eklund, M., K. Leino, and A. Siitonen. 2002. Clinical *Escherichia coli* strains carrying *stx* genes: *stx* variants and *stx*-positive virulence profiles. *J. Clin. Microbiol.* **40**:4585–4593.
- Friedrich, A. W., M. Bielaszewska, W. L. Zhang, M. Pulz, T. Kuczus, A. Ammon, and H. Karch. 2002. *Escherichia coli* harboring Shiga toxin 2 gene variants: frequency and association with clinical symptoms. *J. Infect. Dis.* **185**:74–84.
- Gavin, P. J., L. R. Peterson, A. C. Pasquariello, J. Blackburn, M. G. Hamming, K. J. Kuo, and R. B. Thomson, Jr. 2004. Evaluation of performance and potential clinical impact of ProSpecT Shiga toxin *Escherichia coli* microplate assay for detection of Shiga toxin-producing *E. coli* in stool samples. *J. Clin. Microbiol.* **42**:1652–1656.
- Holland, J. L., L. Louie, A. E. Simor, and M. Louie. 2000. PCR detection of *Escherichia coli* O157:H7 directly from stools: evaluation of commercial extraction methods for purifying fecal DNA. *J. Clin. Microbiol.* **38**:4108–4113.
- Johnson, K. E., C. M. Thorpe, and C. L. Sears. 2006. The emerging clinical importance of non-O157 Shiga toxin-producing *Escherichia coli*. *Clin. Infect. Dis.* **43**:1587–1595.
- Karch, H., P. I. Tarr, and M. Bielaszewska. 2005. Enterohaemorrhagic *Escherichia coli* in human medicine. *Int. J. Med. Microbiol.* **295**:405–418.
- Kehl, S. C. 2002. Role of the laboratory in the diagnosis of enterohaemorrhagic *Escherichia coli* infections. *J. Clin. Microbiol.* **40**:2711–2715.
- Nataro, J. P., and J. B. Kaper. 1998. Diarrheagenic *Escherichia coli*. *Clin. Microbiol. Rev.* **11**:142–201.
- Ostroff, S. M., P. I. Tarr, M. A. Neill, J. H. Lewis, N. Hargrett-Bean, and J. M. Kobayashi. 1989. Toxin genotypes and plasmid profiles as determinants of systemic sequelae in *Escherichia coli* O157:H7 infections. *J. Infect. Dis.* **160**:994–998.
- Park, C. H., H. J. Kim, and D. L. Hixon. 2002. Importance of testing stool specimens for Shiga toxins. *J. Clin. Microbiol.* **40**:3542–3543.
- Paton, A. W., and J. C. Paton. 1998. Detection and characterization of Shiga toxin-producing *Escherichia coli* by using multiplex PCR assays for *stx1*, *stx2*, *eaeA*, enterohaemorrhagic *E. coli hlyA*, *rfbO111*, and *rfbO157*. *J. Clin. Microbiol.* **36**:598–602.
- Perna, N. T., G. Plunkett III, V. Burland, B. Mau, J. D. Glasner, D. J. Rose, G. F. Mayhew, P. S. Evans, J. Gregor, H. A. Kirkpatrick, G. Posfai, J. Hackett, S. Link, A. Boutin, Y. Shao, L. Miller, E. J. Grobeck, N. W. Davis, A. Lim, E. T. Dimalanta, K. D. Potamouis, J. Apodaca, T. S. Anantharaman, J. Lin, G. Yen, D. C. Schwartz, R. A. Welch, and F. R. Blattner. 2001. Genome sequence of enterohaemorrhagic *Escherichia coli* O157:H7. *Nature* **409**:529–533.
- Persson, S., K. E. Olsen, S. Ethelberg, and F. Scheutz. 2007. Subtyping method for *Escherichia coli* Shiga toxin (verocytotoxin) 2 variants and correlations to clinical manifestations. *J. Clin. Microbiol.* **45**:2020–2024.
- Pierard, D., G. Muyldermans, L. Moriau, D. Stevens, and S. Lauwers. 1998. Identification of new verocytotoxin type 2 variant B-subunit genes in human and animal *Escherichia coli* isolates. *J. Clin. Microbiol.* **36**:3317–3322.
- Schmidt, H., J. Scheef, S. Morabito, A. Caprioli, L. H. Wieler, and H. Karch. 2000. A new Shiga toxin 2 variant (Stx2f) from *Escherichia coli* isolated from pigeons. *Appl. Environ. Microbiol.* **66**:1205–1208.
- Schuurman, T., A. Roovers, W. K. van der Zwaluw, A. A. van Zwet, L. J. Sabbe, A. M. Kooistra-Smid, and Y. T. van Duynhoven. 2007. Evaluation of 5'-nuclease and hybridization probe assays for the detection of Shiga toxin-producing *Escherichia coli* in human stools. *J. Microbiol. Methods* **70**:406–415.
- Serna, A. T., and E. C. Boedeker. 2008. Pathogenesis and treatment of Shiga toxin-producing *Escherichia coli* infections. *Curr. Opin. Gastroenterol.* **24**:38–47.
- Sloan, L. M., B. J. Duresko, D. R. Gustafson, and J. E. Rosenblatt. 2008. Comparison of real-time PCR for detection of the *tcdC* gene with four toxin immunoassays and culture in diagnosis of *Clostridium difficile* infection. *J. Clin. Microbiol.* **46**:1996–2001.
- Smith, T. F., M. J. Espy, J. Mandrekar, M. F. Jones, F. R. Cockerill, and R. Patel. 2007. Quantitative real-time polymerase chain reaction for evaluating DNAemia due to cytomegalovirus, Epstein-Barr virus, and BK virus in solid-organ transplant recipients. *Clin. Infect. Dis.* **45**:1056–1061.
- Tarr, P. I., C. A. Gordon, and W. L. Chandler. 2005. Shiga-toxin-producing *Escherichia coli* and haemolytic uraemic syndrome. *Lancet* **365**:1073–1086.
- Tzipori, S., A. Sheoran, D. Akiyoshi, A. Donohue-Rolfe, and H. Trachtman. 2004. Antibody therapy in the management of Shiga toxin-induced hemolytic uremic syndrome. *Clin. Microbiol. Rev.* **17**:926–941.
- Wong, C. S., S. Jelacic, R. L. Habeeb, S. L. Watkins, and P. I. Tarr. 2000. The risk of the hemolytic-uremic syndrome after antibiotic treatment of *Escherichia coli* O157:H7 infections. *N. Engl. J. Med.* **342**:1930–1936.