

**Appendix B: Information about patient safety events available in the MAUDE database.**

<b>Database field</b>	<b>Example</b>
MDR Report Key	<i>1570801</i>
MDR Text Key	<i>7441362</i>
Text	<i>There is a drop down menu for medication dosing frequency with more than 225 options that has an inaccurate and error prone override, requiring scrolling though the options arranged in alphabetical order or counterintuitive terms. This promoted errors in dosing frequency for uncommon programs, such as three days per week. We found dosing errors with Digoxin, a drug of narrow therapeutic ratio, in which the patient received 4X excess.</i>
Brand name	<i>Millenium</i>
Generic name	<i>CPOE</i>
Manufacturer name	<i>Cerner</i>