

SUPPORTING INFORMATION S1

Supplemental Material and Methods

Density gradient ultracentrifugation

For density gradient analysis of detergent-resistant material (DRM), cells were lysed in DRM lysis buffer (150 mmol/L NaCl, 25 mmol/L Tris, 10 mmol/L ATP, 5 mmol/L dithiothreitol, 5 mmol/L EDTA, 1.0% Triton X-100, pH 7.4) supplemented with protease inhibitor cocktail (Roche). After Benzonase treatment to remove DNA (50 units/ml lysis buffer for 30 min on ice), lysates were homogenized by 10 passages through a 22 gauge syringe needle, and then adjusted with Optiprep to 40%. Samples were placed on the bottom of an ultracentrifuge tube and overlaid with a discontinuous Optiprep gradient (30-5% in DRM lysis buffer). After ultracentrifugation at 200,000 x g for 4 hours in a fixed angle rotor, 11-12 fractions were collected from top and prepared for SDS-PAGE.

Figure legend

Co-localization of phospho-ERK1/2 with filamentous actin and ERK1/2 distribution between caveolar and non-caveolar fractions. (A) A7r5 cells were stimulated with either DPBA or FCS for 10 minutes, or left unstimulated, and then fixed and stained for immunofluorescence microscopy with a phospho-ERK1/2 antibody. Cells were co-stained with phalloidin to visualize actin filaments and with DAPI to visualize nuclei. Please note the filamentous ERK1/2 staining after DPBA stimulation (E and H, arrowheads), and the stronger nuclear ERK1/2 staining after FCS stimulation (I and L, arrows) compared to DPBA stimulation (E and H, arrows). Scale bar, 20 μ M. (B) Unstimulated A7r5 lysates were subjected to density gradient ultracentrifugation, and fractions (numbered from top to bottom) were analyzed for total ERK1/2 in western blots. The caveolar fractions were determined by co-staining the membranes with a caveolin-1 antibody. The percentage of ERK1/2 in the caveolar fraction was determined by densitometry (n=3).

A

pERK1/2

phalloidin


DAPI

Merge

unstim

DPBA

FCS



B

