

**Title:** High Frequency Components of Hemodynamic Shear Stress Profiles Are A Major Determinant of Shear-Mediated Platelet Activation In Therapeutic Blood Recirculating Devices

**Authors:** Filippo Consolo<sup>1,2\*</sup>, Jawaad Sheriff<sup>3</sup>, Silvia Gorla<sup>1</sup>, Nicolò Magri<sup>1</sup>, Danny Bluestein<sup>3</sup>, Federico Pappalardo<sup>2</sup>, Marvin J Slepian<sup>3,4</sup>, Gianfranco B Fiore<sup>1</sup>, Alberto Redaelli<sup>1</sup>.

**Affiliations:**

<sup>1</sup> Department of Electronics, Information and Bioengineering, Politecnico di Milano, Italy.

<sup>2</sup> Anesthesia and Cardiothoracic Intensive Care, IRCCS San Raffaele Scientific Institute, Vita Salute University, Italy.

<sup>3</sup> Department of Biomedical Engineering, Stony Brook University, USA.

<sup>4</sup> Department of Medicine and Biomedical Engineering, Sarver Heart Center, The University of Arizona, USA.

**SUPPLEMENTARY INFORMATION**

**The Platelet Activity State (PAS) assay**

The extent of activation of HSD-stimulated platelets was quantified via the PAS assay<sup>27</sup>. The PAS assay utilizes as substrate acetylated prothrombin, which upon exposure to damaged activated platelet membranes, as a result of elevated shear exposure, is converted to thrombin that is measured as a surrogate marker of activation<sup>26</sup>. The PAS assay allows characterizing the dynamic response of platelets by measuring the level of activation at increasing number of stimulation cycles<sup>11,12,26-30</sup>.

To characterize the dynamics of activation, at each time step (0, 2, 5, and 10 min), 25  $\mu$ l of GFP drawn from the HSD stimulation chamber were added to a 100- $\mu$ l tube containing (final concentrations) 5,000 platelets/ $\mu$ l, 200 nM Ac-FII, 5 mM  $\text{Ca}^{2+}$ , and 100 pM FXa and the tube was incubated at 37°C for 10 min; then, a 10- $\mu$ l sample was assayed for thrombin generation in a 96-wells microplate reader (Multiskan GO, Thermo Fisher Scientific Inc., Waltham, MA, USA), using 0.3 mM Chromozym-TH (Tosyl-Gly-Pro-Arg-4-nitranilide acetate, Roche Life Science, Milan, Italy) as the thrombin-specific chromogenic peptide substrate. Kinetic absorbance readings were performed at room temperature at 405-nm wavelength for 8 min<sup>29</sup>. The PAS value was calculated as the slope of the linear fitting of the absorbance-time data points over the 8-min kinetic reading. PAS values were normalized against those obtained by sonicating non-stimulated platelets with a microprobe sonicator (HD 2070 Sonoplus, Bandelin Electronic GmbH & Co. KG, Berlin, Germany). The sonication step is meant to yield platelets with maximal prothrombinase activity; thus, normalized PAS values (expressed in percentage, PAS [%]) represent the bulk activity as a fraction of the thrombin generation rate of sonicated platelets (100%)<sup>28</sup>. Sonication conditions (10 W for 10 s) were optimized for bovine platelets starting from the protocol reported for human platelets<sup>28</sup>.