1 Accuracy of Samsung Smartphone Integrated Pulse Oximetry Meets Full 2 FDA Clearance Standards for Clinical Use. 3 Sara H. Browne^{1,4} *†, Mike Bernstein²*, Philip E. Bickler ³ 4 5 6 Division of Infectious Diseases and Global Public Health, University of California San Diego, La 7 Jolla, CA, USA¹, Physio Monitor, LLC, San Ramon, CA, USA², Department of Anesthesia, 8 University of California San Francisco, San Francisco, CA, USA. ³Specialists in Global Health, 9 Encinitas, CA, USA,4 10 * These Authors contributed equally to this work 11 12 † Corresponding Author 13 14 Funding Sources: This research was funded by a NIH grant supplement R01MH110057-04S to 15 SHB and by the non-profit, Specialists in Global Health (https://sigh.global/). 16 17 18 **Abstract** 19 Background: Pulse oximetry is used as an assessment tool to gauge the severity of COVID-19 infection and identify patients at risk of poor outcomes. 1,2,3,4 The pandemic highlights the need for 20 21 accurate pulse oximetry, particularly at home, as infection rates increase in multiple global regions 22 including the UK, USA and South Africa⁵. Over 100 million Samsung smartphones containing 23 dedicated biosensors (Maxim Integrated Inc. San Jose, CA) and preloaded Apps to perform pulse

25

26

27

28

29

30

31

32

33

34

35

36

37

38

39

40

41

42

43

44

45

46

oximetry, are in use globally. We performed detailed in human hypoxia testing on the Samsung S9 smartphone to determine if this integrated hardware meets full FDA/ISO requirements for clinical pulse oximetry. Methods: The accuracy of integrated pulse oximetry in the Samsung 9 smartphone during stable arterial oxygen saturations (SaO₂) between 70% and 100% was evaluated in 12 healthy subjects. Inspired oxygen, nitrogen, and carbon dioxide partial pressures were monitored and adjusted via a partial rebreathing circuit to achieve stable target SaO₂ plateaus between 70% and 100%. Arterial blood samples were taken at each plateau and saturation measured on each blood sample using ABL-90FLEX blood gas analyzer. Bias, calculated from smartphone readings minus the corresponding arterial blood sample, was reported as root mean square deviation (RMSD). Findings: The RMSD of the over 257 data points based on blood sample analysis obtained from 12 human volunteers tested was 2.6%. *Interpretation:* Evaluation of the smartphone pulse oximeter performance is within requirements of <3.5% RMSD blood oxygen saturation (SpO₂) value for FDA/ISO clearance for clinical pulse oximetry. This is the first report of smartphone derived pulse oximetry measurements that meet full FDA/ISO accuracy certification requirements. Both Samsung S9 and S10 contain the same integrated pulse oximeter, thus over 100 million smartphones in current global circulation could be used to obtain clinically accurate spot SpO₂ measurements to support at home assessment of COVID-19 patients. Introduction Pulse oximetry supports the triage and initial management of symptomatic adults during respiratory infection pandemics^{6,7}. Pulse oximetry is being used during the current pandemic as an assessment

48

49

50

51

52

53

54

55

56

57

58

59

60

61

62

63

64

65

66

67

68

69

tool to identify patients with COVID-19 at risk of poor outcomes ¹ and gauge the severity of infection.^{2,3,4} The pandemic has highlighted the need for accurate pulse oximetry, particularly at home. In conjunction with COVID symptom diaries, at home oximetry is being used to detect early deterioration of patients in primary and community care settings.^{8,9} Primary care assessment pathways within major healthcare systems have developed guidance directives centered around SpO₂ measurements taken at home due to isolation; these readings with respiratory rates (RR) and mental state assessments obtained during virtual visits are being used to generate national early warning scores (NEWS).8 Remote healthcare providers also guide patients on early warning signs of worsening COVID-19 infection centered on deteriorating at home SpO₂ measurements over time, as well as exertion oximetry (exercise induced hypoxia). 8,10,11 The term 'silent hypoxia' associated with COVID-19 infection has appeared in the literature to describe the lack of dyspnea in some patients. ¹⁰ These reports serve as a reminder that the diagnosis of hypoxemia, even severe hypoxemia, is sometimes missed without objective measures such as non-invasive oximetry or arterial blood gases (ABGs).¹¹ The broad availability of accurate pulse oximetry is evidently of considerable significance during the current pandemic; furthermore, inequity in the distribution of accurate oximetry devices globally has been well documented by the WHO¹². FDA/ISO cleared pulse oximeter devices are expensive and it is likely that the majority of at home pulse oximeters currently in use, including those dispensed by large healthcare systems, are inexpensive devices that may be associated with wide variability in accuracy. 13 Samsung Galaxy S9 and 10 smartphones include the same dedicated hardware using high grade integrated biosensors and preloaded proprietary Samsung Apps that process the sensor data,

71

72

73

74

75

76

77

78

79

80

81

82

83

84

85

86

87

88

89

90

91

92

calculate saturation and heartrate and then display the user's pulse oximetry measurement ¹⁴. In 2018, an estimated 45 million S9 smartphones were shipped worldwide. ¹⁵ In addition reports indicate that 71 million S10 smartphones were shipped by the end of 2019. Given the increased global demand for accurate pulse oximetry associated with the COVID-19 pandemic, we performed detailed in human testing to determine if these smartphones meet the FDA/ISO requirements for clinical pulse oximetry. FDA 510K clearance and ISO CE marking require assessment of accuracy during hypoxemia. This type of testing was introduced by John Severinghaus MD of the University of California, San Francisco (UCSF) and refined by Dr. Philip Bicker UCSF^{16,17} During this testing detailed evaluation of bias between smartphone pulse oximeter readings and arterial blood gas measurements are made. Methods A Samsung S9+ smartphone (Samsung, Seoul, South Korea) containing Maxim Integrated biosensors, part number MAX86916 (Maxim Integrated, San Jose, CA) associated with the proprietary Samsung Health App was the test instrument. The test instrument has phone biosensors include a specialized dedicated photodetector and 2 precision wavelength LEDs equivalent to those used in clinical oximeters that in this case are dedicated to the oximetry function (See Figure 1). The associated App processes the sensor data, calculates saturation and heartrate then displays it for the user view and record if desired. Evaluation of test instrument oximeter performance was done using controlled steady-state hypoxia at the UCSF Hypoxia Research Laboratory. The work was sanctioned by the UCSF Committee on Human Research, protocol 10-00437, and conformed to all internationally accepted standards for

94

95

96

97

98

99

100

101

102

103

104

105

106

107

108

109

110

111

112

113

114

115

the protection of human subjects. We enrolled 12 participants, 4 male, 8 female, including 3 participants with darkly pigmented skin. All volunteers provided written informed consent. The left forefinger of each subject was placed over the smartphone sensor system (see Figure 1). During laboratory testing a series of desaturations are performed over 30 minutes. To enable test Participant to hold their finger continuously in position over a 30 minute period of time a silicone boot was used during Laboratory testing. Figure 2 shows the silicone boot attached to a plastic cell phone case that was utilized in this test. A 22-gauge radial artery cannula was placed in each participant after lidocaine local anesthesia. The participants then breathed a mixture of air, nitrogen and carbon dioxide to produce stable levels of arterial saturation between 100% and 70%. When oxygenation was stable, as assessed from endtidal gas and reference pulse oximeters, 2 arterial blood samples were taken, 30 seconds apart. Saturation and other oximetry parameters was measured on each blood sample with an ABL-90FLEX blood gas analyzer (Radiometer, Brea, CA, USA). Statistical Analysis: Bias was computed as smartphone readings minus the corresponding arterial blood sample value. Bias is reported as root mean square deviation (RMSD). A plot was generated following Bland and Altman with adjustments for multiple measurements for each individual according to the "Method Where the True Value Varies." The FDA 2013 guidance prescribes the use of a Bland Altman error plot with specifically calculated limits of agreement.¹⁸ **Results** In this study no plateaus were rejected for lack of stability between the first blood sample and the second. Eighteen readings from the device under test were rejected because the results were

117

118

119

120

121

122

123

124

125

126

127

128

129

130

131

132

133

134

135

136

137

138

delayed more than 15 seconds past the sample time. Twelve subjects completed the study and there were no adverse events. Figure 3 displays a modified Bland-Altman Plot of the entire data set. The RMSD of the over 257 data points based on blood sample analysis obtained from 12 human volunteers tested is 2.6%. **Discussion:** Interpretation of findings: Evaluation of oximeter performance during controlled steady-state hypoxia using the Samsung S9+ smartphone, containing Maxim Integrated biosensors and Samsung Health App, revealed that readings had an RSMD of 2.6% from over 257 simultaneous ABGs. This finding is well within requirements for FDA/ISO clearance for clinical pulse oximetry, which is <3.5% RMSD SpO₂ value. This is the first report of smartphone derived pulse oximetry measurements that meets full FDA/ISO accuracy certification requirements. These findings indicate that an estimated 45 million Samsung S9/9+ smartphones, with existing embedded dedicated hardware and preloaded Apps, in current circulation may be used to take accurate clinical grade pulse oximeter readings. Furthermore, S10 series smartphones contain the same biosensor hardware and preloaded Apps with the same proprietary Samsung algorithm, with reports stating 71 million such smartphones were shipped by the end of 2019. ¹⁹ Consequently, the implications of our study for global access to accurate pulse oximetry during the current COVID-19 pandemic are considerable. Up to this point studies evaluating smartphone pulse oximetry have shown variable accuracy with the minority supporting any clinical use. Jordan et al. compared pulse oximetry Apps downloaded on to iPhones to reference monitors within an Emergency room setting and found the Apps

140

141

142

143

144

145

146

147

148

149

150

151

152

153

154

155

156

157

158

159

160

provided inaccurate measurements with a sensitivity for detection SpO₂ <94% on the reference monitor ranging from 0-69%.²⁰ Two of these Apps utilized the onboard light and camera lens (Pox and Ox) and one used an external device that plugged into the iPhone. Alexander et al., ²¹ evaluated two other smartphone Apps (Pulse Oximeter and Pulse Oximeter Pro) again downloaded onto an iPhone and using onboard light and camera lens to perform SpO₂ measurements, predominantly in a Pre-operative setting. Comparison of these Apps with clinical devices utilized by the Dept of Anesthesia showed similar mean and median values, but wide and highly significant variance in measurements. ²¹ Tayfur et al., compared SpO₂ measurements obtained on 101 hospitalized patients (43% having pulmonary disease) using Samsung Galaxy S8 with that measured by simultaneous ABGs and reported measurements as highly correlated and having a small bias with narrow levels of agreement. 22 Modi et al. conducted 6 minute walk pre- and post- oximetry testing to perform a comparison of a Massimo-radical7 device, a clip Kenek sensor connected to an iPhone, and a Samsung Galaxy 8. 23 They reported a correlation between SpO₂ measurements (r = 0.62 - 0.72, p<0.001) and provided values from 28 out of 47 participants on a Bland Altman analysis that used an average of the SpO₂ obtained by the FDA approved Massimo-radical7 reference device and the smartphone readings as 'true' SpO₂, making this data difficult to evaluate. Most recently, Browne et al reported findings from clinical studies on 320 participants that utilized a repeated measure, nested factorial design to evaluate the accuracy and precision of a smartphone model, containing Maxim Integrated biosensors and App in comparison to Welch Allyn reference units.²⁴ Their results indicated a wide range of adult persons utilized a smartphone model to perform repeated pulse oximeter spot checks with an accuracy and precision equivalent to hospital grade clip sensors. ²⁴ This clinical study had few hypoxemic patients, but did report human steady state hypoxia testing of

162

163

164

165

166

167

168

169

170

171

172

173

174

175

176

177

178

179

180

181

182

the smartphone model versus a portable reference device that indicated FDA/ISO requirements were likely to be met if full testing was conducted.²⁴ These widely varying findings are associated with studies that did not follow uniform protocols, were in different patient populations and none performed as we did full FDA/ISO required human hypoxia testing using ABGs as the oxygen saturation reference. Aside from these limitations, the major determinant of accuracy differences observed across these studies likely reflects differences in the quality of the underlying hardware and proprietary algorithms within the smartphone tested. Whilst traditional pulse oximeters in hospital settings rely on the transmission of light through cutaneous tissues such as the finger or ear lobe, smartphones systems utilize reflected light detected by a sensor on the same surface as the emitter. ²⁵ Within such systems how the reflected light signal is obtained and analyzed is critical. Poorly performing pulse oximetry Apps used the onboard light and camera lens to obtain reflected light to detect PPG signals and conformational changes associated with hemoglobin binding. Signal obtained from camera associated sensors have relatively high signal to noise ratios and may even block near infrared light. ²⁵ In contrast, the smartphone evaluated in our study has dedicated hardware specifically for pulse oximetry function (see Figure 1) with 660nM red and 910nM infrared LEDs used to take measurements, a photodetector connected to an extremely low noise analog channel allowing measurement on a broad range of skin colors, and placement entirely separate from other phone devices such as the camera. Studies testing similar hardware and algorithms reported high SpO₂ measurement accuracy. ^{24, 22} All future studies evaluating smartphone oximetry accuracy should specify the type of embedded hardware including biosensor manufacturer, App proprietor and site of placement within the smartphone, as these features constitute the actual oximetry device

184

185

186

187

188

189

190

191

192

193

194

195

196

197

198

199

200

201

202

203

204

205

Implications for Clinical Assessment: Our introduction described the widespread health system guidance on the use of at home SpO₂ measurement to grade the severity of COVID-19 infection. Currently, infection numbers are increasing rapidly in association with the presence of the UK variant in multiple regions globally, resulting in greater demand for effective, inexpensive at home monitoring to identify patient deterioration^{26,27}. Media reports indicate consumer purchase of oximeters has increased considerably since the start of the current pandemic. ²⁵ The foremost concern of experts regarding home pulse oximetry during the current pandemic is the accuracy of SpO₂ measurements, particularly as saturation falls below 90%.²⁵ Data on the accuracy of inexpensive stand-alone oximeters is limited, as there has been little regulatory oversight. 13,25 The best study available on the accuracy of inexpensive finger oximeters reported wide variability, with the majority of devices tested demonstrating highly inaccurate readings during hypoxia.¹³ Consumers are essentially unaware of these differences and healthcare workers only find reliable accuracy information for more expensive devices (>\$150) on the market. ²⁵ As the current cost of FDA cleared pulse oximeters varies from hundreds to greater than a thousand USD, it is highly likely the majority of at home pulse oximeters in use are not FDA/ISO cleared. In this regard our findings indicate that pulse oximetry measurements obtained by the Samsung Galaxy S9/9+ meet prescribed standards for clinical use and are more accurate than stand-alone oximeters for which no reliable information is available. Pulmonologist practical guidance during COVID-19 at home monitoring states oximetry devices should provide some indication of pulse signal strength and that measurements should be taken and recorded two to three times a day at rest. ²⁵ Using any oximeter, the finger must be positioned accurately relative to the sensor. The Samsung Healthcare App does provide assessment of adequate

207

208

209

210

211

212

213

214

215

216

217

218

219

220

221

222

223

224

225

226

227

228

finger placement and adequate pulse strength and does not take a reading unless both are appropriate. This feature is important as accurately positioning the finger over the sensor on the back of the phone, then turning the hand to view the reading on the front, is more difficult than using a finger clip sensor. The oximetry function on Samsung S9/9+ smartphones is designed for spot check only, where the user holds their finger in place for approximately 30 seconds and gets a single reading. While multiple spot checks during the day can be made, the manner in which these readings are displayed under a label of 'stress test' within the Samsung Health App may be confusing to patients. However, the App does contain a 'share' function which enables approved providers to receive daily SpO₂ and HR data that can be used to support at home assessment; and historical data is stored. Collaboration between software developers and smartphone manufacturers could easily improve this design to include instruction to take 2-3 readings a day at rest and warm extremities before measurement, in line with recommendations²⁵, as well as present graphical data trends over time for transmission to healthcare practitioners. Study Limitations: This study used healthy volunteers under steady state hypoxia testing (SpO₂) approximately 70-100%) and thus provides no data on accuracy below 70%. Data from healthy volunteers provides no information on accuracy of SpO₂ measurement under conditions of reduced extremity pulsatile blood flow, associated with vasoconstriction due to Raynauds, hypotension or peripheral vascular disease, the latter is common in persons with diabetes and coronary heart disease, both risk factors for severe Covid-19 infection ²⁸. Furthermore, measurement accuracy in the presence of diseases associated with hemoglobinopathies, such as sickle cell or thalassemia, or where elevated carboxyhemoglobin (potentially present in heavy smokers) or methemoglobin (associated with use of chloroquine, sulfonamides) my occur, is unknown. This study does not

230

231

232

233

234

235

236

237

238

239

240

241

242

243

244

245

246

247

248

249

250

251

provide any evidence on the use of this system in children. Finally, the system tested is designed to provide SpO₂ spot checks only, there is no evidence from our study to support use for continuous pulse oximetry. Conclusion: Based on our findings and the substantial need for widespread global access to accurate pulse oximetry during the current pandemic we recommend the high-grade biosensors and Samsung proprietary algorithm embedded in S9 and S10 smartphones pursue FDA 510K clearance and ISO certification for use to obtain spot SpO₂ measurements for clinical use. Review could be prioritized and fast tracked by the FDA/ ISO organizations under the 'Emergency Use Authorizations' for 'Remote or Wearable Patient Monitoring Devices associated with the Covid-19 pandemic'. 30 This may substantially support global efforts to respond to the current Covid-19, particularly in LMIC settings, where access to accurate spot pulse oximetry is limited 11. Acknowledgements: Funding was provided by NIH grant supplement R01MH110057-04S to SHB and by the non-profit, Specialists in Global Health (https://sigh.global/). We thank Samsung Electronics, Suwon-si, Korea, for critical appraisal of this manuscript. References 1. Richardson S, Hirsch JS, Narasimhan M, et al. Presenting characteristics, comorbidities, and Outcomes among 5700 Patients Hospitalized with COVID-19 in the New York City Area. JAMA. 2020;323:2052-2059. doi:10.1001/jama.2020.6775 2. Gandhi RT, Lynch JB, del Rio C. Mild or Moderate Covid-19. Solomon CG, ed. New England Journal of Medicine. April 2020. doi:10.1056/nejmcp2009249

252 3. Berlin DA, Gulick RM, Martinez FJ. Severe Covid-19. Solomon CG, ed. New England 253 Journal of Medicine. May 2020. doi:10.1056/nejmcp2009575 254 4. Nicola, M., O'Neill, N., Sohrabi, C., & Khan, M. (2020). Evidence Based Management 255 Guideline for the COVID-19 Pandemic - Review article. *International Journal of Surgery* 256 (London, England). doi:10.1016/j.ijsu.2020.04.001 257 5. The latest global coronavirus statistics, charts and maps. (n.d.). Retrieved from 258 graphics.reuters.com/world-coronavirus-tracker-and-259 maps/ 6. Sala, H., Roca, J. S., Zerbo, G., Garcia, R., Cabral, G., Fernandez, A., ... Rizzo, O. (2011). 260 261 Initial Clinical Management of Symptomatic Adult Patients during Influenza A (H1N1) 262 Epidemics. The Journal of Emergency Medicine, 41(4), 435-440. 263 doi:10.1016/j.jemermed.2010.05.093 264 7. Verhoeven, D., Teijaro, J. R., & Farber, D. L. (2009). Pulse-oximetry accurately predicts 265 lung pathology and the immune response during influenza infection. Virology, 390(2), 151-266 156. doi:10.1016/j.virol.2009.05.004 267 8. Coronavirus: Pulse oximetry to detect early deterioration of patients with COVID-19 in 268 primary and community care settings. (2020). Retrieved from 269 https://www.england.nhs.uk/coronavirus/publication/pulse-oximetry-to-detect-early-270 deterioration-of-patients-with-covid-19-in-primary-and-community-care-settings/ 271 9. NHS England. (2020). Annex 2: Remote Monitoring COVID-19 Diary. Retrieved from NHS

England website: https://www.arstubiedriba.lv/wp-content/uploads/2020/11/NHS-COVID-

272

273

19 oksimetrs.pdf

274 10. Fuglebjerg, N. J. U., Jensen, T. O., Hoyer, N., Ryrsø, C. K., Lindegaard, B., & Harboe, Z. B. 275 (2020). Silent hypoxia in patients with SARS CoV-2 infection before hospital discharge. 276 International Journal of Infectious Diseases, 99, 100-101, doi:10.1016/j.ijid.2020.07.014 277 11. Bickler, P. E., Feiner, J. R., Lipnick, M. S., & Mckleroy, W. (2020). "Silent" Presentation of 278 Hypoxemia and Cardiorespiratory Compensation in COVID-19. Anesthesiology. 279 doi:10.1097/aln.0000000000003578 280 12. WHO. (2008), Global Pulse Oximetry Project Background Document. Retrieved from 281 https://www.who.int/patientsafety/events/08/1st pulse oximetry meeting background doc. 282 pdf 283 13. Lipnick, M. S., Feiner, J. R., Au, P., Bernstein, M., & Bickler, P. E. (2016). The Accuracy of 284 6 Inexpensive Pulse Oximeters Not Cleared by the Food and Drug Administration. 285 Anesthesia & Analgesia, 123(2), 338-345. doi:10.1213/ane.000000000001300 286 14. Personal Communication Samsung Electronics, Suwon-si, Korea. 287 15. O'Dea, S. (2020). Samsung Galaxy S series shipments worldwide 2016-2018. Retrieved 288 from https://www.statista.com/statistics/864691/samsung-galaxy-s-series-smartphone-289 shipments-worldwide/ 290 16. Severinghaus, J. W., Naifeh, K. H., & Koh, S. O. (1989). Errors in 14 pulse oximeters 291 during profound hypoxia. Journal of Clinical Monitoring, 5(2), 72-81. 292 doi:https://doi.org/10.1007/bf01617877

17. Severinghaus, J. W. (2009). Gadgeteering for Health Care. Anesthesiology, 110(4), 721-

293

294

728.

296

297

298

299

300

301

302

303

304

305

306

307

308

309

310

311

312

313

314

315

316

18. Bland, J. M., & Altman, D. G. (2007). Agreement Between Methods of Measurement with Multiple Observations Per Individual. Journal of Biopharmaceutical Statistics, 17(4), 571-582. doi:10.1080/10543400701329422 19. Canalys Newsroom- Global smartphone market Q4 and full year 2019. (2019). Retrieved from Canalys.com website: https://www.canalys.com/newsroom/canalys-globalsmartphone-market-q4-2019 20. Jordan, T. B., Meyers, C. L., & Schrading, W. A. (2020). The utility of iPhone oximetry apps: A comparison with standard pulse oximetry measurement in the emergency department. American Journal of Emergency Medicine, 38(5), 925-928. doi:10.1016/j.ajem.2019.07.020 21. Alexander, J. C., Minhajuddin, A., & Joshi, G. P. (2017). Comparison of smartphone application-based vital sign monitors without external hardware versus those used in clinical practice: a prospective trial. Journal of Clinical Monitoring, 825-831. doi:10.1007/s10877-016-9889-6 22. Tayfur, I., & Afacan, M. A. (2019). Reliability of smartphone measurements of vital parameters: A prospective study using a reference method. American Journal of Emergency Medicine, 37(8), 1527-1530. doi:10.1016/j.ajem.2019.03.021 23. Modi, A. M., Kiourkas, R., Li, J., & Scott, J. B. (2020). Reliability of Smartphone Pulse Oximetry in Subjects at Risk for Hypoxemia. *Respiratory Care*, respcare.07670. doi:10.4187/respcare.07670 24. Browne, S. H., Bernstein, M., Pan, S. C., Garcia, J. G., Easson, C. A., Huang, C. C., . . . Vaida, F. (2020). Maxim Integrated Smartphone Sensor with App Meets FDA/ISO

318

319

320

321

322

323

324

325

326

327

328

329

330

331

332

333

334

335

336

337

338

339

Standards for Clinical Pulse Oximetry and can be Reliably Utilized by a Wide Range of Patients. Chest. doi:10.1016/j.chest.2020.08.2104 25. Luks, A. M., & Swenson, E. R. (2020). Pulse Oximetry for Monitoring Patients with COVID-19 at Home. Potential Pitfalls and Practical Guidance. Ann America Thoracic Society, 17(9), 1040-1046. doi:10.1513/AnnalsATS.202005-418FR 26. Mahase, E. (2020). Covid-19: What have we learnt about the new variant in the UK? BMJ, 371. https://doi.org/10.1136/bmj.m4944 27. Wise, J. (2020). Covid-19: New coronavirus variant is identified in UK. BMJ, 371. https://doi.org/10.1136/bmj.m4857 28. Daily, L. (2020). What is a pulse oximeter, and does the coronavirus pandemic mean you need one? Retrieved from https://www.washingtonpost.com/lifestyle/wellness/pulseoximeter-covid-19-coronavirus/2020/05/18/5b6f8a98-96df-11ea-9f5e-56d8239bf9ad story.html 29. Du Y, Tu L, Zhu P, et al. Clinical Features of 85 Fatal Cases of COVID-19 from Wuhan: A Retrospective Observational Study. Am J Respir Crit Care Med. 2020;201(11):1372-1379. doi:10.1164/rccm.202003-0543OC 30. Center for Devices and Health. (2020). Coronavirus Disease 2019 (COVID-19) Emergency Use Authorizations for Medical Devices. Retrieved from https://www.fda.gov/medicaldevices/emergency-use-authorizations-medical-devices/coronavirus-disease-2019-covid-19emergency-use-authorizations-medical-devices

FIGURES: SMARTPHONE INTEGRATED PULSE OXIMETER

Figure 1: (A) shows the location of the SpO₂ sensor on the back of the smartphone making clear that the sensor is entirely separate from the other phone devices such as the camera. In the expanded view (B) the dedicated photo sensor and LEDs for the pulse oximetry function are shown. The device has a total of 4 LEDs two of which are reserved for future additional functionality. The arrow indicates the 660 nM red and 910nM infrared LEDs used in this measurement. The photodetector is connected to an extremely low noise analog channel allowing measurement on a broad range of skin colors.



Figure 2: Plastic phone case with silicon boot attached used to hold each participant's finger in position continuously for over 30 minutes.



Figure 3: Plot of data collected from 12 human volunteers during this test of the bias in blood oxygen measurements. Each point on the plot represents paired blood oxygen levels recorded by the hemoximeter associated with the Samsung Phone reading. The upper and lower limits of agreement per Bland Altman 2007¹⁷ are shown in red and blue lines respectively.



