

# IF YOU TALK ABOUT IMMUNOCHEMICALS— WE WILL UNDERSTAND YOU!

## 各种免疫试验用品应有尽有

אם אתה מדבר על אימונוכימיקלים - אנחנו נבין אותך!

We speak English, French, Chinese, Hebrew and Swahili, but mostly we speak Immunochemicals, and that is the language we speak best.

Does it seem strange that our goats graze the hills leading to Jerusalem and yet if we don't leave the rabbit houses before 4:30 p.m., we get caught in the five o'clock rush back to the lab in Kankakee? Strange, perhaps, but our animals are certainly productive! Jointly, Miles-Yeda, Rehovot, Israel, and Miles Research Division, Kankakee, Illinois, produce over 150 different antisera, immunodiffusion plates, immunochemical kits and immunoabsorbents.

We find our international environment stimulating, too. Our involvement with the Weizmann Institute of Science (we are right on the campus) plus the technological advances made at our Kankakee labs provide insight and background to produce the products you need.

We invite your questions, comments or suggestions. Again, we don't care what the accent is—if you speak *Immunochemicals*, we will understand you.

Reach us by phone in the U.S.A. at 815 939-4417, or in Israel at 03 95-29-22, or by letter through any of our offices listed.

Ask for P. K. or Tsvi.

Dr. P. K. Chung  
Dr. T. Hirshfeld

RESEARCH DIVISION



Miles Laboratories, Inc.  
Post Office Box 272  
Kankakee, Illinois 60901  
Phone: 815 939-4417

Miles Laboratories Ltd.  
Post Office Box 37  
Stoke Court, Stoke Poges  
Slough, England SL 2 4 LY  
Phone: Farnham Common 2151

Miles-Yeda Ltd.  
Post Office Box 1122  
Kiryat Weizmann  
Rehovot, Israel  
Phone: 03 95-29-22



# ADVERTISING RATES

Cover and Special Positions	£90
Full Page	£70
Half Page	£40
Quarter Page	£25
Agency Commission	10%
Publisher's Discount	10%

All communications regarding advertising matters should be addressed to:

Advertisement Section,  
Biochemical Society,  
7 Warwick Court,  
London, WC1R 5DP

Telephone: 01-242 1076  
(4 lines)

# THE BIOCHEMICAL JOURNAL

# LIPID ANALYSIS

William W Christie

*Hannah Research Institute, Ayr*

352 pages 53 illustrations  
October 1973 £6.00 hard cover

This book critically examines the literature and brings together in a systematic manner the best of the procedures that have been developed for separating, identifying and determining lipid classes and their component parts.

Newcomers will find this book a useful guide through the potential complexities of lipid analysis, and experts in the field will find it a valuable reference work.



Headington Hill Hall, Oxford OX3 0BW

**Pergamon**

## Automated Analysis of Drugs and Other Substances of Pharmaceutical Interest

Edited by C. T. Rhodes, BPharm, PhD, FRIC, MPS  
and R. E. Hone, BPharm, MPS

The editors of this book are convinced that the development and use of automation in the field of quantitative drug analysis is inevitable in the very near future. The twelve contributors are experienced research scientists; their practical advice will consequently be of special value to scientists considering conversion from manual to automated methods.

1973 300 pp., illustrated 0 407 11151 4 cased £5.95  
0 407 11150 6 limp £3.95

## Drugs in Anaesthetic Practice 4th Edition

F. G. Wood-Smith, MA, MB(Cantab), FFARCS.;  
M. D. Vickers, MB, BS(Lond), FFARCS; and  
H. C. Stewart, MA, MD(Cantab), PhD(Lond), FRCP,  
FFARCS

This is the only medium-sized book which covers the whole of this subject and gives not only the theoretical background and science, but also actual practical help on specific drugs and therapeutic problems. It has been thoroughly revised, and some sections have been completely rewritten. The major part of the book deals with various classes and groups of drugs, and a complete grounding in basic physical chemistry and physiological regulation is included.

1973 650 pp., illustrated 0 407 15502 3 £6.50

Obtainable from booksellers or from  
**The Butterworth Group**  
88 Kingsway, London WC2B 6AB  
Showroom: 4-5 Bell Yard, WC2



# Hoppe-Seyler's Zeitschrift für Physiologische Chemie

Editors in Chief

A. BUTENANDT · F. LYNEN · G. WEITZEL

Subscription Rates

For one volume (12 parts) DM 420,—

Vol. 354 No. 6

Contents

June 1973

Regulation of transcription and translation in eukaryotes (Abstracts)

24. MOSBACHER KOLLOQUIUM DER GESELLSCHAFT FÜR BIOLOGISCHE CHEMIE

The chemical structure of the nucleus as the basis of its function (Abstracts)

17. KONFERENZ DER GESELLSCHAFT FÜR BIOLOGISCHE CHEMIE

The effect of a non-peptide interchain crosslink on the reoxidation of reduced insulin

D. BRANDENBURG and A. WOLLMER

Proliferation and energy metabolism of Ehrlich ascites tumor cells in a glucose-free medium

H. P. KRAUSE and FR. SCHNEIDER

On the mechanism of ketogenesis and its control, I: On a possible role of acetylaceto-CoA thiolase in the control of ketone body production

W. HUTH, C. DIERICH, V. v. OEYHAUSEN and W. SEUBERT

Chlorination of the cystine- and tyrosine peptides from the C-terminus of insulin A-chain

S. STOEV, B. ALEKSIEV and A. VALKOVA

Preparative isolation of disc-electrophoretically pure components of hordein

K. KLING

Studies on polypeptides, IV. The synthesis of C-peptide of human proinsulin

V. K. NATHANI

A combination of micro-disc electrophoresis with antigen-antibody crossed electrophoresis (identification and quantitative determination of individual serum proteins)

W. GIEBEL and H. SAECHTLING

Contents of deoxyribonucleoside triphosphates in Yoshida ascites tumor cells

E. LISS, S. BRAMMER and H. SCHMIDT

SHORT COMMUNICATION

Crystallisation of the hemoglobin of carp (*Cyprinus carpio*)

T. KLEINSCHMIDT, B. GRUJIC-INJAC and G. BRAUNITZER

Indexed in Current Contents



Walter de Gruyter · Berlin · New York

# DNA Polymerase I


(EC 2.7.7.7)

Boehringer Mannheim have DNA Polymerase I that you can rely upon in every way. For convenience of use the enzyme is supplied in 50% glycerol in a polypropylene vial that can be kept happily at  $-20^{\circ}\text{C}$  during shipment and storage.


Our Grade I preparation with a specific activity of 2500–5000 U/mg<sup>1</sup> corresponds to stage 7 of the Jovin et al. preparation<sup>2</sup>. Before it leaves our Research Laboratories in Tutzing it undergoes exacting quality control procedures including SDS disc-electrophoresis which shows just one heavy and a few very faint bands. A fragment of the DNA Polymerase I enzyme—Enzyme A according to Klenow<sup>3</sup> is the latest addition to our range. It has a molecular weight of only 75,000, as opposed to 109,000 for the DNA Polymerase I, due to the removal by proteolytic activity of the 5'→3' Exonuclease.

#### References

1. Richardson, J., Biol. Chem. 239, 222 (1964).
2. Jovin et al, J. Biol. Chem. 244, 2996 (1969).
3. Klenow, H. et al, Proc. Nat. Acad. Sci. 65, 168 (1970).



DNA Polymerase is no exception to the Boehringer Mannheim dictum 'that for each enzyme its own substrate, co-factor and inhibitor is supplied'. Examples from the range of complementary biochemicals available for DNA Polymerase I are deoxytriphosphates including 5'-bromo dUTP, Nalidixic acid and Poly d (G-C). Analytical specifications and applications for these and the complete range of our biochemicals for Molecular Biology and Structural Studies are in these illustrated brochures.



Please write or telephone for your own personal copy to:


**The Boehringer Corporation (London) Ltd.**  
Bilton House, Uxbridge Road, London, W5 2TZ  
Tel: 01-579 6944 (5 lines).

# Pharmacia Peristaltic Pump P-3

**A cassette pump  
electronically controlled for  
constant liquid flows in  
chromatography and  
gradient formation**

**You can**

- obtain three different flow rates at the same time through the three individually tensioned pumping channels
- change tubes in seconds with the snap-in cassette
- forget about pulsation thanks to the large driven rollers
- rely on a constant motor speed, electronically compensated for load and temperature
- use any flow rate between 0.6 and 400 ml/h per channel with the same gear box
- reverse the flow instantly
- push for maximum flow and clear pump lines rapidly



Pharmacia (Great Britain) Ltd.  
Paramount House  
75 Uxbridge Road  
EALING LONDON W5 5SS




**Pharmacia  
Fine Chemicals**

# Ideas that cannot wait

Since publishing our first *Handbook of Biochemicals*, we have added many new and interesting biochemicals, some of which are too exciting to await our second handbook.




**85,852-8**  
**3-Guanidino-L-alanine**  
Competitive inhibitor of arginine  
decarboxylase in *E. coli*.  
*J. Biol. Chem.*, **243**, 1671 (1968).  
**250mg \$10.00**




**85,847-1 H**  
**5-Fluorouracil**  
 Antitumor agent.  
*Cancer Res.*, 33, 128 (1973);  
*J. Bacteriol.*, 113, 1311 (1973).  
**1g \$6.00 5g \$24.00**




**85,787-4**  
**(Methylazoxy)methanol acetate**  
**Induces tumors.**  
*Fed. Proc.*, 31, 1485 (1972).  
**250mg \$28.00 1g \$90.00**




**16,431-3**  
**3-O-Methyldopamine hydrochloride**  
**(3-methoxytyramine)**  
 Metabolite of L-DOPA; building  
 block for substituted catecholamines.  
*Nature*, **230**, 533 (1971);  
*J. Med. Chem.*, **15**, 692 (1972).  
 1g \$12.00    5g \$40.00




**85,788-2**  
**3-Deazauridine**  
Exhibits cytotoxic and  
antitumor activity.  
*J. Med. Chem.*, **16**, 294 (1973).  
**250mg \$20.00 1g \$70.00**




**85,755-6**  
**4-(Ethoxymethylene)-2-phenyl-5(4H)-oxazolone**  
Synthetic hapten used to induce cutaneous inflammatory reaction. *Immunology*, 23, 289 (1972); *J. Nat. Cancer Inst.*, 48, 1697 (1972).  
**5g \$18.00 10g \$26.00**



**85,845-5**  
**3-Isobutyl-1-methylxanthine**  
Inhibitor of cyclic AMP  
phosphodiesterase.  
*Biochem J.*, **122**, 115 (1971).  
**1g \$10.00**



**85,843-9**  
**Erucic acid**  
Inhibitor of fatty acid  
oxidation in the heart.  
*Biochim. Biophys. Acta*,  
**280**, 506 (1972).  
**5g \$5.00**



**85,773-4**  
***m*-Aminobenzamidine hydrochloride**  
Used in the affinity chromatography  
of human complement protein.  
*Immunochemistry*, **8**, 539 (1971).  
**1g \$5.00**

**Aldrich Chemical Company, Inc.**  
*Craftsmen in Chemistry*



**Home Office:  
Aldrich Chemical Co., Inc.  
940 W. St. Paul Ave.  
Milwaukee, Wisconsin 53233**

**In Great Britain:  
Ralph N. Emanuel Ltd.  
264 Water Rd., Wembley Middx.,  
HAO 1PY, England**

**In Continental Europe:  
Aldrich-Europe  
B-2340 Beerse,  
Belgium**

**In Germany:  
EGA-Chemie KG  
7924 Steinheim am Albuch  
Germany**

*Printed in Great Britain by  
William Clowes & Sons Limited, London, Colchester and Beccles*