

# In Vitro Alternatives for Ocular Irritation

Rodger D. Curren and John W. Harbell

Institute for In Vitro Sciences, Inc., Gaithersburg, Maryland

The necessity of using animals to test whether new chemicals and products are eye irritants has been questioned with increasing frequency and fervor over the last 20 years. During this time many new nonanimal methods have been proposed as reliable alternatives to the traditional rabbit (Draize) test. To date, however, none of these nonanimal (*in vitro*) tests have become universally accepted as a complete replacement for the Draize test. To understand why a complete replacement has not been found, one has to first understand the reasonably complex structure of the eye, the standard Draize scoring scale—which is based on a qualitative evaluation of three different tissues—the differences between human and rabbit eyes, the intrinsic variability of the animal test, and the details of the different *in vitro* tests that have been proposed as replacements. The *in vitro* tests vary from relatively simple assays using single cells to more sophisticated assays that use discarded animal tissue or artificially constructed human tissue. It is clear that appropriately designed *in vitro* tests will eventually give more useful mechanistic information about ocular injury from which we can more comfortably predict the risk of human eye irritation from new products and ingredients. — *Environ Health Perspect* 106(Suppl 2):485–492 (1998). <http://ehpnet1.niehs.nih.gov/docs/1998/Suppl-2/485-492/curren/abstract.html>

Key words: alternatives, *in vitro*, eye, irritation, toxicology, ocular

## Introduction

The necessity of using animals to test whether new chemicals and products are eye irritants has been questioned with increasing frequency and fervor over the last 20 years. Admittedly the questions are very complex, and strong social, political, ethical, and scientific arguments have been raised on both sides of the issue. During this process, numerous nonanimal methods have been proposed as reliable alternatives to the traditional animal tests. However, before such tests come into common use, they must be carefully evaluated to determine if, in fact, these new methods can replace or reduce the use of whole animals. Such evaluation involves investigating the basic details of ocular irritation, reviewing what type of information is currently obtained from the animal tests, understanding how the design of *in vitro* systems relates to the

animal model, and only then determining what real progress has been made in the search for alternatives to traditional ocular irritation testing.

## Structure of the Eye


In interpreting the results from any toxicologic study, there must be some basic knowledge of the organ system being studied—at the very least an understanding of its morphology, cellular constituents, and normal function—that allows one to determine whether an injury has occurred and what the consequences of that injury are. Because there are many similarities between the animal and the human eye, and because the human eye is the organ we are trying to protect, the human eye will be used here as an example for the discussion of ocular structure. Where differences exist between the eyes of humans and rabbits (the usual target species of ocular irritation testing), they will be noted.

Figure 1 depicts a human eye in both a normal front view and in cross section. The latter view more clearly shows the tissues that are of concern to toxicologists and ophthalmologists. Perhaps the most important tissue is the cornea: the exterior surface that is exposed to the outside environment. The normally transparent cornea allows light to freely enter the eye and eventually be focused

on the retina. If the cornea becomes cloudy (opaque)—as can happen after accidental exposure to strongly irritating chemicals—light can no longer pass easily into the eye and vision becomes impaired or even completely blocked. Although the eyelids offer the cornea some protection, it is still very susceptible to injury.

About 80% of the cornea's structure is the stroma—a regular array of macromolecules through which light can easily pass as a consequence of the stroma's high degree of order and exact level of hydration. Maintenance of this very important hydration level (75–80% water) is the responsibility of two active cell layers, a single-cell-thick endothelium covering the inside surface of the cornea and a much thicker epithelium that covers the outside surface. These cell layers work together to keep additional water from entering the cornea, which would result in swelling and opacity. The epithelium also has a second function of providing a physical barrier against the entry of foreign materials. If the epithelium is injured, corneal opacity can result. However, minor opacities can often be reversed because the epithelium can repair itself either by movement of surrounding cells to cover the wound or by the actual replacement of damaged tissue through new cell division. In contrast, the endothelium is generally not capable of repair. Therefore, if these cells suffer cytotoxic damage there can be significant consequences, e.g., permanent blindness. It is this relationship between the induction of cellular damage and resulting ocular irritation or other injury that is the basis for many of our current *in vitro* ocular irritation screening systems.

Another delicate tissue of the eye is the conjunctiva, the nonkeratinized squamous



**Figure 1.** Tissues of the eye, diagram of the human eye. A cross section through the center of the eye is shown on the left and the normal front view is on the right. The three major tissues of interest for eye irritation studies are illustrated.

This paper was prepared as background for the 13th Meeting of the Scientific Group on Methodologies for the Safety Evaluation of Chemicals (SGOMSEC): Alternative Testing Methodologies held 26–31 January 1997 in Ispra, Italy. Manuscript received at EHP 9 May 1997; accepted 16 January 1998.

Address correspondence to Dr. R.D. Curren, Institute for In Vitro Sciences, Inc., 21 Firstfield Road, Suite 220, Gaithersburg, MD 20878. Telephone: (301) 947-6527. Fax: (301) 947-6538. E-mail: rcurren@iivs.org

Abbreviations used: BCOP, Bovine Cornea Opacity and Permeability Assay; CAM, chloroallantoic membrane; IRAG, U.S. Interagency Regulatory Alternatives Group; TEA, Tissue Equivalent Assay.

epithelium that lines the inner surfaces of the eyelids and much of the external surface of the ocular globe (it is continuous with the cornea). The conjunctiva is highly vascularized and may become quite inflamed after exposure to irritating material. Mildly irritating chemicals or other products often cause conjunctivitis without any associated corneal damage.

A third important ocular tissue is the iris (the colored part of the eye), which, by constricting or dilating, controls the amount of light that enters the eye and is eventually focused on the retina. The iris lies under the cornea within the aqueous humor. In some cases foreign materials penetrate completely through the cornea and interact with the iris. The iris may then become very inflamed and may lose its ability to react to light, seriously damaging the ability to see.

Observations of the degree of injury to each of these tissues in the animal model are incorporated as part of the scoring system of most common eye irritation protocols. The details of these scoring systems will be discussed in "In Vivo Ocular Irritation Testing."

### In Vivo Ocular Irritation Testing

It is important to understand how manufacturers assure themselves that their products will not pose an unacceptable risk to the eyes of their customers. Generally the process consists of several steps. First, the maximum potential hazard of the ingredient or formulation to the ocular tissue is determined. Second, the actual use of the product is considered, estimating the probability that it may inadvertently enter the eye. Third, a final safety assessment takes into account benefits, risks, and the impact of the instructions for use that generally accompany the product. Although the entire process is important, it is the first stage of this process—generally termed hazard identification—and the development of improved *in vitro* systems to detect such hazards that are important to this discussion.

For obvious ethical reasons, tests using animals, rather than volunteer human subjects, have generally been used to assist toxicologists in determining the degree of danger a material poses to the eye. Although it is possible to expose human eyes to dilute forms of materials whose chemical properties are well known and generally regarded as safe, it is obvious because of the risk of severe injury that this cannot be

done with novel materials whose toxic properties are as yet unknown.

The albino rabbit has historically been the animal of choice for testing potential eye irritants, primarily because its large eyes make it easy to observe damage. In addition, it has a large conjunctival sac (accentuated by loose lids) that easily accepts test material and holds it against the eyes. However, because of several striking differences, the rabbit is far from a perfect model for humans. One difference is the presence of a nictitating membrane, or third eyelid, in the rabbit. This membrane moves laterally across the eye, likely causing the kinetics of removal of many test materials to differ from humans. Another difference is that the conjunctival sac of the rabbit is much larger than in humans, meaning that more test material can be placed in a rabbit's eye than would be likely to ever get into the human eye during an accidental exposure. Additionally, the rabbit cornea is somewhat thinner than that of humans and there is less tear production to aid in washing out a foreign material. For these and other reasons, the rabbit is generally considered an overly sensitive model for humans. Although this may be considered a positive aspect of the rabbit model because it adds a margin of safety to the risk assessment, it also presents the problem of inappropriate hazard assessment and suggests that a more predictive model would be beneficial.

The conduct of the animal test now needs to be examined in detail to help us understand the subjective nature of the test and appreciate the difficulties faced in developing and validating *in vitro* models. To test a material for potential ocular irritation, the lower lid of the animal is pulled away from the eyeball, and 100  $\mu$ l of a liquid (100 mg of a solid) is placed in the resulting conjunctival sac. The lids are then held together for a few seconds to ensure contact between the test material and the ocular tissues. The animal's eyes are carefully observed, first at 1-hr and then at 24-hr intervals for up to 14 days. It is important to highlight again a major difference between the structure of a rabbit's eye and that of a human eye. Because 100  $\mu$ l will not fit into the human eye, the animal's eyes are being exposed to much more material than might actually enter the human eye from an accidental exposure. The low-volume eye test (1), which uses one-tenth of the material normally applied to the rabbit eye, is reported to better predict the response of human eyes and to be less hazardous to the animal.

The time that the test material is in contact with the eye is not controlled during the Draize test because the material is only removed by the natural processes of tearing and blinking. Therefore, time of exposure may differ with each test material, which makes it difficult to develop an *in vitro* model.

At various standard time periods after instillation, the three major tissues of the eye (cornea, conjunctiva, and iris) are observed macroscopically for injury. Each tissue is observed for different signs and the degree of injury is recorded according to a standard scale. For the cornea, the degree of opacity and the area of the eye involved are recorded. The iris is examined for inflammation and the conjunctiva—a mucus membrane—is examined for redness, chemosis, and any exudate. Generally mild responses of the conjunctiva alone are not serious unless the test material is designed to be applied to or around the eyes.

### Draize Scoring Scale

The fact that three ocular tissues can be affected by chemical treatment makes simple scoring and evaluation of ocular damage difficult. Draize (2) proposed in 1944 what has become a solution to the problem. He devised an individual numerical scoring system for each of the three ocular tissues of interest and then proposed a special weighting system to combine the scores into a single eye irritation score. Table 1 shows specifically how Draize reduced the evaluation of a very complicated type of injury to a single number.

This awareness of the various ocular tissues and the ways they respond to injury is very important because we need to understand exactly what score an *in vitro* eye irritation assay is supposed to predict. However, the complete Draize scoring system is not generally used to classify materials for regulatory purposes. The European Union classification scheme, for example, uses only discrete categories such as R36 (irritating to eyes), and R41 (risk of serious

**Table 1.** The Draize test<sup>a</sup> for assessing ocular irritation.

Eye tissue	Maximum
Cornea score = opacity (0–4) $\times$ area (0–4) $\times$ 5	80
Iris score = grading value (0–2) $\times$ 5	10
Conjunctiva score = [redness (0–3) + edema (0–4) + discharge (0–3)] $\times$ 2	20
Total score	110

<sup>a</sup>An illustration of how individual subjective observations of injury to three ocular tissues are converted into a single numerical score estimating total eye injury.

damage to eyes) to define the amount of danger that exposure to a material represents (3). However, a continuous scale like that of Draize gives considerably more information about the severity of the hazard than does an abbreviated classification system, and this prediction of many levels of severity is what most people envisage as being supplied by a replacement *in vitro* test.

However, the single number (Draize score) presented as the eye irritation potential of a chemical or formulation is not exact (4). The subjective nature of the gross observations made during the scoring of the test, plus normal animal-to-animal variability, make it virtually impossible to routinely reproduce the final Draize score, especially for midrange irritants. This does not mean that Draize scores are completely meaningless. Repeated trials will generally generate scores within an acceptable range. However, a single Draize score should not be viewed as an exact predictor of eye irritation potential. Bruner et al. (5) presented a summary of the effect of this *in vivo* variability on the evaluation of *in vitro* tests.

With this degree of variability as the norm, it becomes very difficult to have confidence that a Draize test will be able to reveal any real differences between two mild materials. Yet product developers and toxicologists often need to be able to detect with confidence small differences between candidate products. Thus the need for alternative tests that may provide more precise data springs not just from concerns for animal welfare, but also from the impreciseness of the animal test and the desire of toxicologists for better and more accurate tools.

### Development of Alternative Models

The need to develop alternative *in vitro* tests for eye irritation has been apparent for some time. These alternative assays have used a diverse set of human and animal cells, tissues, and even biochemical matrices (Table 2). What is the strategy that researchers have used to develop *in vitro* assays that model the animal assay? First, only a few *in vitro* assays actually attempt to

model the entire eye. In fact, most *in vitro* tests that have been proposed are reductionist, i.e., they tend to model only one small part of the complex process of eye irritation. This has led to a situation where an *in vitro* test that may measure only one specific type of damage is compared to a Draize score that covers several types of damage in several tissues. Success at this type of comparison is more than might reasonably be expected from a single *in vitro* test; thus it makes sense to think of an eventual test battery with several *in vitro* tests, each one capable of detecting a different type of damage. This type of approach has recently been used by the U.S. Interagency Regulatory Alternatives Group (IRAG). This group led an extensive international evaluation of the state-of-the-art of *in vitro* eye irritation tests based on a comparison of *in vitro* scores to individual animal tissue scores, not to the total combined Draize score (6).

### Ex Vivo Models


Figure 2 illustrates the continuum of reductionist relationships between the

**Table 2.** Characteristics of common *in vitro* assays.

Assay, reference	General description	Method of applying test material
Neutral red release (16), neutral red uptake (17)	Target cells (primary or continuous; fibroblasts or epithelial-like) are grown in submerged monolayer culture	Generally increasing dilutions of test material are added to growth medium until a predetermined end point (generally cytotoxicity) is reached
Fluorescein leakage (14,15)	Target cells (primary or continuous; fibroblasts or epithelial-like) are grown in submerged monolayer culture; medium may be removed for dosing	Either increasing dilutions of test material are added to growth medium or cell surface for a set time, or a single concentration is added for varying times. End point is induction of permeability of the monolayer
BCOP (10)	Living bovine corneas are treated with test material and changes in opacity and permeability are measured by instrument	Test materials are applied neat or at in-use concentrations directly to the epithelial surface of the cornea
HetCam (25)	Chorioallantoic membrane of a chicken egg is treated	Test materials are applied neat or at in-use concentrations directly to the membrane and damage to the membrane is recorded
Tissue equivalent assay (11)	Three-dimensional reconstructed tissue (often human) is grown with top surface exposed to air	Test materials are applied neat or at in-use concentrations directly to the tissue construct and cell killing is measured
Enucleated chicken eye (23)	Isolated eye of a chicken is treated and subsequent damage recorded	Test materials are applied neat or at in-use concentrations directly to the tissue construct
Cytosensor microphysiometer <sup>a</sup> (23)	Cells held over or on a coated sensor are treated and changes in cellular metabolism are recorded in real time	Generally increasing dilutions of test material added to growth medium until a predetermined end point (decrease in metabolism) is reached
Irritection (Eytex) <sup>b</sup> (23)	End point is precipitation of protein in a nonviable commercially supplied matrix. Meant to mimic opacity formation in the cornea	Either dilutions or neat test material is added to a membrane bullet over a responding protein matrix
Pollen tube growth (26)	Tobacco pollen is allowed to germinate in the presence of test material	Dilutions of test material used; end point is inhibition of pollen tube elongation
Red blood cell (19)	Red blood cells are exposed to test material	Dilutions of test material used; lysis (release of hemoglobin) and hemoglobin denaturation are monitored
SIRC (27)	Target cells (continuous cell line derived from rabbit cornea) grown in submerged culture at clonal densities	Generally increasing dilutions of test material added to growth medium until a predetermined end point (generally cytotoxicity) is reached

Abbreviations: HetCam, hen's egg test on the chorioallantoic membrane; SIRC, Staatens Seruminstituut Rabbit Corneal—a fibroblastlike cell line derived from a rabbit cornea.

<sup>a</sup>Molecular Devices Corporation, Menlo Park, CA. <sup>b</sup>InVitro International, Irvine, CA.



**Figure 2.** Relationship of *in vitro* models to animal tissue. Diagram illustrating how several specific *in vitro* assays have been developed to model different parts of the animal eye. The eye is essentially dissected into its component parts and an *in vitro* assay is developed to represent each part so that specific types of eye injury can be detected and the mechanism of action understood.

whole animal and the *in vitro* model. Some assays focus on a first stage of reduction, i.e., looking just at the isolated eye without any associated conjunctiva that would be present in the animal. In this type of assay, test material is applied directly to the excised eye and any resulting injury is recorded. Such assays are the enucleated rabbit eye test (7) and the enucleated chicken eye assay (8). These tests have been developed by using the normally discarded eyes from food production species such as chickens or from laboratory animals that have been used for other purposes.

A second group of *in vitro* tests models only the cornea, which is a logical approach because maintenance of an intact transparent cornea is the major concern of ocular safety studies (severe damage to the cornea can lead to permanent blindness), and because damage to the cornea contributes more to the Draize score than does damage to any other ocular tissue. In fact, one recent study indicated that corneal score alone is an excellent predictor of total Draize score (9). An example of an assay that focuses primarily on detecting corneal damage is the bovine cornea opacity and permeability (BCOP) test (10). This assay uses corneas isolated from cattle used for meat production. With this model the amount of corneal opacity that has been induced by a test material can be quantitatively measured with an optical instrument, as opposed to the subjective estimation of opacity made by gross observation in the animal test. Damage to the cornea's barrier function can also be measured with the BCOP assay.

Another *in vitro* model designed to mimic corneal response is the tissue equivalent assay (TEA) (11). This assay uses a reconstructed, nonkeratinized epithelial-like tissue made of human cells,

upon which test materials can be directly placed. Such tissue is commercially available (MatTek Corporation, Ashland, MA and Skin Ethic, Nice, France), although it can also be produced by the individual investigator (12). The structure of this three-dimensional model is meant to simulate the epithelial covering of the cornea. Damage is estimated by measuring the viability of the human cells after treatment. Because this model is only a reconstructed tissue, it is not clear how closely it mimics the response of the epithelial layer of a normal cornea, but the model has been reasonably predictive of ocular irritation (Draize data) in recent studies (13).

### Cell-Based Assays

One step further in the reduction of the animal eye into less complicated *in vitro* models is accomplished using single cell or monolayer culture assays, which generally use epithelial cells similar to those that make up the outer surface of the cornea. If these cells are injured or killed in the animal, chemicals can more easily penetrate into the stroma of the cornea and cause additional damage. This penetration phenomenon is modeled by an *in vitro* assay called the fluorescein leakage test (14,15), in which a single layer of cells acts as a barrier to a common dye, fluorescein. If the cells are damaged or killed by a test chemical they lose their ability to act as a barrier to fluorescein. The subsequent movement of this dye through the cell layer can be measured and is an indicator of the amount of cell damage.

On a simpler level, one can simply observe the amount of cell killing that occurs in a single layer of human or animal cells (a cytotoxicity assay) and use this information to infer damage that might occur to the eye. Examples are the neutral

red release test (16), which uses short exposure times (such as might occur with an accidental splash to the eye followed by a quick rinse), and the neutral red uptake assay (17), which looks at longer exposure periods. An extensive review of these assays has been published as part of the U.S. IRAG evaluation (18).

Eye irritation can be reduced further by looking at damage to only the cell membranes. An example of this type of test is the red blood cell assay, where red blood cells are exposed to test material and membrane lysis is quantitated by measuring the amount of hemoglobin released (19). However, we should remember that the more reductionist the *in vitro* assay is, the more likely that it will only respond to certain classes of chemicals that are likely to cause eye irritation by the same mechanism. Thus only certain types of eye damage will be predicted by each *in vitro* test.

### Performance of the *in Vitro* Tests

Practical experience with the performance of various *in vitro* tests either in validation trials or in everyday use has revealed that considerable care must be taken when using the tests in a routine safety testing program. Each test seems to exhibit a slightly different level of sensitivity to correctly predict only a specific range of chemical classes. To consistently give a correct prediction, an *in vitro* test must do at least two things. First it must appropriately model the exposure kinetics, i.e., it must accept the test material in the same physical form as the animal test; it must be able to be exposed to the same concentration as in the animal test; and it must remain in contact with the test material for the same amount of time. Second, the end points that are developed for the *in vitro* assay must be predictive of the underlying *in vivo* tissue responses and this relationship must be clearly understood.

Currently, not all (or even the majority) of the *in vitro* tests fully meet these criteria. Nonetheless, if the use and interpretation of the tests is approached in an empirical fashion and attention is paid to certain key factors, the results can provide significant information for toxicologic evaluations. These key factors include type of product(s) to be used, physical characteristics of the product(s), expected level of toxicity, resolution required, intended use of the resulting data, and resources available.

An example of how to apply the above considerations is provided by examining the second factor, physical characteristics

of a product. Because there are two general forms of *in vitro* assays—those in which the substrate is completely immersed in growth medium (e.g., the neutral red uptake assay) and those in which the target surface is available for direct application of the test material (e.g., TEA)—the water solubility of a test material should first be considered. If a material is not water soluble, it would be fruitless to attempt to test it in a neutral red uptake assay because the test material will likely never actually come in contact with the target tissue. However, a topical application assay would be the logical choice because in this situation the test material will be applied directly to the surface of the target cells, ensuring exposure similar to what would occur *in vivo*. Examples of which assays are most suitable depending on the water solubility and form of the test material are as follows:

#### Water-soluble formulations

- BCOP assay
- Fluorescein leakage assays
- Neutral red uptake/release assays
- Chorioallantoic membrane (CAM)-based systems

#### Hydrophobic formulations

- Topical application assays
- BCOP assays
- Fluorescein leakage assays
- CAM-based systems

Also to be considered is whether to use dilution-based assays, i.e. assays in which serial dilutions of the test material are applied to the target tissue and the end point is the concentration that causes a certain response, or assays in which test material is only applied undiluted (neat) or at its in-use concentration.

Both types of assays have strengths and weaknesses, as can be seen in Tables 3 and 4.

A third consideration is the expected level of toxicity possessed by the test material. Ocular toxicity ranges from very slight conjunctival redness to full corrosive destruction of the three primary tissues. For a single *in vitro* test to address this full range, with the desired resolution, would be challenging. Most *in vitro* assays are designed to balance resolution with dynamic range. Dilution-based assays rely on the changes in concentration to provide sensitivity and dynamic range. However, they are limited in the types of materials that can be tested. In contrast, the topical application assays, in which the test materials are applied neat, use time of exposure (tissue constructs) or the robustness of the tissue (BCOP) to provide dynamic range.

**Table 3.** Advantages and disadvantages of dilution-based assays.

Advantages	Disadvantages
Rapid to execute	Cannot be used with water-insoluble materials
Most are machine scored	Dilution effects may mask toxicity of neat material
Generally very cost effective; materials are often batched (grouped together)	Change in the physical form, e.g., solids to solutions
Seem to work well with surfactants	Buffering effects of the medium may effect toxicity significantly
Often differentiate well between very mild materials	Possible reaction of the test material with the solvent

In these assays serial dilutions of the test material are applied to the test system and the end points are the concentrations of test material that cause a selected response.

**Table 4.** Advantages and disadvantages of topical application assays.

Advantages	Disadvantages
Material is tested in its native form, i.e., in the same form as an <i>in vivo</i> exposure	Test substrate can often be expensive
Exposure of the target tissue can be assured	Exposure times may be inconveniently long
In some models, exposure time can be selected to match expected <i>in vivo</i> exposure	

In these assays only the neat or in-use concentration test material is applied to the test system.


The tissue constructs provide high resolution for assessing potentially very mild (e.g., eye area cosmetics) to moderately irritating materials by extending the exposure times. Their resolution declines for the more aggressive materials because very short exposures (often a few seconds) are sufficient to kill the tissue. In contrast, the bovine cornea, with many layers of epithelium, does not resolve very mild products without excessively long exposures. However, it has the robustness to discriminate at the medium to high end of toxicity (20) (Figure 3). The double end points of opacity and permeability help the assay span the range from shampoos to industrial cleaners.

### Choosing the Appropriate Assay

The foregoing discussion illustrates that the choice of an *in vitro* ocular irritation assay is not simple. Unless the appropriate test is chosen, the subsequent results may be poor predictors of the actual ocular irritation potential of a test material. The previously mentioned key factors are extremely important and must be considered every time an *in vitro* test is contemplated. Table 5 gives some additional information (such as laboratory resources required and the technical skill needed) that should also be considered before choosing an *in vitro* ocular irritation test.

### Need for a Prediction Model

How will we eventually determine that an *in vitro* test produces valid information



**Figure 3.** Overlapping sensitivity of two types of *in vitro* assays. Diagram illustrating how the sensitivity of the BCOP assay, and assays using tissue constructs, cover different areas of the ocular irritation scale. Tissue constructs are effective at the milder end, whereas the BCOP can correctly characterize differences between more irritating materials.


that is at least as good as that produced by animal tests? Generally the accuracy of an *in vitro* model is assessed by identifying a number of materials that have been tested for ocular irritation in animals and then retesting the same materials in an *in vitro* assay. Both sets of data are then graphed, plotting the *in vitro* data on one axis and the animal data on the other axis, as in Figure 4. Each point on this graph represents the eye irritation score of a single material tested both *in vitro* and *in vivo*.

Knowing that both assays give variable results, we use error bars rather than a single point to represent the range of scores that might be expected if we repeated each test several times. The data in Figure 4 are from the U.S. Cosmetics, Toiletries, and Fragrance Association Phase III evaluation of surfactants and surfactant-based

**Table 5.** Further characterization of common *in vitro* assays

Test, reference	End point	<i>In vivo</i> tissue or irritation scale modeled	Resources needed	Technician skill level
BCOP (10)				
(Total score)	Opacity + permeability	Draize MAS	Specialized equipment, commercially available, spectrophotometer	General laboratory skills
(Opacity)	Opacity	Cornea		
(Permeability)	Permeability	Cornea/conjunctiva?		
Fluorescein leakage (14,15)	Permeability	Conjunctiva?, cornea	Specialized equipment	Tissue culture skills
HetCam (25)	Vascular damage, coagulation	Draize MAS, conjunctiva	General lab equipment	General laboratory skills
CAMVA (28)	Vascular damage, coagulation	Draize MAS, conjunctiva	General lab equipment	General laboratory skills
Tissue equivalent assay (11)	Cytotoxicity	Draize MAS, corneal epithelium damage	General tissue culture lab equipment	Some tissue culture skills
Neutral red uptake (17)	Cytotoxicity	Conjunctiva, corneal epithelium damage	General tissue culture lab equipment, 96-well plate reader	Tissue culture skills
Neutral red release (16)	Cytotoxicity/membrane damage?	Conjunctiva, corneal epithelium damage	General tissue culture lab equipment, 96-well plate reader	Tissue culture skills
Enucleated chicken eye (23)	Opacity, corneal swelling	Corneal damage	Specialized equipment	General laboratory skills
Cytosensor microphysiometer (23)	Cellular metabolism	Conjunctiva, Draize MAS	Specialized equipment (expensive)	Tissue culture skills
Irritation (Eytex) <sup>a</sup> (23)	Precipitation	Draize MAS	Specialized equipment	General laboratory skills
Pollen tube growth (26)	Cytotoxicity	Draize MAS	General lab equipment, spectrophotometer	General laboratory skills
Red blood cell (19)	Membrane lysis	Draize MAS	General lab equipment, spectrophotometer	General laboratory skills
SIRC (27)	Cytotoxicity	Conjunctiva	General tissue culture lab equipment	Tissue culture skills

Abbreviations: CAMVA, chorioallantoic membrane vascular assay; MAS, maximum average score. <sup>a</sup>InVitro International.



**Figure 4.** Surfactant formulations—BCOP permeability score. Graph showing the relationship between the Draize score of a number of surfactant-containing formulations with permeability measurements made in the *in vitro* BCOP assay. The error bars represent  $\pm 1$  standard deviation.  $r = 0.79$ .

products (21) and illustrate how data obtained from the *in vitro* BCOP assay compare with the data obtained using animals. This relationship between *in*

*vitro* and *in vivo* scores, i.e., the algorithm that allows one to predict an *in vivo* score from an *in vitro* score, is currently called the prediction model (22). Without such

a known relationship it is impossible not only to use an *in vitro* test correctly but also to conduct a rigorous validation (22).

Notice the variability of *in vivo* data in Figure 4, especially when compared to the more easily reproducible *in vitro* data. However, if a relationship is found to be good enough, the model could be used in a validation study to test the validity of both the assay and its prediction model. However, the difficult question is how good does this relationship have to be? Bruner et al. (5) provide a discussion of what could be expected if one Draize test is used to predict the results of a second Draize test.

Recently two large international validation trials of *in vitro* ocular irritation assays have been completed (23,24). The results from these studies were mixed, as might be expected from previous discussion in this paper concerning the difficulties involved in validating against the Draize test. Whereas the results from the first study [sponsored by the European Commission and the British Home Office (23)] were not encouraging, the results from the second [sponsored by the European Cosmetic, Toiletry and Perfumery Association (24)] indicated that three of the assays tested satisfied one or

more of the predetermined success criteria. There were major differences in the scope of each study, especially in the range of chemical classes tested; the second study was limited to chemicals commonly used in the cosmetics industry. More validation studies of this type can be expected in the coming years. It is prudent for researchers and toxicologists searching for appropriate *in vitro* assays to keep abreast of the findings of current validation studies—and to take notice of all the key factors involved in assay selection before deciding which test to use with any given set of test materials.

## Summary

Before a perfect, quantitative *in vitro* ocular irritation model is available, much basic work still needs to be done to understand mechanistically how injury happens in the human eye and how to model these

mechanisms *in vitro*. However, several of the appropriate *in vitro* models may already be in development. The upcoming results from current and planned validation studies may tell us just how close any of these assays are to meeting our goals.

Several points need to be emphasized as we evaluate the state of readiness of *in vitro* eye irritation assays:

- The eye is a very intricate organ made up of multiple tissues, each of which responds differently to injury. Current animal tests for ocular irritation use a complex scoring system involving three important ocular tissues.
- *In vitro* tests have generally been designed to model only one—or just a few—ocular tissues, not the whole eye. This is very helpful in obtaining more detailed mechanistic information about the process of eye irritation. However,
- Validation of *in vitro* ocular irritation assays will be difficult because the animal test is not very reproducible and because the animal test scores represent a combination of subjective observations of multiple ocular tissues.
- Different *in vitro* tests are suitable for different types of test materials and different ranges of toxicity. Careful consideration must be given to choosing the correct *in vitro* test for the required purpose.

## REFERENCES AND NOTES

1. Freeberg FE, Nixon GA, Reer PJ, Weaver JE, Bruce RD, Griffith JF, Sanders LW 3d. Human and rabbit eye responses to chemical insult. *Fundam Appl Toxicol* 7:626–634 (1986).
2. Draize JH, Woodard G, Calvery HO. Methods for the study of irritation and toxicity of substances applied topically to the skin and mucous membranes. *J Pharmacol Exp Ther* 82:377–390 (1944).
3. Annexes I, II, III and IV to Commission Directive 93/21/EEC of 27 April 1993 adapting to technical progress for the 18th time Council Directive 67/548/EEC on the approximation of the laws, regulations and administrative provisions relating to the classification, packaging and labeling of dangerous substances. *Off J Europ Commun* L110A:1–86 (1993).
4. Weil CS, Scala RA. Study of intra- and interlaboratory variability in the results of rabbit eye and skin irritation tests. *Toxicol Appl Pharmacol* 19:276–360 (1971).
5. Bruner LH, Carr GJ, Chamberlain M, Curren RD. Validation of alternative methods for toxicity testing. *Toxicol In Vitro* 10:479–501 (1996).
6. Bradlaw J, Gupta K, Green S, Hill R, Wilcox N. Practical application of non-whole animal alternatives: summary of IRAG workshop on eye irritation testing. Interagency Regulatory Alternatives Group. *Food Chem Toxicol* 35(1):175–178 (1997).
7. Burton ABD, York M, Lawrence RS. The *in vitro* assessment of severe eye irritants. *Food Cosmet Toxicol* 19:471–480 (1981).
8. Prinsen MK, Koeter HBWM. Justification of the enucleated eye test with eyes from slaughterhouse animals as an alternative to the Draize eye irritation test with rabbits. *Food Chem Toxicol* 31:69–76 (1993).
9. Lovell DP. Principal component analysis of Draize eye irritation tissue scores from 72 samples of 55 chemicals in the ECE-TOC data bank. *Toxicol In Vitro* 10:609–618 (1996).
10. Gautheron P, Dukic M, Alix D, Sina JF. Bovine corneal opacity and permeability test: an *in vitro* assay of ocular irritancy. *Fundam Appl Toxicol* 18:442–449 (1992).
11. Osborne R, Perkins MA, Roberts DA. Development and intralaboratory evaluation of an *in vitro* human cell-based test to aid ocular irritancy assessments. *Fundam Appl Toxicol* 28:139–153 (1995).
12. Kruszewski FH, Walker TL, Ward SL, Dipasquale LC. Progress in the use of human ocular tissues for *in vitro* alternative methods. *Comments Toxicol* 5:203–224 (1995).
13. De Silva O. Alternatives to eye irritation evaluation: the industry perspectives. In: *Alternatives to Animal Testing* (Lisansky SG, Macmillan R, Dupuis J, eds). Newbury, U.K.:CPL Press, 1996;67–73.
14. Tchao R. Trans-epithelial permeability of fluorescein *in vitro* as an assay to determine eye irritants. In: *Alternative Methods in Toxicology*. Vol 6 (Goldberg AM, ed). New York:Mary Ann Liebert, 1988.
15. Shaw AJ, Clothier RH, Balls M. Loss of trans-epithelial impermeability of a confluent monolayer of Madin-Darby canine kidney (MDCK) cells as a determinant of ocular irritancy potential. *ATLA* 18:145–151 (1990).
16. Reader SJ, Blackwell V, O'Hara R, Clothier RH, Griffin G, Balls M. A vital dye release method for assessing the short term cytotoxic effects of chemicals and formulations. *ATLA* 17:28–37 (1989).
17. Borenfreund E, Puerner JA. Toxicity determined *in vitro* by morphological alterations and neutral red absorption. *Toxicol Lett* 24:119–124 (1985).
18. Harbell JW, Koontz SW, Lewis RW, Lovell D, Acosta D. IRAG working group 4. Cell cytotoxicity assays. Interagency Regulatory Alternatives Group. *Food Chem Toxicol* 35:79–126 (1997).
19. Pape WJW, Pfannenbecker U, Hoppe U. Validation of the red blood cell test system as an *in vitro* assay for the rapid screening of irritation potential of surfactants. *Mol Toxicol* 1:525–536 (1987).
20. Swanson JE, Donnelly TA, Harbell JW, Huggins J. Prediction of ocular irritancy of full-strength cleaners and strippers by tissue equivalent and bovine corneal assays. *J Toxicol Cut Ocular Toxicol* 14:179–195 (1995).
21. Gettings SD, Lordo RA, Hintze KL, Bagley DM, Casterton PL, Chudkowski M, Curren RD, Demetrulias JL, DiPasquale LC, Earl LK et al. The CTFA evaluation of alternatives program: an evaluation of *in vitro* alternatives to the Draize primary eye irritation test. (Phase III). Surfactant-based formulations. *Food Chem Toxicol* 34:79–117 (1996).
22. Bruner LH, Carr G, Chamberlain M, Curren R. No prediction model, no validation study. *ATLA* 24:139–142 (1996).

23. Balls M, Botham PA, Bruner LH, Spielmann H. The EC/HO international validation study on alternatives to the Draize eye irritation test. *Toxicol In Vitro* 9: 871-929 (1995).
24. Brantom PG, Bruner LH, Chamberlain M, De Silva O, Dupuis J, Earl LK, Louell DP, Pape WJW, Uttley M, Bagley DM et al. A summary report of the COLIPA international validation study on alternatives to the Draize rabbit eye irritation test. *Toxicol In Vitro* 11:141-179 (1997).
25. Luepke NP. Hen's egg chorioallantoic membrane test for irritation potential. *Food Chem Toxicol* 23:287-291 (1985).
26. Kristen U, Kappler R. The pollen tube growth test. In: *In Vitro Testing Protocols. Methods in Molecular Biology*. Vol. 43 (O'Hare S, Atterwill CK, eds). Totowa, NJ:Humana Press, 1995.
27. North-Root H, Yackovich F, Demetrulias J, Gacula M, Heinze JE. Evaluation of an *in vitro* cell toxicity test using rabbit corneal cells to predict the eye irritation potential of surfactants. *Toxicol Lett* 14:207-212 (1982).
28. Bagley DM, Waters D, Kong BM. Development of a 10-day chorioallantoic membrane vascular assay as an alternative to the Draize rabbit eye irritation test. *Food Chem Toxicol* 32:1155-1160 (1994).