



SERVALL ANGLE CENTRIFUGES

Ask for Bulletin L-61
with complete data on all models



Large capacity
G type machines


Type G-1, for 6 bottles of 250 cc each
18 tubes of 60 cc each
12 tubes of 15 cc each
78 tubes of 10 cc each
Type G-2, for 4 bottles of 500 cc each
Type G-3, for 10 tubes of 100 cc each



Ask for
Bulletin L-62

STERLING AUTOMATIC AND HAND PIPETTES

Dispense accurately, at a set speed, pre-determined amounts of 0.1 cc to 10.0 cc.
Easy to Operate—Time Saving



Ask for
Bulletin L-63

ORIGINAL—ODHNER

Portable Calculator for the Scientist
Many different models

Efficient—Sturdy—Low Priced

Manufacturers and Distributors:

IVAN SORVALL, Inc.

210 FIFTH AVE.

NEW YORK 10, N. Y.

New and Valuable
Reference Work

Encyclopedia of MEDICAL SOURCES

By

EMERSON CROSBY KELLY

M.D., F.A.C.S.

*Associate Professor of Surgery
Albany Medical College, Editor, Medical Classes*


After several experiences in finding errors or incomplete information in references to original work, Dr. Kelly began to keep a list of references to medical eponyms and original investigations. The present compilation is the outgrowth of this list. In it he gives accurately with title and complete reference, the outstanding contributions to medical science, together with medical eponyms such as culture media, methods, signs, syndromes, etc. In all, the list includes over 5000 names and 150,000 references. More than ninety percent of the papers have been consulted in the original. The index includes disorders, operations, etc., and provides in effect a completely cross-referenced volume of medical sources. Not only will the book be of great use to medical students and practicing physicians, but librarians and research workers will find it invaluable.

481 pages \$7.50

**THE WILLIAMS
& WILKINS COMPANY**

Mt. Royal and Guilford Aves.
Baltimore 2, Maryland

(In writing to advertisers, please mention the journal it helps.)



International Refrigerated Centrifuge

Centrifugation under controlled temperatures at or well below 0°C . offers new possibilities to the biochemist. The cancer program—virus investigations—poliomyelitis research—enzyme research—low temperature solvent fractionation—all have benefited by the use of this new research tool. Countless other applications are waiting the guidance of the research scientist.

The International Portable Refrigerated Centrifuge, Model PR-1, is now equipped

with a larger capacity compressor. Centrifuge speeds and capacities remain unchanged, but much lower controlled temperatures are now possible. For example, with six 100 ml. tubes at 5000 r.p.m., temperatures as low as -15°C . can be maintained; with four 250 ml. bottles at 2600 r.p.m., -8°C .; with four 25 ml. tubes at 18,000 r.p.m., -5°C .

Here is a scientifically designed research instrument definitely worth investigating. Write for descriptive Bulletin RC, 1949.

INTERNATIONAL EQUIPMENT COMPANY
Boston 35, Massachusetts

(In writing to advertisers, please mention the journal it helps.)



for Microbiological Assay
of
VITAMINS and AMINO ACIDS

Bacto dehydrated media containing all the necessary nutrients and growth factors for the microbiological assay of vitamins and amino acids are now available from Difco. These basal media require only the addition of graduated amounts of the substance under assay to obtain linear growth of the test organism for construction of the standard curve. The vitamin or amino acid content of the material under assay is determined by adding appropriate concentrations of the test substance to the basal medium and comparing the growth response obtained with the standard.

BACTO-FOLIC ACID ASSAY MEDIUM
BACTO-NIACIN ASSAY MEDIUM
BACTO-PANTOTHENATE ASSAY MEDIUM
BACTO-PYRIDOXINE ASSAY MEDIUM

BACTO-RIBOFLAVIN ASSAY MEDIUM
BACTO-THIAMIN ASSAY MEDIUM
BACTO-TRYPTOPHANE ASSAY MEDIUM
BACTO-VITAMIN B12 ASSAY MEDIUM

“ ” ”

The method employed in carrying stock cultures of the test organisms and preparing the inoculum for microbiological assay is important. The following media have been developed especially for carrying stock cultures and for preparation of the inoculum.

BACTO-MICRO ASSAY CULTURE AGAR

BACTO-MICRO INOCULUM BROTH

BACTO-NEUROSPORA CULTURE AGAR

“ ” ”

BACTO-VITAMIN FREE CASAMINO ACIDS, dehydrated, is an acid hydrolysate of vitamin free casein prepared especially for laboratories investigating microbiological assay of vitamins.

Descriptive literature is available upon request.

Specify “DIFCO”

The Trade Name of the Pioneers in the Research and Development of Bacto-Peptone and Dehydrated Culture Media

DIFCO LABORATORIES
DETROIT 1, MICHIGAN