

Textbook of **VIROLOGY**

*Essentials of Virus and
Rickettsial Diseases*
SECOND EDITION

By A. J. Rhodes, M.D., F.R.C.P., and
C. E. van Rooyen, M.D., D.Sc., M.R.C.P.

Second edition rewritten from cover to cover,
including all the new developments in the field
of virology.

Covers: fundamental characteristics of virus
infections; exanthemata; skin, respiratory, ve-
nereal, eye diseases; arthropod-borne and
tropical fevers; infectious hepatitis and
serum jaundice; Coxsackie viruses; encephalo-
myocarditis and encephalomyelitis of ani-
mals; neurotropic virus diseases; rickettsial
infections.

560 pp. 76 figs. \$8.00

The Williams & Wilkins C.

Mt. Royal and Guilford Aves.
Baltimore 2, Maryland

BBL

Prepared Culture Media

for
Isolation and Cultivation
of

M. tuberculosis

LOWENSTEIN-JENSEN

DORSET

PETRAGNANI

ATS Medium

Folder PM-3007 Sent on Request

**BALTIMORE BIOLOGICAL
LABORATORY, INC.**

1640 Gorsuch Ave. • Baltimore 18, Md.

CUMULATIVE INDEX TO THE JOURNAL OF BACTERIOLOGY

Volumes 31 to 64

INCREASE THE USEFULNESS of your volumes of the Journal.

Locate the EXACT INFORMATION you want—quickly and easily.

Turn to ANY SUBJECT—any author—from any issue—with a minimum of time and effort.

Find ALL the articles that have appeared in the Journal on the topic of particular interest to you.

THE WILLIAMS & WILKINS COMPANY

Mt. Royal and Guilford Aves., Baltimore 2, Maryland

Please send: CUMULATIVE INDEX to the JOURNAL OF BACTERIOLOGY, Vols. 31-64. I enclose \$7.50.

Name _____

Address _____

In writing to advertisers, please mention the journal—it helps.

A New International Centrifuge



FOR THE *Modern Laboratory*

The new INTERNATIONAL MODEL U is more than just a Centrifuge in a cabinet. International's experienced engineering has produced a Centrifuge with years-ahead features and unmatched performance.

YOU will want to know more about the new aerodynamic guard bowl ventilating system—the new stainless steel interior—the stainless steel auxiliary cover that saves valuable bench space—the electric tachometer—the electric time clock—the unique air-filtering system—the spring-loaded hinge—the new vibration dampening system which results in smooth

operation without the “shaker” action of over-flexible mountings—the versatility of 84 interchangeable accessory combinations.

Let us tell you more about this completely new Centrifuge, handsomely designed to fit in with your wall counters and cabinets. Write for Bulletin J explaining all these features in detail.

INTERNATIONAL EQUIPMENT COMPANY

1284 SOLDIERS FIELD ROAD, BOSTON 35, MASS.

The logo consists of the word "DIFCO" in a bold, serif font, enclosed within a dark, horizontally-oriented oval shape with a slightly irregular, hand-drawn appearance.

Bacto-Dehydrated Media *for Microbiological Assay of* **VITAMINS and AMINO ACIDS**

These media contain all the necessary nutriment for the growth of specified test organisms for the microbiological assay of vitamins and amino acids except for the component under assay. These basal media require only the addition of specified increasing amounts of the vitamin or amino acid being assayed to obtain a growth response which may be measured by acidimetric or turbidimetric methods for the construction of standard curves. The vitamin or amino acid content of the material under assay is determined by adding appropriate concentrations of the test substance to the basal medium and comparing the growth response obtained with that of the standard.

BACTO-RIBOFLAVIN ASSAY MEDIUM	BACTO-CHOLINE ASSAY MEDIUM
BACTO-NIACIN ASSAY MEDIUM	BACTO-CF ASSAY MEDIUM
BACTO-THIAMIN ASSAY MEDIUM	BACTO-TRYPTOPHANE ASSAY MEDIUM
BACTO-PANTOTHENATE ASSAY MEDIUM	BACTO-LEUCINE ASSAY MEDIUM
BACTO-B ₁₂ ASSAY MEDIUM USP	BACTO-METHIONINE ASSAY MEDIUM
BACTO-CS VITAMIN B ₁₂ AGAR	BACTO-LYSINE ASSAY MEDIUM
BACTO-FOLIC ACID ASSAY MEDIUM	BACTO-ISOLEUCINE ASSAY MEDIUM
BACTO-PYRIDOXINE ASSAY MEDIUM	BACTO-ARGININE ASSAY MEDIUM
BACTO-BIOTIN ASSAY MEDIUM	BACTO-TYROSINE ASSAY MEDIUM
	BACTO-CYSTINE ASSAY MEDIUM

The method employed in carrying stock cultures of the test organisms and preparing the inoculum for microbiological assay is important. The following media have been developed especially for carrying stock cultures and for preparation of the inoculum:

BACTO-MICRO ASSAY CULTURE AGAR	BACTO-B ₁₂ CULTURE AGAR USP
BACTO-MICRO INOCULUM BROTH	BACTO-B ₁₂ INOCULUM BROTH USP
BACTO-NEUROSPORA CULTURE AGAR	

BACTO-VITAMIN FREE CASAMINO ACIDS, dehydrated, is an acid hydrolysate of vitamin free casein prepared especially for laboratories investigating microbiological assay of vitamins.

Descriptive literature available upon request

DIFCO LABORATORIES
DETROIT 1, MICHIGAN