Low copy number plasmids for regulated low-level expression of cloned genes in *Escherichia coli* with blue/white insert screening capability

Claude G.Lerner* and Masayori Inouye


Department of Biochemistry, Robert Wood Johnson Medical School – UMDNJ at Rutgers, 675 Hoes Lane, Piscataway, NJ 08854- 5635, USA

Submitted June 13, 1990

We have constructed pCL1920 and pCL1921, a pair of low copy number plasmids which contain a 580 bp BstUI fragment that carries the *lac* promoter/operator, multiple cloning sites and *lacZ* fragment of pUC19 (1) cloned in place of the polylinker region in pGB2 (2), a pSC101 derived plasmid which confers spectinomycin (50 μ g/ml) and streptomycin (100 μ g/ml) resistance in Escherichia coli. All multiple cloning sites indicated are unique except for an additional EcoRI site as shown in the figure. pCL1920 and pCL1921 contain the BstUI fragment in opposite orientations with respect to the pGB2 sequences. In the absence of inducer the pCL1920/21 vectors do not produce detectable levels of β -galactosidase in JM105 (lacI^q lacZ Δ M15) (1) cells (less than 2 Miller units) (3). In the presence of 2 mM IPTG (isopropyl- β -D-thiogalactopyranoside) the β - galactosidase levels of the pCL1920/21 [JM105] transformants rose to 11 units. while the pUC19 [JM105] transformants produced 470 units; a 43 fold increase. These results are consistent with the expected 40 fold difference in plasmid copy number between pCL1920/21 (5 copies per cell) compared to that of the pUC vectors (200 copies per cell). Thus the pCL1920 and pCL1921 vectors allow regulated low-level expression of genes inserted downstream of the lac promoter-operator when transformed into strains containing the lac gene. They should also be useful for cloning genes which may be deleterious at high copy number. Since the pCL1920/21 vectors are compatible with ColE1 derived plasmids they can be used to form stable co-transformants together with pBR322 or pUC derived plasmids. For blue/white screening of inserts competent host cells with the lacZ Δ M15 gene are used, and the transformation mixture is plated on LB, spectinomycin plates pre-spread with 5 μ l of 0.2 M IPTG and 25 μ l 40 mg/ml X-gal (5-bromo-4-chloro-indolyl- β -D-galactopyranoside) per plate. We will provide these plasmids in CL83, a recA-(recA56) derivative of JM83 (1).

REFERENCES

- 1. Yanisch-Perron, C., Vieira, J. and Messing, J. (1985) Gene 33, 103-119.
- 2. Churchward, G., Belin, D. and Nagamine, Y. (1984) Gene 31, 165-171.
- Miller, J.H. (1972) Experiments in Molecular Genetics, pp. 352-355. Cold Spring Harbor Laboratories, Cold Spring Harbor, NY.



^{*} To whom correspondence should be addressed