

Supplementary Methods:

Oligonucleotide probe retrieval assay:

Transfection: Control and WRN-depleted cells were transiently transfected (*via* precipitation of DNA with calcium phosphate) with primer-template DNA probes in which the template strand harbors a single *TaqI* (TCGA) restriction site, and the primer strand contains a 5' terminal biotin moiety. Both strands are protected from exonucleolytic digestion by the presence of dideoxycytidine at their 3' ends (Fig. 9). Twenty-four hrs post-transfection, DNA was extracted from total cell (cyto in Fig. 9) or nuclear lysates. Transfected DNA probes were isolated from total DNA by binding to streptavidin-coated magnetic beads (Invitrogen).

DNA probes/substrates:

- a) Primer-template partial duplex DNA was generated by annealing a single-stranded oligonucleotide to M13mp18 ssDNA. As illustrated in Fig. 9, the oligonucleotide (43 nt) contains sequence that, in addition to being complementary to M13 DNA at its 3' end, has a unique barcoded sequence and a 5' terminal biotin cap. Linear M13mp18 ssDNA (949 nt) was obtained from circular phage DNA by annealing two synthetic oligonucleotides at base positions 2071-2110 and 3021-3060, and digesting the resultant double-stranded segments with the restriction endonucleases, *BaeGI* and *AfeI*, respectively. The 5'-biotinylated oligomer was hybridized to the linear ssDNA fragment (position 2670-2693) generated by restriction digestion.
- b) The D-loop structure was formed by annealing a synthetic oligomer (43 nt) to a restriction fragment derived from digesting plasmid pCEP4 DNA (Invitrogen) with *AatII*. The sequence at the 3' end of the oligomer is complementary to the insertion site of plasmid DNA (position 6272-6295). This oligomer also has a unique barcoded sequence and a biotin residue at the 5'-terminus (Fig. 9). The 5'-biotin tagged oligomer was mixed in 100-fold molar excess over plasmid DNA to increase the probability of

base pairing with plasmid DNA and facilitate displacement of the complementary plasmid DNA strand upon denaturation (Fig. 9).

qPCR: A sensitive quantitative PCR protocol, developed and implemented in our laboratory, was used to measure the mutation frequency at the *TaqI* restriction site (1,2). *In vivo* elongated DNA substrates were first extensively digested with *TaqI* restriction enzyme to eliminate DNA containing wild-type sequence. The DNA was then amplified using primers flanking the *TaqI* restriction site (Table 1). This generated a cycling threshold value for amplification of the test/target site [Ct (test)]. Similarly, using internal control primers that flank a region adjacent to the *TaqI* site (Table 1), we obtained control Ct values. The mutation frequency (MF) was calculated from the measured Ct values; MF (per base) = $1/(2^{\Delta C_t} \times 4)$, where $\Delta C_t = C_t(\text{control}) - C_t(\text{test})$.

References:

1. Bielas, J. H., and Loeb, L. A. (2005) *Nat Methods* **2**, 285-290
2. Wright, J. H., Modjeski, K. L., Bielas, J. H., Preston, B. D., Fausto, N., Loeb, L. A., and Campbell, J. S. (2011) *Nucleic Acids Res* **39**, e73

Supplementary Table 1: Primers Used in Oligonucleotide Retrieval Assay


Primer used to anneal to linear M13 ssDNA:
5'-biotin-GCACGTCAGGCACGGCGTCACATTCAACCGATTGAGGGAGGGA-3'
Oligonucleotides hybridized to M13 ssDNA for restriction cleavage:
Bae G1: 5'-TCAGTGCCTTGAGTAACAGTGCCCGTATAAACAGTTAATG-3'
Afe 1: 5'-AAGTCAGAGGGTAATTGAGCGCTAATATCAGAGAGATAAC-3'
qPCR Primers for M13 ssDNA elongation:
1. Test: 5'-GCACGTCAGGCACGGCGTC-3' (forward)
5'-AACGCGCTACAGTCTGACGCTAAA-3' (reverse)
2. Control: 5'-GCACGTCAGGCACGGCGTC-3' (forward)
5'-TCCGGCCTTGCTAATGGTAATGGT-3' (reverse)
Primer used to form D-loop DNA:
5'-biotin GCACGTCAGGCACGGCGTCTGAAGATCAGTTGGGTGCACGAGT-3'
qPCR Primers for D-loop elongation:
1. Test: 5'-GCACGTCAGGCACGGCGTC-3' (forward)
5'-AATAGTGTATGCGGCGACCGAGTT-3' (reverse)
2. Control: 5'-AACTCGGTGCGCCGCATACACTATT-3' (forward)
5'-TGTCAGAAGTAAGTTGGCCGCACT-3' (reverse)

Supplementary Figure Legends:

Supplementary Fig. 1: Synergistic degradation of bubble DNA requires the simultaneous presence of WRN and Pol δ . WRN or Pol δ (50-100 fmols of each) was incubated with bubble DNA at 37 °C for 20 min. Thereafter, the reaction mixtures were incubated at 63 °C for 10 min to render the enzymes inactive. Control reactions demonstrated that this temperature does not denature bubble DNA, but does abolish the enzymatic activities of WRN and Pol δ . Following inactivation, Pol δ or WRN was added to reactions that contained WRN or Pol δ , respectively. The reactions were incubated for an additional 20 min at 37 °C. As a control, WRN and Pol δ were incubated together with bubble DNA. W, WRN; δ , Pol δ .

Supplementary Fig. 2: Synergism with Pol δ requires more than just the exonuclease activity of WRN. Increasing amounts (0.3, 1.5, 6 pmols, indicated as 1, 2, and 3, respectively) of a minimal exonuclease domain of WRN (aa 1-240) were pre-incubated with Pol δ holoenzyme (60 fmol). The reactions were initiated by the addition of bubble DNA and incubated at 37 °C for 20 min as described in the legend to Fig. 1.

Supplementary Fig. 3: Co-operation of WRN with DNA Pol δ is not observed on all DNA substrates and with other 3'→5' exonucleases. WRN, Pol δ , or WRN + Pol δ (~6 fmol of each) were incubated with a D-loop mimic (**A**), or Trex 1 (0.25 fmols) and/or WRN (15 fmols) were incubated with bubble DNA (**B**) at 37 °C and processed as described in Materials and Methods. Band intensities of hydrolytic products were quantified and plotted using ImageJ software.





Supplementary Figure 1
Kamath-Loeb et al.

NWRN	-	1	2	3	-	1	2	3
δ	-	-	-	-	+	+	+	+



Supplementary Figure 2
Kamath-Loeb et al.



Supplementary Figure 3
Kamath-Loeb et al.