

# An Azobenzene Photoswitch Sheds Light on Turn Nucleation in Amyloid- $\beta$ Self-Assembly

Todd M. Doran,<sup>1</sup> Elizabeth A. Anderson,<sup>1</sup> Sarah E. Latchney,<sup>2</sup> Lisa A. Opanashuk,<sup>2</sup>

Bradley L. Nilsson<sup>\*,1</sup>


Department of Chemistry, University of Rochester, Rochester, NY 14627, USA

## ELECTRONIC SUPPORTING INFORMATION


### Contents:

<b>Figure S1.</b> Analytical HPLC trace (215 nm) of <i>trans cis-2</i> in water .....	S2
<b>Figure S2.</b> Analytical HPLC trace (215 nm) of <i>trans cis-3</i> in water .....	S2
<b>Figure S3.</b> Digital images of aggregate pellets from peptides <b>1</b> and <i>trans-2</i> and <i>cis-2</i> after Congo red staining .....	S3
<b>Figure S4.</b> Congo red absorbance spectra in the presence and absence of <i>trans-2</i> and <i>cis-2</i> .....	S3
<b>Figure S5.</b> FT-IR spectra of aggregates derived from photoisomerized peptides .....	S4
<b>Figure S6.</b> Analytical HPLC trace (215 nm) of peptide <b>1</b> .....	S4
<b>Figure S7.</b> Analytical HPLC trace (215 nm) of <i>trans-2</i> .....	S5
<b>Figure S8.</b> Analytical HPLC trace (215 nm) of <i>cis-2</i> .....	S5
<b>Figure S9.</b> Analytical HPLC trace (215 nm) of <i>trans-3</i> .....	S5
<b>Figure S10 .</b> Analytical HPLC trace (215 nm) of <i>cis-3</i> .....	S6
<b>Table S1.</b> Analytical HPLC conditions for peptides <b>1–3</b> .....	S6
<b>Figure S11.</b> MALDI-TOF mass spectrum of peptide <b>1</b> .....	S7
<b>Figure S12.</b> MALDI-TOF mass spectrum of peptide <b>2</b> .....	S7
<b>Figure S13.</b> MALDI-TOF mass spectrum of peptide <b>3</b> .....	S8
<b>Table S2.</b> Calculated and observed <i>m/z</i> ratios for peptides <b>1–3</b> .....	S8
<b>Figure S14.</b> Concentration curve for peptide <b>1</b> .....	S9
<b>Figure S15.</b> Concentration curve for peptides <b>2–3</b> .....	S9


**Figure S1.** Analytical HPLC trace (215 nm) of *trans*-**2** to *cis*-**2** conversion in water. The peak at 12.1 min is *cis*-**2** and the peak at 12.2 minutes is residual *trans*-**2**.




**Figure S2.** Analytical HPLC trace (215 nm) of *trans*-**3** to *cis*-**3** conversion in water. The peak at 11.9 minutes is *cis*-**3** and the peak at 12.3 minutes is residual *trans*-**3**.




**Figure S3.** Digital images of aggregate pellets from peptides **1** and *trans-2* and *cis-2* after Congo red staining




**Figure S4.** A. Congo red absorbance spectra in the presence and absence of *trans-2* and *cis-2*. B. Congo red absorbance spectra for *trans-2* and *cis-2* with the spectrum for Congo red alone removed.




**Figure S5.** FT-IR spectra of photoswitched aggregates of peptide **2**. *Trans*-**2** was self-assembled into fibrils and photoirradiated to induce switching to *cis*-**2**, resulting in a coincident change in FTIR (significant loss of signal). The *cis*-**2** aggregates were then re-irradiated to *trans*-**2**, providing restoration of the characteristic  $\beta$ -sheet amide I stretch.




**Figure S6.** Analytical HPLC trace (215 nm) of A $\beta$ 42




**Figure S7.** Analytical HPLC trace (215 nm) of *cis-2*




**Figure S8.** Analytical HPLC trace (215 nm) of *trans-2*



**Figure S9.** Analytical HPLC trace (215 nm) of *cis-3*




**Figure S10.** Analytical HPLC trace (215 nm) of *trans*-3




**Table S1.** Analytical HPLC conditions for peptides **1–3**

Peptide	Retention time (min)	Gradient (solution A: water/0.5%TFA; solution B: acetonitrile/0.5% TFA)
<b>1</b>	11.80	Isocratic 5% B over 5 minutes, increase 5-95% B over 15 minutes, maintain at 95% B over 5 minutes
<i>Trans</i> -2	12.27	Isocratic 5% B over 5 minutes, increase 5-95% B over 15 minutes, maintain at 95% B over 5 minutes
<i>Cis</i> -2	12.14	Isocratic 5% B over 5 minutes, increase 5-95% B over 15 minutes, maintain at 95% B over 5 minutes
<i>Trans</i> -3	12.00	Isocratic 5% B over 5 minutes, increase 5-95% B over 15 minutes, maintain at 95% B over 5 minutes
<i>Cis</i> -3	11.9	Isocratic 5% B over 5 minutes, increase 5-95% B over 15 minutes, maintain at 95% B over 5 minutes


**Figure S11.** MALDI-TOF mass spectrum of A $\beta$ 42



**Figure S12.** MALDI-TOF mass spectrum of peptide **2**.



**Figure S13.** MALDI-TOF mass spectrum of peptide **3**




**Table S1.** Calculated and observed  $m/z$  ratios for peptides **1–3**

Peptide	Calculated $m/z$	Observed $m/z$
<b>1</b>	4515.0	4516.2
<b>2</b>	4507.0	4500.5
<b>3</b>	4565.1	4566.5



**Figure S14.** Concentration curve for peptide 1



**Figure S15.** Concentration curve for peptides 2–3

