



Figure S4. AT-FLAG can be immunoprecipitated by either anti-AT or anti-FLAG.

Cells stably expressing Sortilin or Sort.LAYA were transiently transfected with pATM-FLAG (ATM), pATZ-FLAG (ATZ) or the empty vector pRc/RSV (V). After 24 h, cells were labeled for 20 min with ^{35}S -amino acids, followed by a 20 min chase. AT levels in cell lysates were analyzed by immunoprecipitation with either anti-FLAG or anti-AT and subject to SDS-PAGE analysis. Arrowheads denote the positions of non-specific species.