



### Supplementary Figure S1


**Supplemental Figure S1. Examples of peak detection in biological probes.** Examples of detected and quantified BP metabolite peaks in skin models are marked in blue and shown for each of the analyzed BP derivatives.



**7-OH-BP**




**BP**




**Supplementary Figure S2**

**A**



**B**



**Supplementary Figure S2. The metabolic capacity toward cinnamic aldehyde in MatTek *in vitro* 3D skin models is comparable to human skin.** (A) Rates of cinnamic aldehyde conversion into cinnamic alcohol in the MatTek EpiDerm and EpiDerm full-thickness model compared to the rate in human dermatomized skin. (B) Rates of cinnamic aldehyde conversion into cinnamic acid in skin models and skin. All 3 models were topically exposed to cinnamic aldehyde (1  $\mu\text{mol}/\text{cm}^2$ ) diluted in DMSO. Metabolite levels were determined by HPLC after incubation for 1 and 2 h, essentially as described elsewhere (Smith *et al.*, 2001). This figure displays the sum of metabolites detected in both incubation media and LiberaseFT-digested tissue homogenates ( $n=3$ ; obtained with 3 different batches/donors, mean  $\pm$  SD). After cinnamic aldehyde exposure for 1 h the overall recovery was found at 96, 82, and 85% of the initially applied dose for EpiDerm, EpiDermFT and skin *ex vivo*, respectively.