

**Table 6. Patient Infected Status – Female *Chlamydia trachomatis***

PIS	NAAT1		NAAT2		Xpert			Symptom Status		Total
	SW	UR	SW	UR	VS	ES	UR	Symp	Asymp	
NI <sup>b</sup>	-	-	-	-	-	-	-	534	1044	1579
NI	-	-	-	-	IND <sup>c</sup>	-	-	4	5	9
NI	-	-	-	-	-	IND	-	2	8	10
NI	-	-	-	-	-	-	IND	1	1	2
NI	-	-	-	-	+	+	-	0	1	1
NI	-	-	-	-	+	-	-	1	1	2
NI	-	-	-	-	-	+	-	1	0	1
NI	-	-	-	EQ <sup>d</sup>	-	-	-	1	8	9
NI	-	-	EQ	-	-	-	-	2	2	4
NI	-	-	-	+	-	-	-	0	5	5
NI	-	-	+	-	-	-	-	2	0	2
NI <sup>f</sup>	-	+	-	+	-	-	+	1	2	3
NI	-	-	+	-	-	+	-	0	1	1
NI	-	+	-	-	-	-	-	1	0	1
NI	+	-	-	-	-	-	-	1	4	5
NI	+	-	-	-	+	-	-	0	1	1
NI	+	-	-	-	-	+	+	0	2	2
NI	+	-	-	-	-	+	-	0	1	1
NI	+	+	-	-	-	-	+	0	1	1
NI	+	+	-	-	+	+	+	0	1	1
NI	+	+	-	-	+	-	+	1	0	1
NI	-	-	EQ	-	-	-	IND	1	0	1
NI	-	-	-	EQ	-	IND	IND	1	0	1
NI	-	EQ	-	EQ	-	-	-	1	0	1
<b>Total Non-Infected</b>								555	1088	1643
I <sup>e</sup>	+	+	+	+	+	+	+	28	43	71
I <sup>f</sup>	-	+	-	+	-	-	+	1	2	3
I	+	-	+	+	+	+	+	1	0	1
I	+	+	-	+	+	-	+	0	1	1
I	+	+	+	-	+	+	+	1	1	2
I	+	+	+	-	+	+	-	0	1	1
I	+	+	+	+	-	+	+	0	1	1
I	+	+	+	+	+	-	+	0	1	1
I	+	+	+	+	+	IND	+	0	1	1
<b>Total Infected</b>								31	51	82

<sup>a</sup>PIS = Patient Infected Status; SW = Swab; UR = urine; VS = Vaginal Swab; ES = Endocervical Swab

<sup>b</sup>NI = Non-infected

<sup>c</sup>IND = Indeterminate – ERROR, INVALID or NO RESULT by Xpert CT/NG Assay; specimens with IND results by Xpert are not included in the 2x2 tables for that specimen type.

<sup>d</sup>EQ = Equivocal result for this individual specimen type only; PIS status determined based on remaining specimens.

<sup>e</sup>I = Infected

<sup>f</sup>These samples are infected for urine and non-infected for swabs. In this table they appear twice.

**Table 7. Patient Infected Status – Female *Neisseria gonorrhoeae***

PIS	NAAT1		NAAT2		Xpert			Symptom Status		Total
	SW	UR	SW	UR	VS	ES	UR	Symp	Asymp	
NI <sup>b</sup>	-	-	-	-	-	-	-	553	1096	1649
NI	-	-	-	-	IND <sup>c</sup>	-	-	4	5	9
NI	-	-	-	-	-	IND	-	2	9	11
NI	-	-	-	-	-	-	IND	1	1	2
NI	-	-	-	-	+	-	+	0	1	1
NI	-	-	-	-	+	-	-	1	0	1
NI	-	-	EQ <sup>d</sup>	-	-	-	-	2	2	4
NI	-	-	-	+	-	-	-	0	1	1
NI	-	-	+	-	-	-	-	4	2	6
NI <sup>f</sup>	-	+	-	+	-	-	+	1	0	1
NI	-	-	+	+	-	-	-	1	0	1
NI	+	-	-	-	-	-	-	1	0	1
NI	-	-	EQ	-	-	-	IND	1	0	1
NI	-	-	-	EQ	-	-	-	2	8	10
NI	-	-	-	EQ	-	IND	IND	1	0	1
NI	-	EQ	-	EQ	-	-	-	1	0	1
<b>Total Non-Infected</b>								575	1125	1700
I <sup>e</sup>	+	+	+	+	+	+	+	8	7	15
I	+	+	+	-	+	+	+	1	2	3
I <sup>f</sup>	-	+	-	+	-	-	+	1	0	1
I	+	-	+	-	+	+	-	0	1	1
I	+	-	+	+	+	+	+	0	1	1
I	+	+	-	+	+	+	+	1	0	1
I	+	+	+	EQ	+	+	+	0	1	1
<b>Total Infected</b>								11	12	23

<sup>a</sup>PIS = Patient Infected Status; SW = Swab; UR = urine; VS = Vaginal Swab; ES = Endocervical Swab

<sup>b</sup>NI = Non-infected

<sup>c</sup>IND = Indeterminate – ERROR, INVALID or NO RESULT by Xpert CT/NG Assay; specimens with IND results by Xpert are not included in the 2x2 tables for that specimen type.

<sup>d</sup>EQ = Equivocal result for this individual specimen type only; PIS status determined based on remaining specimens.

<sup>e</sup>I = Infected

<sup>f</sup>These samples are infected for urine and non-infected for swabs. In this table they appear twice.