

Analytical Evaluation of GeneXpert CT/NG, the First Genetic Point-of-Care Assay for Simultaneous Detection of *Neisseria gonorrhoeae* and *Chlamydia trachomatis*

Sepehr N. Tabrizi,^{a,b} Magnus Unemo,^c Daniel Golparian,^c Jimmy Twin,^b Athena E. Limnios,^d Monica Lahra,^d Rebecca Guy,^e
on behalf of the TTANGO Investigators

Department of Microbiology and Infectious Diseases, The Royal Women's Hospital, Parkville, Victoria, Australia^a; Department of Obstetrics and Gynaecology, University of Melbourne, The Royal Women's Hospital and Murdoch Children's Research Institute, Parkville, Victoria, Australia^b; WHO Collaborating Centre for Gonorrhoea and other STIs, Department of Laboratory Medicine, Microbiology, Örebro University Hospital, Örebro, Sweden^c; WHO Collaborating Centre for STD, Microbiology Department, The Prince of Wales Hospital, Randwick, Sydney, Australia^d; Department of Microbiology and the Kirby Institute, University of NSW, Sydney, Australia^e

GeneXpert CT/NG was evaluated with 372 characterized bacterial strains. Sensitivity of 10 genome copies/reaction was obtained for both agents. Four *Neisseria mucosa* and two *Neisseria subflava* isolates were positive for one of two gonococcal targets; however, the assay flagged all as negative. The assay was analytically highly sensitive and specific.

Point-of-care (POC) assays which are affordable, sensitive, specific, user friendly, rapid and robust, equipment free, and deliverable (ASSURED) (1) are needed in many settings for the accurate diagnosis of *Neisseria gonorrhoeae* and *Chlamydia trachomatis* infections. In the past, POC tests have generally relied upon lateral flow technology, and most tests have demonstrated suboptimal performance (2, 3). Cepheid GeneXpert CT/NG is a new Food and Drug Administration (FDA)-approved rapid molecular assay for simultaneous detection of *N. gonorrhoeae* and *C. trachomatis*. This assay is the first genetic point-of-care assay that amplifies one chromosomal target (CT1) for the detection of *C. trachomatis*, two chromosomal targets (NG2 and NG4) for detection of *N. gonorrhoeae*, a single-copy human gene which should be present in each specimen to act as a sample adequacy control (SAC), and *Bacillus globigii* DNA added to each cartridge to serve as a sample-processing/internal control (SPC).

Prior to implementing any assay for diagnostics, it is important to comprehensively evaluate its performance. As Cepheid GeneXpert CT/NG is a molecular assay, consistent with evaluations of laboratory-based nucleic acid amplification tests (NAATs), it should initially be evaluated with characterized strains to ensure that the analytical sensitivity and specificity are optimal. Several NAATs, including commercial assays, have experienced problems with false positivity and false negativity due to cross-reaction (4–7), mutations (8, 9), and genetic exchanges with related organisms (10–12). In particular, diagnosis of extragenital specimens, such as pharyngeal and rectal samples, can result in suboptimal specificity due to cross-reaction with commensal *Neisseria* species or *Neisseria meningitidis* (4, 5, 7, 13–15), which result in low positive predictive values, especially in low-prevalence populations.

In order to ensure accuracy of NAAT results, strategies such as utilization of dual targets in the assays (16, 17) or inclusion of a supplemental “confirmation” have been incorporated in the testing protocol (15, 18, 19). These approaches are, in particular, important for detection from the rectum and oropharynx, which is crucial particularly in population groups of men who have sex with men (20) or with history of oral and/or anal sex. Currently, extragenital specimens are not licensed for use by the FDA; however, several studies have conducted an evaluation and shown ap-

propriate sensitivity with these types of samples using the available commercial NAATs (21–28). Nevertheless, most studies evaluating the sensitivity and specificity of the gonococcal and chlamydial NAATs utilize clinical samples whereby the gold standard for confirmation of *N. gonorrhoeae* is a consensus NAAT gold standard (22). In some cases, cross-reactions are not evident, particularly if the gold-standard NAATs also cross-react with the nongonococcal isolate, giving a false-positive result in the assay.

In this study, we evaluated the analytical sensitivity and specificity of this new POC assay, in particular to evaluate cross-reactivity with nongonococcal *Neisseria* species or other closely related bacteria.

N. gonorrhoeae isolates and isolates of nongonococcal *Neisseria* or closely related bacterial species were obtained from two international *Neisseria* reference laboratories to allow more geographically, temporally, and genetically diverse isolates to be tested. There were 89 isolates from the *Neisseria* Reference Laboratory at the World Health Organization Collaborating Centre for STD in Sydney, Australia, and 246 isolates from the WHO Collaborating Centre for Gonorrhoea and other STIs, Swedish Reference Laboratory for Pathogenic *Neisseria*, Örebro, Sweden. Assay performance was evaluated using 111 well-characterized ATCC and clinical isolates of *N. gonorrhoeae* (including 3 mutants with an *N. meningitidis* *porA* gene) (10, 12) from Asia, Asia-Pacific, Australia, North America, South America, Africa, and Europe. In addition, we tested 236 isolates of nongonococcal *Neisseria* species or closely related bacterial species and 25 *C. trachomatis* strains representing all 15 main serovars (Table 1). All gonococcal and chlamydial isolates were prepared and placed at –80°C prior to testing as described previously (5, 29).

The *N. gonorrhoeae* and *C. trachomatis* genome copy numbers were determined for each DNA extract by quantitative PCR

Received 24 March 2013 Accepted 27 March 2013

Published ahead of print 3 April 2013

Address correspondence to Sepehr N. Tabrizi, sepehr.tabrizi@thewomens.org.au.

Copyright © 2013, American Society for Microbiology. All Rights Reserved.

doi:10.1128/JCM.00806-13

TABLE 1 Positivity rate for the detection of characterized isolates of *Neisseria gonorrhoeae*, nongonococcal *Neisseria* and closely related species, and *Chlamydia trachomatis*

Bacterial species tested	No. of isolates tested	No. of GeneXpert CT/NG positive results	
		CT	NG
<i>N. gonorrhoeae</i>	108	0	108
<i>N. gonorrhoeae</i> <i>porA</i> pseudogene mutant ^a	3	0	3
All <i>N. gonorrhoeae</i>	111	0	111
<i>N. animalis</i>	1	0	0
<i>N. caviae</i>	1	0	0
<i>N. cinerea</i>	15	0	0
<i>N. elongata</i>	1	0	0
<i>N. flava</i>	1	0	0
<i>N. flavescens</i>	6	0	0
<i>N. kochii</i> ^b	4	0	0
<i>N. lactamica</i>	30	0	0
<i>N. meningitidis</i> ^c	97	0	0
<i>N. mucosa</i> ^d	11	0	0 ^d
<i>N. perflava</i>	1	0	0
<i>N. pharyngis</i>	1	0	0
<i>N. polysacchara</i>	2	0	0
<i>N. sicca</i>	4	0	0
<i>N. subflava</i> ^e	42	0	0 ^e
<i>N. weaveri</i>	1	0	0
Other commensal <i>Neisseria</i> ^f	5	0	0
<i>Kingella denitrificans</i>	3	0	0
<i>Moraxella catarrhalis</i>	6	0	0
<i>Moraxella osloensis</i>	1	0	0
Other <i>Moraxella</i> species ^f	3	0	0
All commensal isolates	236	0	0
<i>C. trachomatis</i> serovar A	1	1	0
<i>C. trachomatis</i> serovar B	1	1	0
<i>C. trachomatis</i> serovar Ba	1	1	0
<i>C. trachomatis</i> serovar C	1	1	0
<i>C. trachomatis</i> serovar D	1	1	0
<i>C. trachomatis</i> serovar E	1	1	0
<i>C. trachomatis</i> serovar F	1	1	0
<i>C. trachomatis</i> serovar G	1	1	0
<i>C. trachomatis</i> serovar H	1	1	0
<i>C. trachomatis</i> serovar I	1	1	0
<i>C. trachomatis</i> serovar J	1	1	0
<i>C. trachomatis</i> serovar K	1	1	0
<i>C. trachomatis</i> serovar L1	2	2	0
<i>C. trachomatis</i> serovar L2	2	2	0
<i>C. trachomatis</i> serovar L2b	1	1	0
<i>C. trachomatis</i> serovar L3	2	2	0
<i>C. trachomatis</i> serovar nvCT 2006–2011 ^g	6	6	0
All <i>C. trachomatis</i>	25	25	0

^a Representing three isolates that harbor an *N. meningitidis* *porA* gene sequence instead of the gonococcal *porA* pseudogene (10, 12).

^b *N. kochii* (or *N. gonorrhoeae* subspecies *kochii*) is a very rare commensal organism.

^c Representing isolates of all major meningococcal clones spreading worldwide, including serogroups A, B, C, W-135, Y, and 29E.

^d Four out of 11 *N. mucosa* strains amplified 1 (NG4) of the 2 NG targets.

^e Two out of 41 *N. subflava* strains amplified 1 (NG4) of the 2 NG targets.

^f *Neisseria* or *Moraxella* species that were not able to be identified to the species level.

^g Representing the Swedish new variant of *C. trachomatis* (8, 9) isolates obtained over a 6-year period.

(qPCR) in conjunction with a quantified standard (Advanced Biotechnologies Inc. [ABI], Columbia, MD). To each cartridge, 1 ml of GeneXpert CT/NG swab transport reagent, 5 μ l of extracted DNA from quantitated bacterial strains, and 2 μ l of A549 (human alveolar basal epithelial) cell line DNA were added to simulate sample conditions. Approximately 1×10^5 copies of the *N. gonorrhoeae* genome, 2×10^5 copies of other *Neisseria* or related strains, and 5×10^4 copies of *C. trachomatis* were tested in each cartridge.

Strict procedures were followed to avoid specimen contamination and carryover. All strains were deidentified and blind tested.

All 111 *N. gonorrhoeae* isolates were positive for both NG targets (NG2 and NG4) included in the assay and were interpreted as positive. The limit of detection of *N. gonorrhoeae* was 10 genome copies per reaction. Overall, 236 isolates of nongonococcal *Neisseria* species or closely related species were tested. The results of GeneXpert CT/NG were matched with the corresponding organism (Table 1). None of these isolates tested resulted in false-positive results. However, four out of 11 *Neisseria mucosa* isolates and two of 42 *Neisseria subflava* isolates were positive in one (NG4) of two NG targets. Upon retesting, the NG4 target was still positive for all 6 strains. GeneXpert CT/NG assay software interpreted all of these 6 isolates, reactive to one of the two *N. gonorrhoeae* targets, as negative. Accordingly, the NG4 target appears to have some homology to sequences in some of the *N. mucosa* and *N. subflava* strains. As strains of these two species can be present in extragenital samples, the presence of two targets is especially useful when analyzing such samples. The GeneXpert CT/NG assay was recently evaluated for rectal samples and was shown to have high clinical specificity in this sample type (30).

All 15 serovars of *C. trachomatis* were detectable to 10 genome copies per reaction. The GeneXpert CT/NG assay was also able to detect the Swedish new variant of *C. trachomatis* (nvCT) (8, 9) and the L2b strain. The presence of human DNA in all samples was documented by virtue of the sample adequacy control, and valid test results were generated from each analyte. Assay time was 87 min for each sample.

Like commercial and in-house laboratory-based NAAT systems, the GeneXpert CT/NG assay relies on detection of nucleic acid sequences, and thus, any changes in the target may result in loss of detection. There are several examples of such changes in *C. trachomatis* (8, 9) and *N. gonorrhoeae* (10–12). The closed, self-contained, fully integrated and automated GeneXpert CT/NG system was analytically highly sensitive and specific for detection of all types of genetically different *C. trachomatis* and *N. gonorrhoeae* strains. Also, the incorporation of two *N. gonorrhoeae* targets, both of which would need to be positive for reporting a positive result, allows reporting with more confidence.

ACKNOWLEDGMENTS

We thank Cepheid for providing the diagnostic kits for this study.

Test, Treat and Go (TTANGO) Investigators include Rebecca Guy, James Ward, John Kaldor, Basil Donovan, David Wilson, Handan Wand, David Regan, and Louise Causer from the Kirby Institute, UNSW, Belinda Hengel from Apunipima Health Council, Lisa Natoli and David Anderson from the Burnet Institute, David Whitley from the Queensland Pediatric Infectious Diseases (QPID) Laboratory, Sepehr N. Tabrizi from the Royal Women's Hospital, Mark Shephard from Flinders University, and Christopher Fairley from the University of Melbourne and Melbourne Sexual Health Centre.

TTANGO is funded by Australian National Health and Medical Research Council project grant number 1009902.

REFERENCES

- Peeling RW, Holmes KK, Mabey D, Ronald A. 2006. Rapid tests for sexually transmitted infections (STIs): the way forward. *Sex. Transm. Infect.* 82(Suppl 5):v1–v6.
- Hislop J, Quayyum Z, Flett G, Boachie C, Fraser C, Mowatt G. 2010. Systematic review of the clinical effectiveness and cost-effectiveness of rapid point-of-care tests for the detection of genital chlamydia infection in women and men. *Health Technol. Assess.* 14:1–97.
- van Dommelen L, van Tiel FH, Ouburg S, Brouwers EE, Terporten PH, Savelkoul PH, Morre SA, Bruggeman CA, Hoebe CJ. 2010. Alarming poor performance in *Chlamydia trachomatis* point-of-care testing. *Sex. Transm. Infect.* 86:355–359.
- Palmer HM, Mallinson H, Wood RL, Herring AJ. 2003. Evaluation of the specificities of five DNA amplification methods for the detection of *Neisseria gonorrhoeae*. *J. Clin. Microbiol.* 41:835–837.
- Tabrizi SN, Unemo M, Limnios AE, Hogan TR, Hjelmenvoll SO, Garland SM, Tapsall J. 2011. Evaluation of six commercial nucleic acid amplification tests for detection of *Neisseria gonorrhoeae* and other *Neisseria* species. *J. Clin. Microbiol.* 49:3610–3615.
- Whiley DM, Garland SM, Harnett G, Lum G, Smith DW, Tabrizi SN, Sloots TP, Tapsall JW. 2008. Exploring 'best practice' for nucleic acid detection of *Neisseria gonorrhoeae*. *Sex. Health* 5:17–23.
- Upton A, Bromhead C, Whiley DM. 13 March 2013. A *Neisseria gonorrhoeae* false-positive result from a pharyngeal swab using the Roche Cobas 4800 CT/NG assay, New Zealand 2012. *J. Clin. Microbiol.* doi:10.1128/JCM.00485-13.
- Ripa T, Nilsson P. 2006. A variant of *Chlamydia trachomatis* with deletion in cryptic plasmid: implications for use of PCR diagnostic tests. *Euro Surveill.* 11(45):pii=3076. <http://www.eurosurveillance.org/ViewArticle.aspx?ArticleId=3076>.
- Ripa T, Nilsson PA. 2007. A *Chlamydia trachomatis* strain with a 377-bp deletion in the cryptic plasmid causing false-negative nucleic acid amplification tests. *Sex. Transm. Dis.* 34:255–256.
- Golparian D, Johansson E, Unemo M. 2012. Clinical *Neisseria gonorrhoeae* isolate with a *N. meningitidis* *porA* gene and no prolyliminopeptidase activity, Sweden, 2011: danger of false-negative genetic and culture diagnostic results. *Euro Surveill.* 17(9):pii=20102. <http://www.eurosurveillance.org/ViewArticle.aspx?ArticleId=20102>.
- Ison CA, Golparian D, Saunders P, Chisholm S, Unemo M. 13 December 2012. Evolution of *Neisseria gonorrhoeae* is a continuing challenge for molecular detection of gonorrhoea: false negative gonococcal *porA* mutants are spreading internationally. *Sex. Transm. Infect.* doi:10.1136/sextrans-2012-050829.
- Whiley DM, Limnios A, Moon NJ, Gehrig N, Goire N, Hogan T, Lam A, Jacob K, Lambert SB, Nissen MD, Sloots TP. 2011. False-negative results using *Neisseria gonorrhoeae* *porA* pseudogene PCR—a clinical gonococcal isolate with an *N. meningitidis* *porA* sequence, Australia, March 2011. *Euro Surveill.* 16(21):pii=19874. <http://www.eurosurveillance.org/ViewArticle.aspx?ArticleId=19874>.
- Alexander S. 2009. The challenges of detecting gonorrhoea and chlamydia in rectal and pharyngeal sites: could we, should we, be doing more? *Sex. Transm. Infect.* 85:159–160.
- McNally LP, Templeton DJ, Jin F, Grulich AE, Donovan B, Whiley DM, Cunningham PH. 2008. Low positive predictive value of a nucleic acid amplification test for nongenital *Neisseria gonorrhoeae* infection in homosexual men. *Clin. Infect. Dis.* 47:e25–e27.
- Smith DW, Tapsall JW, Lum G. 2005. Guidelines for the use and interpretation of nucleic acid detection tests for *Neisseria gonorrhoeae* in Australia: a position paper on behalf of the Public Health Laboratory Network. *Commun. Dis. Intell. Q. Rep.* 29:358–365.
- Moller JK, Pedersen LN, Persson K. 2010. Comparison of the Abbott RealTime CT new formulation assay with two other commercial assays for detection of wild-type and new variant strains of *Chlamydia trachomatis*. *J. Clin. Microbiol.* 48:440–443.
- Stevens MP, Tan SE, Horvath L, Fairley CK, Garland SM, Tabrizi SN. 2008. Absence of a *Chlamydia trachomatis* variant, harbouring a deletion in the cryptic plasmid, in clients of a sexually transmissible infection clinic and antenatal patients in Melbourne. *Commun. Dis. Intell. Q. Rep.* 32:77–81.
- Johnson RE, Newhall WJ, Papp JR, Knapp JS, Black CM, Gift TL, Steece R, Markowitz LE, Devine OJ, Walsh CM, Wang S, Gunter DC, Irwin KL, DeLisle S, Berman SM. 2002. Screening tests to detect *Chlamydia trachomatis* and *Neisseria gonorrhoeae* infections—2002. *MMWR Recomm. Rep.* 51(RR-15):1–38.
- Workowski KA, Berman S. 2010. Sexually transmitted diseases treatment guidelines, 2010. *MMWR Recomm. Rep.* 59(RR-12):1–110.
- Weinstock H, Workowski KA. 2009. Pharyngeal gonorrhea: an important reservoir of infection? *Clin. Infect. Dis.* 49:1798–1800.
- Bachmann LH, Johnson RE, Cheng H, Markowitz L, Papp JR, Palella FJ, Jr, Hook EW, III. 2010. Nucleic acid amplification tests for diagnosis of *Neisseria gonorrhoeae* and *Chlamydia trachomatis* rectal infections. *J. Clin. Microbiol.* 48:1827–1832.
- Bachmann LH, Johnson RE, Cheng H, Markowitz LE, Papp JR, Hook EW, III. 2009. Nucleic acid amplification tests for diagnosis of *Neisseria gonorrhoeae* oropharyngeal infections. *J. Clin. Microbiol.* 47:902–907.
- Lister NA, Smith A, Tabrizi S, Hayes P, Medland NA, Garland S, Fairley CK. 2003. Screening for *Neisseria gonorrhoeae* and *Chlamydia trachomatis* in men who have sex with men at male-only saunas. *Sex. Transm. Dis.* 30:886–889.
- Lister NA, Tabrizi SN, Fairley CK, Garland S. 2004. Validation of Roche COBAS AmpliCor assay for detection of *Chlamydia trachomatis* in rectal and pharyngeal specimens by an *omp1* PCR assay. *J. Clin. Microbiol.* 42:239–241.
- Moncada J, Schachter J, Liska S, Shayevich C, Klausner JD. 2009. Evaluation of self-collected glans and rectal swabs from men who have sex with men for detection of *Chlamydia trachomatis* and *Neisseria gonorrhoeae* by use of nucleic acid amplification tests. *J. Clin. Microbiol.* 47:1657–1662.
- Nagasawa Z, Ikeda-Dantsuji Y, Niwa T, Miyakoshi H, Nagayama A. 2010. Evaluation of APTIMA Combo 2 for cross-reactivity with oropharyngeal *Neisseria* species and other microorganisms. *Clin. Chim. Acta* 411:776–778.
- Ota KV, Tamari IE, Smieja M, Jamieson F, Jones KE, Towns L, Juzkiw J, Richardson SE. 2009. Detection of *Neisseria gonorrhoeae* and *Chlamydia trachomatis* in pharyngeal and rectal specimens using the BD Probetec ET system, the Gen-Probe Aptima Combo 2 assay and culture. *Sex. Transm. Infect.* 85:182–186.
- Schachter J, Moncada J, Liska S, Shayevich C, Klausner JD. 2008. Nucleic acid amplification tests in the diagnosis of chlamydial and gonococcal infections of the oropharynx and rectum in men who have sex with men. *Sex. Transm. Dis.* 35:637–642.
- Stevens MP, Twin J, Fairley CK, Donovan B, Tan SE, Yu J, Garland SM, Tabrizi SN. 2010. Development and evaluation of an *ompA* quantitative real-time PCR assay for *Chlamydia trachomatis* serovar determination. *J. Clin. Microbiol.* 48:2060–2065.
- Goldenberg SD, Finn J, Sedudzi E, White JA, Tong CY. 2012. Performance of the GeneXpert CT/NG assay compared to that of the Aptima AC2 assay for detection of rectal *Chlamydia trachomatis* and *Neisseria gonorrhoeae* by use of residual Aptima samples. *J. Clin. Microbiol.* 50:3867–3869.