



## A Comparison of the Analytical level of agreement of Nine Treponemal Assays for Syphilis and Possible Implications for Screening Algorithms

Journal:	<i>BMJ Open</i>
Manuscript ID:	bmjopen-2013-003347
Article Type:	Research
Date Submitted by the Author:	05-Jun-2013
Complete List of Authors:	Castro, Arnold; Centers for Disease Control and Prevention, laboratory Reference and Research Branch Jost, Heather; Centers for Disease Control, LRRB Cox, David; Centers for Disease Control and Prevention, Fakile, Yetunde; Centers for Disease Control and Prevention, Kikkert, Susan; Centers for Disease Control and Prevention, Tun, Ye; Centers for Disease Control and Prevention, Zaidi, Akbar; Centers for Disease Control and Prevention, Park, Mahin; Georgia Department of Community Health,
<b>Primary Subject Heading</b>:	Infectious diseases
Secondary Subject Heading:	Diagnostics, Infectious diseases, Immunology (including allergy)
Keywords:	Microbiology < BASIC SCIENCES, Immunology < BASIC SCIENCES, MICROBIOLOGY

SCHOLARONE™  
Manuscripts

Only

**Title page**

**A Comparison of the Analytical level of agreement of Nine Treponemal Assays for Syphilis and Possible Implications for Screening Algorithms**

**Corresponding author:**

Arnold R Castro. Ph.D. Centers for Disease Control and Prevention. 1600 Clifton Rd. MS A-12. Atlanta, Ga. 30333. USA. Phone: 404-639-2874. Fax: 404-639-4664. Email: [acastro@cdc.gov](mailto:acastro@cdc.gov).

**Co-authors:**

Heather A. Jost. B.S. Centers for Disease Control and Prevention. Atlanta, Ga. USA. Tel: 404-639-2462. Fax: 404-639-4664. Email: [hjost@cdc.gov](mailto:hjost@cdc.gov)

David L Cox. Centers for Disease Control and Prevention. Atlanta, Ga. USA. Tel: 770-595-1277. Email: [drdavide@me.com](mailto:drdavide@me.com)

Yetunde F. Fakile. Ph.D. Centers for Disease Control and Prevention. Atlanta, Ga. USA. Tel: 404-639-3784. Fax: 404-639-4664. Email: [yxfo@cdc.gov](mailto:yxfo@cdc.gov)

Susan Kikkert. B.A. Centers for Disease Control and Prevention. Atlanta, Ga. USA. Tel: 404-639-2871. Fax: 404-639-4664. Email: [sjk4@cdc.gov](mailto:sjk4@cdc.gov)

Ye Tun. Ph.D. Centers for Disease Control and Prevention. Atlanta, Ga. USA. Tel: 404-639-7640. Fax: 404-639-xxxx. Email: [gwr5@cdc.gov](mailto:gwr5@cdc.gov)

Akbar A. Zaidi. Ph.D. Centers for Disease Control and Prevention. Atlanta, Ga. USA. Tel: 404-639-8185. Fax: 404-498-3420. Email: [aaz1@cdc.gov](mailto:aaz1@cdc.gov)

Mahin M. Park. Ph.D. Georgia Department of Community Health. Email: [mparkdhr.state.ga.us](mailto:mparkdhr.state.ga.us)

**Address for Reprints:**

Arnold R Castro. Ph.D. Centers for Disease Control and Prevention. 1600 Clifton Rd. MS A-12. Atlanta, Ga. 30333. USA. Phone: 404-639-2874. Fax: 404-639-4664. Email: [acastro@cdc.gov](mailto:acastro@cdc.gov).

**Word counts:**

Summary: 16

Abstract: 240

Text: 1,703

Number of references: 8

Number of tables: 2

Keywords: Comparison analytical agreement of syphilis treponemal assays

## Summary:

This study describes the effect of the comparison of eight qualitative and quantitative treponemal syphilis assays.

## Cover Letter

Dear Sirs,

We are including for your consideration, the manuscript entitled “A Comparison of the Analytical level of agreement of Nine Treponemal Assays for Syphilis and Possible Implications for Screening Algorithms.” This work describes the effect that occurs when eight treponemal syphilis assays are evaluated both qualitatively and quantitatively to demonstrate the variability on analytical sensitivities and to determine the feasibility in selecting the proper assay for a reverse algorithm system.

## Competing Interests

One of the functions of the Laboratory Reference and Research Branch (LRRB) of the Division of STD Prevention of the National Center for HIV/Aids, Viral Hepatitis, STD and TB Prevention of the Centers of Disease Control and Prevention is to evaluate public domain procedures to determine its diagnostic value in proving or disproving a particular claim.

## Sources of support

There are no sources of support for this study that require acknowledgement.

## Duplication

Please note that, there is no related paper or publications in reference to this work that describes the effect of the analytical level of agreement of nine treponemal assays for syphilis serology and its possible implications in selecting the proper assays for the reverse screening algorithm.

**Acknowledgments:** The authors would like to thank the Georgia Department of Health laboratories for supplying the panel of serum samples and the management of the Division of STD Prevention and the Laboratory Reference and Research Branch at CDC for supplying all the necessary diagnostic reagents and materials for this study.

1  
2  
3  
4  
5  
6  
7  
8  
9  
10  
11  
12  
13  
14  
15  
16  
17  
18  
19  
20  
21  
22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

**Funding:**

The Centers for Disease Control and Prevention provided financial support for this study in the form of allowing the principal investigator and CDC coauthors to devote time and effort in conducting this study.

**Conflict of interests:** There is no conflict of interest related to this study.

**Ethics approval:** This study was performed on serum samples supplied by the Georgia Department of Health laboratories from which all patient identifiers had been removed. Protocol for studies using these sera were reviewed by institutional review board of the CDC (protocol no 2018) and determined to be exempt from further review.

**Contributors:**

DLC was the project microbiologist, developed the test protocol, design the test parameters and functional components of the test.

ARC assisted in the developing of the test protocol, designed the test parameters and was the corresponding author.

HAI performed the evaluation of the assays.

YT assisted in the evaluation of the assays.

KK assisted in performing the evaluation.

YFF assisted in performing the evaluation.

AAZ assisted in the statistic evaluation of the data.

All authors contributed to the write up and critically reviewed the manuscript and approved the final draft.

**Provenance and peer review:** Not commissioned; externally peer reviewed.

## Abstract

**Objective:** The serological diagnosis of syphilis requires the detection of two distinct antibodies, the nontreponemal and treponemal. CDC recommends screening first with a nontreponemal test such as (RPR/VDRL), and then confirming those results with one of several treponemal tests (FTA-ABS, EIA, CIA, TP-PA or POC). Due to the high volume of samples processed by some laboratories using automated systems, the screening with treponemal assays and confirming with nontreponemal tests is becoming the established norm. The purpose of this study was to evaluate eight treponemal assays using TP-PA as the predicate assay.

**Methods:** Two hundred ninety stored serum samples were tested qualitatively according to the manufacturer's directions.

**Results:** Concordance with specimens tested as reactive or nonreactive using TP-PA was: FTA-ABS 94.5% - 100%, Trep-Sure 100% - 98.9%, Bioelisa 100% - 98.9%, INNO-LIA 99.1% - 99.4%, Bioline 100% - 98.9%, Captia IgG 100% - 97.2%, Trep-ID 100%-100%, and LIAISON 100% - 99.4%. In order to properly evaluate the performance of these assays, the analytical sensitivity was determined by endpoint titration of serial dilutions of the reactive serum samples in normal sera. The median endpoint titer varied from 1:4 for FTA-ABS to 1:512 for Trep-Sure.

**Conclusions:** The performance of the treponemal serological assays was comparable when using medium and high titer sera. However, the varying performance on specimen dilutions suggest that there may be differences in sensitivity with low titer sera that are more prevalent in primary and late syphilis cases.

## Introduction

Syphilis is a sexually transmitted disease caused by the bacterium *Treponema pallidum*. The clinical diagnosis is difficult due to the complexity of manifestations which requires the aid of serological interpretation. The serological diagnosis of syphilis depends upon the detection of two distinct antibodies, the nontreponemal or heterophile antibodies (reagin) directed against cardiolipin released from damaged host cells and from the treponemes themselves. These antibodies can also be present in other diseases and human conditions such as (Lupus, malaria, HIV, IV drug users, etc.) (1). The presence of nontreponemal antibodies is indicative of active infection, thus a reduction in titer can suggest a successful antibiotic therapy, and a significant increase can indicate a possible relapse or reinfection (2). The treponemal antibodies are primarily directed against specific lipoprotein antigens

of the bacterium such as 15-, 17-, and 47-kDa. Even after treatment or time, these treponemal antibodies usually remain present for life. A positive treponemal test cannot distinguish between active, treated, and old cases of syphilis.

The traditional algorithm is to screen with a nontreponemal assay such as VDRL or RPR and those serum samples found reactive are then confirmed using a treponemal test. With the introduction of automated systems which are appropriate for high output volumes, the reverse algorithm is gaining acceptability and its usefulness is arguable justified (6). However, the adoption of this reverse algorithm has increased the number of discrepant results between the screening and confirmatory tests (4, 5). One reason the traditional algorithm was established was to avoid the detection of previously treated cases, especially in low prevalence settings where the majority of sera will test negative and not require further evaluation, unless a patient has been reinfected or in rare circumstances, their non-treponemal antibody titers remain serofast. With the application of the reverse algorithm even in low prevalence settings all previously treated cases will likely be detected. When the screening treponemal test is positive and the confirmatory nontreponemal test is negative, the CDC guidelines recommend using a second treponemal test to validate the results of the screening test. The selection of a second confirmatory test introduces the possibility of having an analytically less sensitive test which would suggest a false positive screening or a less sensitive second confirmatory test. Then it is important to determine if all available treponemal assays are comparable in analytical sensitivity to avoid uncertainty in their selection (3).

**Methods**

Two hundred ninety stored serum samples randomly selected from our serum bank were used for this study. Originally the serum samples were obtained from the Georgia Public Health Laboratory with all identifiers removed and under a CDC IRB approval. Due to the prevalence of high lipid contents of the sera, the samples were treated with Cleanascite HC (PureBiotech. Middlesex, NJ) (7). The patterns of reactivity of the test panel were determined at CDC by a quantitative RPR (Becton Dickinson, Baltimore, MD) and a qualitative TP-PA (Fujirebio Diagnostics, Inc., Malvern, PA). All serum samples were divided into aliquots and frozen at -20°C so that the number of freeze-thaw cycles was consistent for each sera tested. Eight of the treponemal assays: FTA-ABS; (Zeus

Scientific, Raritan, NJ), LIAISON Treponema Assay (DiaSorin, Inc., Stillwater, MN), SD Bioline Syphilis 3.0 Rapid Point of Care (POC) test (Standard Diagnostics, Inc., Korea), INNO-LIA immunoblot (Innogenetics, Gent, Belgium), Bioelisa (BioKit, Barcelona, Spain) and the Captia IgG, Trep-ID, and Trep-Sure (Trinity Biotech, Jamestown, NY) were tested qualitatively, then quantitatively with reactive samples. In the quantitative test, two-fold serial dilutions (from 1:2 to 1:16,384) were prepared using normal nonreactive human plasma converted to serum (8) as a diluent. All dilutions were treated as neat sera and tested with the nine different assays according to the manufacturer's package insert. The test was performed blind, the EIA results were collated and compared to the TP-PA used as the predicate for the qualitative test.

## Statistical Methods.

The endpoint titer data was converted to dilution data from the exponential scale (titer) to a linear scale (dilution) by taking the logarithm of the endpoint titer, i.e.,  $\log_2(\text{titer}) = \text{dilution}$ . This transformation reduces the variance within the test and between the tests. A generalized linear model was fitted to these data by the Generalized Estimating Equations (GEE) method to test for the differences between the mean dilutions of these tests. Adjustments for multiple comparisons between the mean dilutions of these tests were made by Tukey-Kramer method. A p-value of  $<0.05$  was used as a significance level for the differences. All statistical analyses were conducted using SAS, version 9.3 (SAS Institute, Cary, North Carolina).

## Results

All treponemal syphilis assays were performed according to the manufacturer's direction circulars using two hundred ninety stored serum samples of which 109 were reactive and 181 nonreactive for TP-PA. Concordance with specimens tested as reactive or nonreactive using TP-PA was: FTA-ABS 94.5% - 100%, INNO-LIA 99.1% - 99.4%, LIAISON 100% - 99.4%, Trep-Sure 100% - 98.9%, Bioelisa 100% - 98.9%, Bioline 100% - 98.9%, Captia IgG 100% - 97.2%, Trep-ID 100%-100%. These results indicated that most of the assays were comparable to TP-PA when the serum samples were tested undiluted. However, FTA-ABS failed to detect six TP-PA reactive samples and INNO-LIA one, while the other four assays missed none. One hundred and nine stored serum samples



reactive by TP-PA were serially diluted in nonreactive human plasma converted to serum (8), with the purpose of demonstrating the endpoint titer of each sample with the nine different assays. With a simple comparison the endpoint titers indicated that Trep-Sure had a median value of 1:512, Bioelisa 1:128, LIAISON 1:64, Trep-ID 1:64, Bioline 1:32, INNO-LIA 1:16, TP-PA 1:16, CAPTIA IgG 1:8 and FTA-Abs 1:4 (Table 2).

Of the 109 TP-PA reactive sera, 95(85%) were nonreactive for the confirmatory RPR test. Fourteen samples were reactive for both TP-PA and RPR. All other assays were able to detect treponemal antibody in those undiluted sera that were reactive by the TP-PA. However, five sera found to be reactive by the CAPTIA IgG test and one by the BIOLINE were found to be non-reactive by all other tests suggesting false positives. All other sera found to be non-reactive by the TP-PA but reactive by another test could be considered as true positives since the sera were found to be reactive by at least two other tests (see Table 1). If the TP-PA would have been selected as the screening test for the reverse algorithm then eight samples reactive by the other assays would have been missed resulting in a false serological result.

Mean dilutions and their standard errors of eight quantitative treponemal assays and the TP-PA estimated by the generalized linear model are shown in Table 2. There were significant ( $p<0.0001$ ) differences between the mean dilutions of all tests, except between TP-PA and INNO-LIA ( $p=1.000$ ) and Trep-ID and LIASON ( $p=1.000$ ). This analysis was repeated after subtracting the FTA-Abs titer from the titer of the corresponding test. In this analysis also there was no significant difference between TP-PA and INNO-LIA ( $p=1.000$ ) and Trep-ID and LIASON ( $p=1.000$ ), and the differences among the rest of the pairs were significant ( $p<0.0001$ ). (Table 2)

**Discussion:**

The reverse sequence syphilis screening algorithm is being used because it takes advantage of automated systems. However, the selection of these assays to screen serum samples with suspected cases of syphilis presents problems that were not clearly anticipated (6). In cases where the treponemal test is positive and the confirmatory RPR test is negative, the CDC recommends that a second treponemal test be used to determine the serological validity of the screening test (5). If the second treponemal test is negative, then it is inconclusive whether the first screening test is a false



positive or a more sensitive test. This situation presents a dilemma for physicians and understandably they will be asking for clarification as to whether it is or not a case of syphilis in asymptomatic situations. Consequently, it would be advisable for the laboratorian to consider the selection of more sensitive screening and confirmatory treponemal assays to avoid having discrepant results. This is not an easy task because of the different assay platforms presently used for screening in the reverse algorithm system ranging from FTA-ABS, EIA, CIA, TP-PA and POC. In this study an attempt was made to determine the analytical level of agreement of the assays with the purpose of better understanding the abilities of the assays for detection of low level antibodies. When the nine assays were tested quantitatively, the FTA-ABS was the treponemal test with the lowest median endpoint titer (1:4). In contrast, Trep-Sure had a median endpoint titer of (1:512) (Table2). Since the FTA-ABS was the least reactive, its reevaluation should be considered.

Further studies with documented clinical samples are necessary to demonstrate the correlation between the antibody titer and the stage of the disease. With these comparison results we then could infer that there are marked differences between the analytical sensitivities of these assays especially with low titer sera such as those found in primary or latent cases. The variability of results with the different assays may be due in part to the composition of their systems. Those assays using the double antigen or sandwich configuration such as Trep-Sure, Bioelisa, Trep-ID and LIAISON where the recombinant antigen such as Tp15, 17, 47 KD or combination of them are both attached to the ELISA plate or latex particle and are directly conjugated to horse radish peroxidase (HRP) or chemiluminescence (CIA). This configuration tends to be more sensitive because they are able to detect both IgG and IgM antibodies in the sera of patients with syphilis. In addition, the detection of IgM antibodies amplifies the signal due to its multiple binding sites. On the other hand, normally those assays using the indirect or anti-human system such as Captia (IgG) and FTA-ABS are limited in that they only detect one single immunoglobulin.

When selecting the reverse algorithm system for syphilis diagnostic serology and in cases when the treponemal test is reactive and the RPR nonreactive, it would be advisable to consider the selection of more analytically sensitive screening and confirmatory treponemal assays preferably those based on the double antigen or sandwich method, such as Trep-Sure, Bioelisa, Trep-ID and LIASON to avoid the necessity of using a less sensitive second confirmatory treponemal assay.

REFERENCES

1. Larsen, S. A. (1989). "Syphilis." Clin.Lab.Med. **9**: 545-557.

2. Larsen, S. A., V.Pope, et al. (1998). A Manual of Tests for Syphilis. Washington, DC, American Public Health Association.

3. Zhang WM, Yen-Lieberman, Means K, et al. The Impact of Analytical Sensitivity on Screening Algorithms for Syphilis. Clinical Chemistry. 2012; 58:6.

4. Binnicker MJ, Jespersen DJ, Rollins LO. Direct comparison of the traditional and reverse syphilis screening algorithms in a population with a low prevalence of syphilis. J Clin Microbiol 2012;**50**:148-50.

5. Centers for Disease, C. and Prevention (2011). "Discordant results from reverse sequence syphilis screening--five laboratories, United States, 2006-2010." MMWR Morb Mortal Wkly Rep **60**(5): 133-137.

6. Binnicker MJ. Which algorithm should be used to screen for syphilis?. Curr Opin Infect Dis 2012; 25:79-85.

7. Castro AR, Morrill WE, Pope V. Lipid Removal from Human Serum Samples. Clin Diag Lab Immonol. 2000;7:197-199.

8. Castro AR, Kikkert SE, Fears MB, Pope V. Defibrination of Blood Plasma for use in the Serological Tests for Syphilis. Clin Diag Lab Immonol. 2002;9:1376-1378.

TABLE 1. Percent concordance of eight qualitative syphilis treponemal test compared to the predicate TP-PA

TP-PA	Total no of samples	FTA-ABS		INNOLIA		Liaison		Trep-Sure		Bioelisa		Bioline		Captia (IgG)		Trep-ID	
Predicate		R <sup>a</sup>	NR <sup>b</sup>	R	NR	R	NR	R	NR	R	NR	R	NR	R	NR	R	NR
Reactive	109	103	6(94.4)	108	1(99.1)	109	0(100)	109	0(100)	109	0(100)	109	0(100)	109	0(100)	109	0(100)
Nonreactive	181	0	181(100)	1	180(99.4)	1	180(99.4)	2	179(98.9)	2	179(98.9)	2	179(98.9)	5	176(97.2)	0	181(100)

<sup>a</sup> R, reactive<sup>b</sup> NR, nonreactive

TABLE 2. Mean Dilutions (Quantitative end point titers) for nine treponemal syphilis assays of TP-PA reactive sera

Test	Median Titer	Estimates of Mean Dilutions and 95% Confidence limits				Mean Dilution difference from FTA-ABS	
		Mean Dilution	Std Error	Lower 95% CL	Upper 95% CL	Mean Difference	Std Error
FTA-ABS*	4	1.9266	0.1655	1.6022	2.2510	n/a	n/a
Captia IgG*	8	3.1560	0.1552	2.8517	3.4602	1.229	0.082
TP-PA**	16	4.0734	0.1624	3.7551	4.3917	2.147	0.102
INNO-LIA*	16	4.0917	0.1927	3.7141	4.4694	2.165	0.123
Bioline**	64	5.6606	0.1886	5.2910	6.0301	3.734	0.135
Trep-ID**	64	6.2752	0.1906	5.9016	6.6489	4.349	0.123
LIAISON**	128	6.3119	0.1603	5.9977	6.6262	4.385	0.105
Bioelisa**	128	7.1560	0.1789	6.8053	7.5066	5.229	0.112
Trep-Sure**	512	8.8899	0.1747	8.5475	9.2323	6.963	0.1182

\* IgG only

\*\* IgG/IgM detection



## A Comparison of the Analytical level of agreement of Nine Treponemal Assays for Syphilis and Possible Implications for Screening Algorithms

Journal:	<i>BMJ Open</i>
Manuscript ID:	bmjopen-2013-003347.R1
Article Type:	Research
Date Submitted by the Author:	28-Jul-2013
Complete List of Authors:	Castro, Arnold; Centers for Disease Control and Prevention, laboratory Reference and Research Branch Jost, Heather; Centers for Disease Control, LRRB Cox, David; Centers for Disease Control and Prevention, Fakile, Yetunde; Centers for Disease Control and Prevention, Kikkert, Susan; Centers for Disease Control and Prevention, Tun, Ye; Centers for Disease Control and Prevention, Zaidi, Akbar; Centers for Disease Control and Prevention, Park, Mahin; Georgia Department of Community Health,
<b>Primary Subject Heading</b>:	Infectious diseases
Secondary Subject Heading:	Diagnostics, Infectious diseases, Immunology (including allergy)
Keywords:	Microbiology < BASIC SCIENCES, Immunology < BASIC SCIENCES, MICROBIOLOGY

SCHOLARONE™  
Manuscripts

Only

1  
2  
3  
4  
5  
6  
7  
8  
9  
10  
11  
12  
13  
14  
15  
16  
17  
18  
19  
20  
21  
22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

**Abstract**

**Objective:** The serological diagnosis of syphilis requires the detection of two distinct antibodies, the nontreponemal and treponemal. CDC recommends screening first with a nontreponemal test such as (RPR/VDRL), and then confirming those results with one of several treponemal tests (FTA-ABS, EIA, CIA, TP-PA or POC). Due to the high volume of samples processed by some laboratories using automated systems, the screening with treponemal assays and confirming with nontreponemal tests is becoming the established norm. The purpose of this study was to evaluate eight treponemal assays using TP-PA as the predicate assay.

**Methods:** Two hundred ninety stored serum samples were tested qualitatively according to the manufacturer's directions.

**Results:** Concordance with specimens tested as reactive or nonreactive using TP-PA was: FTA-ABS 94.5% and 100%, Trep-Sure 100% and 98.9%, Bioelisa 100% and 98.9%, INNO-LIA 99.1% and 99.4%, Bioline 100% and 98.9%, Captia IgG 100% and 97.2%, Trcp-ID 100% and 100%, LIAISON 100% and 99.4% respectively. In order to properly evaluate the performance of these assays, the analytical sensitivity was determined by endpoint titration of serial dilutions of the reactive serum samples in normal sera. The median endpoint titer varied from 1:4 for FTA-ABS to 1:512 for Trep-Sure.

**Conclusions:** The performance of the treponemal serological assays was comparable when using medium and high titer sera. However, the varying performance on specimen dilutions suggest that there may be differences in sensitivity with low titer sera that are more prevalent in primary and late syphilis cases.

**Introduction**

Syphilis is a sexually transmitted disease caused by the bacterium *Treponema pallidum*. The clinical diagnosis is difficult due to the complexity of manifestations which requires the aid of serological interpretation. The serological diagnosis of syphilis depends upon the detection of two distinct antibodies, the nontreponemal or heterophile antibodies (reagin) directed against cardiolipin released from damaged host cells and from the treponemes themselves. These antibodies can also be present in other diseases and human conditions such as (Lupus, malaria, HIV, IV drug users, etc.) (1). The presence of nontreponemal antibodies is indicative of active infection, thus a reduction in titer can suggest a successful antibiotic therapy, and a significant increase can indicate a possible relapse or reinfection (2). The treponemal antibodies are primarily directed against specific lipoprotein antigens

of the bacterium such as 15,17, and 47-kDa. Even after treatment or time, these treponemal antibodies usually remain present for life. A positive treponemal test cannot distinguish between active, treated, and old cases of syphilis.

The traditional algorithm is to screen with a nontreponemal assay such as VDRL or RPR and those serum samples found reactive are then confirmed using a treponemal test. With the introduction of automated systems which are appropriate for high output volumes, the reverse algorithm is gaining acceptability and its usefulness is arguable justified (6). However, the adoption of this reverse algorithm has increased the number of discrepant results between the screening and confirmatory tests (4, 5). One reason the traditional algorithm was established was to avoid the detection of previously treated cases, especially in *low* prevalence settings where the majority of sera will test negative and not require further evaluation, unless a patient has been reinfected or in rare circumstances, their non-treponemal antibody titers remain serofast. With the application of the reverse algorithm even in low prevalence settings all previously treated cases will likely be detected. When the screening treponemal test is positive and the confirmatory nontreponemal test is negative, the CDC guidelines recommend using a second treponemal test to validate the results of the screening test. The selection of a second confirmatory test introduces the possibility of having an analytically less sensitive test which would suggest a false positive screening or a less sensitive second confirmatory test. Then it is important to determine if all available treponemal assays are comparable in analytical sensitivity to avoid uncertainty in their selection (3).

## Methods

Two hundred ninety stored serum samples randomly selected from our serum bank were used for this study. Originally the serum samples were obtained from the Georgia Public Health Laboratory with all identifiers removed and under a CDC IRB approval. Due to the prevalence of high lipid contents of the sera, the samples were treated with Cleanascite HC (PureRiotech, Middlesex, NJ) (7). The patterns of reactivity of the test panel were determined at CDC by a quantitative RPR (Becton Dickinson, Baltimore, MD) and a qualitative TP-PA (Fujirebio Diagnostics, Inc., Malvern, PA). All serum samples were divided into aliquots and frozen at -20°C so that the number of freeze-thaw cycles was consistent for each sera tested. Eight of the treponemal assays: FTA-ABS; (Zeus

Scientific, Raritan, NJ), LIAISON Treponema Assay (DiaSorin, Inc., Stillwater, MN), SD Bioline Syphilis 3.0 Rapid Point of Care (POC) test (Standard Diagnostics, Inc., Korea), INNO-LIA immunoblot (Innogenetics, Gent, Belgium), Bioelisa (BioKit, Barcelona, Spain) and the Captia IgG, Trep-ID. and Trep-Surc (Trinity Biotech, Jamestown, NY) were tested qualitatively, then quantitatively with reactive samples. In the quantitative test, two-fold serial dilutions (from 1:2 to 1:16,384) were prepared using normal nonreactive human plasma converted to serum (8) as a diluent. All dilutions were treated as neat sera and tested with the nine different assays according to the manufacturer's package insert. The test was performed blind, the EIA results were collated and compared to the TP-PA used as the predicate for the qualitative test.

**Statistical Methods.**

The endpoint titer data was converted to dilution data from the exponential scale (titer) to a linear scale (dilution) by taking the logarithm of the endpoint titer, i.e.,  $\log_2(\text{titer}) = \text{dilution}$ . This transformation reduces the variance within the test and between the tests. A generalized linear model was fitted to these data by the Generalized Estimating Equations (GEE) method to test for the differences between the mean dilutions of these tests. Adjustments for multiple comparisons between the mean dilutions of these tests were made by Tukey-Kramer method. A p-value of  $<0.05$  was used as a significance level for the differences. All statistical analyses were conducted using SAS, version 9.3 (SAS Institute, Cary, North Carolina).

**Results**

All treponemal syphilis assays were performed according to the manufacturer's direction circulars using two hundred ninety stored serum samples of which 109 were reactive and 181 nonreactive for TP-PA. Concordance with specimens tested as reactive or nonreactive using TP-PA was: FTA-ABS 94.5% and 100%, INNO-LIA 99.1% and 99.4%, LIAISON 100% and 99.4%. Trep-Sure 100% and 98.9%, Bioelisa 100% and 98.9%, Bioline 100% and 98.9%, Captia IgG 100% and 97.2%, Trep-ID 100% and 100% respectively.

These results indicated that most of the assays were comparable to TP-PA when the serum samples were tested undiluted. However, FTA-ABS failed to detect six TP-PA reactive samples and INNO-



88 LIA one, while the other four assays missed none. Of the 109 TP-PA reactive sera, 95(85%) were  
 1 nonreactive for the confirmatory RPR test. Fourteen samples were reactive for both TP-PA and RPR. All  
 2 other assays were able to detect treponemal antibody in those undiluted sera that were reactive by the TP-  
 3 PA. However, five sera found to be reactive by the Captia IgG test and one by the BIOLINE were found to  
 4 be non-reactive by all other tests suggesting false positives. All other sera found to be non-reactive by the  
 5 TP-PA but reactive by another test could be considered as true positives since the sera were found to be  
 6 reactive by at least two other tests (see Table 1). If the TP-PA would have been selected as the screening test  
 7 for the reverse algorithm then eight samples reactive by the other assays would have been missed  
 8 resulting in a false serological result.  
 9  
 10  
 11  
 12

One hundred and nine stored serum samples reactive by TP-PA were serially diluted in nonreactive  
 human plasma converted to serum (8), with the purpose of demonstrating the endpoint titer of each sample  
 with the nine different assays. With a simple comparison the endpoint titers indicated that Trep-Sure had a  
 median value of 1:512, Bioelisa 1:128, LIAISON 1:64, Trep-ID 1:64, Bioline 1:32, INNO-LIA 1:16, TP-PA  
 1:16, Captia IgG 1:8 and FTA-ABS 1:4 (Table 2).

103 Mean dilutions and their standard errors of eight quantitative treponemal assays and the TP-PA  
 104 estimated by the generalized linear model are shown in Table 2. There were significant ( $p < 0.0001$ )  
 105 differences between the mean dilutions of all tests, except between TP-PA and INNO-LIA ( $p = 1.000$ )  
 106 and Trep-ID and LIASON ( $p = 1.000$ ). This analysis was repeated after subtracting the FTA-ABS  
 107 titer from the titer of the corresponding test. In this analysis also there was no significant difference  
 108 between TP-PA and INNO-LIA ( $p = 1.000$ ) and Trep-ID and LIASON ( $p = 1.000$ ), and the differences  
 109 among the rest of the pairs were significant ( $p < 0.0001$ ). (Table 2)

# 111 Discussion:

112 The reverse sequence syphilis screening algorithm is being used because it takes advantage of  
 113 automated systems. However the selection of these assays to screen serum samples with suspected  
 114 cases of syphilis presents problems that were not clearly anticipated (6). In cases where the  
 115 treponemal test is positive and the confirmatory RPR test is negative, the CDC recommends that a  
 116 second treponemal test be used to determine the serological validity of the screening test (5). If the  
 117 second treponemal test is negative, then it is inconclusive whether the first screening test is a false

positive or a more sensitive test. This situation presents a dilemma for physicians and understandably they will be asking for clarification as to whether it is or not a case of syphilis in asymptomatic situations. Consequently, it would be advisable for the laboratorian to consider the selection of more sensitive screening and confirmatory treponemal assays to avoid having discrepant results. This is not an easy task because of the different assay platforms presently used for screening in the reverse algorithm system ranging from FTA-ABS, EIA, CIA, TP-PA and POC. In this study an attempt was made to determine the analytical level of agreement of the assays with the purpose of better understanding the abilities of the assays for detection of low level antibodies. When the nine assays were tested quantitatively, the FTA-ABS was the treponemal test with the lowest median endpoint titer (1:4). In contrast, Trep-Sure had a median endpoint titer of (1:512) (Table 2). Since the FTA-ABS was the least reactive, its reevaluation should be considered.

Further studies with documented clinical samples are necessary to demonstrate the correlation between the antibody titer and the stage of the disease. With these comparison results we then could infer that there are marked differences between the analytical sensitivities of these assays especially with low titer sera such as those found in primary or latent cases. The variability of results with the different assays may be due in part to the composition of their systems. Those assays using the double antigen or sandwich configuration such as Trep-Sure, Bioelisa, Trep-ID and LIAISON where the recombinant antigen such as Tp15, 17, 47 KD or combination of them are both attached to the ELISA plate or latex particle and are directly conjugated to horse radish peroxidase (HRP) or chemiluminescence (CIA). This configuration tends to be more sensitive because they are able to detect both IgG and IgM antibodies in the sera of patients with syphilis. In addition, the detection of IgM antibodies amplifies the signal due to its multiple binding sites. On the other hand, normally those assays using the indirect or anti-human system such as Captia (IgG) and FTA-ABS are limited in that they only detect one single immunoglobulin.

When selecting the reverse algorithm system for syphilis diagnostic serology and in cases when the treponemal test is reactive and the RPR nonreactive, it would be advisable to consider the selection of more analytically sensitive screening and confirmatory treponemal assays preferably those based on the double antigen or sandwich method, such as Trep-Sure, Bioelisa, Trep-ID and LIASON to avoid the necessity of using a less sensitive second confirmatory treponemal assay.

## REFERENCES

1. Larsen, S. A. (1989). "Syphilis." Clin.Lab.Med. 9: 545-557.
2. Larsen, S. A., V.Pope. et al. (1998). A Manual of Tests for Syphilis. Washington, DC, American Public Health Association.
3. Zhang WM, Yen-Lieberman, Means K,et al. The Impact of Analytical Sensitivity on Screening Algorithms for Syphilis. Clinical Chemistry. 2012; 58:6.
4. Binnicker MJ, Jespersen DJ, Rollins LO. Direct comparison of the traditional and reverse syphilis screening algorithms in a population with a low prevalence of syphilis. J Clin Microbiol 2012;50:148-50.
5. Centers for Disease, C. and Prevention (2011). "Discordant results from reverse sequence syphilis screening--five laboratories, United States, 2006-2010." MMWR Morb Mortal Wkly Rep 60(5): 133-137.
6. Binnicker MJ. Which algorithm should be used to screen for syphilis?. Curr Opin Infect Dis 2012; 25:79-85.
7. Castro AR, Morrill WE, Pope V. Lipid Removal from Human Serum Samples. Clin Diag Lab Immunol. 2000;7:197-199.
8. Castro AR, Kikkert SE, Fears MB, Pope V. Defibrination of Blood Plasma for use in the Serological Tests for Syphilis. Clin Diag Lab Immunol. 2002;9:1376-1378.

TABLE 1. Percent concordance of eight qualitative syphilis treponemal test compared to the predicate TP-PA

TP-PA Predicate	Total no of samples	FTA-ABS		INNOLIA		Liaison		Trep-Sure		Bioelisa		Bioline		Captia (IgG)		Trep-ID	
		R <sup>a</sup>	NR <sup>b</sup>	R	NR	R	NR	R	NR	R	NR	R	NR	R	NR	R	NR
Reactive	109	103	6(94.4%)	108	1(99.1%)	109	0(100%)	109	0(100%)	109	0(100%)	109	0(100%)	109	0(100%)	109	0(100%)
Nonreactive	181	0	181(100%)	1	180(99.4%)	1	180(99.4%)	2	179(98.9%)	2	179(98.9%)	2	179(98.9%)	5	176(97.2%)	0	181(100%)

<sup>a</sup> R, reactive

<sup>b</sup> NR, nonreactive

TABLE 2. Mean Dilutions (Quantitative end point titers) for nine treponemal syphilis assays of TP-PA reactive sera

Test	Median Titer	Estimates of Mean Dilutions and 95% Confidence limits				Mean Dilution difference from FTA-ABS	
		Mean Dilution	Std Error	Lower 95% CL	Upper 95% CL	Mean Difference	Std Error
FTA-ABS*	4	1.9266	0.1655	1.6022	2.2510	n/a	n/a
Captia IgG*	8	3.1560	0.1552	2.8517	3.4602	1.229	0.082
TP-PA**	16	4.0734	0.1624	3.7551	4.3917	2.147	0.102
INNO-LIA*	16	4.0917	0.1927	3.7141	4.4694	2.165	0.123
Bioline**	64	5.6606	0.1886	5.2910	6.0301	3.734	0.135
Trep-ID**	64	6.2752	0.1906	5.9016	6.6489	4.349	0.123
LIAISON**	128	6.3119	0.1603	5.9977	6.6262	4.385	0.105
Bioelisa**	128	7.1560	0.1789	6.8053	7.5066	5.229	0.112
Trep-Sure**	512	8.8899	0.1747	8.5475	9.2323	6.963	0.1182

\* IgG only

\*\* IgG/IgM detection

## Abstract

**Objective:** The serological diagnosis of syphilis requires the detection of two distinct antibodies, the nontreponemal and treponemal. CDC recommends screening first with a nontreponemal test such as (RPR/VDRL), and then confirming those results with one of several treponemal tests (FTA-ABS, EIA, CIA, TP-PA or POC). Due to the high volume of samples processed by some laboratories using automated systems, the screening with treponemal assays and confirming with nontreponemal tests is becoming the established norm. The purpose of this study was to evaluate eight treponemal assays using TP-PA as the predicate assay.

**Methods:** Two hundred ninety stored serum samples were tested qualitatively according to the manufacturer's directions.

**Results:** Concordance with specimens tested as reactive or nonreactive using TP-PA was: FTA-ABS 94.5% and 100%, Trep-Sure 100% and 98.9%, Bioelisa 100% and 98.9%, INNO-LIA 99.1% and 99.4%, Bioline 100% and 98.9%, Captia IgG 100% and 97.2%, Trcp-ID 100% and 100%, LIAISON 100% and 99.4% respectively. In order to properly evaluate the performance of these assays, the analytical sensitivity was determined by endpoint titration of serial dilutions of the reactive serum samples in normal sera. The median endpoint titer varied from 1:4 for FTA-ABS to 1:512 for Trep-Sure.

**Conclusions:** The performance of the treponemal serological assays was comparable when using medium and high titer sera. However, the varying performance on specimen dilutions suggest that there may be differences in sensitivity with low titer sera that are more prevalent in primary and late syphilis cases.

## Introduction

Syphilis is a sexually transmitted disease caused by the bacterium *Treponema pallidum*. The clinical diagnosis is difficult due to the complexity of manifestations which requires the aid of serological interpretation. The serological diagnosis of syphilis depends upon the detection of two distinct antibodies, the nontreponemal or heterophile antibodies (reagin) directed against cardiolipin released from damaged host cells and from the treponemes themselves. These antibodies can also be present in other diseases and human conditions such as (Lupus, malaria, HIV, IV drug users, etc.) (1). The presence of nontreponemal antibodies is indicative of active infection, thus a reduction in titer can suggest a successful antibiotic therapy, and a significant increase can indicate a possible relapse or reinfection (2). The treponemal antibodies are primarily directed against specific lipoprotein antigens

of the bacterium such as 15,17, and 47-kDa. Even after treatment or time, these treponemal antibodies usually remain present for life. A positive treponemal test cannot distinguish between active, treated, and old cases of syphilis.

The traditional algorithm is to screen with a nontreponemal assay such as VDRL or RPR and those serum samples found reactive are then confirmed using a treponemal test. With the introduction of automated systems which are appropriate for high output volumes, the reverse algorithm is gaining acceptability and its usefulness is arguable justified (6). However, the adoption of this reverse algorithm has increased the number of discrepant results between the screening and confirmatory tests (4, 5). One reason the traditional algorithm was established was to avoid the detection of previously treated cases, especially in *low* prevalence settings where the majority of sera will test negative and not require further evaluation, unless a patient has been reinfected or in rare circumstances, their non-treponemal antibody titers remain serofast. With the application of the reverse algorithm even in low prevalence settings all previously treated cases will likely be detected. When the screening treponemal test is positive and the confirmatory nontreponemal test is negative, the CDC guidelines recommend using a second treponemal test to validate the results of the screening test. The selection of a second confirmatory test introduces the possibility of having an analytically less sensitive test which would suggest a false positive screening or a less sensitive second confirmatory test. Then it is important to determine if all available treponemal assays are comparable in analytical sensitivity to avoid uncertainty in their selection (3).

## Methods

Two hundred ninety stored serum samples randomly selected from our serum bank were used for this study. Originally the serum samples were obtained from the Georgia Public Health Laboratory with all identifiers removed and under a CDC IRB approval. Due to the prevalence of high lipid contents of the sera, the samples were treated with Cleanascite HC (PureRiotech, Middlesex, NJ) (7). The patterns of reactivity of the test panel were determined at CDC by a quantitative RPR (Becton Dickinson, Baltimore, MD) and a qualitative TP-PA (Fujirebio Diagnostics, Inc., Malvern, PA). All serum samples were divided into aliquots and frozen at -20°C so that the number of freeze-thaw cycles was consistent for each sera tested. Eight of the treponemal assays: FTA-ABS; (Zeus



Scientific, Raritan, NJ), LIAISON Treponema Assay (DiaSorin, Inc., Stillwater, MN), SD Bioline Syphilis 3.0 Rapid Point of Care (POC) test (Standard Diagnostics, Inc., Korea), INNO-LIA immunoblot (Innogenetics, Gent, Belgium), Bioelisa (BioKit, Barcelona, Spain) and the Captia IgG, Trep-ID. and Trep-Sure (Trinity Biotech, Jamestown, NY) were tested qualitatively, then quantitatively with reactive samples. In the quantitative test, two-fold serial dilutions (from 1:2 to 1:16,384) were prepared using normal nonreactive human plasma converted to serum (8) as a diluent. All dilutions were treated as neat sera and tested with the nine different assays according to the manufacturer's package insert. The test was performed blind, the EIA results were collated and compared to the TP-PA used as the predicate for the qualitative test.

### Statistical Methods.

The endpoint titer data was converted to dilution data from the exponential scale (titer) to a linear scale (dilution) by taking the logarithm of the endpoint titer, i.e.,  $\log_2(\text{titer}) = \text{dilution}$ . This transformation reduces the variance within the test and between the tests. A generalized linear model was fitted to these data by the Generalized Estimating Equations (GEE) method to test for the differences between the mean dilutions of these tests. Adjustments for multiple comparisons between the mean dilutions of these tests were made by Tukey-Kramer method. A p-value of  $<0.05$  was used as a significance level for the differences. All statistical analyses were conducted using SAS, version 9.3 (SAS Institute, Cary, North Carolina).

### Results

All treponemal syphilis assays were performed according to the manufacturer's direction circulars using two hundred ninety stored serum samples of which 109 were reactive and 181 nonreactive for TP-PA. Concordance with specimens tested as reactive or nonreactive using TP-PA was: FTA-ABS 94.5% and 100%, INNO-LIA 99.1% and 99.4%, LIAISON 100% and 99.4%. Trep-Sure 100% and 98.9%, Bioelisa 100% and 98.9%, Bioline 100% and 98.9%, Captia IgG 100% and 97.2%, Trep-ID 100% and 100% respectively.

These results indicated that most of the assays were comparable to TP-PA when the serum samples were tested undiluted. However, FTA-ABS failed to detect six TP-PA reactive samples and INNO-



LIA one, while the other four assays missed none. Of the 109 TP-PA reactive sera, 95(85%) were nonreactive for the confirmatory RPR test. Fourteen samples were reactive for both TP-PA and RPR. All other assays were able to detect treponemal antibody in those undiluted sera that were reactive by the TP-PA. However, five sera found to be reactive by the Captia IgG test and one by the BIOLINE were found to be non-reactive by all other tests suggesting false positives. All other sera found to be non-reactive by the TP-PA but reactive by another test could be considered as true positives since the sera were found to be reactive by at least two other tests (see Table 1). If the TP-PA would have been selected as the screening test for the reverse algorithm then eight samples reactive by the other assays would have been missed resulting in a false serological result.

One hundred and nine stored serum samples reactive by TP-PA were serially diluted in nonreactive human plasma converted to serum (8), with the purpose of demonstrating the endpoint titer of each sample with the nine different assays. With a simple comparison the endpoint titers indicated that Trep-Sure had a median value of 1:512, Bioelisa 1:128, LIAISON 1:64, Trep-ID 1:64, Bioline 1:32, INNO-LIA 1:16, TP-PA 1:16, Captia IgG 1:8 and FTA-ABS 1:4 (Table 2).

Mean dilutions and their standard errors of eight quantitative treponemal assays and the TP-PA estimated by the generalized linear model are shown in Table 2. There were significant ( $p < 0.0001$ ) differences between the mean dilutions of all tests, except between TP-PA and INNO-LIA ( $p = 1.000$ ) and Trep-ID and LIASON ( $p = 1.000$ ). This analysis was repeated after subtracting the FTA-ABS titer from the titer of the corresponding test. In this analysis also there was no significant difference between TP-PA and INNO-LIA ( $p = 1.000$ ) and Trep-ID and LIASON ( $p = 1.000$ ), and the differences among the rest of the pairs were significant ( $p < 0.0001$ ). (Table 2)

## Discussion:

The reverse sequence syphilis screening algorithm is being used because it takes advantage of automated systems. However the selection of these assays to screen serum samples with suspected cases of syphilis presents problems that were not clearly anticipated (6). In cases where the treponemal test is positive and the confirmatory RPR test is negative, the CDC recommends that a second treponemal test be used to determine the serological validity of the screening test (5). If the second treponemal test is negative, then it is inconclusive whether the first screening test is a false

positive or a more sensitive test. This situation presents a dilemma for physicians and understandably they will be asking for clarification as to whether it is or not a case of syphilis in asymptomatic situations. Consequently, it would be advisable for the laboratorian to consider the selection of more sensitive screening and confirmatory treponemal assays to avoid having discrepant results. This is not an easy task because of the different assay platforms presently used for screening in the reverse algorithm system ranging from FTA-ABS, EIA, CIA, TP-PA and POC. In this study an attempt was made to determine the analytical level of agreement of the assays with the purpose of better understanding the abilities of the assays for detection of low level antibodies. When the nine assays were tested quantitatively, the FTA-ABS was the treponemal test with the lowest median endpoint titer (1:4). In contrast, Trep-Sure had a median endpoint titer of (1:512) (Table 2). Since the FTA-ABS was the least reactive, its reevaluation should be considered.

Further studies with documented clinical samples are necessary to demonstrate the correlation between the antibody titer and the stage of the disease. With these comparison results we then could infer that there are marked differences between the analytical sensitivities of these assays especially with low titer sera such as those found in primary or latent cases. The variability of results with the different assays may be due in part to the composition of their systems. Those assays using the double antigen or sandwich configuration such as Trep-Sure, Bioelisa, Trep-ID and LIAISON where the recombinant antigen such as Tp15, 17, 47 KD or combination of them are both attached to the ELISA plate or latex particle and are directly conjugated to horse radish peroxidase (HRP) or chemiluminescence (CIA). This configuration tends to be more sensitive because they are able to detect both IgG and IgM antibodies in the sera of patients with syphilis. In addition, the detection of IgM antibodies amplifies the signal due to its multiple binding sites. On the other hand, normally those assays using the indirect or anti-human system such as Captia (IgG) and FTA-ABS are limited in that they only detect one single immunoglobulin.

When selecting the reverse algorithm system for syphilis diagnostic serology and in cases when the treponemal test is reactive and the RPR nonreactive, it would be advisable to consider the selection of more analytically sensitive screening and confirmatory treponemal assays preferably those based on the double antigen or sandwich method, such as Trep-Sure, Bioelisa, Trep-ID and LIASON to avoid the necessity of using a less sensitive second confirmatory treponemal assay.

147

1  
2  
3 148  
4 149  
5 150  
6 151  
7  
8 152  
9 153  
10  
11 154  
12  
13 155  
14 156  
15 157  
16 158  
17 159  
18  
19 160  
20 161  
21 162  
22  
23 163  
24 164  
25 165  
26 166  
27 167  
28 168  
29 169  
30  
31 170  
32  
33  
34 171  
35  
36  
37  
38  
39  
40  
41  
42 172  
43  
44  
45  
46  
47  
48  
49 173  
50  
51  
52  
53 174  
54  
55  
56  
57 175  
58  
59  
60

REFERENCES

1. Larsen, S. A. (1989). "Syphilis." Clin.Lab.Med. 9: 545-557.

2. Larsen, S. A., V.Pope. et al. (1998). A Manual of Tests for Syphilis. Washington, DC, American Public Health Association.

3. Zhang WM, Yen-Lieberman, Means K,et al. The Impact of Analytical Sensitivity on Screening Algorithms for Syphilis. Clinical Chemistry. 2012; 58:6.

4. Binnicker MJ, Jespersen DJ, Rollins LO. Direct comparison of the traditional and reverse syphilis screening algorithms in a population with a low prevalence of syphilis. J Clin Microbiol 2012;50:148-50.

5. Centers for Disease, C. and Prevention (2011). "Discordant results from reverse sequence syphilis screening--five laboratories, United States, 2006-2010." MMWR Morb Mortal Wkly Rep 60(5): 133-137.

6. Binnicker MJ. Which algorithm should be used to screen for syphilis?. Curr Opin Infect Dis 2012; 25:79-85.

7. Castro AR, Morrill WE, Pope V. Lipid Removal from Human Serum Samples. Clin Diag Lab Immunol. 2000;7:197-199.

8. Castro AR, Kikkert SE, Fears MB, Pope V. Defibrination of Blood Plasma for use in the Serological Tests for Syphilis. Clin Diag Lab Immunol. 2002;9:1376-1378.

TABLE 1. Percent concordance of eight qualitative syphilis treponemal test compared to the predicate TP-PA

TP-PA Predicate	Total no of samples	FTA-ABS		INNOLIA		Liaison		Trep-Sure		Bioelisa		Bioline		Captia (IgG)		Trep-ID	
		R <sup>a</sup>	NR <sup>b</sup>	R	NR	R	NR	R	NR	R	NR	R	NR	R	NR	R	NR
Reactive	109	103	6(94.4%)	108	1(99.1%)	109	0(100%)	109	0(100%)	109	0(100%)	109	0(100%)	109	0(100%)	109	0(100%)
Nonreactive	181	0	181(100%)	1	180(99.4%)	1	180(99.4%)	2	179(98.9%)	2	179(98.9%)	2	179(98.9%)	5	176(97.2%)	0	181(100%)

<sup>a</sup> R, reactive<sup>b</sup> NR, nonreactive

TABLE 2. Mean Dilutions (Quantitative end point titers) for nine treponemal syphilis assays of TP-PA reactive sera

Test	Median Titer	Estimates of Mean Dilutions and 95% Confidence limits				Mean Dilution difference from FTA-ABS	
		Mean Dilution	Std Error	Lower 95% CL	Upper 95% CL	Mean Difference	Std Error
FTA-ABS*	4	1.9266	0.1655	1.6022	2.2510	n/a	n/a
Captia IgG*	8	3.1560	0.1552	2.8517	3.4602	1.229	0.082
TP-PA**	16	4.0734	0.1624	3.7551	4.3917	2.147	0.102
INNO-LIA*	16	4.0917	0.1927	3.7141	4.4694	2.165	0.123
Bioline**	64	5.6606	0.1886	5.2910	6.0301	3.734	0.135
Trep-ID**	64	6.2752	0.1906	5.9016	6.6489	4.349	0.123
LIAISON**	128	6.3119	0.1603	5.9977	6.6262	4.385	0.105
Bioelisa**	128	7.1560	0.1789	6.8053	7.5066	5.229	0.112
Trep-Sure**	512	8.8899	0.1747	8.5475	9.2323	6.963	0.1182

\* IgG only

\*\* IgG/IgM detection