

## **Extended Experimental Procedures**

### **Cells, viruses and reagents**

Human colon cancer cell line HT-29, Vero, 293T, and 3T3-SA were maintained at 37° C in 5% CO<sub>2</sub> using Dulbecco's modified Eagle's medium (DMEM) containing 10% FBS (Atlanta Biologicals), 4.5 g/mL glucose, 2 mM L-glutamine, 100 U/mL penicillin and 100 U/mL streptomycin (Invitrogen). U937 cells were maintained in complete RPMI medium containing 10% FBS, 2mM L-glutamine, 100 U/mL penicillin and 100 U/mL streptomycin. The R1-null mutant HSV1  $\Delta$ ICP6 and its parental strain KOS were kindly provided by Dr. Sandra Weller (University of Connecticut Health Center, Farmington, CT). HSV viral stocks were prepared by infecting HT-29 cells stably expressing ICP10 or ICP6 at an MOI of 0.01 and were titrated by plaque assay on monolayer cultures of Vero cells. MCMV K181 and MCMV-M45*mut*RHIM viruses were generated as described (Upton et al., 2010). z-VAD-fmk was from Enzo Life Sciences, RIP1 inhibitor GSK'963, RIP3 inhibitor GSK'840 and IAP antagonist SMAC007 (Mandal et al., 2014), as well as the pRIP1 S166-specific antibody, were provided by GlaxoSmithKline (contact Peter Gough, peter.j.gough@gsk.com), IAP antagonist BV6 was provided by Domagoj Vucic (Genentech), and recombinant human TNF was from R&D.

### **Plasmids, transfection and transduction**

The ICP6 ORF was PCR-amplified from pKHF plasmid (a gift from Dr. Sandra Weller) and then inserted into pCMV10-3XFLAG vector (Sigma) in-frame with an amino-terminal epitope tag. Truncation mutants were generated from pCMV10-3XFLAG-ICP6 as a template using PCR to amplify nucleotides 1-1467, 1-1887, 1-2505, 1-3318, 1-3348, 730-1887 or 730-3414, to produce ICP6(1-489), ICP6(1-629), ICP6(1-835), ICP6(1-1106), ICP6(1-1116), ICP6(244-629) and ICP6( $\Delta$ 1-243), respectively. ICP6*mut*RHIM was generated by overlap extension PCR to change ICP6 aa73-76 from VQCG to AAAA. ICP6(G865/867/870A) was generated by overlap extension PCR to change ICP6 aa865, aa867 and aa870 from G to A. All plasmids were verified by DNA

sequencing. Expression vectors encoding Myc-tagged or FLAG-tagged RIP1, RIP2, RIP3, RIP4 have been previously described (Kaiser and Offermann, 2005; Kaiser et al., 2008; Meylan et al., 2004; Rebsamen et al., 2009). The RIP1 intermediate domain (aa301-558) was PCR amplified and cloned into the expression vectors pcDNA3-6XMyC and pCMV10-3XFLAG. To mutate the RHIM in RIP3 or in RIP1, the amino acids VQVG or IQIG, which correspond to RIP3aa458-461 or RIP1aa539-542, were changed to AAAA by overlap extension PCR. ICP10 ORF was PCR-amplified from plasmid AdCMV5-GFP-R1 (a gift from Dr. Yves Langelier) and then inserted into pCMV10-3XFLAG vector (Sigma) resulting in an in-frame amino-terminal 3XFLAG epitope tag. This FLAG-tagged ICP10 construct was used as a template using PCR to amplify nucleotides 748-3435 (aa250-1144) and produce ICP10( $\Delta$ 1-249). Epitope-tagged ICP6 and ICP10 and mutant forms were subcloned into pLV-EF1 $\alpha$ -MCS-IRES-Puro (Biosettia) and HT-29 cells stably expressing these proteins were generated by lentiviral transduction (Upton et al., 2010). Transfections were performed with Lipofectamine 2000 (Invitrogen) and DNA at a 2:1 ratio in Opti-MEM (Invitrogen).

## Supplemental Figures

**Figure S1, related to Figure 1.** (A) Viability of HT-29 cells infected with KOS or  $\Delta$ ICP6 (MOI=10) and treated from 2 to 14 hpi with T, T+S, T+V, T+S+V or anti-Fas antibody (F; 5 ng/ml), F+S, F+V, F+S+V. (B) Viability of HT-29 cells infected with  $\Delta$ ICP6 (MOI=10) and treated from 4 to 19 hpi with T, T+V, T+S7, T+S7+V, in the absence or presence GSK'963 (1  $\mu$ M) or GSK'840 (3  $\mu$ M). (C) Viability of HT-29 cells infected with mock, KOS, or  $\Delta$ ICP6 and treated with T+S+V from 2 to 24 hpi in the presence or absence of phosphonoformate (PFA; 300  $\mu$ g/ml). (D-F) Replication levels of KOS (left) and  $\Delta$ ICP6 (right) virus in HT-29 cells (D, MOI=1; E, MOI=0.5; F, MOI=0.1) in the absence (DMSO control) or presence of T+S7, T+S7+V or T+S7+V with RIP1 kinase inhibitor GSK'963(1  $\mu$ M) or RIP3 kinase inhibitor GSK'840 (3  $\mu$ M), as shown. (G) IB to detect cFLIP, cIAP1, cIAP2 and  $\beta$ -actin in KOS or  $\Delta$ ICP6-infected HT-29 cell lysates collected at 2, 4, 6, 8, 10 hpi (MOI=5).


**Figure S2, related to Figure 2.** (A) Photomicrographs of HT-29-EV, HT-29-ICP6 and HT-29-ICP10 cell cultures at 0, 4, 8, 12 h post-treatment with T+S+V in the presence of Sytox Green (62.5 nM), a nucleic acid detecting fluorescent dye, assessed on an IncuCyte instrument. Original magnification, 200x. (B) Viability of HT-29-EV or HT-29-M45 cells treated with T+S+V for 24h. (C) Viability of HT-29-EV, HT-29-ICP6 or HT-29-ICP10 cells treated with F+S+V for 24 h. (D) Viability of HT-29-EV or HT-29-ICP10 cells treated with T+S+V or F+S+V as indicated for 12, 24 or 48 h. (E) Viability of U937-EV or U937-ICP10 cells treated with T, T+S, T+V or T+S+V as indicated for 18 h. (F) Viability of 3T3-SA-EV or 3T3-SA-ICP10 cell lines treated for 18 h with T (25 ng/mL) and/or V (25  $\mu$ M). (G) Viability of 3T3-SA-EV or 3T3-SA-ICP6 cells treated for 18 h with T and/or V. (H) Viability of 3T3-SA-EV or 3T3-SA-M45 cells treated for 18 h with T+V (left) or infected with parental K181 and M45mutRHIM virus (MOI=10) (right). Insets depict IB detection of FLAG-ICP10 or FLAG-ICP6 in the transduced whole cell lysates. (I and J) HSV2 ICP10 interaction with RIP1 and RIP3. IB/IP to detect FLAG-ICP10 interaction with Myc-tagged versions of RIP1, RIP1(301-558) intermediate domain, RIP2, RIP3, and RIP4. The top panel shows 293T cell lysates subjected to anti-Myc IP followed by anti-FLAG and anti-Myc IB analysis. (K) Reciprocal IP to detect ICP10 interaction with RIP1 and RIP3. Lysates of cells transfected with Myc-RIP1 or Myc-RIP3, with or without FLAG-ICP10 subjected to IP with anti-FLAG antibody, followed by IB with anti-Myc or anti-FLAG antibody. (L) IP/IB to detect Myc-ICP10 interaction with FLAG-tagged versions of RIP1, RIP1(301-558) and RIP1(301-558, *mutRHIM*). Lysates were subjected to IP with anti-FLAG antibody followed by IB with anti-Myc or anti-FLAG antibody. (M) IP/IB to detect HA-ICP10 with FLAG-tagged RIP3 or RIP3(*mutRHIM*). Lysates were subjected to IP with anti-FLAG antibody followed by IB with anti-HA or anti-FLAG antibodies. The lower panels depict IB of input whole cell lysates.




**Figure S3, related to Figure 3.** (A) IB to detect HT-29 cells stably expressing the indicated FLAG-tagged ICP6 constructs. (B) IP/IB to detect interaction of Myc-tagged ICP6, ICP10, UL45 or M45, or EV, with V5His-Casp8 in transfected 293T cells. Lysates were subjected to IP with anti-Myc antibody followed by IB with anti-His or anti-Myc antibody. The vertical line shows where lanes from the original gel were brought adjacent.

**Figure S4, related to Figure 4.** (A) Viability of HT-29-EV, HT-29-ICP6, HT-29-ICP6*mut*RHIM and HT-29-ICP6(1-629) cells infected with  $\Delta$ ICP6 (MOI=5) and treated from 1 to 24 hpi with T, T+V, T+S7, in the absence or presence GSK'963 (1  $\mu$ M) or GSK'840 (3  $\mu$ M) as shown. (B) Viability of HT-29 cells infected with HSV1 (F strain) or HSV1*mut*RHIM (MOI=10) or mock-infected and then treated from 1 to 24 hpi with T+S7, T+S7+GSK'963, T+S7+GSK'840. (C-E) Viability of HT-29-EV, HT-29-ICP10, HT-29-ICP10( $\Delta$ 1-249) cells, either mock-infected (C) or infected with  $\Delta$ ICP6 (D) or KOS (E) virus (MOI=5) and treated from 2 to 20 hpi with T, S and/or V. (F) IB to detect HT-29 cells stably expressing FLAG-tagged ICP10, ICP10( $\Delta$ 1-249). (G) IB to detect Casp3 cleavage (Cl-Casp3) and MLKL phosphorylation (p-MLKL) in HT-29 cells stably expressing EV, ICP10( $\Delta$ 1-249), ICP6( $\Delta$ 1-243) and ICP6*mut*RHIM treated with T+S for 8h. (H) Viral plaque numbers of HT-29 cells infected with HSV1 and HSV1*mut*RHIM virus for 72 hs in the absence or presence of T, T+GSK'963 (1  $\mu$ M), or PFA (150  $\mu$ g/ml).


Supplemental Figure 1



Supplemental Figure 2



Supplemental Figure 3



Supplemental Figure 4

