



SERVALL**ANGLE CENTRIFUGES**

*Ask for Bulletin L-121
with complete data on all models*



Large capacity
G-type machines


Type G/1, for 6 bottles of 250 cc each
18 tubes of 60 cc each
42 tubes of 15 cc each
78 tubes of 10 cc each
Type G/2, for 4 bottles of 500 cc each
Type G/3, for 10 tubes of 100 cc each



*Ask for
Bulletin L-122*

**STERLING AUTOMATIC
AND
HAND PIPETTES**

Dispense accurately, at a set speed, pre-determined amounts of 0.1 cc to 10.0 cc.
Easy to Operate—Time Saving



*Ask for
Bulletin L-123*

ORIGINAL—ODHNER

Portable Calculator for the Scientist
Many different models

Efficient—Sturdy—Low Priced

Manufacturers and Distributors:

IVAN SORVALL, Inc.

210 FIFTH AVE.

NEW YORK 10, N. Y.

PERCIVAL'S MEDICAL ETHICS

Edited by Chauncey D. Leake

University of Texas, Medical Branch

Thomas Percival in 1792 gave to the medical profession a code of ethics upon which present day codes have been largely based. In this printing, interesting additions have been made in the form of several essays by Professor Leake on the ethics and moral obligations of physicians, and the influence of Percival's code today.

I. Of Professional Conduct, Relative to Hospitals, or Other Medical Charities

II. Of Professional Conduct in Private or General Practice

III. Of the Conduct of Physicians Towards Apothecaries

IV. Of Professional Duties, in Certain Cases Which Require a Knowledge of Law

Selection from the "Notes and Illustrations" Appended to Percival's Code

305 pages

8 illus.

\$2.00




**THE WILLIAMS
& WILKINS CO.**
BALTIMORE 2, MARYLAND

New Lusteroid Tubes

FOR HIGH SPEED CENTRIFUGING

For Angle Heads on Multispeed Attachment


When used in conjunction with the Multispeed Attachment on the Reinforced Size 1, Type SB, Reinforced Size 2 and Model PR-1 Centrifuges, Head No. 295 swings six 10 ml. tubes at 18,000 r.p.m.; Head No. 296 swings four 25 ml. tubes at 18,000 r.p.m. At these speeds the problem of glassware breakage can be overcome by the use of Lusteroid (celluloid) tubes. As illustrated here, they are designed as inserts for aluminum or stainless steel shields which fit directly into the holes in the head. Plastic caps are available with tapered flanges which not only effectively seal the Lusteroid tubes but also prevent the tubes collapsing at the top under the high centrifugal force (25,000 x gravity).




300 658-671



HEAD NO. 295



298 662-664




HEAD NO. 296


295 Head, 6-place, Duralumin for No. 300 aluminum tubes and No. 658 Lusteroid tubes	\$ 75.00
300 Tube, Aluminum, with flanged cover. Capacity 10 ml. No rubber cushion required. For No. 295 head only	7.50
658 Lusteroid Tube, 10 ml., for No. 300 tube	.10
671 Plastic Cap for No. 658 Lusteroid tube	.80
296 Head, 4-place, Duralumin, for No. 298 stainless steel tubes and No. 664 Lusteroid tubes	83.00
298 Tube, Stainless Steel and Cover. Capacity 25 ml., for No. 296 head only	8.00
664 Lusteroid Tube, 25 ml., for No. 298 tube	.15
662 Plastic Cap for No. 664 Lusteroid tube	.85

For Large Capacity Angle Heads

Head No. 840 swings six 125 ml. tubes at speeds up to 5000 r.p.m. in the Size 1, Type SB, Size 2 and Model PR-1 Centrifuges. Head No. 845 swings eight 125 ml. tubes at speeds up to 4400 r.p.m. At these speeds Lusteroid tubes eliminate the problem of glassware breakage. When using the Lusteroid tubes, special rubber sleeve cushions No. 690 are also required. Plastic Caps No. 691 are available to seal the Lusteroid tubes and prevent collapsing at the top.



690 660 691



HEAD NO. 840-845

840 Head, 6-place	\$185.00
845 Head, 8-place	200.00
660 Lusteroid Tube, 125 ml., round bottom for use in No. 690 cushion	.30
690 Rubber Sleeve Cushion for No. 660 Lusteroid tube	2.30
691 Plastic Cap for No. 660 Lusteroid tube	1.00

Prompt delivery is available from your Laboratory Apparatus Supply Dealer.

INTERNATIONAL EQUIPMENT COMPANY

1284 SOLDIERS FIELD ROAD, BOSTON 35, MASS.



DIFCO

PEPTONES AND HYDROLYSATES

for Microbiological Culture Media

Difco peptones and hydrolysates are prepared to meet the diversified nutritional and biochemical requirements of microorganisms. A single peptone cannot provide the essential nitrogenous nutriment for all cultural and metabolic processes. Difco peptones and hydrolysates have been developed to meet the various requirements for microbiological procedures including growth, preparation of toxin and vaccines, and for the study of microbial metabolism.

BACTO-PEPTONE is most widely used in the preparation of routine culture media. It is rich in readily available forms of nitrogen, and in a one per cent solution it is sparklingly clear with a reaction of pH 7.0. Bacto-Peptone is specified in the formulae recommended in "Standard Methods of Water Analysis" of the American Public Health Association.

BACTO-TRYPTOSE is a peptone which was originally developed for cultivation of the *Brucella*. A two per cent solution, as the sole source of nitrogen, is an excellent substitute for the meat infusion generally employed for propagation of the *Streptococci*, *Pneumococci*, *Meningococci* and other discriminative bacteria.

PROTEOSE PEPTONE is universally used in the preparation of media for elaboration of diphtheria toxin. It is also an excellent nutriment for use in media designed for the production of other bacterial toxins such as those of scarlet fever and botulinus.

PROTEOSE PEPTONE NO. 3 is particularly suitable for use in culture media employed for isolation and propagation of *Neisseria gonorrhoeae* and *Corynebacterium diphtheriae*.

NEOPEPTONE is especially recommended as an ingredient of culture media for isolation and study of *Streptococci* and *Pneumococci*. Media prepared with Neopeptone are most satisfactory for cultivation of *Streptococci* in the smooth or mucoid phase. Growth of organisms from relatively small inocula is readily initiated in media containing this peptone.

BACTO-TRYPTONE is especially adapted for the elaboration of Indol and for use in media to detect Hydrogen Sulfide production. This peptone is recommended in "Standard Methods for the Examination of Dairy Products" of the American Public Health Association for the medium used in milk counts.

BACTO-CASITONE is a pancreatic digest of casein recommended for use in the preparation of media for sterility testing according to standard procedures, and for the preparation of media requiring an enzymatic hydrolyzed casein. Bacto-Casitone can be used for the detection of Indol production by microorganisms, since it has a high tryptophane content.

BACTO-CASAMINO ACIDS is an acid hydrolyzed casein originally prepared for the production of diphtheria toxin in peptone free media, and since adapted for other toxins and vaccines. Iron, copper, and other heavy metals have been removed from Bacto-Casamino Acids. This hydrolysate is excellently suited for nutritional studies of microorganisms.

BACTO-CASAMINO ACIDS TECHNICAL is an acid hydrolyzed casein recommended for use in media where amino acids mixture are required as a source of nitrogen, and the sodium chloride content is not a limiting factor. Bacto-Casamino Acids Technical is of particular value in nutritional studies of bacteria.

BACTO-VITAMIN FREE CASAMINO ACIDS is acid hydrolyzed casein, free from vitamins. It is recommended for the preparation of media for microbiological assay of vitamins and tryptophane. This hydrolysate is a valuable aid in studies of growth requirements of microorganisms.

Specify "DIFCO"

DIFCO LABORATORIES

DETROIT 1, MICHIGAN