

World-Wide Abstract Service **EXCERPTA MEDICA**

*An International Medical
Abstract Service published
in English and covering the
whole field of Clinical
and Experimental Medicine*

I. Anatomy, Anthropology, Embryology and Histology	\$22.50
II. Physiology, Biochemistry and Pharmacology	45.00
III. Endocrinology	15.00
IV. Medical Microbiology and Hygiene	22.50
V. General Pathology and Pathological Anatomy	37.50
VI. Internal Medicine	37.50
VII. Pediatrics	15.00
VIII. Neurology and Psychiatry	22.50
IX. Surgery	25.00
X. Obstetrics and Gynecology	15.00
XI. Oto-rhino-laryngology	15.00
XII. Ophthalmology	15.00
XIII. Dermatology and Venereology ...	25.00
XIV. Radiology	15.00
XV. Tuberculosis	15.00

Special rate to libraries ordering
all 15 sections \$250.00

*Published monthly, one volume
a year. Subscription by volume
only, payable in advance.
Send for a sample copy of
the section which covers
your specialty*


**THE WILLIAMS &
WILKINS COMPANY**

Mt. Royal and Guilford Aves.
Baltimore 2, Md.

Sole distributor for the United States,
Canada and Central America

SERVALL ANGLE CENTRIFUGES


*Ask for Bulletin L-121
with complete data on all models*



Large capacity
G-type machines

Type G/1, for 6 bottles of 250 cc each
18 tubes of 60 cc each
42 tubes of 15 cc each
78 tubes of 10 cc each


Type G/2, for 4 bottles of 500 cc each
Type G/3, for 10 tubes of 100 cc each



*Ask for
Bulletin L-122*

STERLING AUTOMATIC AND HAND PIPETTES

Dispense accurately, at a set speed, pre-
determined amounts of 0.1 cc to 10.0 cc.
Easy to Operate—Time Saving



*Ask for
Bulletin L-123*

ORIGINAL—ODHNER

Portable Calculator for the Scientist
Many different models

Efficient—Sturdy—Low Priced

Manufacturers and Distributors:


IVAN SORVALL, Inc.

210 FIFTH AVE.


NEW YORK 10, N. Y.

International Accessories

Only an International Centrifuge can offer you such a wide range of interchangeable accessories. No other Centrifuge is so adaptable to your requirements present and future.

Only genuine International heads, metal shields, trunnion rings, cups, and rubber cushions can give you peak performance from your International Centrifuge. Fifty years of engineering experience insure safety, reliability and smooth operation. Each accessory is stamped with the  monogram and each individual piece is test run under overload for your protection.

Illustrated here are a few of the accessories for the International Size 2. Send us the serial number of your Centrifuge and we will gladly send you a descriptive Bulletin listing all of the accessories usable on your machine.



INTERNATIONAL EQUIPMENT CO.

Boston 35, Mass., U. S. A.



DIFCO

Penicillin Inactivation

with

BACTO-PENASE

Bacto-Penase is a sterile thermolabile biological extract possessing the capacity to neutralize the antibiotic properties of all known types of *penicillin*. It is non-toxic to microorganisms and is therefore particularly satisfactory for use in the sterility testing of *penicillin* or products containing *penicillin*, and in determining the microbial counts on *penicillin*-containing products. BACTO-PENASE is recommended as a penicillinase for use in official procedures for sterility testing of *penicillin*. It is also recommended for use in media for culturing organisms from blood and other body fluids containing *penicillin* and for estimation of *penicillin* levels in such fluids.

PENICILLIN STERILITY TESTS

BACTO-PENASE in a concentration of 1 ml. per 15 ml. of Bacto-Fluid Thioglycollate Medium inactivates 50,000 units of *penicillin* in 40-48 hours incubation or 10,000 units of *penicillin* in 18 hours incubation.

BLOOD CULTURE

BACTO-PENASE is recommended for use in media for culturing blood or other body fluids containing *penicillin*. It is suggested that 1 ml. of BACTO-PENASE be added to 100 ml. of culture media such as Bacto-Brain Heart Infusion with PAB or Bacto-Brain Heart Infusion with PAB and Agar for this purpose. Media so prepared will inactivate the maximum amount of *penicillin* and sulfonamides present in 10 ml. or less of blood or other body fluids, and will also inactivate up to 1000 units of *streptomycin*.

Stability. BACTO-PENASE and media containing this penicillinase are stable when stored at refrigerator temperatures (2°-6°C.).

Bacto-Penase Concentrate. BACTO-PENASE is available in concentrated form as BACTO-PENASE CONCENTRATE, 1 ml. of which neutralizes 500,000 units of *penicillin*. BACTO-PENASE CONCENTRATE is recommended in the sterility testing of *penicillin*.

BACTO-PENASE and BACTO-PENASE CONCENTRATE are supplied in a standard package of one-half dozen screw cap tubes of 20 ml. each.

Descriptive literature is available upon request.

Specify "DIFCO"

THE TRADE NAME OF THE PIONEERS

In the Research and Development of Bacto-Peptone and Dehydrated Culture Media

DIFCO LABORATORIES

DETROIT 1, MICHIGAN