


## SUPPLEMENTARY DATA

### Supplementary appendix 1

FASTQ files were further processed using the CLC Genomics Server v5.5 software. Reads were trimmed to remove bases with a phred-score  $> 0.01$ , all ambiguous reads, the last 3' base, reads with a length below 30 and finally all duplicate reads were merged into one. Mapping was performed using default settings except length fraction and similarity fraction were set to 0.95. Variants were called using Probabilistic Variant Detection (minimum coverage 10, variant probability  $> 95.0$ ). To identify variants specific for eosinophils, the eosinophil variants were filtered against the lymphocytes mapped reads. The following filter criteria were applied: Minimum read count  $> 3$ , frequency  $> 30\%$ , forward reverse balance  $> 0.2$ , control frequency  $< 10\%$  and control coverage  $> 9$ .

We allowed a control frequency  $< 10\%$  since eosinophil sample preparation can contain small contaminations from lymphocytes. Subsequently, the lists of variants were imported into Ingenuity Variant Analysis for further filtering. Here, the variants were filtered to remove variants in the top 0.2% most exonically variable 100base windows in healthy public genomes and variants in the top 1% most exonically variable genes in healthy public genomes (1000 Genomes). Next, only variants predicted to be potentially deleterious were included, being variants with frameshift, in-frame indel, nonsense mutations, missense or likely splice site loss (up to 10 bases into intronic region). Finally, variants passing all these criteria were manually inspected in CLC Genomic Workbench selecting only true *de novo* mutations.



**Supplementary Figure S1: Validation of the DNA methylation data obtained using Illumina's 450K BeadChip array.** The DNA methylation level of two CpG sites was verified using pyrosequencing and for each CpG site the verified DNA methylation level correlated highly between pyrosequencing and Illuminas 450K BeadChip array with a coefficient of determination higher than 0.91.

**Supplementary Table S1: Assays used for Sanger sequencing for confirmation of mutations in diagnostic samples and validating the genome-wide DNA methylation analysis, respectively**

Gene	Primer sequence		Genomic location	
	Location of probes	Primer sequence*	Pyrosequencing primer	No. of CpG sites analyzed
<i>CDH17</i>		F: 5'TCATGCCCGGACTGTCTACAG R: 5'AGAAACAATGCCCTTCCAAGGGTG		
<i>PUF60</i>		F: 5'CTGGGTTGACCGGTCTTTCC R: 5'TCTATGATGGGCTGGGCCTG		
<i>LMLN</i>		F: 5'AAGCAAATTACAGCATGGCTGAG R: 5'CTTTTGTCTCTGCTGATCAATCCA		
<i>AQP12A</i>		F: 5'TCCCTGCAGAAAGTTCCTCATGGC R: 5'TGCAGGTCAGTACTGAGCTCCCA		
<i>PCSK1</i>		F: 5'AGACCGAAAGCGCTTCACTGA R: 5'AGCAAGATAGGAGAGAAAGCCAGA		
<b>Probe no.</b>	<b>Location of probes</b>	<b>Primer sequence*</b>	<b>Pyrosequencing primer</b>	<b>No. of CpG sites analyzed</b>
cg00620824( <i>HLA-C</i> )	1500TSS	F: 5'-gaagTaggggttgtaTtgtTattgTaa R: 5'-AtccaaataaataAcacactActtaAat	F: 5'-ggtttgtaTtgtTattgTaaT	1
cg05030953( <i>HLA-C</i> )	1500TSS	F: 5'aaggagTagaggaagaattTaaagTagt R: 5'-aAcccattaAttttaaAcaAtcacaca	F: 5'-agTTtgTagggggt	1

\* In the DNA methylation assays, bisulfite converted non-CpG cytosines are indicated as T/A on the sense/antisense strand, respectively

**Supplementary Table S2: The 285 probes corresponding to 128 unique genes that were differentially methylated in samples from patients with known and suspected clonal eosinophilia (S samples) and patients with reactive eosinophilia (R samples). The table is sorted by gene name.**

**Supplementary Table S3: Differentially methylated oncogenic signature genes in our dataset**

Gene name	Genomic annotation*	CpG neighborhood	Comment
<i>TSKS</i>	Body	Island	hypermethylation in S
<i>SPTLC2</i>	Body	Open	hypermethylation in S
<i>C1orf109</i>	Promoter	Shore	hypermethylation in S
<i>HCG9</i>	Promoter	Shore	hypermethylation in S
<i>CREB3L2</i>	3'UTR	Open	hypermethylation in S
<i>PM20D1</i>	Promoter	Island	hypermethylation in S
<i>GSTM1</i>	Promoter	Shore	hypermethylation in S
<i>RPS6KA2</i>	Body	Open	hypermethylation in S
<i>PF4</i>	Promoter	Island	hypermethylation in R
<i>TRIM41</i>	Promoter	Shore	hypermethylation in R
<i>PTH1R</i>	Body	Shelf	hypermethylation in R
<i>LRRC61</i>	Promoter	Island	hypermethylation in R
<i>ANKRD53</i>	Promoter	Island	hypermethylation in R
<i>TP53I13</i>	3'UTR	Island	hypermethylation in R
<i>SLC17A3</i>	Promoter	Open	hypermethylation in R
<i>EPS8L1</i>	Body	Island	hypermethylation in R
<i>SPTBN1</i>	Body	Shore	hypermethylation in R
<i>C3orf32</i>	Promoter	Open	hypermethylation in R
<i>KIAA1274</i>	Promoter	Open	hypermethylation in R
<i>SERHL</i>	Promoter	Island	hypermethylation in R
<i>KCNK3</i>	Body	Island	hypermethylation in R
<i>OXT</i>	Promoter	Island	hypermethylation in R
<i>TCF7L2</i>	Body	Shore	hypermethylation in R
<i>RAI1</i>	Promoter	Island	hypermethylation in R
<i>HSD17B1</i>	Body	Island	hypermethylation in R
<i>BTG2</i>	Promoter	Shore	hypermethylation in R
<i>FN3K</i>	Body	Island	hypermethylation in R
<i>TNXB</i>	Body	Island	hypermethylation in R
<i>CRIP2</i>	Body	Island	hypermethylation in R
<i>SLC39A4</i>	Body	Island	hypermethylation in R
<i>C3</i>	Promoter	Open	hypermethylation in R

\* As per UCSC genome Table. Promoter is 2000bp across TSS