

Supplementary information

Deciphering the relationship among phosphate dynamics, electron-dense body and lipid accumulation in the green alga *Parachlorella kessleri*

Shuhei Ota^{1,2}, Mai Yoshihara¹, Tomokazu Yamazaki^{1,2}, Tsuyoshi Takeshita¹, Aiko Hirata³, Mami Konomi⁴, Kenshiro Oshima⁵, Masahira Hattori⁵, Kateřina Bišová⁶, Vilém Zachleder⁶, Shigeyuki Kawano^{1,2*}

¹Department of Integrated Biosciences, Graduate School of Frontier Sciences, University of Tokyo, Kashiwa, Chiba, 277-8562, Japan. ²CREST, Japan Science and Technology Agency, Tokyo, Japan. ³Bioimaging Center, Graduate School of Frontier Science, University of Tokyo, Kashiwa, Chiba, 277-8562, Japan. ⁴Hitachi High-Technologies Corporation, Science & Medical Systems Business Group, Nishi-shinbashi, Tokyo, 105-8717, Japan. ⁵Center for Omics and Bioinformatics, Graduate School of Frontier Sciences, University of Tokyo, Kashiwa, Chiba, 277-8561, Japan. ⁶Institute of Microbiology, CAS, Centre Algatech, Laboratory of Cell Cycles of Algae, Třeboň, Czech Republic.




Figure S1. PolyP visualization of *P. kessleri* cells with high-concentration DAPI staining. (a) DAPI staining of cells growing on TAP medium. (b) DAPI staining of cells growing on dSTAP medium. (c) DAPI staining of cells growing on dNTAP medium (N-depleted medium). (d) DAPI staining of cells growing on dPTAP medium. Poly-P is stained yellow.




Figure S2. Rescaled representation of a growth curve in S-deficient medium. Values are means \pm standard deviation (S.D.) of four independent assays from the same batch culture.




Figure S3. Remaining and stored phosphate dynamics. (a) The amount of phosphate remaining in the medium. (b) Percent of stored poly-P in each medium.




Figure S4. EDX analysis of DBs in TAP and the resin background. EDX spectrometry of DBs in cells cultured in TAP (a) and the resin background (b). P accumulation in cells cultured in TAP medium was analyzed linearly by EDX (c). (Insets) STEM images showing the site (red lines) of the line analysis in c. Asterisks in c show the starting points of the EDX linear analysis.




Figure S5. Row count data from the EDX analysis. (a) Non-normalized EDX counts corresponding to Fig. 4a (dSTAP). (b) Non-normalized EDX counts corresponding to Fig. 4b (dPTAP).




Figure S6. Ultrastructure of DBs in P-depleted cells. TEM observations of cells in 6-day-old P-depleted culture (a and b). Cp: chloroplast, Db: electron-dense body, St: starch.

Table S1. Recipe of TAP medium

Component	Stock Solution (mg/mL)	Quantity Used*
NH ₄ Cl	200	2.0 mL
CaCl ₂ · 2H ₂ O	100	0.5 mL
MgSO ₄ · 7H ₂ O	300	1.0 mL
K ₂ HPO ₄	100	1.0 mL
KH ₂ PO ₄	100	1.0 mL
Hutner's trace elements	See following recipe	1.0 mL
Acetic acid	-	1.0 mL
Tris (hydroxymethyl) aminomethane	-	2.42 g

*To prepare: make up to 1 liter with distilled water (DW). Autoclave at 121°C for 20 minutes.

Recipe of Hutner's trace elements

Component	Quantity Used (per 100 mL DW)**
Na ₂ EDTA · 2H ₂ O	5.000 g
ZnSO ₄ · 7H ₂ O	2.200 g
H ₃ BO ₃	1.140 g
MnCl ₂ · 4H ₂ O	0.506 g
FeSO ₄ · 7H ₂ O	0.499 g
CoCl ₂ · 6H ₂ O	0.161 g
CuSO ₄ · 5H ₂ O	0.157 g
(NH ₄) ₆ Mo ₇ O ₂₄ · 4H ₂ O	0.110 g
KOH	~1.6 g

**Dissolve each of the above components, and then bring the final volume up to 100 mL with DW.

Table S2. Recipe of dSTAP medium

Component	Stock Solution (mg/mL)	Quantity Used*
NH ₄ Cl	200	2.0 mL
CaCl ₂ · 2H ₂ O	100	0.5 mL
MgCl ₂ · 6H ₂ O	254	1.0 mL
K ₂ HPO ₄	100	1.0 mL
KH ₂ PO ₄	100	1.0 mL
Hutner's trace elements for dSTAP	See following recipe	1.0 mL
Acetic acid	-	1.0 mL
Tris (hydroxymethyl) aminomethane	-	2.42 g

*To prepare: make up to 1 liter with DW. Autoclave at 121°C for 20 minutes.

Recipe of Hutner's trace elements for dSTAP

Component	Quantity Used (per 100 mL DW)**
Na ₂ EDTA · 2H ₂ O	5.000 g
ZnCl ₂	1.040 g
H ₃ BO ₃	1.140 g
MnCl ₂ · 4H ₂ O	0.506 g
FeCl ₃ · 7H ₂ O	0.485 g
CoCl ₂ · 6H ₂ O	0.161 g
CuCl ₂ · 2H ₂ O	0.107 g
(NH ₄) ₆ Mo ₇ O ₂₄ · 4H ₂ O	0.110 g
KOH	~1.6 g

**Dissolve each of the above components, and then bring the final volume up to 100 mL with DW.

Table S3. Recipe of dPTAP medium

Component	Stock Solution (mg/mL DW)	Quantity Used*
NH ₄ Cl	200	2.0 mL
CaCl ₂ · 2H ₂ O	100	0.5 mL
MgCl ₂ · 6H ₂ O	300	1.0 mL
KCl	140	1.0 mL
Hutner's trace elements	See the recipe in Table S1	1.0 mL
Acetic acid	-	1.0 mL
Tris (hydroxymethyl) aminomethane	-	2.42 g

* To prepare: make up to 1 liter with DW. Autoclave at 121°C for 20 minutes.

Movies S1-S10. 3D-representation movies corresponding to Figure 3. Color legends are as shown in Figure 3. Movies S1–S2: 3D-representation movies of a zero-control cell. Movie S3: 3D-representation movie of DBs in a zero-control cell. Movies S4–S5: 3D-representation movies of a starch cell. Movie S6: 3D-representation movie of DBs in a starch cell. Movies S7–S9: 3D-representation movies of a lipid cell. Movie S10: 3D-representation movie of DBs in the lipid cell.