

Spectral Data Processing Algorithm

Igor Pro Procedure:

This simple procedure is designed to take columns of spectral data and take the area under the curve as a function of the change in wavelength. This takes advantage of a built-in Igor Pro 6.3 function called 'areaXY'.

Function BIR5(W, B1, B2, B3, B4, B5, NSB, SP, index) //input function – designed to intergrate input waves B1, B2, B3, B4, B5, NSB, SP where ΔX is determined by the ΔW – usually about 0.3 nm
//The number of input background variables should be altered if using a data set that has more or less responses recorded.

//Variables

Wave NSB, SP //non-specific binding and specific binding spectra

Wave B1, B2, B3, B4, B5 //background spectra

Wave W //wavelength range scanned

Variable index

Make /O/N = 100 NSB_int, SP_int, B1_int, B2_int, B3_int, B4_int, B5_int, ratio

//generates the empty waves to allocate the resulting values into

//if processing larger data sets (greater than 100 values expected), the size of these waves should be changed accordingly

NSB_int[index] = areaXY(W, NSB) //calculates the area under the curve of NSB and using the wavelength as the limits of integration...trapezoidal method

//index is using an assigned value to direct results to a specific location in the wave NSB_int – following commands are engineered in an identical fashion, thus assigning values from one experiment to the same row in each wave

SP_int[index] = areaXY(W, SP) //calculates the area under the curve of SP and using the wavelength as the limits of integration...trapezoidal method

B1_int[index] = areaXY (W,B1)

B2_int[index] = areaXY (W,B2)

B3_int[index] = areaXY (W,B3)

B4_int[index] = areaXY (W,B4)

B5_int[index] = areaXY (W,B5)

ratio[index] = SP_int[index]/NSB_int[index]

End