



## Supporting Information

© 2016 The Authors. Published by Wiley-VCH Verlag GmbH & Co. KGaA, Weinheim

### **Novel Natural Oximes and Oxime Esters with a Vibralactone Backbone from the Basidiomycete *Boreostereum vibrans***

He-Ping Chen,<sup>[a, c]</sup> Zhen-Zhu Zhao,<sup>[a, c]</sup> Zheng-Hui Li,<sup>[a, b]</sup> Ze-Jun Dong,<sup>[a]</sup> Kun Wei,<sup>[a]</sup> Xue Bai,<sup>[a]</sup> Ling Zhang,<sup>[a]</sup> Chun-Nan Wen,<sup>[a, c]</sup> Tao Feng,<sup>\*[a, b]</sup> and Ji-Kai Liu<sup>\*[a, b]</sup>

open\_201500198\_sm\_miscellaneous\_information.pdf

# Supporting Information

## Table of Contents

1. Bioassays .....	6
1.1 Assessment of the inhibitory activity against pancreatic lipase. ....	6
1.2 Results of the inhibitory activity against pancreatic lipase of compounds <b>1-16</b> . ....	7
1.3 Cytotoxicity assay. ....	8
2. Semisynthesis of compound <b>1</b> from vibralactone.....	9
2.1 From vibralactone to (1 <i>R</i> ,5 <i>S</i> )-1-(3-methylbut-2-en-1-yl)-7-oxo-6-oxabicyclo[3.2.0]hept-2-ene-3-carbaldehyde ( <b>17</b> ).....	9
2.2 From (1 <i>R</i> ,5 <i>S</i> )-1-(3-methylbut-2-en-1-yl)-7-oxo-6-oxabicyclo[3.2.0]hept-2-ene-3-carbaldehyde ( <b>17</b> ) to vibralactoxime A ( <b>1a</b> ).....	9
Figure A. $^1\text{H}$ NMR spectrum of <b>17</b> (600 MHz, acetone-d <sub>6</sub> ). ....	10
Figure B. $^1\text{H}$ NMR spectrum of <b>1a</b> (600 MHz, acetone-d <sub>6</sub> ) .....	10
Figure C. The CD and UV spectra of vibralactone (solvent: MeCN). ....	11
Figure D. The CD and UV spectra of <b>1</b> (solvent: MeCN).....	12
Figure E. The CD and UV spectra of synthesized <b>1a</b> (solvent: MeCN).....	13
3. Representative HPLC spectra of some final compounds. ....	14
Figure F. HPLC spectra of compound <b>1</b> ( $t_R$ = 8.71 min). ....	14
Figure G. HPLC spectra of compound <b>2</b> ( $t_R$ = 8.07 min). ....	14
Figure H. HPLC spectra of compound <b>4</b> ( $t_R$ = 9.82 min). ....	15
Figure I. HPLC spectra of compound <b>6</b> ( $t_R$ = 8.92 min).....	15
Figure J. HPLC spectra of compound <b>8</b> ( $t_R$ = 9.19 min).....	16
Figure K. HPLC spectra of compound <b>13</b> ( $t_R$ = 8.73 min). ....	16
Table 1S. $^1\text{H}$ (600 MHz) and $^{13}\text{C}$ NMR (150 MHz) spectral data of <b>1</b> and <b>2</b> in acetone-d <sub>6</sub> .	17
Table 2S. $^1\text{H}$ (600 MHz) and $^{13}\text{C}$ NMR (150 MHz) spectral data of <b>3</b> in acetone-d <sub>6</sub> .....	18
Table 3S. $^1\text{H}$ (600 MHz) and $^{13}\text{C}$ NMR (150 MHz) spectral data of <b>4</b> and <b>5</b> in acetone-d <sub>6</sub> .	19
Table 4S. $^1\text{H}$ (600 MHz) and $^{13}\text{C}$ NMR (150 MHz) spectral data of <b>6</b> and <b>7</b> in acetone-d <sub>6</sub> .	20
Table 5S. $^1\text{H}$ (600 MHz) and $^{13}\text{C}$ NMR (150 MHz) spectral data of <b>8</b> and <b>13</b> in acetone-d <sub>6</sub> .	21

Table 6S. $^1\text{H}$ (600 MHz) and $^{13}\text{C}$ NMR (150 MHz) spectral data of <b>9</b> and <b>14</b> in acetone- $d_6$ .	22
Table 7S. $^1\text{H}$ (600 MHz) and $^{13}\text{C}$ NMR (150 MHz) spectral data of <b>10</b> in acetone- $d_6$ .....	24
Table 8S. Characteristic $^1\text{H}$ NMR (600 MHz) spectral data of <b>11</b> , <b>12</b> , <b>15</b> and <b>16</b> in acetone- $d_6$ .....	25
4. Supplementary NMR spectra. ....	26
Figure 1S: $^1\text{H}$ NMR spectrum of <b>1</b> (acetone- $d_6$ , 600MHz).....	26
Figure 2S: $^{13}\text{C}$ NMR and DEPT spectra of <b>1</b> (acetone- $d_6$ , 150MHz) .....	27
Figure 3S: HSQC spectrum of <b>1</b> .....	28
Figure 4S: $^1\text{H}$ - $^1\text{H}$ COSYspectrum of <b>1</b> .....	29
Figure 5S: HMBC spectrum of <b>1</b> .....	30
Figure 6S: ROESYspectrum of <b>1</b> .....	31
Figure 7S: HREIMSspectrum of <b>1</b> .....	32
Figure 8S: Optical rotation of compound <b>1</b> .....	33
Figure 9S: $^1\text{H}$ NMR spectrum of <b>2</b> (acetone- $d_6$ , 600MHz) .....	34
Figure 10S: $^{13}\text{C}$ NMR and DEPT spectra of <b>2</b> (acetone- $d_6$ , 150MHz) .....	35
Figure 11S: HSQC spectrum of <b>2</b> .....	36
Figure 12S: $^1\text{H}$ - $^1\text{H}$ COSY spectrum of <b>2</b> .....	37
Figure 13S: HMBC spectrum of <b>2</b> .....	38
Figure 14S: ROESYspectrum of <b>2</b> .....	39
Figure 15S: HREIMSSpectrum of <b>2</b> .....	40
Figure 16S: $^1\text{H}$ NMR spectrum of <b>3</b> (acetone- $d_6$ , 600MHz).....	41
Figure 17S: Enlarged view of $^1\text{H}$ NMR spectra of <b>3</b> (acetone- $d_6$ , 600MHz) .....	42
Figure 18S: $^{13}\text{C}$ NMR and DEPT spectra of <b>3</b> (acetone- $d_6$ , 150MHz) .....	43
Figure 19S: Enlarged view of $^{13}\text{C}$ NMR spectrum of <b>3</b> (acetone- $d_6$ , 150MHz) .....	44
Figure 20S: HSQC spectrum of <b>3</b> .....	45
Figure 21S: $^1\text{H}$ - $^1\text{H}$ COSYspectrum of <b>3</b> .....	46
Figure 22S: HMBC spectrum of <b>3</b> .....	47
Figure 23S: Enlarged view of HMBC spectrum of <b>3</b> .....	48
Figure 24S: ROESYspectrum of <b>3</b> .....	49
Figure 25S: HRESIMS spectrum of <b>3</b> .....	50

Figure 26S: $^1\text{H}$ NMR spectrum of <b>4</b> (acetone- $d_6$ , 600MHz).....	51
Figure 27S: $^{13}\text{C}$ NMR and DEPT spectra of <b>4</b> (acetone- $d_6$ , 150MHz) .....	52
Figure 28S: Enlarged view of $^{13}\text{C}$ NMR spectrum of <b>4</b> (acetone- $d_6$ , 150MHz) .....	53
Figure 29S: HSQC spectrum of <b>4</b> .....	54
Figure 30S: $^1\text{H}$ - $^1\text{H}$ COSYspectrum of <b>4</b> .....	55
Figure 31S: HMBC spectrum of <b>4</b> .....	56
Figure 32S: Enlarged view of HMBC spectra of <b>4</b> . ....	57
Figure 33S: ROESYspectrum of <b>4</b> .....	58
Figure 34S: HREIMSspectrum of <b>4</b> .....	59
Figure 35S: IR spectrum of <b>4</b> .....	60
Figure 36S: $^1\text{H}$ NMR spectrum of <b>5</b> (acetone- $d_6$ , 600MHz).....	61
Figure 37S: $^{13}\text{C}$ NMR and DEPT spectra of <b>5</b> (acetone- $d_6$ , 150MHz) .....	62
Figure 38S: Enlarged view of $^{13}\text{C}$ NMR spectrum of <b>5</b> (acetone- $d_6$ , 150MHz) .....	63
Figure 39S: HSQC spectrum of <b>5</b> .....	64
Figure 40S: $^1\text{H}$ - $^1\text{H}$ COSYspectrum of <b>5</b> .....	65
Figure 41S: HMBC spectrum of <b>5</b> .....	66
Figure 42S: Enlarged view of HMBC spectra of <b>5</b> . ....	67
Figure 43S: ROESYspectrum of <b>5</b> .....	68
Figure 44S: IR spectrum of <b>5</b> .....	69
Figure 45S: HRESIMS spectrum of <b>5</b> .....	70
Figure 46S: $^1\text{H}$ NMR spectrum of <b>6</b> (acetone- $d_6$ , 600MHz).....	71
Figure 47S: $^{13}\text{C}$ NMR and DEPT spectra of <b>6</b> (acetone- $d_6$ , 150MHz) .....	72
Figure 48S: HSQC spectrum of <b>6</b> .....	73
Figure 49S: $^1\text{H}$ - $^1\text{H}$ COSYspectrum of <b>6</b> .....	74
Figure 50S: HMBC spectrum of <b>6</b> .....	75
Figure 51S: Enlarged view of HMBC spectra of <b>6</b> . ....	76
Figure 52S: ROESYspectrum of <b>6</b> .....	77
Figure 53S: HRESIMSSpectrum of <b>6</b> .....	78
Figure 54S: $^1\text{H}$ NMR spectrum of <b>7</b> (acetone- $d_6$ , 600MHz).....	80
Figure 55S: $^{13}\text{C}$ NMR and DEPT spectra of <b>7</b> (acetone- $d_6$ , 150MHz) .....	81
Figure 56S: HSQC spectrum of <b>7</b> .....	82

Figure 57S: $^1\text{H}$ - $^1\text{H}$ COSY spectrum of <b>7</b> .....	83
Figure 58S: Enlarged view of $^1\text{H}$ - $^1\text{H}$ COSY spectrum.....	84
Figure 59S: HMBC spectrum of <b>7</b> .....	85
Figure 60S: ROESY spectrum of <b>7</b> .....	86
Figure 61S: HRESIMSspectrum of <b>7</b> .....	87
Figure 62S: $^1\text{H}$ NMR spectrum of <b>8</b> (acetone- $d_6$ , 600MHz).....	88
Figure 63S: $^{13}\text{C}$ NMR and DEPT spectra of <b>8</b> (acetone- $d_6$ , 150MHz) .....	89
Figure 64S: HSQC spectrum of <b>8</b> .....	90
Figure 65S: $^1\text{H}$ - $^1\text{H}$ COSYspectrum of <b>8</b> .....	91
Figure 66S: HMBC spectrum of <b>8</b> .....	92
Figure 67S: ROESYspectrum of <b>8</b> .....	93
Figure 68S: HREIMSSpectrum of <b>8</b> .....	94
Figure 69S: $^1\text{H}$ NMR spectrum of <b>9</b> (acetone- $d_6$ , 600MHz).....	95
Figure 70S: $^{13}\text{C}$ NMR and DEPT spectra of <b>9</b> (acetone- $d_6$ , 150MHz) .....	96
Figure 71S: HSQC spectrum of <b>9</b> .....	97
Figure 72S: $^1\text{H}$ - $^1\text{H}$ COSYspectrum of <b>9</b> .....	98
Figure 73S: HMBC spectrum of <b>9</b> .....	99
Figure 74S: ROESYspectrum of <b>9</b> .....	100
Figure 75S: HREIMSSpectrum of <b>9</b> .....	101
Figure 76S: $^1\text{H}$ NMR spectrum of <b>10</b> (acetone- $d_6$ , 600MHz).....	102
Figure 77S: Enlarged view of $^1\text{H}$ NMR spectra of <b>10</b> (acetone- $d_6$ , 600MHz) .....	103
Figure 78S: $^{13}\text{C}$ NMR and DEPT spectra of <b>10</b> (acetone- $d_6$ , 150MHz) .....	104
Figure 79S: Enlarged view of $^{13}\text{C}$ NMR and DEPT spectra of <b>10</b> (acetone- $d_6$ , 150MHz) 105	105
Figure 80S: HSQC spectrum of <b>10</b> .....	106
Figure 81S: Enlarged view of HSQC spectrum of <b>10</b> .....	107
Figure 82S: $^1\text{H}$ - $^1\text{H}$ COSYspectrum of <b>10</b> .....	108
Figure 83S: HMBC spectrum of <b>10</b> .....	109
Figure 84S: Enlarged view of HMBC spectrum of <b>10</b> .....	110
Figure 85S: ROESYspectrum of <b>10</b> .....	111
Figure 86S: HRESIMSspectrum of <b>10</b> .....	112
Figure 87S: $^1\text{H}$ NMR spectrum of <b>11</b> (acetone- $d_6$ , 600MHz).....	113

Figure 88S: Enlarged view of $^1\text{H}$ NMR spectrum of <b>11</b> (acetone- $d_6$ , 600MHz) .....	114
Figure 89S: HRESIMSspectrum of <b>11</b> .....	115
Figure 90S: $^1\text{H}$ NMR spectrum of <b>12</b> (acetone- $d_6$ , 600MHz).....	116
Figure 91S: Enlarged view of $^1\text{H}$ NMR spectrum of <b>12</b> (acetone- $d_6$ , 600MHz) .....	117
Figure 92S: HRESIMSspectrum of <b>12</b> .....	118
Figure 93S: $^1\text{H}$ NMR spectrum of <b>13</b> (acetone- $d_6$ , 600MHz).....	119
Figure 94S: $^{13}\text{C}$ NMR and DEPT spectra of <b>13</b> (acetone- $d_6$ , 150MHz) .....	120
Figure 95S: Enlarged view of $^{13}\text{C}$ NMR and DEPT spectra of <b>13</b> (acetone- $d_6$ , 150MHz) 121	
Figure 96S: HSQC spectrum of <b>13</b> .....	122
Figure 97S: $^1\text{H}$ - $^1\text{H}$ COSY spectrum of <b>13</b> .....	123
Figure 98S: HMBC spectrum of <b>13</b> .....	124
Figure 99S: ROESY spectrum of <b>13</b> .....	125
Figure 100S: HREIMS spectrum of <b>13</b> .....	126
Figure 101S: $^1\text{H}$ NMR spectrum of <b>14</b> (acetone- $d_6$ , 600MHz).....	127
Figure 102S: $^{13}\text{C}$ NMR and DEPT spectra of <b>14</b> (acetone- $d_6$ , 150MHz) .....	128
Figure 103S: HSQC spectrum of <b>14</b> .....	129
Figure 104S: $^1\text{H}$ - $^1\text{H}$ COSYspectrum of <b>14</b> .....	130
Figure 105S: HMBC spectrum of <b>14</b> .....	131
Figure 106S: Enlarged view of HMBC spectrum of <b>14</b> .....	132
Figure 107S: ROESYspectrum of <b>14</b> .....	133
Figure 108S: HRESIMSspectrum of <b>14</b> .....	134
Figure 109S: $^1\text{H}$ NMR spectrum of <b>15</b> (acetone- $d_6$ , 600MHz).....	135
Figure 110S: Enlarged view of $^1\text{H}$ NMR spectra of <b>15</b> (acetone- $d_6$ , 600MHz) .....	136
Figure 111S: HRESIMSspectrum of <b>15</b> .....	137
Figure 112S: $^1\text{H}$ NMR spectrum of <b>16</b> (acetone- $d_6$ , 600MHz).....	138
Figure 113S: Enlarged view of $^1\text{H}$ NMR spectrum of <b>16</b> (acetone- $d_6$ , 600MHz) .....	139
Figure 114S: HRESIMSspectrum of <b>16</b> .....	140

## **1. Bioassays**

### **1.1 Assessment of the inhibitory activity against pancreatic lipase.**

Porcine pancreatic lipase was used to evaluate the inhibitory activity of the monascus pigment derivatives. *p*-Nitrophenyl butyrate (*p*-NPS) was used as the substrate. Firstly, an enzyme-buffer solution was prepared by adding 30 $\mu$ L of lipase solution (in 10 mM morpholine propane sulfonic acid (MOPS) and 1 mM EDTA, pH 6.8) to 850  $\mu$ L of Tris buffer (100 mM Tris–HCl and 5 mM CaCl<sub>2</sub>, pH 7.0). Pigment solutions were prepared by dissolving each pigment in a mixture (1:1) of ethanol and distilled water. Subsequently, 100  $\mu$ L of the pigment solution was mixed with 880  $\mu$ L of the enzyme-buffer solution. Following incubation of the enzyme-pigment mixture for 15 min at 37 °C, 20  $\mu$ L of the substrate solution (10 mM *p*-NPB in dimethyl formamide) was added. Enzymatic reactions were conducted for 15 min at 37 °C. The hydrolysis of *p*-NPB to *p*-nitrophenol was monitored at 400 nm using a spectrophotometer. One unit of the enzyme was defined as the amount required to liberate 1 mol of *p*-nitrophenol under standard assay conditions.

**Reference:** K. Wei, G. Q. Wang, X. Bai, Y. F. Niu, H. P. Chen, C. N. Wen, Z. H. Li, Z. J. Dong, Z. L. Zuo, W. Y. Xiong, J. K. Liu, *Nat. Prod. Bioprospect.* **2015**, 5, 129-157

## 1.2 Results of the inhibitory activity against pancreatic lipase of compounds 1-16.

Compound	Concentration ( $\mu\text{M}$ )	PPL inhibitory rate (%)	$\text{IC}_{50}$ ( $\mu\text{M}$ )
Orlistat	0.005	86.024	0.0018
	0.002	51.526	
	0.001	28.284	
Vibralactone	10000	91.800	48.668
	1000	87.038	
	100	57.678	
	10	33.130	
<b>1</b>	50	71.318	23.061
	10	26.847	
	2	5.197	
<b>2</b>	50	-13.600	
<b>3</b>	50	-14.573	
<b>4</b>	50	84.188	20.580
	10	22.026	
	2	-27.792	
<b>5</b>	50	78.439	28.619
	10	-3.936	
	2	-22.889	
<b>6</b>	50	19.483	
<b>7</b>	50	89.319	16.813
	10	31.203	
	2	-24.384	
<b>8</b>	50	45.843	
<b>9</b>	50	79.213	23.189
	10	17.979	
	2	-16.907	
<b>10</b>	50	67.824	19.812
	10	36.842	
	2	-35.413	
<b>11</b>	50	60.084	11.108
	10	49.160	
	2	8.359	
<b>12</b>	50	18.841	
<b>13</b>	50	-35.847	
<b>14</b>	50	-34.719	
<b>15</b>	50	-2.875	
<b>16</b>	50	-8.293	

Note: The compounds which PPL inhibitory rate less than 50% on the concentration 50  $\mu\text{M}$  were not to test their  $\text{IC}_{50}$  values.


### **1.3 Cytotoxicity assay.**

The following human tumor cell lines were used: HL-60, SMMC-7721, A-549, MCF-7, and SW-480. All cells were cultured in RPMI-1640 medium (Hyclone, Logan, UT, USA), supplemented with 10% fetal bovine serum (FBS, Hyclone) at 37 °C in a humidified atmosphere with 5% CO<sub>2</sub>. Cell viability was assessed by conducting colorimetric measurements of the amount of insoluble formazan formed in living cells based on the reduction of 3-(4,5-dimethylthiazol-2-yl)-2,5-diphenyltetrazolium bromide (MTT) (Sigma, St. Louis, MO, USA). Briefly, 100 µL of adherent cells was seeded into each well of a 96-well cell culture plate and allowed to adhere for 12 h before test compound addition, while suspended cells were seeded just before test compound addition, both with an initial density of 1 ×10<sup>5</sup> cells/mL in 100 µL of medium. Each tumor cell line was exposed to the test compound at various concentrations in triplicate for 48 h, with cis-platin and paclitaxel (Sigma) as positive control. After the incubation, MTT (100 µg) was added to each well, and the incubation continued for 4 h at 37 °C. The cells were lysed with 100 µL of 20% SDS-50% DMF after removal of 100 µL of medium. The optical density of the lysate was measured at 595 nm in a 96-well microtiter plate reader (Bio-Rad 680). The IC<sub>50</sub> value of each compound was calculated by Reed and Muench's method.

**References:** a) M. C. Alley, D. A. Scudiero, A. Monks, M. L. Hursey, M. J. Czerwinski, D. L. Fine, B. J. Abbott, J. G. Mayo, R. H. Shoemaker, M. R. Boyd, *Cancer Res.* **1988**, *48*, 589-601; b) J. J. Zhang, X. W. Yang, J. Z. Ma, X. Liu, L. X. Yang, S. C. Yang, G. Xu, *Nat. Prod. Bioprospect.* **2014**, *4*, 73-79.

## 2. Semisynthesis of compound 1 from vibralactone.


**2.1** From vibralactone to  
**(1*R*,5*S*)-1-(3-methylbut-2-en-1-yl)-7-oxo-6-oxabicyclo[3.2.0]hept-2-ene-3-carbaldehyde**  
**(17)**



To a solution of vibralactone (16 mg, 0.077 mmol) in  $\text{CH}_2\text{Cl}_2$  (1 mL) was added PCC (17.2 mg, 0.080 mmol). The mixture was stirred at 0 °C for 30 min, during which the progress of the reaction was monitored by TLC. The reaction mixture was passed through a short silica gel pad (eluent: ethyl acetate) and then purified by prep-HPLC [method: 30-50% (acetonitrile/water), 20 min, 8 mL/min, detected UV wavelength: 205 nm, 247 nm,  $t_R = 9.85$  min] to afford **17** (12.4 mg, 78% yield):  $^1\text{H}$  NMR (600 MHz, acetone- $d_6$ ):  $\delta_{\text{H}}$  6.93 (s, H-2), 2.88 (overlapped, H-4a), 2.94 (dd,  $J = 19.0, 13.0$  Hz, H-4b), 5.08 (d,  $J = 5.5$  Hz, H-5), 2.64 (dd,  $J = 15.0, 7.4$  Hz, H-8a), 2.78 (dd,  $J = 15.0, 7.4$  Hz, H-8b), 5.21 (br. t,  $J = 7.4$  Hz, H-9), 1.68 (s, 3H, H-11), 1.71 (s, 3H, H-12), 9.86 (s, H-13).

**Reference:** M. Y. Jiang, F. Wang, X. L. Yang, L. Z. Fang, Z. J. Dong, H. J. Zhu, J. K. Liu, *Chem. Pharm. Bull.* **2008**, 56, 1286-1288.


**2.2** From  
**(1*R*,5*S*)-1-(3-methylbut-2-en-1-yl)-7-oxo-6-oxabicyclo[3.2.0]hept-2-ene-3-carbaldehyde**  
**(17) to vibralactoxime A (1a)**




A solution of  $\text{NH}_2\text{OH} \cdot \text{HCl}$  (4.6 mg, 0.067 mmol) and sodium acetate (5.5 mg, 0.067 mmol) in  $\text{EtOH}/\text{H}_2\text{O}$  (1:1, 0.6 mL) was added dropwise to a solution aldehyde **17** (6.8 mg, 0.033 mmol) in  $\text{EtOH}$  (0.3 mL) at 0 °C. The mixture was stirred at 0 °C for 20 min, during which the progress of the reaction was monitored by TLC. The mixture was concentrated under reduced pressure then purified by prep-HPLC [method: 30-50% (acetonitrile/water), 20 min, 8 mL/min, detected UV wavelength: 247 nm,  $t_R = 15.6$  min] to afford **1a** (6.56 mg, 90% yield).

**Reference:** R. A. Aungst, Jr, C. Chan, R. L. Funk, *Org. Lett.* **2001**, 3, 2611-2613.

**Figure A.**  $^1\text{H}$  NMR spectrum of **17** (600 MHz, acetone- $d_6$ ).




**Figure B.**  $^1\text{H}$  NMR spectrum of **1a** (600 MHz, acetone- $d_6$ )



**Figure C.** The CD and UV spectra of vibralactone (solvent: MeCN).

LCVIB



File: CD LCVIB-1mm(195-400)15112612.dsx

ProBinaryX

Attributes :

- Time Stamp :Thu Nov 26 17:09:58 2015

- File ID : {257931F1-1A3F-4bdb-9C09-1407AE93E278}

- Is CFR Compliant : false

- Original unaltered data

Remarks:

- HV (CDCC channel): 0 v
- Time per point: 1 s
- Description: Sample 1
- Concentration: 0.0960mg/mL Acetonitrile
- Pathlength: 1 mm

Settings:

- Time-per-point: 1s (25us x 40000)


- Wavelength: 195nm - 400nm

- Step Size: 1nm

- Bandwidth: 1nm

**Figure D.** The CD and UV spectra of 1 (solvent: MeCN).

LCBVB-3A1



File: CD LCBVB-3A1-1mm(195-400)15112614.dsx

ProBinaryX

Attributes :

- Time Stamp :Thu Nov 26 17:42:25 2015

- File ID : {C095BDC9-0152-4636-892A-6F1935FC6F54}

- Is CFR Compliant : false

- Original unaltered data

Remarks:


- HV (CDDC channel): 0 v
- Time per point: 1 s
- Description: Sample 1
- Concentration: 0.0300mg/mL Acetonitrile
- Pathlength: 1 mm

Settings:

- Time-per-point: 1s (25us x 40000)
- Wavelength: 195nm - 400nm
- Step Size: 1nm
- Bandwidth: 1nm

**Figure E.** The CD and UV spectra of synthesized 1a (solvent: MeCN).

LCOX-1



File: CD LCOX-1-1mm(195-400)15112613.dsx

ProBinaryX

Attributes :

- Time Stamp :Thu Nov 26 17:23:50 2015

- File ID : {6FFFC8A3-7792-4b78-A7CA-A33391B5F007}

- Is CFR Compliant : false

- Original unaltered data

Remarks:

- HV (CDDC channel): 0 v

- Time per point: 1 s

- Description: Sample 1

- Concentration: 0.1700mg/mL Acetonitrile

- Pathlength: 1 mm

Settings:

- Time-per-point: 1s (25us x 40000)

- Wavelength: 195nm - 400nm

- Step Size: 1nm

- Bandwidth: 1nm


### 3. Representative HPLC spectra of some final compounds.

The HPLC profile were carried out on an agilent 1100 liquid chromatography system equipped with a Zorbax SB-C18 column ( $4.9 \times 150$  mm, particle size  $5\mu\text{m}$ ) and a DAD detector. Acetonitrile and water were used as the mobile phase. The method used to analyze these compounds was as follows:


Time	B% (Acetonitrile)	Flowrate (mL/min)
0	0	1
12	100	1
15	100	1
18	0	1

The detected wavelength of these compounds were at 247 and 205 nm.


**Figure F.** HPLC spectra of compound **1** ( $t_R = 8.71$  min).




**Figure G.** HPLC spectra of compound **2** ( $t_R = 8.07$  min).




**Figure H.** HPLC spectra of compound **4** ( $t_R = 9.82$  min).




**Figure I.** HPLC spectra of compound **6** ( $t_R = 8.92$  min).




**Figure J.** HPLC spectra of compound **8** ( $t_R = 9.19$  min).




**Figure K.** HPLC spectra of compound **13** ( $t_R = 8.73$  min).





**Table 1S.**  $^1\text{H}$  (600 MHz) and  $^{13}\text{C}$  NMR (150 MHz) spectral data of **1** and **2** in acetone- $d_6$ .


No.	<b>1</b>		<b>2</b>	
	$\delta_{\text{C}}$ , type	$\delta_{\text{H}}$ , multi.	$\delta_{\text{C}}$ , type	$\delta_{\text{H}}$ , multi.
1	76.5, s		66.7, s	
2	131.3, d	6.04, s	136.4, d	6.06, s
3	141.1, s		138.5, s	
4	36.7, t	2.88, d (19.2) 2.94, dd (19.2, 5.4)	39.9, t	2.54, dd (17.0, 1.8) 2.84, overlapped
5	78.9, d	4.97, d (5.4)	78.5, d	4.27, ddd (6.0, 5.9, 1.8)
7	172.4, s		173.5, s	
8	28.0, t	2.52, dd (15.0, 7.4) 2.67, dd (15.0, 7.4)	35.7, t	2.20, dd (13.9, 7.6) 2.62, dd (13.9, 7.6)
9	118.5, d	5.19, t (7.4)	120.4, d	5.06, t (7.6)
10	136.4, s		134.9, s	
11	18.0, q	1.66, s, 3H	17.9, q	1.57, s, 3H
12	25.9, q	1.70, s, 3H	26.0, q	1.66, s, 3H
13	146.5, d	7.99, s	147.3, d	7.93, s
7-OCH <sub>3</sub>			51.6, s	3.63, s, 3H
N-OH		10.53, s,		10.20, s
5-OH				4.18, d (5.9)



**Table 2S.**  $^1\text{H}$  (600 MHz) and  $^{13}\text{C}$  NMR (150 MHz) spectral data of **3** in acetone- $d_6$ .

<b>3</b>					
No.	$\delta_{\text{C}}$ , type	$\delta_{\text{H}}$ , multi.	No.	$\delta_{\text{C}}$ , type	$\delta_{\text{H}}$ , multi.
1	77.1, s		1'	61.4, s	
2	138.9, d	6.48, s	2'	64.9, d	3.77, s
3	139.3, s		3'	65.3, s	
4	36.6, t	3.08, dd (18.4, 5.8) 2.99, d (18.4)	4'	37.2, t	2.44, dd (15.2, 6.1) 2.10, d (15.2)
5	79.1, d	5.08, d (5.8)	5'	76.8, d	3.95, dd (7.5, 6.1 )
7	171.7, s		7'	168.7, s	
8	28.0, t	2.74, dd (14.8, 7.5) 2.60, overlapped	8'	32.8, t	2.60, overlapped 2.23, dd (14.3, 7.5)
9	118.4, d	5.20, br.t (7.5)	9'	118.8, d	5.16, br. t (7.5)
10	136.8, s		10'	136.4, s	
11	18.1, q	1.67, s, 3H	11'	18.1, q	1.63, s, 3H
12	26.1, q	1.71, s, 3H	12'	26.0, q	1.69, s, 3H
13	154.2, d	8.44, s	13'	64.0, t	4.49, d, (12.5) 4.20, d (12.5)
CH <sub>3</sub> CO-	20.6, q	2.07, overlapped*			
CH <sub>3</sub> CO-	170.8, s				
5'-OH		2.98, d (7.5)			

Note: \* overlapped in solvent residual peaks.




**4** Vibralactoxime D

**5** Vibralactoxime E


**Table 3S.**  $^1\text{H}$  (600 MHz) and  $^{13}\text{C}$  NMR (150 MHz) spectral data of **4** and **5** in acetone- $d_6$ .

No.	<b>4</b>		<b>5</b>	
	$\delta_{\text{C}}$ , type	$\delta_{\text{H}}$ , multi.	$\delta_{\text{C}}$ , type	$\delta_{\text{H}}$ , multi.
1	77.0, s		77.0,s	
2	138.9, d	6.48, s	138.6,d	6.46,s
3	139.2, s		139.3,s	
4	36.5, t	3.08, dd (19.0, 5.7) 2.97, d (19.0)	36.5,t	3.08,br.dd(18.7,5.8) 2.97,d (18.7)
5	79.0, d	5.07, d (5.6)	79.1,d	5.07,d (5.8)
7	171.6, s		171.6,s	
8	27.9, t	2.73, overlapped 2.59, overlapped	27.9,t	2.72,overlapped 2.58,dd (15.5,7.5)
9	118.3, d	5.20, br. t (7.4)	118.3,d	5.20,br.t (7.5)
10	136.7, s		136.7,s	
11	18.0, q	1.67, s, 3H	18.1,q	1.67,s, 3H <sup>a</sup>
12	25.9, q	1.71, s, 3H	25.9,q	1.71,s, 3H <sup>b</sup>
13	154.3, d	8.45, s	154.0,d	8.42,s
1'	67.5, s		66.9,s	
2'	149.4, d	6.82, s	143.8,d	6.56,s
3'	115.0, s		137.2,s	
4'	43.3, t	3.11, dd (17.6, 5.8) 2.61, dd (17.6, 1.7)	39.8,t	3.00,dd(17.0,6.0) 2.68,d(17.0)
5'	77.8, d	4.44, ddd (5.8, 5.7, 1.7)	78.6,d	4.45,br.t (6.0)
7'	168.2, s		168.3,s	
8'	34.8, t	2.74, overlapped 2.33, dd (14.2, 7.8)	34.8,t	2.75,overlapped 2.34,dd (14.5, 7.5)
9'	119.0, d	5.11, br. t (7.8)	119.2,d	5.12,br.t (7.5)
10'	136.3, s		136.3,s	
11'	18.0, q	1.60, s, 3H	18.0,q	1.59,s, 3H <sup>a</sup>
12'	26.0, q	1.68, s, 3H	26.0,q	1.69,s, 3H <sup>b</sup>
13'	116.8, s		155.0,d	8.40,s
1"			67.4,s	
2"			149.5,d	6.83,s
3"			114.9,s	
4"			43.3,t	3.11,br.dd(16.4,6.0) 2.62,d (16.4)
5"			77.9,d	4.43,br.t (6.0)
7"			169.1,s	
8"			35.3,t	2.76,dd(14.5,7.5) 2.33,dd (14.5,7.5)
9"			119.7,d	5.12,br.t (7.5)
10"			135.7,s	
11"			18.0,q	1.61,s, 3H <sup>a</sup>
12"			26.0,q	1.67,s, 3H <sup>b</sup>
13"			116.8,s	
5'-OH		4.83, d (5.7)		4.61,d (6.0)
5"-OH				4.82,d (6.0)

Note: <sup>a</sup> Interchangeable; <sup>b</sup> Interchangable




**6** Vibralactoxime F



**7** Vibralactoxime G


**Table 4S.**  $^1\text{H}$  (600 MHz) and  $^{13}\text{C}$  NMR (150 MHz) spectral data of **6** and **7** in acetone- $d_6$ .

No.	<b>6</b>		<b>7</b>	
	$\delta_{\text{C}}$ , type	$\delta_{\text{H}}$ , multi.	$\delta_{\text{C}}$ , type	$\delta_{\text{H}}$ , multi.
1	77.0, s		77.0, s	
2	138.2, d	6.44, s	138.6, d	6.46, s
3	139.5, s		139.3, s	
4	36.5, t	3.07, br. dd (19.0, 5.6) 2.97, d (19.0)	36.5, t	3.07, ddd (18.8, 5.5, 2.2) 2.97, d (18.2)
5	79.1, d	5.06, d (5.6)	79.1, d	5.07, d (5.6)
7	171.7, s		171.6, s	
8	27.9, t	2.72, overlapped 2.58, dd (15.0, 7.5)	27.9, t	2.72, overlapped 2.58, dd (15.0, 7.3)
9	118.3, d	5.20, br. t (7.5)	118.3, d	5.20, br. t (7.5)
10	136.7, s		136.7, s	
11	18.0, q	1.66, s, 3H	18.1, q	1.67, s, 3H
12	25.9, q	1.71, s, 3H	26.1, q	1.71, s, 3H
13	153.7, d	8.38, s	154.0, d	8.42, s
1'	65.5, s		66.8, s	
2'	125.1, d	5.67, s	143.0, d	6.52, s
3'	146.1, s		137.5, s	
4'	41.6, t	2.72, overlapped 2.35, br. d (16.5)	39.9, t	3.00, ddd (17.7, 6.5, 2.3) 2.68, d (17.7)
5'	79.1, d	4.29, br. t (6.5)	78.7, d	4.43, ddd (6.5, 6.0, 2.3)
7'	170.5, s		169.1, s	
8'	35.8, t	2.65, overlapped 2.21, dd (14.3, 8.0)	35.3, t	2.75, dd (15.0, 7.2) 2.32, dd (15.0, 8.2)
9'	120.6, d	5.12, br. t (8.0)	119.8, d	5.15, br. t (7.5)
10'	134.7, s		134.6, s	
11'	18.1, q	1.58, s, 3H	18.1, q	1.59, s, 3H
12'	26.1, q	1.67, s, 3H	25.9, q	1.67, s, 3H
13'	61.7, t	4.14, dd (14.0, 5.5) 4.10, dd (14.0, 5.5)	154.3, d	8.34, s
1''			65.4, s	
2''			125.2, d	5.67, s
3''			146.0, s	
4''			41.6, t	2.72, overlapped 2.37, d (16.0)
5''			79.2, d	4.29, ddd (6.2, 6.2, 3.0)
7''			170.6, s	
8''			35.7, t	2.65, dd (14.0, 7.0) 2.22, dd (14.0, 8.0)
9''			120.7, d	5.14, br. t (7.5)
10''			135.6, s	
11''			18.0, q	1.59, s, 3H
12''			26.0, q	1.67, s, 3H
13''			61.7, t	4.12, dd (14.0, 5.5) 4.13, dd (14.0, 5.5)
5'-OH		4.24, d (6.5)		4.59, d (6.0)
5"-OH				4.24, d (6.0)
13'-OH		3.91, t (5.5)		
13"-OH				3.91, t (5.5)



**Table 5S.**  $^1\text{H}$  (600 MHz) and  $^{13}\text{C}$  NMR (150 MHz) spectral data of **8** and **13** in acetone- $d_6$ .


No.	<b>8</b>		<b>13</b>	
	$\delta_{\text{C}}$ , type	$\delta_{\text{H}}$ , multi.	$\delta_{\text{C}}$ , type	$\delta_{\text{H}}$ , multi.
1	77.0, s		67.2, s	
2	138.4, d	6.46, s	144.1, d	6.49, s
3	139.4, s		136.9, s	
4	36.5, t	3.07, br. dd (19.0, 5.6) 2.97, d (19.0)	39.7, t	2.94, overlapped 2.61, overlapped
5	79.0, d	5.07, d (5.6)	78.5, d	4.34, overlapped
7	171.7, s		172.9, s	
8	27.9, t	2.72, overlapped 2.58, overlapped	35.4, t	2.66,dd (14.0, 7.7) 2.25, dd (14.0, 7.7)
9	118.3, d	5.20, br. t (7.5)	120.0, d	5.06, br. t (7.5)
10	136.7, s		135.2, s	
11	18.1, q	1.67, s, 3H	18.1, q	1.58, s, 3H
12	26.0, q	1.71, s, 3H	26.0, q	1.66, s, 3H
13	153.8, d	8.40., s	154.6, d	8.33, s
1'	66.3, s		66.2, s	
2'	135.6, d	6.11, s	135.8, d	6.11, s
3'	139.1, s		139.0, s	
4'	40.1, t	2.90, overlapped 2.60, overlapped	40.1, t	2.89, overlapped 2.61, overlapped
5'	78.7, d	4.36,ddd (6.0, 6.0, 2.7)	78.7, d	4.35, overlapped
7'	169.6, s		169.7, s	
8'	35.6, t	2.72, overlapped 2.27, dd (14.0, 8.0)	35.6, t	2.71, dd (14.0, 7.6) 2.27, dd (14.0, 7.6)
9'	120.1, d	5.12, br. t (7.5)	120.2, d	5.13, br. t (7.6)
10'	135.3, s		135.2, s	
11'	18.0, q	1.58, s, 3H	17.9, q	1.59, s, 3H
12'	25.9, q	1.66, s, 3H	26.0, q	1.67, s, 3H
13'	147.2, d	7.95, s	147.2, d	7.96, s
N-OH		10.26, s		10.26, s
5-OH				4.36, overlapped
5'-OH		4.42, d (6.0)		4.41, d (6.0)
7-OCH <sub>3</sub>			51.8, q	3.64, s 3H



**Table 6S.**  $^1\text{H}$  (600 MHz) and  $^{13}\text{C}$  NMR (150 MHz) spectral data of **9** and **14** in acetone- $d_6$ .

No.	<b>9</b>		<b>14</b>	
	$\delta_{\text{C}}$ , type	$\delta_{\text{H}}$ , multi.	$\delta_{\text{C}}$ , type	$\delta_{\text{H}}$ , multi.
1	77.0, s		67.3, s	
2	138.6, d	6.46, s	144.3, d	6.49, s
3	139.3, s		136.8, s	
4	36.5, t	3.07, dd (19.1, 5.6) 2.97, d (19.2)	39.7, t	2.94, dd (17.0, 5.2) 2.62, d (17.0)
5	79.1, d	5.07, d (5.6)	78.4, d	4.34, overlapped
7	171.6, s		172.9, s	
8	27.9, t	2.58, overlapped 2.72, overlapped	35.4, t	2.66, overlapped 2.25, dd (14, 8.0)
9	118.3, d	5.20, t (7.5)	120.0, d	5.06, t (7.5)
10	136.7, s		135.2, s	
11	18.1, q	1.66, s, 3H	18.0, q	1.57, s, 3H
12	26.0, q	1.71, s, 3H	26.1, q	1.65, s, 3H
13	154.0, d	8.42, s	154.6, d	8.35, s
1'	66.8, s		66.7, s	
2'	143.3, d	6.54, s	143.5, d	6.54, s
3'	137.4, s		137.3, s	
4'	39.8, t	3.00, dd (16.0, 6.0) 2.68, overlapped	39.8, t,	2.99, dd (16.7, 6.0) 2.67, overlapped
5'	78.7, d	4.43, ddd (6.0, 6.0, 2.3)	78.6, d	4.42, overlapped
7'	169.1, s		169.2, s	
8'	35.3, t	2.75, overlapped 2.32, dd (14.2, 8.0)	35.3, t	2.75, dd (14.0, 7.0) 2.30, overlapped
9'	119.8, d	5.12, t (7.5)	120.2, d	5.13, t (7.5)
10'	135.6, s		135.5, s	
11'	18.0, q	1.59, s, 3H	18.1, q	1.59, s, 3H
12'	26.0, q	1.67, s, 3H	26.1, q	1.67, s, 3H
13'	154.5, d	8.35, s	154.8, d	8.36, s
1''	66.2, s		66.2, s	
2''	135.7, d	6.12, s	135.7, d	6.12, s
3''	139.0, s		139.0, s	
4''	40.1, t	2.90, overlapped 2.61, overlapped	40.1, t	2.91, overlapped 2.62, overlapped
5''	78.7, d	4.36, ddd (6.0, 6.0, 2.8)	78.7, d	4.35, overlapped

7"	169.7, s		169.7, s	
8"	35.6, t	2.72, overlapped 2.28, dd (14.2, 8.1)	35.6, t	2.73, overlapped 2.27, overlapped
9"	120.2, d	5.13, t (7.5)	119.9, d	5.13, t (7.5)
10"	135.2, s		135.2, s	
11"	18.1, q	1.59, s, 3H	18.1, q	1.59, s, 3H
12"	25.9, q	1.66, s, 3H	26.0, q	1.67, s, 3H
13"	147.2, d	7.96, s	147.1, d	7.96, s
5-OH				4.42, d (6.0)
5'-OH		4.60, d (6.0)		4.57, d (6.0)
5"-OH		4.42, d (6.0)		unobserved
N-OH		10.25, s		10.26, s
-OCH <sub>3</sub>			51.8, q	3.64, s, 3H



**10** Vibralactoxime J

**Table 7S.**  $^1\text{H}$  (600 MHz) and  $^{13}\text{C}$  NMR (150 MHz) spectral data of **10** in acetone- $d_6$ .

No.	$\delta_{\text{C}}$ , type	$\delta_{\text{H}}$ , multi.	No.	$\delta_{\text{C}}$ , type	$\delta_{\text{H}}$ , multi.
1	77.0, s		1''	66.7, s	
2	138.6, d	6.46, s	2''	143.5, d	6.54, s
3	139.3, s		3''	137.4, s	
4	36.5, t	3.07, dd (19.0, 6.0) 2.97, d (19.0)	4''	39.8, t	3.00, overlapped 2.69, dd (17.0, 10.0)
5	79.0, d	5.07, d (5.6)	5''	78.7, d	4.43, overlapped
7	171.6, s		7''	169.2, s	
8	27.9, t	2.72, overlapped 2.58, dd (15.0, 7.4)	8''	35.3, t	2.76, overlapped 2.33, dd (14.5, 7.8)
9	118.3, d	5.20, t (7.5)	9''	119.9, d	5.13, t (7.5)
10	136.7, s		10''	135.6, s	
11	18.0, q	1.59, s, 3H	11''	18.1, q	1.60, s, 3H
12	25.9, q	1.71, s, 3H	12''	26.1, q	1.66, s, 3H
13	154.0, d	8.42, s	13''	154.6, d	8.36, s
1'	66.8, s		1'''	66.2, s	
2'	143.4, d	6.54, s	2'''	135.7, d	6.12, s
3'	137.4, s		3'''	139.0, s	
4'	39.8, t	3.00, overlapped 2.69, dd (17.0, 10.0)	4'''	40.1, t	2.90, overlapped 2.61, d (16.2)
5'	78.7, d	4.43, overlapped	5'''	78.6, d	4.36, ddd (6.0, 6.0, 2.5)
7'	169.1, s		7'''	169.7, s	
8'	35.3, t	2.76, overlapped 2.33, dd (14.5, 7.8)	8'''	35.6, t	2.74, overlapped 2.28, dd (14.8, 8.2)
9'	119.8, d	5.13, t (7.5)	9'''	120.2, d	5.12, t (7.5)
10'	135.6, s		10'''	135.2, s	
11'	18.1, q	1.60, s, 3H	11'''	18.1, q	1.67, s, 3H
12'	26.1, q	26.1, s, 3H	12'''	26.1, q	1.67, s, 3H
13'	154.7, d	8.37, s	13'''	147.2, d	7.96, s
5'-OH		4.60, d (5.5)	5''-OH		4.60, d (5.5)
5'''-OH		4.43, overlapped	N-OH		10.27, s


**Table 8S.** Characteristic  $^1\text{H}$  NMR (600 MHz) spectral data of **11**, **12**, **15** and **16** in acetone- $d_6$ .

No.	<b>11</b>	<b>12</b>	<b>15</b>	<b>16</b>
H-2 of head unit	6.46, s	6.47, s	6.49, s	6.49, s
H-9 of head unit	5.20, t (7.5)	5.19, t (7.5)	5.06, t (7.5)	5.06, t (7.5)
H-13 of head unit	8.41, s	8.42, s	8.35, s	8.35, s
H-2 of tail unit	6.12, s	6.12, s	6.12, s	6.12, s
H-9 of tail unit	5.12, t (7.5)	5.12, t (7.5)	5.13, t (7.5)	5.12, t (7.5)
H-13 of tail unit	7.96, s	7.96, s	7.96, s	7.96, s
H-2 of the other unit	6.54, each s, totally 3H, overlapped	6.56, each s, totally 4H, overlapped	6.55, each s, totally 2H, overlapped	6.55, each s, totally 3H, overlapped
H-9 of the other unit	5.13, t (7.5), totally 3H, overlapped	5.13, t (7.5), totally 4H, overlapped	5.14, t (7.5), totally 2H, overlapped	5.13, t (7.5), totally 3H, overlapped
H-13 of the other unit	8.36, each s, totally 3H, overlapped	8.36, s 8.37, s	8.36, s 8.37, s	8.36, s 8.37, s
		8.38, each s, totally 2H, overlapped		8.38, s
N-OH	10.23, s	10.28, s	10.23, s	10.27, s
OCH <sub>3</sub>			3.85, s, 3H	3.84, s, 3H


#### 4. Supplementary NMR spectra.

Figure 1S:  $^1\text{H}$  NMR spectrum of **1** (acetone- $d_6$ , 600MHz)


1cbb3a1  
h, acetone av600




**Figure 2S:**  $^{13}\text{C}$  NMR and DEPT spectra of **1** (acetone- $d_6$ , 150MHz)




**Figure 3S:** HSQC spectrum of **1**.




**Figure 4S:**  $^1\text{H}$ - $^1\text{H}$  COSY spectrum of **1**.



**Figure 5S:** HMBC spectrum of **1**.



**Figure 6S:** ROESY spectrum of **1**.



**Figure 7S:** HREIMS spectrum of **1**.

**Elemental Composition Report**

**Page 1**

**Single Mass Analysis**

Tolerance = 10.0 PPM / DBE: min = -10.0, max = 120.0

Selected filters: None

Monoisotopic Mass, Odd and Even Electron Ions

12 formula(e) evaluated with 1 results within limits (up to 51 closest results for each mass)

Elements Used:

C: 0-200 H: 0-400 N: 1-1 O: 2-4

lcbvb-3a1

09:38:12 20-Nov-2013

Voltage EI+

100

%

0

220.900 220.950 221.000 221.050 221.100 221.150 221.200 221.250 221.300 m/z

KIB  
M131120EA-03AFAMM 15 (1.378)  
221.1049

Autospec Premier  
P776  
10.1


Minimum: 200.0 10.0 -10.0  
Maximum: 220.0 120.0 10.0

Mass	Calc. Mass	mDa	PPM	DBE	i-FIT	Formula
221.1049	221.1052	-0.3	-1.4	6.0	5546025.5	C12 H15 N O3


**Figure 8S:** Optical rotation of compound 1.

Optical rotation measurement										
			Model : P-1020 (A060460638)							
No.	Sample	Mode	Data	Monitor Blank	Temp. Cell	Date Comment Sample Name	Light Filter	Cycle Time	Integ Time	Operator
No.1	2 (1/3)	Sp.Rot	-192.6320	-0.1830 0.0000	22.9 50.00	Thu Mar 12 14:38:40 2015 0.00190g/mL Acetonitrile Cell LCBVB-3A1	Na 589nm	2 sec	10 sec	
No.2	2 (2/3)	Sp.Rot	-193.3680	-0.1837 0.0000	22.8 50.00	Thu Mar 12 14:38:53 2015 0.00190g/mL Acetonitrile Cell LCBVB-3A1	Na 589nm	2 sec	10 sec	-192.580
No.3	2 (3/3)	Sp.Rot	-190.8420	-0.1813 0.0000	22.8 50.00	Thu Mar 12 14:39:07 2015 0.00190g/mL Acetonitrile Cell LCBVB-3A1	Na 589nm	2 sec	10 sec	


**Figure 9S:**  $^1\text{H}$  NMR spectrum of **2** (acetone- $d_6$ , 600MHz)




**Figure 10S:**  $^{13}\text{C}$  NMR and DEPT spectra of **2** (acetone- $d_6$ , 150MHz).




**Figure 11S:** HSQC spectrum of **2**.




**Figure 12S:**  $^1\text{H}$ - $^1\text{H}$  COSY spectrum of **2**.



**Figure 13S:** HMBC spectrum of **2**.



**Figure 14S:** ROESY spectrum of **2**.



**Figure 15S:** HREIMS spectrum of **2**.

**Elemental Composition Report**

**Page 1**

**Single Mass Analysis**

Tolerance = 10.0 PPM / DBE: min = -10.0, max = 120.0

Selected filters: None

Monoisotopic Mass, Odd and Even Electron Ions

13 formula(e) evaluated with 1 results within limits (up to 51 closest results for each mass)

Elements Used:

C: 0-200 H: 0-400 N: 1-1 O: 3-5

lcbvb-3a2  
14:49:25 31-Dec-2013

Voltage El+

100

%

0

KIB  
M131231EA-02AFAMM 16 (1.470)  
253.1312


Autospec Premier  
P776  
1.78

252.900 252.950 253.000 253.050 253.100 253.150 253.200 253.250 253.300 253.350 m/z


Minimum: -10.0  
Maximum: 200.0 10.0 120.0

Mass	Calc. Mass	mDa	PPM	DBE	i-FIT	Formula
253.1312	253.1314	-0.2	-0.8	5.0	5546025.5	C13 H19 N O4


**Figure 16S:**  $^1\text{H}$  NMR spectrum of **3** (acetone- $d_6$ , 600MHz).




**Figure 17S:** Enlarged view of  $^1\text{H}$  NMR spectra of **3** (acetone- $d_6$ , 600MHz)




**Figure 18S:**  $^{13}\text{C}$  NMR and DEPT spectra of **3** (acetone- $d_6$ , 150MHz)



**Figure 19S:** Enlarged view of  $^{13}\text{C}$  NMR spectrum of **3** (acetone- $d_6$ , 150MHz)



**Figure 20S:** HSQC spectrum of **3**.



**Figure 21S:**  $^1\text{H}$ - $^1\text{H}$  COSY spectrum of **3**.








Figure 22S: HMBC spectrum of 3.



**Figure 23S:** Enlarged view of HMBC spectrum of **3**.



**Figure 24S:** ROESY spectrum of **3**.




**Figure 25S:** HRESIMS spectrum of 3.

## Qualitative Analysis Report

Data Filename	Icox-27.d	Sample Name	Icox-27
Sample Type	Sample	Position	P1-B3
Instrument Name	Instrument 1	User Name	
Acq Method	SIBU.m	Acquired Time	3/10/2015 3:53:53 PM
IRM Calibration Status	Success	DA Method	Default.m
Comment			
Sample Group	Info.		
Acquisition SW	6200 series TOF/6500 series		
Version	Q-TOF B.05.01 (B5125.2)		

### User Spectra



### Peak List

m/z	z	Abund	Formula	Ion
274.2741	1	32309.41		
318.3002	1	19627.41		
510.2102	1	138483.17	C26 H33 N O8	$(M+Na)^+$
511.2137	1	38055.41	C26 H33 N O8	$(M+Na)^+$
526.1843	1	230987.25		
527.1877	1	65481.15		
528.1861	1	24762.72		
542.2361	1	20034.74		
558.2104	1	43224.76		
997.4312	1	30187.51		

### Formula Calculator Element Limits


Element	Min	Max
C	3	60
H	0	120
O	0	30
N	0	10

### Formula Calculator Results


Formula	CalculatedMass	CalculatedMz	Mz	Diff. (mDa)	Diff. (ppm)	DBE
C26 H33 N O8	487.2206	510.2098	510.2102	-0.4	-0.8	11.0000

--- End Of Report ---


**Figure 26S:**  $^1\text{H}$  NMR spectrum of **4** (acetone- $d_6$ , 600MHz).




**Figure 27S:**  $^{13}\text{C}$  NMR and DEPT spectra of **4** (acetone- $d_6$ , 150MHz).




**Figure 28S:** Enlarged view of  $^{13}\text{C}$  NMR spectrum of **4** (acetone- $d_6$ , 150MHz).



**Figure 29S:** HSQC spectrum of **4**.



**Figure 30S:**  $^1\text{H}$ - $^1\text{H}$  COSY spectrum of **4**.



**Figure 31S:** HMBC spectrum of **4**.






Figure 32S: Enlarged view of HMBC spectra of **4**.



**Figure 33S:** ROESY spectrum of **4**.



**Figure 34S:** HREIMS spectrum of **4**.

**Elemental Composition Report**

**Page 1**

**Single Mass Analysis**

Tolerance = 10.0 PPM / DBE: min = -10.0, max = 120.0

Selected filters: None

Monoisotopic Mass, Odd and Even Electron Ions

18 formula(e) evaluated with 1 results within limits (up to 51 closest results for each mass)

Elements Used:

C: 0-200 H: 0-400 N: 2-2 O: 4-6

Icbvb-3a3

14:55:39 31-Dec-2013

Voltage El+

100

%

0

423.800 423.900 424.000 424.100 424.200 424.300 424.400 424.500 424.600 m/z


KIB  
M131231EA-03AFAMM 14 (1.285)  
424.1987

Autospec Premier  
P776  
3.33

Minimum: 200.0 10.0 -10.0  
Maximum: 400.0 40.0 120.0

Mass	Calc. Mass	mDa	PPM	DBE	i-FIT	Formula
424.1987	424.1998	-1.1	-2.6	12.0	5546025.5	C24 H28 N2 O5


Figure 35S: IR spectrum of 4.




Sample : lcbvb-3a3	Frequency Range : 399.246 - 3996.32	Measured on : 17/03/2015	
Technique : KBr压片	Resolution : 4	Instrument : Tensor27	Sample Scans : 16
Customer : 150317IR6	Zerofilling : 2	Acquisition : Double Sided,For	

**Figure 36S:**  $^1\text{H}$  NMR spectrum of **5** (acetone- $d_6$ , 600MHz).


lcox68, 45, 1, 1r  
h Acetone av600




**Figure 37S:**  $^{13}\text{C}$  NMR and DEPT spectra of **5** (acetone- $d_6$ , 150MHz).




**Figure 38S:** Enlarged view of  $^{13}\text{C}$  NMR spectrum of **5** (acetone- $d_6$ , 150MHz).




**Figure 39S:** HSQC spectrum of **5**.




**Figure 40S:**  $^1\text{H}$ - $^1\text{H}$  COSY spectrum of **5**.



**Figure 41S:** HMBC spectrum of **5**.



**Figure 42S:** Enlarged view of HMBC spectra of **5**.



**Figure 43S:** ROESY spectrum of **5**.





Figure 44S: IR spectrum of 5.




**Figure 45S:** HRESIMS spectrum of **5**.

## Qualitative Analysis Report

Data Filename	150209ESIA10.d	Sample Name	Icox-68
Sample Type	Sample	Position	
Instrument Name	Agilent G6230 TOF MS	User Name	KIB
Acq Method	ESI.m	Acquired Time	2/9/2015 4:40:06 PM
IRM Calibration Status	Success	DA Method	ESI.m
Comment			
Sample Group	Info.		
Acquisition SW	6200 series TOF/6500 series		
Version	Q-TOF B.05.01 (B5125.2)		

### User Spectra



### Peak List

m/z	z	Abund
121.0509	1	32167.5
230.2472	1	29344.04
274.2741	1	203815.23
275.2771	1	31088.2
318.3006	1	312160.5
319.3034	1	53387.34
352.1542	1	55843.86
362.3267	1	121212.16
406.3523	1	30843.82
922.0098	1	29313.56

### Formula Calculator Element Limits


Element	Min	Max
C	0	200
H	0	400
O	5	10
Na	1	1
N	3	3

### Formula Calculator Results


Formula	CalculatedMass	CalculatedMz	Mz	Diff. (mDa)	Diff. (ppm)	DBE
C36 H43 N3 Na O8	668.2948	668.2942	668.2938	0.4	0.5	16.5000

--- End Of Report ---


**Figure 46S:**  $^1\text{H}$  NMR spectrum of **6** (acetone- $d_6$ , 600MHz).




**Figure 47S:**  $^{13}\text{C}$  NMR and DEPT spectra of **6** (acetone- $d_6$ , 150MHz).




**Figure 48S:** HSQC spectrum of **6**.




**Figure 49S:**  $^1\text{H}$ - $^1\text{H}$  COSY spectrum of **6**.




**Figure 50S:** HMBC spectrum of **6**.



**Figure 51S:** Enlarged view of HMBC spectra of **6**.



**Figure 52S:** ROESY spectrum of **6**.




**Figure 53S:** HRESIMS spectrum of **6**.

## Qualitative Analysis Report

<b>Data Filename</b>	Icox-74.d	<b>Sample Name</b>	Icox-74
<b>Sample Type</b>	Sample	<b>Position</b>	P1-B7
<b>Instrument Name</b>	Instrument 1	<b>User Name</b>	
<b>Acq Method</b>	SIBU.m	<b>Acquired Time</b>	3/9/2015 10:34:24 AM
<b>IRM Calibration Status</b>	Success	<b>DA Method</b>	Default.m
<b>Comment</b>			
<b>Sample Group</b>		<b>Info.</b>	
<b>Acquisition SW Version</b>	6200 series TOF/6500 series Q-TOF B.05.01 (B5125.2)		

### User Spectra



### Peak List

m/z	z	Abund	Formula	Ion
102.1279	1	93262.92		
103.1311	1	6175.42		
274.2742	1	15380.4		
301.1413	1	6937.4		
318.3003	1	21869.02		
340.2824	1	10073.84		
430.2231	1	11604.17		
452.2046	1	36882.11	C24 H31 N O6	$(M+Na)^+$
453.2076	1	10524.28	C24 H31 N O6	$(M+Na)^+$
468.1786	1	100357.7		
469.1818	1	28506.38		
470.1801	1	10598.65		
500.2042	1	18659.6		
881.4201	1	9003		
897.395	1	10103.26		

### Formula Calculator Element Limits

Element	Min	Max
C	3	60
H	0	120
O	0	30
N	0	5

### Formula Calculator Results

Formula	CalculatedMass	CalculatedMz	Mz	Diff. (mDa)	Diff. (ppm)	DBE
C24 H31 N O6	429.2151	452.2044	452.2046	-0.2	-0.4	10.0000



Agilent Technologies

Page 1 of 2

Printed at: 8:19 AM on: 3/10/2015

## **Qualitative Analysis Report**

---

--- End Of Report ---




**Agilent Technologies**


Page 2 of 2

Printed at: 8:19 AM on: 3/10/2015


**Figure 54S:**  $^1\text{H}$  NMR spectrum of **7** (acetone- $d_6$ , 600MHz).




**Figure 55S:**  $^{13}\text{C}$  NMR and DEPT spectra of **7** (acetone- $d_6$ , 150MHz).




**Figure 56S:** HSQC spectrum of 7.




**Figure 57S:**  $^1\text{H}$ - $^1\text{H}$  COSY spectrum of **7**.




**Figure 58S:** Enlarged view of  $^1\text{H}$ - $^1\text{H}$  COSY spectrum.



**Figure 59S:** HMBC spectrum of **7**.



**Figure 60S:** ROESY spectrum of **7**.




**Figure 61S:** HRESIMS spectrum of 7.

## Qualitative Analysis Report

Data Filename	Icox-67.d	Sample Name	Icox-67
Sample Type	Sample	Position	P1-B5
Instrument Name	Instrument 1	User Name	
Acq Method	SIBU.m	Acquired Time	3/9/2015 10:29:39 AM
IRM Calibration Status	Success	DA Method	Default.m
Comment			
Sample Group	Info.		
Acquisition SW	6200 series TOF/6500 series		
Version	Q-TOF B.05.01 (B5125.2)		

### User Spectra



### Peak List

m/z	z	Abund	Formula	Ion
274.274	1	7256.5		
275.2774	1	1445.26		
296.2561		2391.28		
318.2999	1	6337.43		
319.3036	1	1533.73		
340.2821		2261.19		
384.3075		1078.58		
411.0943		1458.01		
427.0686	1	1445.7		
452.205		1404.92		
468.1782	1	5010.77		
469.1817	1	1631.39		
673.3104	1	1046.02	C36 H46 N2 O9	(M+Na)+
689.2841	1	2269.47		
690.2864	1	1062.86		

### Formula Calculator Element Limits


Element	Min	Max
C	3	60
H	0	120
O	0	30
N	0	5

### Formula Calculator Results


Formula	CalculatedMass	CalculatedMz	Mz	Diff. (mDa)	Diff. (ppm)	DBE
C36 H46 N2 O9	650.3203	673.3096	673.3104	-1.6	-2.4	15.0000

--- End Of Report ---


**Figure 62S:**  $^1\text{H}$  NMR spectrum of **8** (acetone- $d_6$ , 600MHz).



**Figure 63S:**  $^{13}\text{C}$  NMR and DEPT spectra of **8** (acetone- $d_6$ , 150MHz).



**Figure 64S:** HSQC spectrum of **8**.



**Figure 65S:**  $^1\text{H}$ - $^1\text{H}$  COSY spectrum of **8**.




Figure 66S: HMBC spectrum of **8**.





Figure 67S: ROESY spectrum of **8**.



**Figure 68S:** HREIMS spectrum of **8**.

**Elemental Composition Report**

**Page 1**

**Single Mass Analysis**

Tolerance = 10.0 PPM / DBE: min = -10.0, max = 120.0

Selected filters: None

Monoisotopic Mass, Odd and Even Electron Ions

18 formula(e) evaluated with 1 results within limits (up to 51 closest results for each mass)

Elements Used:

C: 0-200 H: 0-400 N: 2-2 O: 5-7

Icbva-5  
15:17:29 31-Dec-2013

Voltage El<sup>+</sup>

100

%

0

441.800      441.900      442.000      442.100      442.200      442.300      442.400      442.500      442.600      m/z


KIB  
M131231EA-05AFAMM 27 (2.479)  
442.2098

Autospec Premier  
P776  
31.3


Minimum: -10.0  
Maximum: 200.0 10.0 120.0

Mass	Calc. Mass	mDa	PPM	DBE	i-FIT	Formula
442.2098	442.2104	-0.6	-1.4	11.0	5546034.0	C24 H30 N2 O6


**Figure 69S:**  $^1\text{H}$  NMR spectrum of **9** (acetone- $d_6$ , 600MHz).




**Figure 70S:**  $^{13}\text{C}$  NMR and DEPT spectra of **9** (acetone- $d_6$ , 150MHz).



**Figure 71S:** HSQC spectrum of **9**.



**Figure 72S:**  $^1\text{H}$ - $^1\text{H}$  COSY spectrum of **9**.



**Figure 73S:** HMBC spectrum of **9**.





Figure 74S: ROESY spectrum of **9**.



**Figure 75S:** HREIMS spectrum of **9**.

**Elemental Composition Report**

**Page 1**

**Single Mass Analysis**

Tolerance = 10.0 PPM / DBE: min = -10.0, max = 120.0

Selected filters: None

Monoisotopic Mass, Odd and Even Electron Ions

25 formula(e) evaluated with 1 results within limits (up to 51 closest results for each mass)

Elements Used:

C: 0-200 H: 0-400 N: 3-3 O: 8-10

Icbva-6  
15:43:23 31-Dec-2013

Voltage El+

100

%

0

KIB  
M131231EA-06AFAMM 82 (7.527)  
663.3160


Autospec Premier  
P776  
1

662.70 662.80 662.90 663.00 663.10 663.20 663.30 663.40 663.50 663.60 663.70 663.80 663.90 m/z


Minimum: 200.0 10.0 -10.0  
Maximum: 663.3160 0.4 0.6 120.0

Mass	Calc. Mass	mDa	PPM	DBE	i-FIT	Formula
663.3160	663.3156	0.4	0.6	16.0	5546025.5	C <sub>36</sub> H <sub>45</sub> N <sub>3</sub> O <sub>9</sub>


**Figure 76S:**  $^1\text{H}$  NMR spectrum of **10** (acetone- $d_6$ , 600MHz).



**Figure 77S:** Enlarged view of  $^1\text{H}$  NMR spectra of **10** (acetone- $d_6$ , 600MHz)



**Figure 78S:**  $^{13}\text{C}$  NMR and DEPT spectra of **10** (acetone- $d_6$ , 150MHz).



**Figure 79S:** Enlarged view of  $^{13}\text{C}$  NMR and DEPT spectra of **10** (acetone- $d_6$ , 150MHz).





Figure 80S: HSQC spectrum of **10**.



**Figure 81S:** Enlarged view of HSQC spectrum of **10**.




Figure 82S:  $^1\text{H}$ - $^1\text{H}$  COSY spectrum of **10**.







Figure 83S: HMBC spectrum of **10**.



**Figure 84S:** Enlarged view of HMBC spectrum of **10**.



**Figure 85S:** ROESY spectrum of **10**.




**Figure 86S:** HRESIMS spectrum of **10**.

## Qualitative Analysis Report

Data Filename	Icox-72.d	Sample Name	Icox-72
Sample Type	Sample	Position	P1-B6
Instrument Name	Instrument 1	User Name	
Acq Method	SIBU.m	Acquired Time	3/10/2015 8:42:51 AM
IRM Calibration Status	Success	DA Method	Default.m
Comment			
Sample Group		Info.	
Acquisition SW	6200 series TOF/6500 series		
Version	Q-TOF B.05.01 (B5125.2)		

### User Spectra



### Peak List

m/z	z	Abund
102.1281	1	335324.28
103.1312	1	23138.34
163.1231	1	32261.13
274.2743	1	93798.95
318.3006	1	94393.51
319.3035	1	18947.31
362.3267	1	19927.98
432.2388	1	18251.21
437.1937	1	37663.06
453.1679	1	57979.79

### Formula Calculator Element Limits


Element	Min	Max
C	3	60
H	0	120
O	0	30
N	0	10

### Formula Calculator Results


Formula	CalculatedMass	CalculatedMz	Mz	Diff. (mDa)	Diff. (ppm)	DBE
C48 H60 N4 O12	884.4208	907.4100	907.4100	0.8	0.9	21.0000

--- End Of Report ---

**Figure 87S:**  $^1\text{H}$  NMR spectrum of **11** (acetone- $d_6$ , 600MHz).



**Figure 88S:** Enlarged view of  $^1\text{H}$  NMR spectrum of **11** (acetone- $d_6$ , 600MHz)




**Figure 89S:** HRESIMS spectrum of **11**.

## Qualitative Analysis Report

<b>Data Filename</b>	Icox-48.d	<b>Sample Name</b>	Icox-48
<b>Sample Type</b>	Sample	<b>Position</b>	P1-B1
<b>Instrument Name</b>	Instrument 1	<b>User Name</b>	
<b>Acq Method</b>	SIBU.m	<b>Acquired Time</b>	3/10/2015 3:49:12 PM
<b>IRM Calibration Status</b>	Success	<b>DA Method</b>	Default.m
<b>Comment</b>			
<b>Sample Group</b>	<b>Info.</b>		
<b>Acquisition SW</b>	6200 series TOF/6500 series		
<b>Version</b>	Q-TOF B.05.01 (B5125.2)		

### User Spectra



### Peak List

m/z	z	Abund	Formula	Ion
102.1277	1	101863.34		
103.1313	1	8485.26		
144.0805		10690.44		
187.1232		9680.63		
230.248		8409.56		
274.274	1	29855.7		
275.2776	1	6881.42		
318.3003	1	36899.34		
362.3268		10662.17		
1128.5158	1	5871.05	C <sub>60</sub> H <sub>75</sub> N <sub>5</sub> O <sub>15</sub>	(M+Na)+

### Formula Calculator Element Limits


Element	Min	Max
C	3	60
H	0	120
O	0	30
N	0	10

### Formula Calculator Results


Formula	CalculatedMass	CalculatedMz	Mz	Diff. (mDa)	Diff. (ppm)	DBE
C <sub>60</sub> H <sub>75</sub> N <sub>5</sub> O <sub>15</sub>	1105.5260	1128.5152	1128.5158	-0.2	-0.2	26.0000

--- End Of Report ---

**Figure 90S:**  $^1\text{H}$  NMR spectrum of **12** (acetone- $d_6$ , 600MHz).



**Figure 91S:** Enlarged view of  $^1\text{H}$  NMR spectrum of **12** (acetone- $d_6$ , 600MHz)




**Figure 92S:** HRESIMS spectrum of **12**.

## Qualitative Analysis Report

<b>Data Filename</b>	Icox-25.d	<b>Sample Name</b>	Icox-25
<b>Sample Type</b>	Sample	<b>Position</b>	P1-A6
<b>Instrument Name</b>	Instrument 1	<b>User Name</b>	
<b>Acq Method</b>	SIBU.m	<b>Acquired Time</b>	3/23/2015 4:30:40 PM
<b>IRM Calibration Status</b>	Success	<b>DA Method</b>	Default.m
<b>Comment</b>			
<b>Sample Group</b>	<b>Info.</b>		
<b>Acquisition SW</b>	6200 series TOF/6500 series		
<b>Version</b>	Q-TOF B.05.01 (B5125.2)		

### User Spectra



### Peak List

m/z	z	Abund	Formula	Ion
675.3102	2	13176.74		
675.8108	2	13227.66		
683.2972	2	39890.23		
683.7988	2	30519.08		
1327.6368	1	18248.74	C72 H90 N6 O18	(M+H)+
1328.6397	1	15092.47	C72 H90 N6 O18	(M+H)+

### Formula Calculator Element Limits


Element	Min	Max
C	3	100
H	0	200
O	0	30
N	0	10

### Formula Calculator Results


Formula	CalculatedMass	CalculatedMz	Mz	Diff. (mDa)	Diff. (ppm)	DBE
C72 H90 N6 O18	1326.6312	1327.6384	1327.6368	1.7	1.3	31.0000

--- End Of Report ---


**Figure 93S:**  $^1\text{H}$  NMR spectrum of **13** (acetone- $d_6$ , 600MHz).



**Figure 94S:**  $^{13}\text{C}$  NMR and DEPT spectra of **13** (acetone- $d_6$ , 150MHz).



**Figure 95S:** Enlarged view of  $^{13}\text{C}$  NMR and DEPT spectra of **13** (acetone- $d_6$ , 150MHz).



**Figure 96S:** HSQC spectrum of **13**.




Figure 97S:  $^1\text{H}$ - $^1\text{H}$  COSY spectrum of **13**.




Figure 98S: HMBC spectrum of **13**.





Figure 99S: ROESY spectrum of **13**.



**Figure 100S:** HREIMS spectrum of **13**.

**Elemental Composition Report**

**Page 1**

**Single Mass Analysis**

Tolerance = 10.0 PPM / DBE: min = -10.0, max = 120.0

Selected filters: None

Monoisotopic Mass, Odd and Even Electron Ions

19 formula(e) evaluated with 1 results within limits (up to 51 closest results for each mass)

Elements Used:

C: 0-200 H: 0-400 N: 2-2 O: 6-8

Icbla-4

15:07:59 31-Dec-2013

Voltage El+

100

90

80

70

60

50

40

30

20

10

0

%

0

100

90

80

70

60

50

40

30

20

10

0

%

0

100

90

80

70

60

50

40

30

20

10

0

%

0

100

90

80

70

60

50

40

30

20

10

0

%

KIB  
M131231EA-04AFAMM 43 (3.947)  
474.2369


Autospec Premier  
P776  
632

473.800 473.900 474.000 474.100 474.200 474.300 474.400 474.500 474.600 474.700 m/z

Minimum: -10.0  
Maximum: 200.0 10.0 120.0

Mass	Calc. Mass	mDa	PPM	DBE	i-FIT	Formula
474.2369	474.2366	0.3	0.6	10.0	5546332.0	C25 H34 N2 O7

**Figure 101S:**  $^1\text{H}$  NMR spectrum of **14** (acetone- $d_6$ , 600MHz).



**Figure 102S:**  $^{13}\text{C}$  NMR and DEPT spectra of **14** (acetone- $d_6$ , 150MHz).







Figure 103S: HSQC spectrum of **14**.




**Figure 104S:**  $^1\text{H}$ - $^1\text{H}$  COSY spectrum of **14**.




**Figure 105S:** HMBC spectrum of **14**.



**Figure 106S:** Enlarged view of HMBC spectrum of **14**.



**Figure 107S:** ROESY spectrum of **14**.




**Figure 108S:** HRESIMS spectrum of 14.

## Qualitative Analysis Report

Data Filename	Icbva-72630.d	Sample Name	Icbva-72630
Sample Type	Sample	Position	P1-B8
Instrument Name	Instrument 1	User Name	
Acq Method	SIBU.m	Acquired Time	3/10/2015 8:57:06 AM
IRM Calibration Status	Success	DA Method	Default.m
Comment			
Sample Group		Info.	
Acquisition SW	6200 series TOF/6500 series		
Version	Q-TOF B.05.01 (B5125.2)		

### User Spectra



### Peak List

m/z	z	Abund	Formula	Ion
102.128	1	37223.71		
163.1235	1	9297.78		
274.2742	1	17934.31		
318.3004	1	19673.47		
319.3035	1	4342.35		
340.2822	1	5525.16		
513.2003	1	12429.02		
718.3311	1	5901.76	C37 H49 N3 O10	(M+Na)+
734.3056	1	14048.74		
735.3088	1	5943.43		

### Formula Calculator Element Limits


Element	Min	Max
C	3	60
H	0	120
O	0	30
N	0	10

### Formula Calculator Results


Formula	CalculatedMass	CalculatedMz	Mz	Diff. (mDa)	Diff. (ppm)	DBE
C37 H49 N3 O10	695.3418	718.3310	718.3311	-1.1	-1.6	15.0000

--- End Of Report ---

**Figure 109S:**  $^1\text{H}$  NMR spectrum of **15** (acetone- $d_6$ , 600MHz).



**Figure 110S:** Enlarged view of  $^1\text{H}$  NMR spectra of **15** (acetone- $d_6$ , 600MHz).




**Figure 111S: HRESIMS spectrum of 15.**

## Qualitative Analysis Report

<b>Data Filename</b>	Icox-47.d	<b>Sample Name</b>	Icox-47
<b>Sample Type</b>	Sample	<b>Position</b>	P1-B2
<b>Instrument Name</b>	Instrument 1	<b>User Name</b>	
<b>Acq Method</b>	SIBU.m	<b>Acquired Time</b>	3/10/2015 3:51:32 PM
<b>IRM Calibration Status</b>	Success	<b>DA Method</b>	Default.m
<b>Comment</b>			
<b>Sample Group</b>	<b>Info.</b>		
<b>Acquisition SW</b>	6200 series TOF/6500 series		
<b>Version</b>	Q-TOF B.05.01 (B5125.2)		

### User Spectra



### Peak List

m/z	z	Abund	Formula	Ion
274.2745	1	19914.81		
318.3001	1	22309.2		
479.2151	2	21511.53		
495.1891	1	37214		
513.1998	1	27252.69		
939.436	1	32073.08	C49 H64 N4 O13	(M+Na)+
940.438	1	16767.01	C49 H64 N4 O13	(M+Na)+
955.4105	1	79761.27		
956.4138	1	45534.26		
957.4137	1	19558.31		

### Formula Calculator Element Limits


Element	Min	Max
C	3	60
H	0	120
O	0	30
N	0	10

### Formula Calculator Results


Formula	CalculatedMass	CalculatedMz	Mz	Diff. (mDa)	Diff. (ppm)	DBE
C49 H64 N4 O13	916.4470	939.4362	939.4360	0.7	0.8	20.0000

--- End Of Report ---

**Figure 112S:**  $^1\text{H}$  NMR spectrum of **16** (acetone- $d_6$ , 600MHz).



**Figure 113S:** Enlarged view of  $^1\text{H}$  NMR spectrum of **16** (acetone- $d_6$ , 600MHz).




**Figure 114S:** HRESIMS spectrum of **16**.

## Qualitative Analysis Report

Data Filename	Icox-34.d	Sample Name	Icox-34
Sample Type	Sample	Position	P1-A5
Instrument Name	Instrument 1	User Name	
Acq Method	SIBU.m	Acquired Time	3/23/2015 4:13:22 PM
IRM Calibration Status	Success	DA Method	Default.m
Comment			
Sample Group		Info.	
Acquisition SW	6200 series TOF/6500 series		
Version	Q-TOF B.05.01 (B5125.2)		

### User Spectra



### Peak List

m/z	z	Abund	Formula	Ion
274.274	1	15805.6		
318.3006	1	15193.98		
588.7575	2	31665.42		
589.2593	2	26469.95		
1138.5593	1	15225.89	C61 H79 N5 O16	(M+H)+
1160.541	1	15176.26		

### Formula Calculator Element Limits

Element	Min	Max
C	3	100
H	0	200
O	0	30
N	0	10

### Formula Calculator Results

Formula	CalculatedMass	CalculatedMz	Mz	Diff. (mDa)	Diff. (ppm)	DBE
C61 H79 N5 O16	1137.5522	1138.5595	1138.5593	0.2	0.2	25.0000

--- End Of Report ---