



Supplementary Figure 1. ROC curve analysis for RT-QulC results



Arrow indicates a cutoff of 2.0

- AUC = 0.98 (0.94 – 1.02)
- P value <0.0001

Supplementary Figure 2. ROC curve analysis for EP-QuIC



Arrow indicates a cutoff of 4.0

- AUC = 0.98 (0.95 – 1.01)
- P value <0.0001

Supplemental Table 1 Clinical Classification of Samples.

		Diagnosis	Number of cases
Breakdown of diagnoses:	CJD+	+ CJD	39
		+ CJD w/Alzheimers	4
		+ fCJD	2
	CJD-	Alzheimer's Disease	3
		Hydrocephalus	2
		Delirium	2
		Alzheimer's Disease / Dementia	1
		Probable Alzheimer's Disease	1
		Not CJD	1
		Dementia NOS	1
		Encephalitis	1
		Sub clinical Epilepsy	1
		Vascular Dementia	1
		Whipples Disease	1
		Transverse myelitis	1
		Status epilepticus with Hypothyroidism	1
		Lymphoma	1
		Probable Mitochondrial Disease	1
		MS	1
		Fronto-temporal dementia	1
		Myoclonus NYD	1
		Bipolar Disorder	1
		Septicemia	1
		Atypical brain stem stroke	1
		Psychosis	1
		aspiration pneumonia	1
		Hypoxic encephalopathy	1
		NYD	1
		Non Specific Dementia	1
		Chronic Hypoxia	1
		Coma rigidity, NYD	1
		Delusional-Psychiatric Disease NYD	1
		CVA	1
		HIV	1
		Parkinsons. Plus Syndrome	1
		Grade III Astrocytoma	1
		Not available	1
		Motor Neuro Disease	1
		Progressive cognitive decline	1
		Acute Glysomatosis	1
		Degenerative Brain Disorder, Unknown eitology	1
		Motorneuropathy	1
		Alcoholism	1
		Thiamine deficiency	1
		Paraneoplastic Syndrome	1

Supplementary Table 2. Clinical Diagnoses of False positive, False negative and Indeterminate samples

	Diagnosis	# False Results		# Indeterminate	
		EP-QuIC	RT-QuIC	EP-QuIC	RT-QuIC
CJD+	(A) + CJD		1	1 <sup>a</sup>	1 <sup>a</sup>
CJD-	Alzheimer's Disease / Dementia			1	
	Probable Alzheimer's Disease			1	
	Status epilepticus with Hypothyroidism	1 <sup>b</sup>			1 <sup>b</sup>
	Acute Glycomatosis	1			

<sup>a</sup> same CSF sample

<sup>b</sup> same CSF sample

## Supplementary Table 3. Inter-observer agreement analysis (Kappa test)\*

A

		Test #1 Result		
		Positive	Negative	
Test #2 Result	Positive	<b>a</b>	b	$m_1$
	Negative	c	<b>d</b>	$m_0$
		$n_1$	$n_0$	$n$

(a) and (d) are the number of times that the two tests gave the same result. (b) and (c) are the number of times that the two tests did not agree.

$$\text{Kappa} = \frac{(p_o - p_e)}{(1 - p_e)}$$

Where the observed agreement  $p_o = (a+d)/n$

and the expected agreement  $p_e = [(n_1/n)*(m_1/n)] + [(n_0/n)*(m_0/n)]$

B

		RT-QuIC		
		Positive	Negative	
EP-QuIC	Positive	<b>43</b>	2	45
	Negative	0	<b>42</b>	42
		43	44	87

$$K = (0.98 - 0.50) / (1 - 0.50) = 0.95$$

\* Indeterminate samples (two in RT-QuIC and three in EP-QuIC) were excluded from the analysis. A, Kappa calculation data layout; B, Kappa test for RT-QuIC and EP-QuIC reactions on 87 samples.