

# Human Organotypic Cultured Cardiac Slices: New Platform For High Throughput Preclinical Human Trials


Kang C, Qiao Y, Li G, Baechle K, Camelliti P, Rentschler S, Efimov IR

## Supplementary Table and Figures

Preparation	Pro	Con
Isolated Primary Cells	<ul style="list-style-type: none"><li>- Only model for ion channels currents</li></ul>	<ul style="list-style-type: none"><li>- Absence of cell to cell coupling</li><li>- Chunk isolation process alters electrophysiology</li></ul>
iPSC- Derived Cardiomyocytes (Monolayers and Tissue Constructs)	<ul style="list-style-type: none"><li>- Highest throughput for therapy screening</li><li>- Genetic manipulation</li><li>- Patient specific sampling</li><li>- Unlimited tissue supply</li></ul>	<ul style="list-style-type: none"><li>- Unable to fully replicate adult tissue characteristics currently</li><li>- Cell alignment and cell coupling is disorganized</li></ul>
Coronary Perfused Intact Tissue	<ul style="list-style-type: none"><li>- Large Scale electrophysiology</li><li>- Tissue conduction properties</li><li>- Full scaled arrhythmia and defibrillation study</li></ul>	<ul style="list-style-type: none"><li>- Requires Intact coronary system</li><li>- Few preparations per heart</li><li>- Unable to be cultured</li></ul>
Tissue Slices	<ul style="list-style-type: none"><li>- Tissue level electrophysiology from anywhere on the heart</li><li>- Viable in organotypic culture for chronic studies</li><li>- Genetic manipulation</li></ul>	<ul style="list-style-type: none"><li>- Lower optical Signal-to-noise ratio compared to wedge preparations</li><li>- Cannot study true ionic currents</li><li>- Non-patient specific</li></ul>


## Supplemental Table 1. Comparison between *ex-vivo* human cardiac models

Pros and cons of each established and currently developing modeling for human *ex-vivo* cardiac study.



### Supplemental Figure 1. Quiescent Left Ventricular Slices

Ventricular slices from donor tissue do not exhibit automaticity. Action potential can only be observed and recorded during electrical stimulation.



**Supplemental Figure 2. Activation and AP Duration maps of  $\alpha_1$ -AR stimulation**

Ventricular slices from donor hearts are treated either acutely (fresh and after control culture) or chronically. (a) Activation map of slices at each condition under 1 Hz pacing. Activation appears slower under acute  $\alpha_1$ -AR stimulation and faster under chronic  $\alpha_1$ -AR stimulation. (b) Representative AP duration maps of slices at each condition under 1 Hz pacing. Acute  $\alpha_1$ -AR stimulation clearly increased AP duration, while chronic  $\alpha_1$ -AR stimulation had the reverse effect.