

Supplementary Methods

Surface sterilization of eggs

Axenic larvae were produced by placing eggs into sterile Petri dishes containing 70% EtOH for 5 min, transferring to a solution of 3% bleach and 0.1% ROCCAL-D (Pfizer) for 3 min, transferring again to 70% EtOH for 5 min, and rinsing 3x in sterile water. Sterile water was produced by autoclaving.

Verification of sterility of axenic larvae and diet

Sterility of larvae was confirmed by culture-based and PCR analysis. A pool of 30 first instars from surface sterilized eggs were rinsed in sterile water and transferred to 100 μ l of sterile PBS and homogenized. Larval homogenates were plated on Luria broth (LB) and brain-heart infusion (BHI) agar plates at 28° C for 24-72 h and subsequently checked for colonies. DNA was also isolated from three pools of 10 larvae using the Gentra Puregene Yeast/Bacteria Kit (Qiagen). Following DNA extraction, the 16S rDNA gene was amplified using either universal bacterial primers 27F (5'-GAGAGTTTGGATCCTGGCTCAG-3') and 1492R (5'-GGTTACCTTGTTACGACTT-3') or *Wolbachia*-specific primers WspF (5'TGGTCCAATAAGTGATGAAGAACTAGCTA-3') and WspR (5'AAAAATTAAACGCTACTCCAGCTTCTGCAC—3'). For each mosquito species, three replicate 10 μ L reactions were conducted, each containing 1 X HotMaster Taq Buffer, 200 μ M of each dNTP, 0.2 μ M of each primer, 1 U of Hotmaster Taq polymerase (5 Prime) and 0.8 μ L of DNA template. Reactions were run on a MJ Mini thermocycler (Bio-Rad) with amplification cycle conditions as follows: denaturation at 94° C for 2 min, followed by 30 cycles of 94° C for 20 s, 50° C for 20 s, and 65° C for 1 min, and final extension at 65° C for 5 min. Products from the three replicate amplifications were pooled and 5 μ L of the combined PCR product was electrophoresed on a 1% (wt/vol) agarose gel. The same methods were used to verify sterility of irradiated diet.