

Maintenance therapy with histamine plus IL-2 induces a striking expansion of two CD56^{bright} NK cell subpopulations in patients with acute myeloid leukemia and supports their activation

SUPPLEMENTARY TABLES AND FIGURES

Supplementary Table S1: List of AML patients and control healthy donors. Compilation of the 11 untreated AML patients (average 62 years, range 23-83 years), 9 AML patients after chemotherapy (average 49 years, range 28-72 years), 9 AML patients treated in addition with HDC plus IL-2 (average 47, range 21-69 years) and the 48 healthy donors (average 43 years, range 23-83 years) from whom blood samples were obtained.


See Supplementary File 1

Supplementary Table S2: List of antibodies used for identification of NK cells by flow cytometry


mAb anti-	Conjugated Fluorochrome	Clone	Company
CD107a	FITC	H4A3	BD Pharmingen
CD16	APC-Cy7	3G8	BD Pharmingen
CD20	PE	2H7	Biolegend
CD3	PerCP-Cy5.5	UCTH1	BD Pharmingen
CD3	V450	SK7	BD Biosciences
CD45	V500	HI30	BD Biosciences
CD56	PE-Cy7	B159	BD Pharmingen
CD57	APC	NK-1	BD Pharmingen
CD94	FITC	HP-3D9	BD Pharmingen
IFN- γ	APC-Cy7	4S.B3	Biolegend
IgG1	APC	IS5-21F5	Miltenyi Biotec
IgG1	FITC	IS5-21F5	Miltenyi Biotec
IgG1	PE	IS5-21F5	Miltenyi Biotec
IL-10	PE	JES3-19F1	BD Pharmingen
IL-13	PerCP-Cy5.5	JES10-5A2	Biolegend
NKG2A	APC	Z199	BD Pharmingen
NKG2C	PE	134591	R&D Systems
NKG2D	APC	1D11	BD Pharmingen
NKp30	APC	p30-15	BD Pharmingen
NKp46	APC	9E2/NKp46	BD Pharmingen
pan-KIR	FITC	180704	R&D Systems
TNF α	PECY7	MAB111	Biolegend


Supplementary Table S3: Combination of markers used for flow cytometry

Cell surface staining					Intracellular staining	
1	2	3	4	5	6	7
CD56	CD56	CD56	CD56	CD56	CD56	CD56
CD3	CD3	CD3	CD3	CD3	CD3	CD3
CD16	CD16	CD16	CD16	CD16	CD45	CD16
CD45	CD94	KIR	CD158a	NKG2D	IFN- γ	IFN- γ
CD14	NKG2A	NKP30	CD57	NKP46	TNF	CD107a
CD20	NKG2C	NKP44	CD158b		IL-10	eFluor-670
					IL-13	
Live/dead	Live/dead	Live/dead	Live/dead	Live/dead	CD107	Live/dead
					Live/dead	




Supplementary Figure S1: Absolute numbers of NK cells in AML patients and healthy donors. Absolute cell numbers of total NK cells and of single NK cell subsets from 12 healthy donors (H), 9 patients after chemotherapy (Cx) and patients undergoing additional HDC plus IL-2 therapy after 1, 2 and 3 weeks (w) of treatment with HDC plus IL-2 (1w: n=9; 2w: n=5; 3w: n=6) are shown. The calculation is based on the total leukocyte count and the percentages obtained by flow cytometry for the proportion of CD3^{neg}CD56^{pos} cells (total NK cells) and the CD56^{bright}CD16^{neg}, CD56^{bright}CD16^{low} and CD56^{dim}CD16^{high} cells. Mean values \pm SD are shown. * $p < 0.05$, ** $p < 0.01$.





Supplementary Figure S3: *In vitro* effects of HDC and IL-2 on the expression of NK cell receptors in single NK cell subsets. PBMC from 4 healthy donors were stimulated with HDC (10^{-5} M), IL-2 (500 U/ml) or a combination of both for 6 days. Percentages of cells expressing the receptors NKG2D, NKp30, NKp46, KIR and NKG2A. were analyzed on the CD56^{bright}CD16^{neg}, CD56^{bright}CD16^{low} and CD56^{dim}CD16^{high} subsets. The percentage of cells positive for the specific markers was determined as the number of positively stained cells minus the number of cells stained with an isotype-matched negative control antibody.



Supplementary Figure S4: Polychromatic and multifunctional tracing method. NK cells were isolated from PBMC of healthy donors by negative selection, stained with anti-CD56-PE and anti-CD16-FITC and sorted by preparative flow cytometry in a FACS Aria. The three subsets were then individually labelled with either of the three cell tracer dyes CFSE (CD56^{bright}CD16^{neg}), eFluor670 (CD56^{bright}CD16^{low}) or CellTrace Violet (CD56^{dim}CD16^{high}) as described in Methods. The differentially labelled subsets were then recombined and a complete PBMC fraction reconstituted by adding the NK cell depleted fraction obtained during negative selection of the NK cells. The reconstituted PBMC fraction was then cultured for the respective experiments.