

Figure S1





Figure S2



Supplemental Figure Legends.

Figure S1. ZIKV MR766 infects primary differentiating CTB from midgestation placentas. CTB were plated on Matrigel and immunostained for NS3 protein at 96 hpi. Nuclei in blue (DAPI). **Related to Figure 1.**

Figure S2. ZIKV MR766 and Nica1-16 strains infect chorionic villus explants. Chorionic villus explants were plated on Matrigel and either infected or mock-infected with MR766 at 2 or 24 h after attachment to the Matrigel then fixed and immunostained at 3 days after plating (i.e., 2 or 3 dpi). (A) Uninfected control explant (8.1 wk) immunostained for TIM1 (red) and CD68 (green), a marker of Hofbauer cells. TIM1 was detected in CD68-positive cells. (B-G) A 7.5 wk explant infected with MR766 2 h after attachment and immunostained for cytokeratin (CK; MAb 7D3), ZIKV proteins E (B, E) and NS3 (C, F) and CD68 (D, G) in adjacent sections. Both E and NS3 were detected in CD68-positive cells, circumscribed by dotted lines. Chorionic villus explants (8.4 wk) were infected with MR766 2 h (H-M) or 24 h (N-S) after attachment to Matrigel, and adjacent sections were immunostained at 3 days after attachment for E protein and cytokeratin (CK; 7D3 MAb) (H, inset K, N, inset Q), NS3 and CK (I, inset L, O, inset R) and NS3 and Ki67 (J, inset M, P, inset S). In explants infected 2 h after attachment, ZIKV proteins were detected in proliferating CTB (K and L, arrowheads). In explants infected 24 h after attachment, ZIKV protein-positive cells were detected in underlying CTB (Q and R, arrowheads). NS3-positive cells (M and S, arrowheads) were in a zone of Ki67-positive proliferating cells (arrows). Nuclei in blue (DAPI). VC: villus core. **Related to Figure 1.**

Table S1. Titers of prototype and Nicaragua ZIKV strains in different cell types of the human placenta. Related to Figure 1.

ZIKV strain	dpi ^a	AmEpC (gestation week) ^b				TBPC ^c	HPF ^d	HUVEC ^d	CTB D206 ^d	CTB D207 ^e
		21.2	22.6	38.6	40.2					
Uganda MR766	3	2x10 ⁶	4x10 ⁶	4.5x10 ⁵	4x10 ⁵	7.5x10 ⁵	3x10 ⁶	2.5x10 ⁶	2.7x10 ³	1.0x10 ⁵
	5	3.5x10 ⁶	2.7x10 ⁷		3.5x10 ⁵	1.5x10 ⁶	2.7x10 ⁶	1.5x10 ⁶	ND	ND
	7	3.75x10 ⁶	2x10 ⁶		5.2x10 ⁴	1.75x10 ⁵	3x10 ⁶	5x10 ⁶	ND	ND
Nica 1-16	3	1.00x10 ⁵	8.75x10 ⁵	3.75x10 ⁴	2x10 ⁴	1x10 ⁴				
	5	1.75x10 ⁵	1.63x10 ⁵	3.38x10 ⁴	2.75x10 ⁴	1.5x10 ⁴	ND	ND	ND	ND
	7	2.0x10 ⁵	2.88x10 ⁵	4.35x10 ⁴	1.75x10 ⁴	1.75x10 ⁵				

^adpi, days postinfection. ^bInfected at MOI 0.001; ^cInfected at MOI 1 (MR766) and MOI 0.004 (Nica1-16);

^dInfected at MOI 1; ^eInfected at MOI 2. Titers were determined by FFA (FFU/ml) on BHK cells after 3 days.

Table S2. Summary of Axl, Tyro3 and TIM1 immunohistochemical staining in biopsies of human placenta and fetal membranes. Related to Figure 5.

Donor	GA (wk)	Axl				Tyro3				TIM1			
		AmEpC	Ch	Villi	Decidua	AmEpC	Ch	Villi	Decidua	AmEpC	Ch	Villi	Decidua
AM5	22	—				++				++			
AM6	20	+++				+				++			
AM7	22	++				+				+++			
AM8	23	++				+							
AM12	17	+++				+++				+++			
AM13	21	+++											
AM18	27	+++				+				+++			
AM19	18	+++				+				+++			
SF147	23	++	+	+	++								
SF148	17	++	+	+	++							++	++
SF149	23	++	+	+	++							++	++
UCSF7	19			+	++								
AM17	40	—				++				++			
SF34	39	—	++	+	+++								
SF38	39	—	++	+	+++					++	++	++ core STB	++ CTB
SF66	32	—	++	+	+++								
SF70	32	—	++	+	+++								
SF88	39								++ CTB	++	++	++ core STB	++ CTB
SF92	41	—	++	++									
SF95A	36	—	++	+									
SF102	39	—	++	+									
SF105	32	—								++	+	++ core STB	+++ CTB
SF107	36	—			++	++	++	—	+++ CTB	++	+	++ core STB	++ CTB
SF111	40	—			++	+++			+++ CTB	+++			+++ CTB
SF114	37	—	++	+	+++	++			+++ CTB				
SF126A	37	—			+++	++	+++	—	+++ CTB	++	+	++ core STB	++ CTB
SF131	36									+++			++ CTB
SF132	36									+++			++ CTB
SF133	38									+++			++ CTB
SF107	36					++	++	—	+++ CTB	+	+	++ core STB	++ CTB

STB, syncytiotrophoblasts. CTB, cytotrophoblasts. Ch, chorion. The donor number is shown in the first column.

**Table S3. Summary of cofactor expression by tissue/cell type.
Related to Figures 3 and 5.**

	AmEpC	Chorion	Villi	Decidua	CTBs
<u>Biopsies</u>					
AXL	10/24	11/11	12/12	12/12	3/3
TYRO3	13/13	3/3	0/3	6/6	6/6
TIM1	17/17	6/6	8/8	12/12	10/10
<u>Primary cells</u>					
AXL	11/13	NA	NA	NA	4/7
TYRO3	9/13	NA	NA	NA	7/7
TIM1	13/13	NA	NA	NA	7/7

Supplemental Experimental Procedures

Antibodies and reagents

Commercial antibodies used were: goat polyclonal anti-Axl and -GATA4 (R&D Systems); rabbit polyclonal anti-Tyro3, -TIM1, -vimentin, and -Ki67 (Abcam), anti-von Willebrand factor (Thermo Scientific), and anti-cytokeratin 19 (Proteintech); and mouse anti-TIM1 (Abcam), anti-CD209 (BD Pharmingen), anti-CD68 and -cytokeratin 7 (Dako), anti-CD31 (Fitzgerald) and anti-actin (Sigma-Aldrich). Rat anti-human cytokeratin MAb (clone 7D3) and mouse anti-MHC class I HLA-G MAb (clone 4H84) were gifts from Susan Fisher (UCSF) (Damsky et al., 1992) and Michael McMaster (UCSF).

Immunofluorescence staining

Cells grown on cover slips were fixed with 4% PFA and permeabilized with 0.1% Triton X-100. Frozen tissues (5 cases) were cut into 5- μ m sections. After blocking with bovine serum albumin (BSA) or 3-5% normal serum matching the secondary antibody source, primary antibodies were incubated overnight followed by secondary antibodies (fluorescein isothiocyanate or rhodamine red-X conjugates; Jackson ImmunoResearch). Images were obtained using a Nikon Eclipse 50i microscope with a SPOT 7.4 Slider camera (Diagnostic Instruments). Images were processed and analyzed using Adobe Photoshop and NIH Image J software.

Immunoblotting

Whole cell lysates were obtained using RIPA buffer (50 mM Tris [pH 7.4], 150 mM NaCl, 1% [v/v] Nonidet-P40, 1 mM EDTA, 0.1% [w/v] SDS, and 1% Na-deoxycholate) supplemented with complete protease inhibitor cocktail (Roche), and immunoblot analysis was performed as described previously (Tabata et al., 2015).

Statistical analysis

An unpaired t-test was used for pair-wise statistical comparisons and a one-way ANOVA with Dunnett's post-test for multiple comparisons. The mean +/- standard deviation (SD) is shown in all bar- and dot plots. All statistical analyses were performed in GraphPad Prism v6.0.

