S3 Table. Spike experiment to assess reduction in ZIKV RNA levels after tissue homogenization. The same ZIKV stock, 'spike input', was added to MEM alone or tissues from non-ZIKV inoculated macaques and homogenized. Mean vRNA levels were compared between the samples with and without the tissue and showed that spiked vRNA levels were reduced in homogenized tubes containing macaque tissues.

		mean log ₁₀ ZIKV RNA copies/gram measured after homogenization	log ₁₀ ZIKV RNA copies/gram post- homogenization loss
	non-homogenized ZIKV spike input	9.6	n/a
	bladder	8.2	1.4
non-ZIKV infected macaque tissues spiked with ZIKV	brain	6.0	3.6
	liver	8.3	1.3
	muscle	8.8	0.8