


Fig. S2 PCR detection of *VmDCLs* deletion mutants using four pairs of primers.



The protoplast of wild type 03-8 was used for *VmDCL1* and *VmDCL2* single gene deletion mutants ($\Delta VmDCL1$ and $\Delta VmDCL2$). M: DM2000 maker. For $\Delta VmDCL1$ detection, Lane 1: product amplified by VmDCL1-5F/VmDCL1-6R for detecting *VmDCL1*. Lane 2: product amplified by H852/H850 for detection of the *HYG* insertion. Lane 3 and lane 4: products amplified by VmDCL1-7F/H855R and H856F/VmDCL1-8R for confirming homologous recombination. For $\Delta VmDCL2$ detection, Lane 1: product amplified by VmDCL2-5F/VmDCL2-6R for detecting *VmDCL2*. Lane 2: product amplified by H852/H850 for detection of the *HYG* insertion. Lane 3 and lane 4: products amplified by VmDCL2-7F/H855R and H856F/VmDCL2-8R for confirming homologous recombination. The protoplast of $\Delta VmDCL1$ was used for *VmDCL1* and *VmDCL2* double genes deletion mutant ($\Delta VmDCL1DCL2$). For $\Delta VmDCL1DCL2$ detection, Lane 1: product amplified by VmDCL2-5F/VmDCL2-6R for detecting *VmDCL2*. Lane 2: product amplified by G850/G852 for detection of the *NEO* insertion. Lane 3 and lane 4: products amplified by VmDCL2-7F/G855R and G856F/VmDCL2-8R for confirming homologous recombination.