

Supporting Information for:

**Carbon nanodots as dual-mode nanosensors for selective
detection of hydrogen peroxide**

Cheng-Long Shen¹, Li-Xia Su¹, Jin-Hao Zang¹, Xin-Jian Li¹, Qing Lou^{1*}, Chong-Xin
Shan^{1,2*}

¹*School of Physics and Engineering, Zhengzhou University, Zhengzhou 450052, China.*

²*State Key Laboratory of Luminescence and Applications, Changchun Institute of Optics, Fine
Mechanics and Physics, Chinese Academy of Sciences, No. 3888 Dongnanhu Road, Changchun,
130033, China.*

List of Figures and Tables




Figure S1. (a) TEM image of CDs after adding 0.5 M H_2O_2 . (b,c) Histogram of distribution from TEM before (b) and (c) after adding 0.5 M H_2O_2 .

From the histogram and TEM image, the diameters of the CDs after adding H_2O_2 enlarge about ten times, which are considered to be the aggregation of CDs from H_2O_2 .




Figure S2. Excited-emission matrix of the CDs with 2500 rpm (a), 5000 rpm(b), 7500 rpm(c) centrifugation and with the addition in 0.5 M H₂O₂ (d).

The excited-emission matrix of the CDs with 2500 rpm, 5000 rpm and 7500 rpm centrifugation can prove the dual-emission fluorescence is from one kind of CDs with the same size.




Figure S3. Absorption spectrum of the CDs normalized at 500 nm before and after adding 0.5 M H_2O_2 .




Figure S4. XRD pattern of the CDs after adding 0.5 M H₂O₂.

The XRD spectra of the CDs before and after adding H₂O₂ alter little.




Figure S5. Raman spectrum of the CDs after adding 0.5 M H₂O₂.




Figure S6. (a) Fluorescence stability of the CDs before and after adding 0.5 M H₂O₂ with the emission at 450 nm and 500 nm, respectively. (b) Photographic images of the CDs and the CDs under 1M H₂O₂ after 3 day.

The fluorescence stability of the CDs before and after adding H₂O₂ changes little.




Figure S7. FTIR spectra of the CDs before and after adding 0.5 M H_2O_2 .




Figure S8. (a,b) XPS (full survey) of the CDs before (a) and after (b) adding 0.5 M H_2O_2 . (c,d) XPS (C1s) of CDs before (c) and after (d) adding 0.5 M H_2O_2 . (e,f) XPS (O1s) of CDs before (e) and after (f) adding 0.5 M H_2O_2 .




Figure S9. The fluorescence intensity of the $112.5 \mu\text{g mL}^{-1}$ CDs diluted for 10, 20, 30 and 40 times.

The high concentration may induce the fluorescence quenching of the CDs. So, the CDs should be diluted with a proper concentration (30x, $3.75 \mu\text{g mL}^{-1}$) for increasing the accuracy of the nanosensor based on the CDs.




Figure S10. The normalized fluorescence intensity of the CDs under different concentrations of H_2O_2 (0, 0.05, 0.1, 0.15, 0.25, 0.5, 1 and 2 M, from left to right).

With increasing the concentrations of H_2O_2 , the fluorescence intensity of the CDs at 500 nm obviously is enhanced.




Figure S11. The liner range of the absorption for sensing H_2O_2 .




Figure S12. Fluorescence spectra of CDs in the presence of different oxidants.




Figure S13. UV-vis spectra of CDs in the presence of different oxidants.




Figure S14. The fluorescence intensity varied with time after injecting H₂O₂.

Table S1. Limit of detection (LOD) of CDs sensor for H₂O₂

1	0.57062	11	0.57339
2	0.57068	12	0.57407
3	0.57217	13	0.5737
4	0.57029	14	0.57469
5	0.57011	15	0.57456
6	0.57134	16	0.5737
7	0.56819	17	0.573
8	0.57313	18	0.57362
9	0.57228	19	0.5729
10	0.56779	20	0.57937

The LOD rule of IUPAC: $LOD = (k \cdot b) / S$

k=3, S was the slope of calibration curve, b was the standard deviation of blank.

Calibration: S = 0.547, b = 0.00256, so LOD = 0.014M.

Table S2. Comparison of analytical performance of different nanosensors for H₂O₂ determination.

Sensing materials	Method	Linear range	LOD ^{a)}	Remark	Ref.
Au-Ag/C nanoclusters	Colorimetry	0.8-90, 90-500 μ M	0.3 μ M	Noble metals	1
SiO ₂ NPs-HRP ^{b)}	Colorimetry	1.2-72 μ M	1.3 μ M	Complicated operation with enzyme reaction	2
CNT/AgNPs	Electrochemistry	0.05-17mM	0.5 μ M	Complicated operation and noble metals	3
Gr-CCS-AgNPs ^{c)}	Electrochemistry	20 μ M- 5.02mM	2.49 μ M	Complicated operation and noble metals	4
Au NC	Fluorescence	0.161-19.32 mM	20 μ M	Difficulty for naked-eye detection and noble metals	5
dLys-AgNCs ^{d)}	Fluorescence	0.8-200 μ M	0.2 μ M	Difficulty for naked-eye detection and noble metals	7
polymer nanoparticles	Fluorescence	6-1000 μ M	2.0 μ M	Difficulty for naked-eye detection	9
CQD-HRP ^{e)}	Fluorescence	0.5-50 μ M	0.2 μ M	Difficulty for naked-eye detection and operation with enzyme reaction	10
Carbon dots	Fluorescence and colorimetry	0.05-0.5M	50 mM	Simple and visualization by naked eye	This work

^{a)} Limit of detection. ^{b)} SiO₂ nanoparticles and horse radish peroxidase. ^{c)} Silver nanoparticles selectively deposited on graphene-colloidal carbon sphere composite. ^{d)} Lysozyme-silver nanoclusters. ^{e)} Carbon quantum dots and horse radish peroxidase.

References:

- 1 L. Zhang, W. Hou, Q. Lu, M. Liu, C. Chen, Y. Zhang and S. Yao, *Anal. Chim. Acta*, 2016, **947**, 23-31.
- 2 J. Pla-Tolós, Y. Moliner-Martinez, C. Molins-Legua and P. Campins-Falcó, *Sensor. Actuat. B-Chem.*, 2016, **231**, 837-846.
- 3 W. Zhao, H.Wang, X. Qin, Z.Wang, Z. Miao, L. Chen, M. Shan, Y. Fang and Q. Chen, *Talanta*, 2009, **80**, 1029-1033.

- 4 H. Wang, H. Wang, T. Li, J. Ma, K. Li and X. Zuo, *Sensor. Actuat. B-Chem.*, 2017, **239**, 1205-1212.
- 5 P. Zhang, Y. Wang and Y. Yin, *Sensors*, 2016, **16**, 1124.
- 6 F. Liu, T. Bing, D. Shangguan, M. Zhao, and N. Shao, *Anal. Chem.*, 2016, **88**, 10631-10638.
- 7 H. Tan, C. Ma, Q. Li, L. Wang, F. Xu, S. Chen and Y. Song, *Analyst*, 2014, **6139**, 5516-5522.
- 8 Y. Zhang, X. Yang and Z. Gao, *RSC Adv.*, 2015, **5**, 21675-21680.