

## **Electronic Supplementary Information (ESI)**

### **A single-chain derivative of the relaxin hormone is a functionally selective agonist of the G protein-coupled receptor, RXFP1**

Mohammed Akhter Hossain,<sup>\*ab</sup> Martina Kocan,<sup>c</sup> Song T. Yao,<sup>a</sup> Simon G. Royce,<sup>d</sup> Vinojini B. Nair,<sup>ab</sup> Christopher Siwek,<sup>c</sup> Nitin A. Patil,<sup>a,b</sup> Ian P. Harrison,<sup>d</sup> K. Johan Rosengren,<sup>e</sup> Stavros Selemidis,<sup>d</sup> Roger J. Summers,<sup>c</sup> John D. Wade,<sup>\*ab</sup> Ross A. D. Bathgate,<sup>\*af</sup> and Chrishan S. Samuel<sup>\*d</sup>

<sup>a</sup>*Florey Institute of Neuroscience and Mental Health, The University of Melbourne, Victoria 3010, Australia. E-mail: [akhter.hossain@florey.edu.au](mailto:akhter.hossain@florey.edu.au), [john.wade@florey.edu.au](mailto:john.wade@florey.edu.au), [bathgate@florey.edu.au](mailto:bathgate@florey.edu.au),*

<sup>b</sup>*School of Chemistry, The University of Melbourne, Victoria 3010, Australia.*


<sup>c</sup>*Monash Institute of Pharmaceutical Sciences, Monash University, Victoria, Australia.*

<sup>d</sup>*Cardiovascular Disease Program, Biomedicine Discovery Institute and Department of Pharmacology, Monash University, Victoria, Australia. E-mail: [chrishan.samuel@monash.edu](mailto:chrishan.samuel@monash.edu)*

<sup>e</sup>*The University of Queensland, School of Biomedical Sciences, Brisbane, QLD 4072, Australia.*

<sup>f</sup>*Department of Biochemistry, The University of Melbourne, Victoria 3010, Australia*

**Figure S1**



**Fig. S1** Schematic representation of the mechanism of RXFP1 activation by H2 relaxin and B7-33. While (A) H2 relaxin was shown to stimulate both cAMP and ERK pathways, (B) B7-33 selectively activated ERK pathway in naturally expressing RXFP1 cells leading to its strong H2-relaxin-like anti-fibrotic action *in vitro* ((myo)fibroblasts) and *in vivo*. Accumulation of cAMP is believed to be related to prostate cancer promotion and activation of ERK and subsequent promotion of MMPs is related to beneficial effects (e.g. anti-fibrosis and vasodilation) of relaxin. B7-33 is thus a novel functionally selective RXFP1 agonist.

**Figure S2**

**A**

		13      17      20
AcB7-33	Ac-VIKLSGRELVRAQ I AISGMSTWS <b>KRSL</b> -NH2	
R13A(AcB7-33)	Ac-VIKLSG <b>A</b> ELVRAQ I AISGMSTWS <b>KRSL</b> -NH2	
R17A(AcB7-33)	Ac-VIKLSGRELV <b>AA</b> Q I AISGMSTWS <b>KRSL</b> -NH2	
I20A(AcB7-33)	Ac-VIKLSGRELVRA <b>QA</b> AISGMSTWS <b>KRSL</b> -NH2	
R13/17A. I20A(AcB7-33)	Ac-VIKLSG <b>A</b> ELV <b>AA</b> Q AISGMSTWS <b>KRSL</b> -NH2	

**B**

Peptides	Human RXFP1 (HEK-293T)	
	Eu-H2 pKi	cAMP pEC50
H2 relaxin	8.96 ± 0.03 (5)	10.49 ± 0.13 (4)
<b>B7-33</b>	5.54 ± 0.13 (5) <sup>#</sup>	5.12 ± 0.06 (3) <sup>#</sup>
AcB7-33	6.00 ± 0.0 (1) <sup>#</sup>	5.40 ± 0.04 (3) <sup>#</sup>
R13A(AcB7-33)	No binding	No activity
R17A(AcB7-33)	No binding	No activity
I20A(AcB7-33)	No binding	No activity
R13/17A. I20A(AcB7-33)	No binding	No activity

#P<0.001 vs H2 relaxin

**Fig. S2** (A) Amino acid sequences of **B7-33** variants with key residues mutations. (B) Pooled binding affinity (pKi) and cAMP activity (pEC50) data. The dramatic reduction in the binding affinity and cAMP activity of mutated analogues for RXFP1 indicated that **B7-33** used the same residues as H2 relaxin to interact with the LRRs of RXFP1.


**Figure S3**

Peptides	Purity (%)
B7-33	97
AcB7-33	95
R13A(AcB7-33)	95
R17A(AcB7-33)	94
I20A(AcB7-33)	98
R13/R17A.I20A(AcB7-33)	97

**Fig. S3** Purity of B7-33 and its variants

**Figure S4**

**1a. HPLC profile of B7-33**




1b. ESI trace of B7-33

Mass Spectrometry

Sample : 2533613.F


H-VIKLSGRELVRAQIAISGMSTWSKRLS-NH<sub>2</sub>



Theoretical MW of parent is 2986.6 dalton  
3.9% of ion count is below range in figure  
MS operated in the positive ion mode

Instrument : Perkin-Elmer Sciex API 100  
Eluent : 0.1% Acetic acid in 60% acetonitrile  
Ion Source : Ionspray  
Detection : Ion counting

*2a. HPLC profile of AcB7-33*




2b. ESI trace of AcB7-33

Mass Spectrometry

Sample : 2533622.F

Ac-VIKLSGRELVRAQIAISGMSTWSKRSL-NH<sub>2</sub>



Theoretical MW of parent is 3028.6 dalton

0.5% of ion count is below range in figure


0.1% of ion count is above range in figure

MS operated in the positive ion mode

Instrument : Perkin-Elmer Sciex API 100  
Eluent : 0.1% Acetic acid in 60% acetonitrile  
Ion Source : Ionspray  
Detection : Ion counting

*3a. HPLC profile of R13A(AcB7-33)*

Analytical RP-HPLC  
Sample : 2533617F  
Result Id : 9664




	RT	PeakArea	Peak-Height	% Area
1	10.27	6762	1336	0.15
2	10.49	19420	1583	0.43
3	11.15	12202	2637	0.27
4	11.27	98987	17621	2.21
5	11.37	56644	18589	1.27
6	11.57	4255815	467032	95.10
7	11.95	25420	4709	0.57

Column : 150 x 4.6 mm, Monitor C18, Column Engineering  
Solvent A : 0.1% trifluoroacetic acid in 100% water  
Solvent B : 0.1% trifluoroacetic acid in 90% acetonitrile (aqueous)  
Gradient : 10%B for 1.0 min; 10-66.6%B (Linear) over 15.0 min; then column reequilibration  
Flow rate : 1.5 mL/minute  
Wavelength : PDA 214.0 nm

Date Acquired: 25-Sep-2014 20:44:33 EST  
Acq Method Set S6\_15cm\_25min\_TFA  
Analysis using Empower 2 Software Build 2154


*3b. ESI trace of R13A(AcB7-33)*



Theoretical MW of parent is 2943.5 dalton  
13.0% of ion count is below range in figure  
0.5% of ion count is above range in figure  
MS operated in the positive ion mode

Instrument : Perkin-Elmer Sciex API 100  
Eluent : 0.1% Acetic acid in 60% acetonitrile  
Ion Source : Ionspray  
Detection : Ion counting

*4a. HPLC profile of R17A(AcB7-33)*




4b. ESI trace of R17A(AcB7-33)

Mass Spectrometry

Sample : 2533616.J

Ac-VIKLSGRELVAAQIAISGMSTWSKRSL-NH<sub>2</sub>



Theoretical MW of parent is 2943.5 dalton


7.1% of ion count is below range in figure

MS operated in the positive ion mode

Instrument : Perkin-Elmer Sciex API 100  
Eluent : 0.1% Acetic acid in 60% acetonitrile  
Ion Source : Ionspray  
Detection : Ion counting

*5a. HPLC profile of I20A(AcB7-33)*

Analytical RP-HPLC  
Sample : 2533618Fr1  
Result Id : 10150




	RT	PeakArea	PeakHeight	% Area
1	9.54	22668	6232	1.09
2	10.36	3659	1479	0.18
3	10.46	2035561	365152	98.17
4	10.71	11582	1986	0.56

Column : 150 x 4.6 mm, Monitor C18, Column Engineering  
Solvent A : 0.1% trifluoroacetic acid in 100% water  
Solvent B : 0.1% trifluoroacetic acid in 90% acetonitrile (aqueous)  
Gradient : 10% B for 1.0 min; 10-66.6% B (Linear) over 15.0 min; then column reequilibration  
Flowrate : 1.5 mL/minute  
Wavelength : PDA 214.0 nm

Date Acquired: 26-Sep-2014 11:33:11 EST  
Acq Method Set S8\_15cm\_25min\_TFA  
Analysis using Empower 2 Software Build 2154


5b. ESI trace of I20A(AcB7-33)



Theoretical MW of parent is 2986.6 dalton  
0.8% of ion count is below range in figure  
0.1% of ion count is above range in figure  
MS operated in the positive ion mode

Instrument : Perkin-Elmer Sciex API 100  
Eluent : 0.1% Acetic acid in 60% acetonitrile  
Ion Source : Ionspray  
Detection : Ion counting

*6a. HPLC profile of R13/17A.I20(AcB7-33)*




6b. ESI trace of R13/17A.I20(AcB7-33)

Mass Spectrometry

Sample : 2533619.E

Ac-VIKLSGAELVAAQAAISGMSTWSKRSL-NH<sub>2</sub>



Theoretical MW of parent is 2816.3 dalton

2.0% of ion count is below range in figure

0.2% of ion count is above range in figure

MS operated in the positive ion mode


Instrument : Perkin-Elmer Sciex API 100  
Eluent : 0.1% Acetic acid in 60% acetonitrile  
Ion Source : Ionspray  
Detection : Ion counting

**Figure S5**

H2 relaxin	B7-33
2 chain peptide, 3 disulfide bridges, 53 amino acids	Single chain peptide, no disulfide bridge, 27 amino acids
6 HPLC purification steps	1 HPLC purification step
Manufactured in 12 days	Manufactured in 2 days

**Fig. S5** Comparison between H2 relaxin and B7-33 and steps involved in their chemical assembly. B7-33 is far easier and cheaper to manufacture by chemical means and thus more cost-effective as a drug.

**Figure S6**



**Fig. S6** cAMP stimulating activity of B7-33 compared with H2 relaxin and the native ligand of RXFP2 (INSL3) in cells over expressing RXFP2. B7-33 did not activate cAMP responses in these cells. The data are expressed as mean  $\pm$  SEM of n=3 independent experiments.