


1 **Supplementary Information “Microglia Innately develop within cerebral organoids”**

2 Ormel et al.

3



Supplementary Figure 1. iPSC5 characterization and neuronal identity in organoids

a to c- Stem cell markers NANOG and TRA1-81 (a), SOX2 and TRA1-60 (b), and OCT4 and SSEA4 (c) expression in iPSC5 line. Scale bar 200 μ m.

d- Karyogram of iPSC 5.

e to g- Pluripotency potential was evaluated by spontaneous differentiation assay followed by immunohistochemistry for smooth muscle actin (SMA, e), β -III Tubulin (TUJ1, f), and α -fetoprotein (AFP, g). Scale bar 200 μ m.

h- mRNA expression of stem cell markers by qRT-PCR (normalized to *ACTB*) with a commercial embryonic stem cell line (HUES6). mRNA expression of iPSC 5 relative to HUES6.


i- Dorsal forebrain identity of organoids assessed by FOXG1 and PAX6 immunostainings at day 31. Representative pictures of organoids from iPSC 1 are shown. (right panel is a close-up of the left panel) Scale bars 200 μ m.

j- Neuronal features of the organoids was evaluated by immunohistochemistry for progenitor marker SOX2 and pan-neuronal marker TUJ1 at day 31. Representative pictures of organoids derived from iPSC 1 are shown. (right panel is a close-up of the left panel) Scale bars 200 μ m.

k, l and m- Presence of TBR2⁺ progenitors (k), cortical deep layer neurons expressing CTIP2 (l) and TBR1 (m), and mature neurons MAP2 (m) in cerebral organoids at day 31. Representative pictures of organoids from iPSC 1 are shown. Scale bars 200 μ m.

n, o and p- Superficial layers were assessed by SATB2 expression at day 31 (n) and presence of astrocytes evaluated by expression of GFAP (o) and S100 β (p) at day 52. Representative pictures of cerebral organoids from iPSC 1 are shown. Scale bars 200 μ m.

q, r, s- Glutamatergic spontaneous excitatory postsynaptic potentials (sEPSCs) detected at a holding potential of -65 mV. Voltage dependent sodium currents, necessary for the generations of spikes, could be elicited (q). Example of a single sodium current evoked at -10 mV is shown (r). IV plot depicting the typical course of the voltage dependency of the sodium currents is shown (s). Cerebral organoids from iPSC 1 were used for the electrophysiology.



Supplementary Figure 2. The quantity of microglia is similar between batches, donors, and timepoints in culture and co-mature with neurons

a- Mask images created by an automated macro in FIJI to quantify the fraction of nuclei, IBA-1, and TUJ1 positive area of a tiled image of a fluorescent staining. The fluorescent channels were split and a separate threshold was applied to enable further analyses.



Representative pictures of cerebral organoids from iPSC 1 are shown. Scale bar 500 μm .

b- Quantification of IBA-1/nuclei ratio from tiled fluorescent images of 2 sections per organoid. The variation between batches is similar to the variation between organoids of the same batch ranging from 0.05-0.2.

c- Percentage of CD11b⁺ cells in organoid single cell suspension when sorted with flow cytometry. Organoids from three donors were used for this experiment (iPSC 1, 3, and 5) at two timepoints (38 and 52 days *in vitro*). The mean percentage of CD11b⁺/CD45⁺ cells of donors iPSC 1, 3, and 5 was 0.83% \pm 0.3 (SD) at both time points for oMG (n = 6)

d- The increase of IBA1⁺ is positively correlated with TUJ1. Each data-point reflects the IBA-1 and TUJ1 fraction, normalized to nuclei, of one tiled image of organoids from three separate batches.

e- Neuronal identity and cyto-architecture is maintained in organoids containing microglia as shown in tiled pictures from sections of one organoid after 66 days in culture. Co-staining for: microglia (IBA-1) and mature neurons (NEUN) (* indicates NEUN⁺ cells in close proximity with IBA-1⁺ cells, left panel); radial glia (PAX6) and a pan-neuronal marker (TUJ1, middle panel); and for a deep cortical layer marker (CTIP2) and a post-mitotic projection neuron marker (TBR1, right panel). Scale bar 500 μ m.





Supplementary Figure 3. oMG express microglia-specific genes that could be measured in the whole organoid but also when the microglia population is enriched

a- Graphs depicting mRNA expression levels of microglia-specific genes in organoids. mRNA levels were determined by qRT-PCR and normalized to the geomean of the reference genes *SDHA2* and *ACTB*. Data represent the mean of four batches consisting of two

68 organoids per batch per time-point. All batches consisted of organoids derived from iPSC 1.
69 Error bars represent the standard error of the mean (SEM).
70 b- Magnetic automated cell sorting validation in oMG (upper panel) and adult MG (lower
71 panel) by qRT-PCR by using a panel of classical microglia genes but also genes that should
72 not be expressed by microglia (*GFAP*, *TUBB3*, *MAP2*, and *OLIG2*). mRNA levels were
73 compared with the flow through fraction (oMG flow through and adult MG flow through).
74 Data was log transformed and scaled for each sample to visualize the expression pattern.
75 mRNA levels were determined by qRT-PCR and normalized to the geomean of the reference
76 genes *SDHA2* and *ACTB*. n = 3 separate experiments in which oMG were isolated with
77 CD11b-magnetic cell sorting from 8 organoids per experiment and also an n = 3 separate
78 experiments to enrich adult MG with CD11b-magnetic cell sorting from fresh human brain
79 tissue. (*p < 0.05)

80



Supplementary Figure 4. Correlation among microglia-like cells

a- Spearman correlation analysis between oMG day 38, oMG day 52, adult MG 1, iPSC and fibroblast samples. DESeq2 rlog transformed raw gene counts of all genes annotated after removal of common genes ($\text{FDR} > 0.05$, sum of raw read counts > 0) between the samples were used as input. Size and color of circles indicate the strength and direction of the correlation, respectively.

b- Unsupervised hierarchical cluster analysis on DESeq2 rlog transformed raw counts of oMG day 38, oMG day 52, and adult MG1 based on all genes after removal of common genes ($\text{FDR} > 0.05$, sum of raw read counts > 0) between samples.

c- Principal component analysis on DESeq2 rlog transformed raw counts of oMG day 38, 52, adult MG, iPSC and fibroblasts.


d- Spearman correlation analysis between oMG day 38, 52, adult MG1, fetal MG, iPSC MG and adult MG2 scaled log transformed FPKM values of genes used in figure 3f. Median log transformed FPKM values for biological replicates were used as input for the correlation analysis.

e- Spearman correlation analysis between oMG day 52, adult MG, fetal MG, and iPSC MG scaled log₂ FPKM values of a panel of transcription factor families that recognize microglia-related motifs. Median log transformed FPKM values for biological replicates were used as input for the correlation analysis.

f, g, and h- Volcano plots show differentially expressed genes ($\text{FDR} < 0.05$ in red) between day 38 oMG vs adult MG1 (f), day 52 vs day 38 oMG (g), and day 52 oMG vs adult MG1 (h).

i, j, and k- Venn diagrams show common expressed ($\text{FDR} > 0.5$, sum of raw read counts > 0) and differentially expressed genes ($\text{FDR} < 0.05$, enriched genes in sample/total amount of differentially expressed genes) between day 52 oMG vs adult MG1 (i), day 38 oMG vs adult

109 MG1 (j), and day 38 vs day 52 oMG (k) after shrinkage correction of the log2Fold change
110 and removal of identified common genes between iPSC, fibroblasts, oMG and adult MG1.
111



Supplementary Figure 5. Flow cytometry gating strategies used for cell sorting and protein expression quantification

a- Gating strategy to sort CD11b+CD45+ single living cells from the organoid/brain single cell suspensions. Representative gates of oMG enrichment of iPSC 3 are shown.

b- Gating strategy to quantify protein expression of CD11b+ cells gated from organoid/brain single cell suspension. Isotype controls of respective antibodies were used to accurately determine the protein expression. Representative gates of CD11b+ stained cells of iPSC 1 are shown.

Abbreviations: SSC-A = sideward scatter area; FSC-A = forward scatter area; FSC-W = forward scatter width.

125 **Supplementary Table 1.** Comparative overview of adaptations in the organoid

126 differentiation protocol used in this study and the original protocol of Lancaster et al. 2014

Description		Ormel et al. 2018	Lancaster et al. 2014
Embryoid bodies generation	medium	hES4 (-P/S; +FBS); Y27 (1:100)	hES4 (-P/S; +FBS); Y27 (1:100)
	plates	AggreWell 800, (300 microwells/well)	V- bottom ULS 96 well plate
	Volume	2 mL (1.75x10 ⁶ cells per mL)	150 µL (6x10 ⁵ cells/mL)
	cells per embryoid body (EB)	Approx. 11500 cells	9000 cells
Germ layer differentiation I	timing	up to day 4	up to day 4
	medium	hESC4 ; Y27 (1:100)	hESC4 ; Y27 (1:100)
	plates	flat bottom ULA 96 well plate	V- bottom ULA 96 well plate
	EB size	320 (+/- 40) µm	> 350-400 µm
Germ layer differentiation II	timing	Days 4 to 6	Days 4 to 6
	medium	hES0	hES0
	plates	flat bottom ULA 96 well plate	V- bottom ULA 96 well plate
	EB size	330 (+/- 37) µm	350-600 µm
Induction of neural ectoderm	timing	days 6 to 12	days 6 to 9
	medium	NIM with 0.1 µg/mL Heparin	NIM with 1 µg/mL Heparin
	plates	flat bottom ULA 96 well plate	flat bottom ULA 24 well plate
	EB size	330-570 µm	500-600 µm
Transfer to Matrigel	timing	Day 13	Day 11
	medium	Differentiation medium without RA	Differentiation medium without RA
	plates	60 mm petri dish	60 mm tissue culture dish
Transfer bioreactor	timing	Day 17	Day 15
	medium	Differentiation medium with RA	Differentiation medium with RA
	platform	Spinning bioreactor	Spinning bioreactor
	speed	27.5 rpm	25 rpm

127

128

129 **Supplementary Table 2.** Gene panels consisting of transcription factor genes important for
130 microglia functioning *in vivo*

Transcription family	Transcription factors important for microglia functioning <i>in vivo</i>
PU.1	<i>SPI1</i>
CTCF	<i>CTCF</i>
IRF	<i>IRF1</i> <i>IRF2</i> <i>IRF3</i> <i>IRF8</i> <i>IRF9</i>
RUNX	<i>RUNX1</i> <i>RUNX2</i>
AP-1	<i>JUN</i> <i>JUNB</i> <i>JUND</i> <i>FOS</i> <i>FOSB</i> <i>FOSL2</i> <i>ATF4</i> <i>BATF</i> <i>BATF2</i> <i>BATF3</i>
C/EBP	<i>CEBPA</i> <i>CEBPB</i> <i>CEBPG</i>
MEF2	<i>MEF2A</i> <i>MEF2B</i> <i>MEF2C</i> <i>MEF2D</i>
SMAD	<i>SMAD3</i>
MAF	<i>MAF</i> <i>MAF1</i> <i>MAFB</i> <i>MAFF</i> <i>MAFG</i> <i>MAFK</i>

131

132

133 **Supplementary Table 3.** Primer sequences used for qRT-PCR experiments

Gene	5'-Forward primer-3'	5'-Reverse primer-3'
<i>ACTB</i>	GCTCCTCCTGAGCGCAAG	CATCTGCTGGAAGGTGGACA
<i>GAPDH</i>	TGTTTCGACAGTCAGCCGCATCTTC	CAGAGTTAAAAGCAGCCCTGGTGA
<i>SOX2 endo</i>	CGAGGGAAATGGGAGGGGTGC	TGCAGCTGTCATTTGCTGTGGGT
<i>SOX2 viral</i>	GCATGACCAGCAGCCAGACCTA	TCTTGACCACGCTGCCCATGCT
<i>NANOG endo</i>	GCCTGTGATTTGTGGGCCTGA	GTGGAAGAATCAGGGCTGTCCTG
<i>OCT4 endo</i>	TGTCTCCGTCAACACTCTGGGC	CCCAAAAACCCTGGCACAACCTCC
<i>OCT4 viral</i>	AACCCCGAGGAAAGCCAGGACA	ACAGCACGCCCAGTGTCACT
<i>C-MYC endo</i>	GCGGGCACTTTGCACTGGAAGT	TTTCAGAGAAGCGGGTCTCTGGCA
<i>C-MYC viral</i>	TACGCCCTGTTGAAGCTGGCTG	TGCACCGAGTCGTAGTCGAGGT
<i>C-MYC total</i>	ACCGAAAATGCACCAGCCCCA	CGATCTGGTCACGCAGGGCAAA
<i>KLF4 endo</i>	TCCCGCCGCTCCATTACCAA	TTTTGCCGCGAGCCCGCGTAA
<i>KLF4 viral</i>	TGGAAGTTCGCCAGAAGCGACG	TTCATGTGCAGAGCCAGGTGGT
<i>dTOMATO</i>	TGAAGATGCGCGGCACCAACT	TGGTGGATCTCGCCCTTCAGCA
<i>SDHA2</i>	GAAGCCCTTTGAGGAGCACT	GTTTTGTGTCATCACGGGTCT
<i>AIF1</i>	AGACGTTTCAGCTACCCTGACTT	GGCCTGTGGCTTTTCTCTTTCTC
<i>PTPRC</i>	GCAGCTAGCAAGTGGTTTGTTC	AAACAGCATGCGTCCTTTCTC
<i>CD68</i>	CTTCTCTCATTCCCCTATGGACA	GAAGGACACATTGTACTCCACC
<i>ITGAM</i>	TGCTTCCTGTTTGGATCCAACCTA	AGAAGGCAATGTCACTATCCTCTGA
<i>CX3CR1</i>	CTTACGATGGCAGCCAGTGA	CAAGGCAGTCCAGGAGAGTT
<i>TREM2</i>	TCAGGAAGGTCCTGGTGGA	GGGTGGGAAGGGGATTCTC
<i>TYROBP</i>	TACGGCCTCTGTGTGTTGAG	CGGAAACAGCGTATCACTGAG
<i>P2RY12</i>	TTTGTGTGTCAAGTTACCTCCG	CTGGTGGTCTTCTGGTAGCG
<i>GPR34</i>	CCTGATGTCCAGTAACATTTCGC	CATGCAGGGAGTATCCTGGT
<i>PROS1</i>	TTGCACTTGTAACCAGGTTGG	CAGGAACAGTGGTAACTTCCAG
<i>MERTK</i>	CTCTGGCGTAGAGCTATCACT	AGGCTGGGTGGTGAAAACA
<i>GAS6</i>	CTCTCTCTGTGGCACTGGTA	CCTTGATCTCCATTAGGGCCAA
<i>HLADRA</i>	CCCAGGGGAAGACCACCTTT	CACCCTGCAGTCGTAAACGT
<i>TLR2</i>	ATCCTCCAATCAGGCTTCTCT	GGACAGGTCAAGGCTTTTACA
<i>TLR3</i>	CAAACACAAGCATTTCGGAATCTG	AAGGAATCGTTACCAACCACATT
<i>TLR4</i>	AGTTGATCTACCAAGCCTTGAGT	GCTGGTTGTCCCAAAATCACTTT
<i>TMEM119</i>	CTTCCTGGATGGGATGTTGGAC	GCACAGACGATGAACATCAGC
<i>GFAP</i>	AGGTCCTATGTGGAGCTTGAC	GCCATTGCCTCATACATCGCT
<i>TUBB3</i>	GGCCTTTGGACATCTCTTC	CTCCGTGTAGTGACCTTG
<i>MAP2</i>	CTCAGCACCGCTAACAGAGG	CATTGGCGCTTCGGACAAG
<i>OLIG2</i>	AGGACAAGAAGCAAATGACAG	TCCATGGCGATGTTGAGG
<i>IL6</i>	TGCAATAACCAACCCCTGACC	TGCGCAGAATGAGATGAGTTG
<i>IL1B</i>	TTTGAGTCTGCCCAGTTCCC	TCAGTTATATCCTGGCCGCC
<i>CD163</i>	TTTGTCAACTTGAGTCCCTTCAC	TCCCGCTACACTTGTTTTAC
<i>MRC1</i>	TGCAGAAGCAAACCAAACCTGTAA	CAGGCCTTAAGCCAACGAAACT
<i>TNF</i>	TGGAGAAGGGTGACCGACTC	TCACAGGGCAATGATCCCAA
<i>RUNX1</i>	AAGACCCTGCCCATCGCTTT	CATCATTGCCAGCCATCACAG
<i>SPI1</i>	GTGCAAAATGGAAGGGTTTCCC	TACTCGTGCCTTTGGCGTTG
<i>CSF1R</i>	ATCAGCATCCGGCTGAAAGT	CTCGAATCCGCACCAGCTCT
<i>IL34</i>	TGCACTGTCACGGGTTTTCT	CCCTCGTAAGGCACACTGAT
<i>CSF1</i>	GCAGGAGTATCACCGAGGAG	CACGAGGTCTCCATCTGACTG
<i>TGFB1</i>	CAATTCTTGCGGATACCTCAG	GCACAACTCCGTGACATCAA
<i>IRF8</i>	ATCAAAAGGAGCCCTTCCCC	ATCAAAAGGAGCCCTTCCCC
<i>TGFB1</i>	TCCAAACCACAGAGTGGGAAC	TCCAAACCACAGAGTGGGAAC
<i>TGFB2</i>	GTATCGCCAGCACGATCCCA	GAAACTTGACTGCACCGTTGTT
<i>CIQA</i>	GAGCACCAGACGGGAAGAAA	TAAGGCCCTTGATGCCTGTC

134

135

136 **Supplementary Table 4.** Antibodies used in this study for immune histo/cytochemistry

Antigen/target	Host species	Dilutions	Provider, article number
SMA	Mouse	1:100	Sigma, A2547
AFP	Rabbit	1:50	Quartett, 2011200530
Brachyury	Goat	1:1000	R&D systems, AF2085-SP
PAX6	Mouse	1:200	DSHB, Pax6-s
NEUN	Mouse	1:300	Abcam, AB104224
TUJ1	Rabbit	1:1000	Sigma, T2200
TUJ1	Mouse	1:1000	Covance, MMS-435P
IBA-1	Rabbit	1:1000	Wakko, 019-19741
IBA-1	Goat	1:1000	Abcam, AB5076
FOXG1	Rabbit	1:100	Abcam, AB18259
CD68	Rabbit	1:100	Invitrogen, MA5-13324
CTIP2	Rat	1:100	Abcam, AB18465
TBR1	Rabbit	1:100	Gift from Robert Hevner
S100 β	Rabbit	1:600	Dako, Z0311
GFAP-pan	Rabbit	1:1000	Dako, Z0334
PSD-95	Mouse	1:300	NeuroMab, 75028
MAP2	Mouse	1:300	Biologend, SMI-52p
PU.1	Rabbit	1:100	Invitrogen, A13971
SATB2	Rabbit	1:300	Abcam, AB34735
NANOG	Rabbit	1:200	
OCT4	Rabbit	1:200	
SSEA4	Mouse	1:200	STEMLight iPSC
TRA1-60	Mouse	1:200	characterization Kit, Cell
TRA1-81	Mouse	1:200	Signalling, 9656S
SOX2	Rabbit	1:200	
Mouse, 568	Donkey	1:1000	ThermoFisher, A10037
Mouse, 555	Donkey	1:1000	ThermoFisher, A31570
Mouse, 488	Donkey	1:1000	ThermoFisher, A21202
Rabbit, 568	Donkey	1:1000	ThermoFisher, A10042
Rabbit, 488	Donkey	1:1000	ThermoFisher, A21206
Goat, 488	Donkey	1:1000	ThermoFisher, A11055
Rat, 488	Donkey	1:1000	ThermoFisher, A21208
Rabbit, Atto647N	Goat	1:200	Sigma Aldrich, 40839
Rabbit, 488	Goat	1:200	ThermoFisher, A11034
Mouse, 594	Goat	1:200	ThermoFisher, A11032

137

138

139