



# Prospective Evaluation of Molecular Assays for Diagnosis of Vaginitis

Sandra S. Richter,<sup>a\*</sup> Joshua Otiso,<sup>a</sup> Oluwatosin J. Goje,<sup>b</sup> Sherilynn Vogel,<sup>a</sup> Jory Aebly,<sup>a</sup> George Keller,<sup>a</sup> Hillary Van Heule,<sup>a</sup> Danielle Wehn,<sup>b</sup> Amy L. Stephens,<sup>b</sup> Salena Zanotti,<sup>b</sup> Tricia Johnson,<sup>a</sup> Sixto M. Leal,<sup>c</sup> Gary W. Procop<sup>a</sup>

<sup>a</sup>Department of Laboratory Medicine, Cleveland Clinic, Cleveland, Ohio, USA

<sup>b</sup>Ob/Gyn & Women's Health Institute, Cleveland Clinic, Cleveland, Ohio, USA

<sup>c</sup>Department of Pathology and Laboratory Medicine, University of Alabama, Birmingham, Alabama, USA

**ABSTRACT** Molecular tests to diagnose conditions involving the disruption of normal microbiota are difficult to optimize. Using Nugent-scored Gram stain (NS) as the reference standard, we evaluated the performance of 3 molecular assays for the diagnosis of bacterial vaginosis (BV) and examined the impact of an incremental increase in bacterial targets. The BD Affirm assay includes a DNA probe for *Gardnerella vaginalis*, the Hologic transcription-mediated amplification (TMA) analyte-specific reagent (ASR) assay adds a second *Lactobacillus* sp. target, and the recently cleared *in vitro* diagnostic use (IVD) Aptima BV assay includes a third target (*Atopobium vaginae*). The diagnosis of vulvovaginal candidiasis (VVC) by the Affirm and *Candida* vaginitis Hologic TMA ASR assays was assessed using microscopy for yeast as the reference standard. From May to December 2018, 111 women with vaginitis symptoms prompting the clinician to order an Affirm test were enrolled with informed consent for the collection of additional specimens. Clinicians accurately predicted BV as the most likely diagnosis for 71% of the 45 patients with BV. Coinfection occurred in 13.5% of patients. For BV, the specificity of the Aptima IVD assay (86.3%) was higher than the Affirm assay (60.6%,  $P = 0.0002$ ), but sensitivities were not significantly different. For VVC, the sensitivity of the ASR assay (100%) was higher than Affirm (75.9%;  $P = 0.023$ ) and the specificity of the Affirm assay (98.8%) was higher than the ASR assay (86.6%;  $P = 0.004$ ).

**KEYWORDS** molecular, vaginitis

Molecular tests to diagnose conditions that represent a disruption of normal microbiota are difficult to optimize. Bacterial vaginosis (BV) is a common cause of vaginitis associated with alteration of the normal microbiome, namely, a decrease in *Lactobacillus* spp. accompanied by an increase in *Gardnerella vaginalis* and other bacteria (1). The reference standard for diagnosing BV is Nugent scoring of a Gram stain prepared with vaginal secretions (2). Research using broad-range molecular techniques has uncovered more bacteria associated with BV (*Atopobium vaginae*, *Megasphaera* spp., *Leptotrichia amnionii*, *Sneathia sanguinegens*, and *Porphyromonas asaccharolytica*), including 3 new species designated "bacterial vaginosis-associated bacteria" (BVAB1, BVAB2, and BVAB3) (3). An evaluation of a quantitative PCR targeting 7 bacteria associated with BV concluded *G. vaginalis* ( $\geq 10^9$  copies/ml) and *A. vaginae* ( $\geq 10^8$  copies/ml) had the best predictive value for diagnosing BV with Nugent score as the reference method (4).

Current options to diagnose BV in our institution include Gram stain with Nugent score or DNA probe (BD Affirm vaginal pathogens III) targeting *G. vaginalis*. The Affirm assay targets 2 additional infectious etiologies of vaginitis, namely, *Candida* spp. and *Trichomonas vaginalis*. The Affirm assay is ergonomically difficult for technologists to

**Citation** Richter SS, Otiso J, Goje OJ, Vogel S, Aebly J, Keller G, Van Heule H, Wehn D, Stephens AL, Zanotti S, Johnson T, Leal SM, Procop GW. 2020. Prospective evaluation of molecular assays for diagnosis of vaginitis. J Clin Microbiol 58:e01264-19. <https://doi.org/10.1128/JCM.01264-19>.

**Editor** Michael J. Loeffelholz, Cepheid

**Copyright** © 2019 American Society for Microbiology. All Rights Reserved.

Address correspondence to Gary W. Procop, [procopg@ccf.org](mailto:procopg@ccf.org).

\* Present address: Sandra S. Richter, bioMérieux, Durham, North Carolina, USA.

**Received** 31 July 2019

**Returned for modification** 22 August 2019

**Accepted** 24 October 2019

**Accepted manuscript posted online** 6 November 2019

**Published** 23 December 2019

perform, and recent microbiome studies suggest that using only a single target to diagnose BV is probably suboptimal.

Using the Nugent score as the gold standard, we evaluated the performance of 3 molecular assays for the diagnosis of BV and examined the impact of an incremental increase in bacterial targets; the Affirm assay targets *G. vaginalis*, the Hologic transcription-mediated amplification (TMA) analyte-specific reagent (ASR) assay adds a second *Lactobacillus* sp. target, and the recently cleared *in vitro* diagnostic use (IVD) Aptima BV assay contains a third target (*A. vaginae*). Additionally, a TMA ASR assay targeting *Candida* spp. was compared to the Affirm assay for the detection of vulvo-vaginal candidiasis (VVC) using microscopy for yeast as the reference standard. The prevalence of coinfections and the ability of clinicians to predict the etiology of vaginitis were also examined.

(The study was presented, in part, as a poster presentation at the ASM Microbe on June 23, 2019 in San Francisco, CA.)

## MATERIALS AND METHODS

**Clinical specimens.** From May to December 2018, Cleveland Clinic outpatients with vaginitis symptoms that would normally prompt their provider to order the BD Affirm assay were invited to participate in this institutional review board (IRB)-approved study. After informed consent was obtained, 2 vaginal swabs (Copan swab in Amies transport for Gram stain and Aptima vaginal swab for BV and VVC testing) were collected in addition to a swab for the Affirm assay. On a questionnaire submitted with patient samples to the main campus microbiology laboratory, providers were asked to record a description of vaginal discharge, whiff test results, and vaginal pH; and use clinical judgment to indicate the most likely syndrome (rank order BV, VVC, and trichomoniasis).

**BD Affirm.** The Affirm VPIII microbial identification test was performed according to the manufacturer's instructions and reported in the patient record according to laboratory protocol. The Affirm assay incorporates single-stranded DNA probes immobilized on a bead embedded in a probe analysis card, which contains a separate bead for each target organism and positive and negative controls. After hybridization, an enzyme conjugate-binding step occurs to enable detection. Observing visible blue color on the target organism bead is a positive result. Negative results indicate <10,000 CFUs of *Candida* cells, <200,000 CFUs of *G. vaginalis*, and <5,000 trichomonads.

**Microscopy.** Every Gram stain was read by 2 technologists who were blind to all of the molecular test results. Bacteria and yeast observed were quantified based on the number present per oil immersion field (i.e., rare, few, moderate, and many) (5). Any quantity of yeast was considered a positive result. The Nugent score system guidelines were applied to interpret the significance of bacteria observed and a third technologist (blind to the initial 2 Nugent score results) reviewed the slide if the initial 2 readings were not concordant (positive, negative, and intermediate). The presence of clue cells was noted but not factored into the assigned score. A predominance of medium-to-large Gram-positive bacilli suggestive of lactobacilli (score, 0 to 3) is consistent with normal vaginal flora. A predominance of Gram-negative or Gram-variable bacilli suggestive of *G. vaginalis*, *Mobiluncus* species, and *Bacteroides* species (score, 7 to 10) is consistent with bacterial vaginosis. An intermediate Nugent score (4–6) suggests a transition from normal vaginal flora and was considered positive for BV if  $\geq 2$  Amsel's criteria (positive whiff test, clue cells on Gram stain, or vaginal pH of >4.5) were present (6).

**Hologic ASR BV assay.** The ASR BV assay for the Panther instrument (Hologic) utilizes transcription-mediated amplification (TMA) technology. After release of RNA in specimens, hybridization occurs to target organisms. After target amplification, results are reported as T-times, utilizing a scale of 0 to 100 (lower values indicate less time to detection; T-time = 100 if no target detected). The ASR BV algorithm uses the *Lactobacillus* results only if a medium organism burden of *G. vaginalis* is detected. There were 2 scenarios when the assay was considered positive for BV, namely, (i) if a high burden of *G. vaginalis* was detected (T-time sample < T-time high control) or (ii) if a medium amount of *G. vaginalis* (T-sample < T-time low control) was detected along with no *Lactobacillus* sp.

**Aptima IVD BV assay.** In May 2019, the remnant specimens (stored at  $-70^{\circ}\text{C}$ ) were tested using a new investigational use only (subsequently cleared as IVD) Aptima BV assay that contains *A. vaginae* along with the *G. vaginalis* and *Lactobacillus* sp. targets in the ASR assay initially tested. Software was installed on the Panther instrument that provided a revised algorithm for interpreting results as negative or positive for BV and allowed an additional (fourth) test to be run on a specimen.

**Aptima ASR *Candida* assay.** The VVC assay utilizes TMA targeting 5 *Candida* species (*Candida albicans*, *Candida parapsilosis*, *Candida tropicalis*, *Candida dubliniensis*, and *Candida glabrata*). Probes provide detection of *C. glabrata* or other *Candida* spp. The *Candida* vaginitis (CV) assay was an ASR when utilized in our study and received FDA approval as a combination test with *T. vaginalis* on May 29, 2019.

**Data analysis.** Using microscopy as the reference standard, we determined the sensitivity, specificity, positive predictive value, and negative predictive value for all assays. Confidence intervals were calculated using the modified Wald method; McNemar's test was applied to analyze differences between test performance (sensitivity and specificity), with *P* values of  $\leq 0.05$  considered statistically significant (<https://www.graphpad.com>). The rate of coinfection was determined using the *T. vaginalis* result from the Affirm assay and microscopy for BV and VVC.

**TABLE 1** Comparison of BV molecular assays to NS Gram stain for 111 specimens<sup>a,b</sup>

Assay	TP	FP	TN	FN	Sensitivity (%)	Specificity <sup>c</sup> (%)	PPV (%)	NPV (%)
Aptima IVD	38	9	57	7	84.4 (70.9–92.6)	86.3 (75.9–92.9)	80.9 (67.2–89.8)	89.1 (78.8–94.9)
BD Affirm	39	26	40	6	86.7 (73.5–94.1)	60.6 (48.5–71.5)	60.0 (47.8–71.0)	87.0 (74.0–94.3)
Hologic ASR	34	12	54	11	75.6 (63.6–87.6)	81.8 (70.7–89.4)	73.9 (59.6–84.5)	83.1 (72.0–90.5)

<sup>a</sup>BV, bacterial vaginosis; NS, Nugent score; TP, true positive; FP, false positive; TN, true negative; FN, false negative; PPV, positive predictive value; NPV, negative predictive value. Values in parentheses are 95% confidence intervals.

<sup>b</sup>Intermediate NSs (4–6) were considered positive for BV if  $\geq 2$  of Amsel's criteria (positive whiff test, clue cells on Gram stain, and vaginal pH of  $>4.5$ ) were present.

<sup>c</sup>The specificity of Aptima IVD was higher than BD Affirm ( $P = 0.0002$ ).

## RESULTS

A total of 111 women were enrolled in the study. To achieve concordance for Nugent score interpretation (positive, negative, or intermediate), a third reader was required for 24% of specimens. There were 40 patients with positive Nugent scores (7–10) consistent with BV and 52 patients with low scores (0 to 3) interpreted as normal vaginal microbiota. Only 5 of the 19 specimens with intermediate Nugent scores were positive by Amsel's criteria and classified as positive for BV; the other 14 specimens were considered negative.

For BV, the sensitivities and specificities were 75.6% (34/45) and 81.8% (54/66) for the Hologic ASR assay, 86.7% (39/45) and 60.6% (40/66) for the Affirm assay, and 84.4% (38/45) and 86.3% (57/66) for the Aptima BV IVD assay, respectively (Table 1). Although the differences in sensitivity did not reach statistical significance, the specificity of the Aptima BV IVD assay was higher than the BD Affirm assay ( $P = 0.0002$ ). Excluding intermediate Nugent score, the sensitivities and specificities for BV were 82.5% (33/40) and 84.6% (44/52) for the Hologic ASR assay, 87.5% (35/40) and 63.4% (33/52) for the Affirm assay, 90% (36/40) and 92.3% (48/52) for the Aptima BV IVD assay, respectively (Table 2). This subset analysis also demonstrated higher specificity of the Aptima BV IVD assay than the BD Affirm assay ( $P = 0.0003$ ).

The Aptima ASR assay detected *Candida* spp. in all 29 specimens with yeast observed by microscopy, while only 22 (75.9%;  $P = 0.023$ ) were detected by the Affirm DNA probe (Table 3). There were more false positives ( $n = 11$ ) and lower specificity with the Aptima ASR assay than Affirm assay (86.6% versus 98.8%;  $P = 0.004$ ) for VVC.

Coinfection occurred in 15 (13.5%) patients (12 BV and VVC, 2 BV and trichomoniasis, and 1 trichomoniasis and VVC). Clinicians accurately chose BV as the most likely diagnosis for 71% of the 45 patients with BV.

## DISCUSSION

Despite the higher cost, there are a number of reasons why molecular assays to diagnose vaginitis have been embraced by clinicians and laboratorians. The high percentage of women that may have BV by Nugent criteria and, yet, are without symptoms (up to 75%) explains why clinicians appreciate having the option of ordering syndromic testing covering multiple etiologies (7–9). The importance of diagnosing and treating BV is underscored by the association with preterm delivery and acquisition of HIV and other infections (10–12).

The gold standard for BV is labor-intensive, and many laboratories have difficulty finding technologists competent and confident in reading Gram stains. The subjective nature of Nugent scoring is reflected by 24% of specimens in the current study needing

**TABLE 2** Comparison of BV molecular assays to positive and negative NS Gram stain for 92 specimens<sup>a</sup>

Assay	TP	FP	TN	FN	Sensitivity (%)	Specificity <sup>b</sup> (%)	PPV (%)	NPV (%)
Aptima IVD	36	4	48	4	90.0 (76.4–96.6)	92.3 (81.3–97.5)	90.0 (76.4–96.6)	92.3 (81.2–97.5)
BD Affirm	35	19	33	5	87.5 (73.4–95.0)	63.4 (49.8–75.2)	64.8 (51.5–76.2)	86.8 (72.2–94.7)
Hologic ASR	33	8	44	7	82.5 (67.7–91.6)	84.6 (72.2–92.3)	80.5 (65.7–90.0)	86.3 (74.0–93.5)

<sup>a</sup>The 19 specimens with intermediate NSs (4–6) were excluded. BV, bacterial vaginosis; NS, Nugent score; TP, true positive; FP, false positive; TN, true negative; FN, false negative; PPV, positive predictive value; NPV, negative predictive value. Values in parentheses are 95% confidence intervals.

<sup>b</sup>The specificity of Aptima IVD was higher than BD Affirm ( $P = 0.0003$ ).

**TABLE 3** Comparison of molecular assays to microscopy for detection of vulvovaginal candidiasis<sup>a</sup>

Assay	TP	FP	TN	FN	Sensitivity <sup>b</sup> (%)	Specificity <sup>c</sup> (%)	PPV (%)	NPV (%)
Aptima ASR	29	11	71	0	100 (86.1 to 100)	86.6 (77.4 to 92.5)	72.5 (57.0 to 84.0)	100 (93.0 to 100)
BD Affirm	22	1	81	7 <sup>d</sup>	75.9 (57.6 to 88.1)	98.8 (92.8 to >99.9)	95.7 (77.3 to >99.9)	92.1 (84.2 to 96.3)

<sup>a</sup>BV, bacterial vaginosis; NS, Nugent score; TP, true positive; FP, false positive; TN, true negative; FN, false negative; PPV, positive predictive value; NPV, negative predictive value. Values in parentheses are 95% confidence intervals.

<sup>b</sup>The sensitivity of Aptima ASR was higher than BD Affirm ( $P = 0.023$ ).

<sup>c</sup>The specificity of BD Affirm was higher than Aptima ASR ( $P = 0.004$ ).

<sup>d</sup>Quantities of yeast observed by microscopy but not detected by Affirm were rare ( $n = 4$ ), few ( $n = 1$ ), moderate ( $n = 1$ ), and many ( $n = 1$ ).

a third reader to achieve a concordant category result. It can be difficult to get agreement for specimens with Nugent scores adjacent to a categorical cutoff value (6 and 7 or 3 and 4). The quality of research studies without multiple readers for each Gram stain should be questioned.

There was little change in Affirm assay performance with the elimination of intermediate Nugent specimens (from 86.7% to 87.5% sensitivity and 60.6% to 63.4% specificity). Although the Affirm assay is FDA approved and performed by many laboratories, there is limited published data on the accuracy of this method to compare with our results. Sensitivity and specificity of the Affirm assay compared with the Nugent criteria determined in a study of 195 Korean women was 75.0% and 88.9% for *G. vaginalis* (13). The Affirm package insert reports 83.8% sensitivity and 100% specificity for *G. vaginalis* compared with Gram stain and 89.0% sensitivity and 99.1% specificity compared with culture of *G. vaginalis*. Culture is not recommended or performed routinely as a test for diagnosing BV since *G. vaginalis* can be recovered from cultures collected from healthy women.

Of the 3 molecular assays we evaluated for BV, Aptima BV IVD demonstrated the highest specificity, which may reflect value for the *A. vaginae* target unique to that assay. It is difficult to predict whether the difference observed in the Aptima IVD and Affirm assay positive predictive values would be seen in other settings (85% of our patients were enrolled from general gynecology clinics and 15% from a vulvovaginitis specialty clinic). The sensitivity and specificity of the Aptima IVD reached 90% only with elimination of intermediate Nugent-scored specimens (Table 2). This approach (eliminating intermediate Nugent-scored specimens) was used for an evaluation of the BD Max assay containing 2 additional BV targets (BVAB2 and Megasphaera-1) and reported 92.7% sensitivity and 91.5% specificity compared with Nugent score (14). Although our study had fewer subjects, the Aptima IVD results (90.0% sensitivity, 92.3% specificity) showed similar performance to that reported for BD Max. Although assays that incorporate more bacterial targets are attractive since they reflect the bacterial diversity that has been reported in BV, it is uncertain whether they will provide better diagnostic accuracy to offset the higher cost usually charged for additional targets. Third party reimbursement of multitarget nucleic acid amplification tests (NAAT) for BV may be denied by insurers who want to see evidence of better patient management or improved clinical outcomes to support NAAT as an alternative to traditional approaches.

For VVC, our Affirm assay sensitivity compared with microscopy was 75.9%; the sensitivity stated in the Affirm package insert is 80.6% compared with culture. A comparison of Affirm assay to KOH microscopy reported similar performance, with a sensitivity of 82% and specificity of 100% (15). Our Affirm assay specificity compared with microscopy (98.8%) was similar to the Affirm package insert comparison to culture (98.2%) and higher than the TMA assay (86.6%); false positives are not surprising when NAAT is compared with microscopy as the reference standard.

A major limitation of our study is the uncertainty regarding whether microscopy is the optimal gold standard for BV and VVC. Hopefully future studies will be performed that are designed to determine whether clinical outcomes are improved in populations where a highly sensitive NAAT method rather than microscopy is used to diagnose vaginitis. More subjects and clinical outcome data would enrich the study but were prevented by limited clinician time to enroll patients and the amount of funding.

Strengths of this evaluation include a relevant, prospective study design with a comparison of a new commercially available molecular assay (Aptima IVD) to an established, commonly ordered assay (Affirm) with little published performance data. The reference standard of microscopy was performed in a rigorous manner (minimum of 2 technologists blind to the molecular results). Comparisons of the performance of the Affirm assay with one target (*G. vaginalis*), Hologic BV ASR with 2 targets (*Lactobacillus* sp. and *G. vaginalis*), and Aptima BV IVD with a third target (*A. vaginae*) to a reference standard of microscopy using current clinical samples are a novel demonstration of the utility that may be associated with individual targets for diagnosing BV. Many laboratories use the Affirm assay but few have the resources needed to conduct a prospective study requiring informed consent to collect additional specimens. The ergonomic complaints reported by technologists performing the Affirm assay illustrates the relevance of comparing those results to more automated NAAT solutions.

Although coinfections in our study population were uncommon, ordering patterns at our institution for Affirm demonstrate clinicians' preference for syndromic tests that cover multiple etiologies. Although workflow and some performance parameters for the IVD TMA assays were better than the DNA probe, studies demonstrating improved clinical outcomes with NAAT testing for vaginitis may be needed to ensure reimbursement.

## ACKNOWLEDGMENTS

Funding of the study was provided by Hologic, Inc.

S.S.R. has been an employee of bioMérieux since August 12, 2019. When employed at Cleveland Clinic, she received research funding from bioMérieux, BD Diagnostics, Roche, Diasorin, Accelerate, and Affinity Biosensors. All other authors have no conflict of interest.

## REFERENCES

- Hill GB. 1993. The microbiology of bacterial vaginosis. *Am J Obstet Gynecol* 69:450. [https://doi.org/10.1016/0002-9378\(93\)90339-K](https://doi.org/10.1016/0002-9378(93)90339-K).
- Nugent RP, Krohn MA, Hillier SL. 1991. Reliability of diagnosing bacterial vaginosis is improved by a standardized method of gram stain interpretation. *J Clin Microbiol* 29:297–301.
- Fredricks DN, Fiedler TL, Marrazzo JM. 2005. Molecular identification of bacteria associated with bacterial vaginosis. *N Engl J Med* 353:1899–1911. <https://doi.org/10.1056/NEJMoa043802>.
- Menard J-P, Fenollar F, Henry M, Bretelle F, Raoult D. 2008. Molecular quantification of *Gardnerella vaginalis* and *Atopobium vaginae* loads to predict bacterial vaginosis. *Clin Infect Dis* 47:33–43. <https://doi.org/10.1086/588661>.
- Chan WW. 2016. Gram stain, p 3.2.1.1–3.2.1.16. In Leber A (ed), *Clinical microbiology procedures handbook*, 4th ed, ASM Press, Washington, DC.
- Amsel R, Totten PA, Spiegel CA, Chen KC, Eschenbach D, Holmes KK. 1983. Non-specific vaginitis: diagnostic criteria and microbial and epidemiologic associations. *Am J Med* 74:14–22. [https://doi.org/10.1016/0002-9343\(83\)91112-9](https://doi.org/10.1016/0002-9343(83)91112-9).
- Klebanoff MA, Schwebke JR, Zhang J, Nansel TR, Yu KF, Andrews WW. 2004. Vulvovaginal symptoms in women with bacterial vaginosis. *Obstet Gynecol* 104:267–272. <https://doi.org/10.1097/01.AOG.0000134783.98382.b0>.
- Yen S, Shafer M-A, Moncada J, Campbell CJ, Flinn SD, Boyer CB. 2003. Bacterial vaginosis in sexually experienced and non-sexually experienced young women entering the military. *Obstet Gynecol* 102:927–933. [https://doi.org/10.1016/s0029-7844\(03\)00858-5](https://doi.org/10.1016/s0029-7844(03)00858-5).
- Schwebke JR, Desmond R. 2007. Natural history of asymptomatic bacterial vaginosis in a high-risk group of women. *Sex Transm Dis* 34:876–877. <https://doi.org/10.1097/OLQ.0b013e318073bd82>.
- Klebanoff MA, Hillier SL, Nugent RP, MacPherson CA, Hauth JC, Carey JC, Harper M, Wapner RJ, Trout W, Moawad A, Leveno KJ, Miodovnik M, Sibai BM, Vandersten JP, Dombrowski MP, O'Sullivan MJ, Varner M, Langer O, National Institute of Child Health and Human Development Maternal-Fetal Medicine Units Network. 2005. Is bacterial vaginosis a stronger risk factor for preterm birth when it is diagnosed earlier in gestation? *Am J Obstet Gynecol* 192:470–477. <https://doi.org/10.1016/j.jog.2004.07.017>.
- Cohen CR, Lingappa JR, Baeten JM, Ngayo MO, Spiegel CA, Hong T, Donnell D, Celum C, Kapiga S, Delany S, Bukusi EA. 2012. Bacterial vaginosis associated with increased risk of female-to-male HIV-1 transmission: a prospective cohort analysis among African couples. *PLoS Med* 9:e1001251. <https://doi.org/10.1371/journal.pmed.1001251>.
- Wiesenfeld HC, Hillier SL, Krohn MA, Landers DV, Sweet RL. 2003. Bacterial vaginosis is a strong predictor of *Neisseria gonorrhoeae* and *Chlamydia trachomatis* infection. *Clin Infect Dis* 36:663–668. <https://doi.org/10.1086/367658>.
- Byun SW, Park YJ, Hur SY. 2016. Affirm VPIII microbial identification test can be used to detect *Gardnerella vaginalis*, *Candida albicans*, and *Trichomonas vaginalis* microbial infections in Korean women. *J Obstet Gynaecol Res* 42:422–426. <https://doi.org/10.1111/jog.12913>.
- Schwebke JR, Gaydos CA, Nyirjesy P, Paradis S, Kotsi S, Cooper CK. 2018. Diagnostic performance of a molecular test versus clinician assessment of vaginitis. *J Clin Microbiol* 56:e00252-18. <https://doi.org/10.1128/JCM.00252-18>.
- Petrikkos G, Makrilakis K, Pappas S. 2007. AFFIRM VPIII in the detection and identification of *Candida* species in vaginitis. *Int J Gynaecol Obstet* 96:39–40. <https://doi.org/10.1016/j.ijgo.2006.09.008>.