

## **Supporting Information**

Parthenolide inhibits ubiquitin-specific peptidase 7 (USP7), Wnt signaling,  
and colorectal cancer cell growth

**Xue Li<sup>1,2</sup>, Lingmei Kong<sup>1,3</sup>, Qihong Yang<sup>1,2</sup>, Aizhu Duan<sup>1,2</sup>, Xiaoman Ju<sup>1,2</sup>,  
Bicheng Cai<sup>1</sup>, Lin Chen<sup>1,2</sup>, Tao An<sup>1,3,\*</sup>, Yan Li<sup>1,3,\*</sup>**

### **Inventory of Supporting Information:**

**Figure S1:** Identification of PTL as a USP7 inhibitor.

**Figure S2:** Determination of PTL binding sites in USP7 using LC-MS/MS.

Figure S1





Figure. S1 Identification of PTL as a USP7 inhibitor. (A) PTL inhibited USP7 enzymatic activity in a concentration dependent manner in the Ub-AMC hydrolysis assay. (B) PTL inhibited USP7 enzymatic activity in a concentration dependent manner in the Ub-R110 hydrolysis assay. (C) Selectivity profile of PTL against a panel of DUBs. Screening was performed at a fixed concentration of 100  $\mu\text{M}$  PTL by Ubiquigent. Data were presented as mean  $\pm$  SD (n=2). (D) HCT116 cells were treated with indicated doses of PTL for 2 h and then labeled with Ub-PA. Individual DUBs were detected using specific antibodies.

## Figure S2

A



B

Sites of PTL modification on recombinant USP7	
adduction sites on USP7	peptide sequence
Cys90	LSESVLSPPC*FVR
Cys315	MKGTC*VEGTIPK
Cys334	MVSYIQC*K
Cys478	WCKFDDDWVSR
Cys488	C*TKEEAIEHNYGGHDDDLVSR
Cys510	HCTNAYM <sup>o</sup> LVYIR
Cys702	SLNYC*GHIYTPISCK
Cys721	DLLPVM* <sup>o</sup> C*DR
Cys799	VDVIFC*DK
Cys896	C*IWLNSQFR
Cys917	EEEEITLYPDKHGC*VR
Cys925	DLLEEC*KK
Cys961	IIGVHQEDELLEC*LSPATSR

\*Indicates modified residue

<sup>o</sup>Oxidation of some amino acids

Figure. S2 Determination of PTL binding sites in USP7 using LC-MS/MS. 3  $\mu$ g recombinant USP7 protein were incubated with 50  $\mu$ M PTL at room temperature for 2 h, and then separated by SDS-PAGE and visualized by coomassie brilliant blue. The protein-containing bands in gel were excised, followed by in-gel digestion and analyzed by LC-MS/MS on Nano LC-LTQ-Orbitrap XL (Thermo Finnigan, San Jose,

CA). Data processing was done using Proteome Discoverer 1.4 software. (A) Sequence coverage achieved in MS peptide mapping experiments: covered sequence (84%) is green-highlighted and ligand-bound cysteines are red-highlighted. (B) The cysteines and corresponding peptides covalently targeted by PTL in MS experiments.