

## Supplement 2

### Analysis of DSS samples for designer stimulants and opioids.

#### *Chemicals and reagents*

Propan-2-ol (Rotisol<sup>®</sup>, ≥ 99.95%, LC-MS grade) and formic acid (HCOOH, > 98%, p.a.) were obtained from Carl Roth (Karlsruhe, Germany), while acetonitrile (ACN, HPLC-Super gradient grade) was purchased from VWR International (Fontenay-sous-Bois, France). Ammonium formate (10 M, 99,995%) and hydrochloric acid (HCl; ≥ 37%, p.a.) were obtained from Sigma Aldrich (Steinheim, Germany). Methanol (MeOH, Chromasolv<sup>™</sup>, LC-MS, ≥ 99.9%) was purchased from Honeywell Riedel-de Haën (Seelze, Germany), while deionized water (H<sub>2</sub>O) was prepared using a Medica<sup>®</sup> Pro single high flow purification system from ELGA LabWater (Celle, Germany). Reference material of the opioids and designer stimulants included in the study (Table 1 and 2) were obtained from certified reference material retailers or obtained online as 'research chemical'. 'Research chemicals' were evaluated for identification and purity using GC-MS, LC-QTOF-MS, and NMR analysis before use.

#### *Solutions*

Working solutions (0.05 ng/mL and 0.5 ng/mL) containing all designer stimulants were prepared in MeOH. The methanolic working solution for opioids contained the opioids at two concentration levels. Level 1 was set to a concentration of 1 µg/mL, while level 2 was set to a concentration of 0.1 µg/mL (Table 2). The internal standard solution (IS) for designer stimulants (concentrations adjusted to yield medium peak intensities with respect to the calibrated concentration range) contained ethylone-D5 (0.1 µg/mL), methamphetamine-D5, methylone-D3, MDE-D5, MDPV-D8, PCP-D5 (0.25 µg/mL each), 2C-C-D6, 2C-I-D6, mephedrone-D3, MDMA-D5 (0.5 µg/mL each), amphetamine-D5, benzyloxy-piperazine-D7, mCPP-D8 (1.0 µg/mL each), norephedrine-D3 (2.5 µg/mL), and methcathinone-D3 (5.0 µg/mL). The IS for opioids contained nortilidine-D3, tapentadol-D3, tramadol-D3, fentanyl-D5, norfentanyl-D5, O-desmethyltramadol-D6, 7-hydroxymitragynine-D3, buprenorphine-D4, dihydrocodeine-D6, EDDP-D3 (perchlorate), hydromorphone-D3, methadone-D9, mitragynine-D3, morphine-D3, norbuprenorphine-D3, oxycodone-D3, oxycodone-D3, and sufentanil-D5 at 100 ng/mL each. Mobile phase A was deionized water with 1% acetonitrile (ACN), 0.1% formic acid, and 2 mM ammonium formate. Mobile phase B consisted of 0.1% formic acid and 2 mM ammonium formate in ACN.

#### *Instrumentation and methods*

The Shimadzu liquid chromatography (LC) system combined a CMB-20AC communications BUS module, a Prominence CTO-20AC column oven (set to 30 °C/40 °C for designer stimulants/opioids), a Nexera X2 SIL-30AC autosampler (10 µl/5 µl injection volume for designer stimulants/opioids, set temperature: 10 °C), a DGU-20A5R degasser unit, two Nexera X2 LC-30AD pumps, and a Prominence LC-2AD pump. Coupled to the LC system was a QTRAP<sup>®</sup> 5500 mass spectrometer (Sciex, Darmstadt, Germany). Ionization was performed in positive electrospray ionization (ESI+), while scheduled multi reaction monitoring (sMRM) was used for data acquisition. MS source parameters were set to the following values:

- a) Designer stimulants: Curtain gas 30 psi, ion source gas (1) 40 psi, ion source gas (2) 65 psi, ion spray voltage +4000 V, temperature 550 °C.
- b) Opioids: Curtain gas 30 psi, ion source gas (1) 60 psi, ion source gas (2) 70 psi, ion spray voltage +4500 V, temperature 500 °C.

Chromatographic separation for stimulants and opioids was achieved on a biphenyl column (100 x 2.1 mm, 2.6 µm particle size, Phenomenex, Aschaffenburg, Germany) with a corresponding guard column (SecurityGuard™ ULTRA Cartridges UHPLC Biphenyl for 2.1 mm ID columns, Phenomenex, Aschaffenburg, Germany). The chromatographic conditions were as follows:

- a) Designer stimulants: The gradient was set to a total flow rate of 0.3 mL/min and started at 5% mobile phase B. After 3 min at 5% mobile phase B was increased to 15% and kept at 15% for 3 min. Within 7 min the percentage of mobile phase B was further increased to 80%. For column cleaning 95% mobile phase B were reached after a total of 18.5 min. After 2.5 min at 95% starting conditions were restored within 0.5 min and kept for additional 4.5 min. For signal enhancement propan-2-ol was added post-column at 0.1 mL/min.
- b) Opioids: With a total flow rate of 0.4 mL/min the gradient started at 1% mobile phase B and these conditions were kept for 0.25 min. Within 1.25 min the gradient was increased to 3% mobile phase B. The gradient reached 22.5% mobile phase B after 3 min. Percentage of mobile phase B was further increased to 27.5% within 0.75 min. At the same time the total flow rate was set to 0.45 mL/min. After 8.5 min the gradient reached 32.5% mobile phase B which was further increased to 95% within 1 min. At that time the total flow rate was set to 0.6 mL/min. These conditions were kept for 1.5 min. Finally, starting conditions were restored within 0.1 min and kept for 1.4 min.

### *Sample preparation*

For calibration, 7 calibration samples consisting of 100 µL serum each were spiked with designer stimulants and opioids. The following concentration levels were obtained: 1.0, 5.0, 10, 20, 50, 70, 100 ng/mL (for designer stimulants and opioids level 1) and 10, 50, 100, 200, 500, 700, 1,000 ng/mL (for opioids level 2). The serum was then pipetted on empty filter paper spots (Whatman 903). After drying, 10 µL of both IS solutions were added and the serum spots were transferred in 5 mL Eppendorf tubes. For extraction, 3 mL of MeOH were added and the tubes were ultrasonicated for 15 min. The methanolic extract was transferred into glass vials and subsequently constricted to 100 µL at 40 °C using a gentle stream of nitrogen. After addition of 100 µL of a mixture (v/v; 3/1) of propan-2-ol and HCl (≥ 37%) the extract was evaporated to dryness and reconstituted in 100 µL mobile phase A/B (v/v; 99/1). For each calibration a blank serum sample (containing only IS) was prepared accordingly. Samples from the subtype C cluster as well as matching non-cluster control samples were spiked with 10 µL of each IS solution and treated as mentioned above.

### *Evaluation of process efficiency*

For evaluation of process efficiency two sets of spiked samples were prepared. Set 1 consisted of 5 blank serum samples (obtained from 5 individuals, 100 µL each) spiked with stimulants and opioids yielding a concentration of 10 ng/mL (100 ng/mL for opioids at level 1). Set 2 consisted of 100 µL MeOH spiked with opioids and designer stimulants and subsequently evaporated to dryness. Process efficiency was evaluated by comparison of absolute peak areas of set 1 and set 2.

**Table 1: Included designer stimulants and the mean process efficiency calculated for each analyte.**

Designer stimulants	Process efficiency [%]	Designer stimulants	Process efficiency [%]
2,5-DMA	72	Butylone	70
2/3/4-Fluoroamphetamine	52	Cathine	57
2/3/4-Fluoromethamphetamine	46	Cathinone	43
2/3/4-Methylmethcathinone	37	Cl-Pseudoephedrine	2.0
2-AI	52	DB-MDBP	74
2-FMC	n.e.	DBZP	83
2-MAPB	60	Desoxyipradrol	50
2-MeO-Ketamine	56	Dimethocaine	0
3,4-CTMP	2.0	Dimethylcathinone	40
3,4-DMA	88	Diphenidine	50
3,4-DMMC	50	Ephedrine	63
3,4-MeO-a-PHP	31	Etaqualone	51
3/4-Me-Buphedrone	49	Ethcathinone	41
3-Cl-Methcathinone	19	Ethylamphetamine	52
3-FMC	20	EthylNaphtidate	50
3-F-Phenetrazine	60	Ethylone	72
3-FPM	61	Ethylphenidate	47
3-MeO-MC	39	Isopentdrone	19
4-APDB	n.e.	Isopenmetrazine	64
4-CAB	61	Isopropylphendiate	54
4-Cl-Methamphetamine	57	Ketamine	72
4-Cl-Methcathinone	26	mCPP	72
4-Cl-PVP	10	MDA	65
4-Ethylethcathinone	56	MDAI	66
4-Ethylmethcathinone	45	MDAT	91
4-F-a-PBP	33	MDMA	70
4-F-a-PVP	33	MDPBP	71
4-F-BF	40	MDPHP	38
4-F-Buphedrone	45	MDPPP	72
4-F-Ethylphenidate	49	MDPV	55
4-F-IPV	40	MEAI	77
4-FMC	33	Mebroqualone	58
4-F-Methylphenidate	73	Mephtramine	28
4-F-PV8	24	Methamnetamine	40
4-F-PV9	10	Methamphetamine	53
4-Me-Methylphenidate	0	Methcathinone	38
4-Me-N-ethylnorpentdrone	32	Methiopropamine	49
4-MeO-BF	39	Methoxetamine	73
4-MeO-PV9	28	Methoxphenidine	41
4-MeO-PVP	43	MethylNaphtidate	60
4-Me-Pentdrone	56	Methylone	65
4-Me-Phenmetrazine	65	Methylphenidate	70
4-Me-PHP	41	MDE	41
4-MMA	54	NEB-Indane-analog	42
4-MTA	815	N-Ethylpentylone	75
5-APB	57	N-Ethylphenmetrazole	52
5-APDI	63	Nitracaine	0
5-BPDI	36	N-Me-2AI	54
5-DBFPV	33	N-Me-bk-MMDA-2	57
5-IAI	77	Norephedrine	46
5-IT	0	NRG-3	38
5-MAPDB	70	PCP	42
5-MBPB	57	Pentdrone	43
5-PPDI	48	Pentylone	70
6/5/4-EAPB	64	Phenetrazine	60
6/5/4-MAPB	59	PMA	66
6-APB	61	PMMA	64
6-APDB	75	Propylphendiate	59
7-APDB	66	PV9	24
a-Naphyrone	29	Pyrovalerone	27
a-ET	67	Ritalinic acid	78
AH-792	43	TFMPP	66
ALEPH-2	415	TMA (3,4,5)	101
ALEPH-4	224	TMA-2 (2,4,5)	97
a-Me-AHP	47	TMA-6 (2,4,6)	89
Amphetamine	54		
a-PAVP	47		
a-PHP	26		
a-PNP	21		
a-PVP	48		
Benzylpiperazine	27		
bk-MDDMA	73		
Buphedrone	45		

**Table 2: Included opioids, their concentration level, and the mean process efficiency calculated for each analyte.**

Opioids	Opioid level	Process efficiency [%]
2-F-Isobutyrfentanyl	2	64
3,4-Methylenedioxy-U-47700	1	59
3-Methylfuranylfentanyl	2	68
4-ANPP	2	38
4-Cl-Isobutyrfentanyl	2	72
4-F-Butyrfentanyl	2	66
4-MeO-Butyrfentanyl	2	70
7-Hydroxymitragynine	1	0
Acetylfentanyl	2	51
Acryloylfentanyl	2	53
AH-792	1	66
Alfentanil	2	53
Benzodioxolfentanyl	2	76
Benzylfentanyl	2	54
Buprenorphine	2	0
Butyrfentanyl	2	67
Carfentanil	2	66
Codeine	1	85
Cyclopentylfentanyl	2	79
Cyclopropylfentanyl	2	60
Desomorphine	1	63
Despropionyl-o-F-fentanyl	2	24
Dextrometorphan	1	57
Dihydrocodeine	1	77
Dihydromorphine	1	62
EDDP	1	41
Fentanyl	2	65
Furanylethylfentanyl	2	40
Furanylfentanyl	2	68
Furanylfentanyl 3-CA-isomer	2	58
Hydrocodone	1	49
Hydromorphone	1	46
Loperamide	1	76
MeO-Acetylfentanyl	2	49
Meptazinol	1	78
Methadone	1	74
Mitragynine	1	30
Morphine	1	74
MT-45	1	69
M-U-47700	1	64
N,N-Bidesmethyl-U-47700	1	70
Nalbuphine	1	61
Naloxone	1	39
Naltrexone	1	45
N-Desmethyltapentadol	1	57
N-Desmethyl-U-47700	1	52
Norbuprenorphine	2	0
Norcodeine	1	92
Norfentanyl	2	50
Normorphine	1	64
Noroxycodone	1	36
Nortilidine	1	66
Noscapine	1	62
Ocfentanil	2	0
O-Desmethyltramadol	1	65
Oxycodone	1	47
Oxymorphone	1	44
Papaverine	1	58
Pentazocine	1	57
Pethidine	1	55
Pholcodine	1	102
Remifentanil	2	67
Sufentanil	2	73
Tapentadol	1	63
THF-fentanyl	2	60
Tilidine	1	56
Tramadol	1	69
U-47700	1	66
U-48800	1	57
U-49900	1	68
Valeryl fentanyl	2	61
W-18	2	0