

## **Generation of Small RNA-Modulated Exosome Mimetics for Bone Regeneration**

*Jiabing Fan<sup>1,2#</sup>, Chung-Sung Lee<sup>1#</sup>, Soyon Kim<sup>1</sup>, Chen Chen<sup>1</sup>, Tara L. Aghaloo<sup>3</sup>, Min Lee<sup>1,2,4\*</sup>*

*1. Division of Advanced Prosthodontics, School of Dentistry, University of California, Los Angeles, California, 90095, USA.*

*2. Weintraub Center for Reconstructive Biotechnology, School of Dentistry, University of California, Los Angeles, California, 90095, USA.*

*3. Division of Diagnostic and Surgical Sciences, School of Dentistry, University of California, Los Angeles, California, 90095, USA.*

*4. Department of Bioengineering, University of California, Los Angeles, California, 90095, USA.*

*# These authors contributed equally to this work*

*\* To whom correspondence should be addressed:*

*Min Lee, PhD*

*Professor*

*Division of Advanced Prosthodontics*

*Weintraub Center for Reconstructive Biotechnology*

*Department of Bioengineering*


*University of California, Los Angeles*

*Phone: 310-825-6674*


*Fax: (310) 825-6345*

*Email: leemin@ucla.edu*


## Supplemental Materials:



**Figure S1: The yields of EMs and EXOs from hMSCs with OM or OMN.** (A): The yields between EM-OM and EXO-OM were measured by the protein concentration. (B): The yields between EM-OMN and EXO-OMN were measured by the protein concentration. \*\*\*p<0.001. EMs, exosome mimetics; EXOs, exosomes; OM, osteogenic medium; OMN, osteogenic medium and noggin suppression; EXO-OM, EXOs from hMSCs with osteogenic medium; EM-OM, EMs from hMSCs with osteogenic medium; EXO-OMN; EXOs from noggin-suppressed hMSCs with osteogenic medium; EM-OMN, EMs from noggin-suppressed hMSCs with osteogenic medium.



**Figure S2: Western-blot analysis of CD63 and cellular protein in EMs and hMSCs.** A western-blot assay was conducted to measure expression of CD63 and cellular protein including Calnexin, Cytochrome c and GM-130 in EMs and parent hMSCs (whole cell lysate). EMs, exosome mimetics.



**Figure S3: Biocompatibility of EM-OMN or EM-OM-laden MeGC hydrogel.** (A,B): Cell viability and proliferation was respectively examined by Live/Dead (A) and AlamarBlue (B) assay at indicated times. Scale bar = 100  $\mu$ m. Ctl, Control; L-EM-OM, hydrogel encapsulated with low dose EM-OM; H-EM-OM, hydrogel encapsulated with high dose EM-OM; L-EM-OMN, hydrogel encapsulated with low dose EM-OMN; H-EM-OMN, hydrogel encapsulated with high dose EM-OMN.

**1. Table S1: Sequence of primers for real-time PCR assay.**

<b>Gene</b>	<b>Forward</b>	<b>Reverse</b>
<i>Runx2</i>	TGGTTACTGTCATGGCGGGTA	TCTCAGATCGTTGAACCTTGCTA
<i>Osterix</i>	CCTCTGCGGGACTCAACAAC	AGCCCATTAGTGCTTGTAAGG
<i>ALP</i>	ACTGGTACTCAGACAACGAGAT	ACGTCAATGTCCCTGATGTTATG
<i>Osteocalcin</i>	CACTCCTCGCCCTATTGGC	CCCTCCTGCTTGGACACAAAG
<i>BMPRII</i>	AGATGACCAGGGAGAAACCAC	CAACATTCTATTGTCCGGCGTA
<i>ID1</i>	ACGGCTGTTACTCACGCCTC	GCTGGAGAATCTCCACCTTGC
<i>GAPDH</i>	TGTGGGCATCAATGGATTTGG	ACACCATGTATTCCGGGTCAAT