

Supplemental Protocols

Protocol 1: Recipes for chemical solutions

Longmire's Solution (Longmire, et al., 1991, 1997; Renshaw, et al., 2015; Williams, et al., 2016)			
Recipe is for total volume: 500 mL			
Stock Chemical Name	[Stock]	[Final]	Volume/Amount
Tris Base pH 8	1M	0.1M	50 mL
Disodium-EDTA pH 8	0.5M	0.1M	100 mL
Sodium chloride (NaCl)	5M	0.01M	1 mL
Sodium dodecyl sulfate (SDS)	10%	0.5%	25 mL
Water			324 mL
Recipe notes:			
<ul style="list-style-type: none">Sodium azide can be added to prevent microbial growth (Williams, et al., 2016)Keep all waste in a properly labelled Hazardous Waste container.			
Lysis Buffer I (Lever, et al., 2015)			
Recipe is for total volume: 250 mL			
Stock Chemical Name	[Stock]	[Final]	Volume/Amount
Tris Base pH 8	1M	0.03M	7.5 mL
Disodium-EDTA pH 8	0.5M	0.03M	15 mL
Guanidium hydrochloride	1M	0.8M	200 mL
Triton X-100	10%	0.5%	12.5 mL
Water			15 mL
Recipe notes:			
<ul style="list-style-type: none">Adjust final solution to pH [10.0] using NaOH pellets ($n = 2$) + ~1mL [1M] NaOHKeep all waste in a properly labelled Hazardous Waste container.			
PEG8000-NaCl Precipitation Solution (Cseke, et al., 2011 "Handbook of Molecular and Cellular Methods in Biology and Medicine, Third Edition"; Lever, et al., 2015)			
Recipe is for total volume: 100 mL			
Stock Chemical Name	[Stock]	[Final]	Volume/Amount
Polyethylene glycol 8000	100%	30%	30 g
Sodium chloride (NaCl)	5M	1.6M	32 mL
Water			68 mL
Recipe notes:			
<ul style="list-style-type: none">No pH adjustment required (unadjusted solution is acidic).Keep all waste in a properly labelled Hazardous Waste container.			

Protocol 2: Preparation of DNase-treated glycogen (200 mg/mL)

Step 1: Prepare DNase reaction buffer reagents

- 1) Make 25mL of [1M] Tris (MW 121.14; CAS no 77-86-1)
 - a. 3.0275g into 20 mL DNA-free water then adjust up to 25 mL final
 - b. NOTE: If solution has a yellow color get better quality Tris
- 2) Make 25mL of [0.25M] MgCl₂ (MW = 95.21;CAS no 7786-30-3)
 - a. 0.5951g into 20 mL DNA-free water then adjust up to 25 mL final
- 3) Make 25mL of [0.05M] CaCl₂ (MW = 110.98; CAS no 10043-52-4)
 - a. 0.138725g into 20 mL DNA-free water then adjust up to 25 mL final

Step 2: Prepare reaction buffer (100mL):

- 1) Add 80mL DNA-free water to 100 mL DNA-free bottle
- 2) Add 1000 μ L [1M] Tris (10 mM final)
- 3) Add 1000 μ L [0.25M] MgCl₂ (2.5 mM final)
- 4) Add 1000 μ L [0.05M] CaCl₂ (0.5 mM final)
- 5) Adjust solution to pH 7.6
 - a. Add 5-10 DROPS of 5% HCl then confirm pH
 - b. Add 3-5 DROPS [1M] NaOH then confirm pH
- 6) Adjust final volume up to 100 mL with MilliQ water
- 7) Adjust final solution to pH 7.6
 - a. Add 3-4 DROPS [1M] NaOH then confirm pH

Step 3: Add 5g oyster glycogen (Sigma-Aldrich cat no. G8751-5G) to empty 50mL tube

- 1) Add reaction buffer up to 25 mL (approx. 20mL)
- 2) Agitate gently until glycogen fully dissolved then adjust volume to 25mL

Step 4: Add 25 mg DNase (1 mg/mL final; Sigma-Aldrich cat no. 10104159001) to solution

- 1) Mix by inversion or gentle agitation (note: do not vortex solution containing DNase)

Step 5: Start DNase reaction by incubating (37°C, 15 min) with gentle mixing

Step 6: Stop DNase reaction by adding 250 μ L [0.5M] EDTA and incubating (75°C, 10min)

Step 7: Divide DNase-treated glycogen into 1mL aliquots and freeze ($\leq 20^{\circ}\text{C}$)

•Note: Remaining reaction buffer can be aliquoted ($\approx 20\text{mL}$) and frozen ($\leq 20^{\circ}\text{C}$)

Reagent	Price (AUD)	Unit	\$ per 25mL (200mg/mL)	\$ per sample [Minimum]	\$ per sample [Maximum]	# Samples @ [Minimum]	# Samples @ [Maximum]
Glycogen	\$131.00	5g	\$131.00	\$0.008	\$0.031	16667	4167
DNase	\$95.50	100mg	\$23.88	\$0.001	\$0.006		
CaCl	\$86.00	100g	\$0.02	\$0.000	\$0.000		
MgCl	\$78.50	100g	\$0.02	\$0.000	\$0.000		
[1M] Tris-HCl	\$2.25	25mL	\$0.09	\$0.000	\$0.000		
Labor (hr)	\$100	2	\$200	\$0.012	\$0.048		
TOTAL			\$355.00	\$0.02	\$0.09		

Protocol 3: Silica purification

Silica matrix preparation:

- 1) Mix 5g Silicon dioxide (Sigma #5631) with 50mL MilliQ water
- 2) Allow to settle overnight¹
- 3) Remove supernatant (containing fine silica particles) and add 50mL MilliQ water
- 4) Allow to settle overnight¹
- 5) Remove supernatant (containing fine silica particles) and add 50mL MilliQ water
 - a. Silica matrix ready @ [100mg/mL]²
 - b. Store at room temp for ≤ 1 year²

Silica-based purification of eDNA samples:

- 1) Add 1 volume of [3M] NaCl and 20 μ L² silica matrix
- 2) Vortex and incubate at room temperature for 15min
- 3) Centrifuge 16,000 x g for 10 sec² and discard supernatant
- 4) Wash twice with 500 μ L washing solution (add and vigorously vortex)
 - a. 50% ethanol, 10 mM Tris-HCl pH 7.5, 100 mM NaCl, 1 mM EDTA^{1,2}
- 5) Elute DNA by resuspending silica matrix with MilliQ water equal to starting volume
- 6) Incubate at room temperature ($\geq 20^{\circ}\text{C}$) for ≥ 30 min then vortex briefly to resuspend
- 7) Transfer slurry onto Econospin™ column (Epoch Life Sciences, Texas USA; #1920)[^]
- 8) Centrifuge 20,000 x g for 10 sec at room temp (see Note in Step 3a)
- 9) Purified samples suitable for qPCR-based detection of target species*

[^]eDNA samples expected to contain trace copies (e.g., < 10) or intended for quantification (e.g., biomass estimation) should utilize Zymo Inhibitor Removal columns instead of Econospin™ columns (“Silica-Zymo”)

*“Silica-Epoch” purified samples can exhibit 1-2 qPCR cycle shift due to trace silica carry-over (species-dependent).

¹Huanca-Mamani, W., Rivera-Cabello, D., & Maita-Maita, J. (2015). A simple, fast, and inexpensive CTAB-PVP-Silica based method for genomic DNA isolation from single, small insect larvae and pupae. *Genetics and Molecular Research*, 14(3), 8001-8007).

²Li, J. F., Li, L., & Sheen, J. (2010). Protocol: a rapid and economical procedure for purification of plasmid or plant DNA with diverse applications in plant biology. *Plant Methods*, 6(1), 1.