




American Society for Microbiology Provides 2020 Guidelines for Detection and Identification of Group B *Streptococcus*

 Laura Filkins,^a Jocelyn R. Hauser,^b Barbara Robinson-Dunn,^c Robert Tibbetts,^d Bobby L. Boyanton, Jr.,^e  Paula Revell^{f,g} on behalf of the American Society for Microbiology Clinical and Public Health Microbiology Committee, Subcommittee on Laboratory Practices

^aDepartment of Pathology, University of Texas Southwestern Medical Center, Dallas, Texas, USA

^bPublic Health Laboratory, DC Department of Forensic Sciences, Washington, DC, USA

^cBeaumont Health System, Royal Oak, Michigan, USA

^dDepartment of Pathology, Henry Ford Health System, Detroit, Michigan, USA

^eDepartment of Pathology, Arkansas Children's Hospital, Little Rock, Arkansas, USA

^fDepartment of Pediatrics, Baylor College of Medicine/Texas Children's Hospital, Houston, Texas, USA

^gDepartment Pathology and Immunology, Baylor College of Medicine/Texas Children's Hospital, Houston, Texas, USA

KEYWORDS EOD, GBS, early-onset disease, group B *Streptococcus*, intrapartum prophylaxis, screening

Maternal colonization with group B *Streptococcus* (GBS) is a primary risk factor for early-onset-disease (EOD) GBS infection in infants, and intrapartum prophylaxis reduces neonatal infection. The guidelines for GBS screening and intrapartum prophylaxis were last revised in 2010. In 2019, the CDC transferred stewardship of the GBS laboratory testing recommendations to the American Society for Microbiology. For the complete laboratory guidelines for GBS screening specimen collection and handling, organism detection and identification, and antimicrobial susceptibility testing (AST), we direct readers to the 2020 *Guidelines for the Detection and Identification of Group B Streptococcus* (1). Here, we highlight updates found in the 2020 GBS laboratory guidelines.

In 2019, the American College of Obstetricians and Gynecologists recommended antepartum screening for GBS during 36 0/7 to 37 6/7 weeks of gestation, a change from 35 to 37 weeks (2). Seven percent of U.S. births occur at or after 41 weeks (3). Shifting specimen collection 1 week later lengthens the 5-week predictive window of GBS screening tests to 41 0/7 weeks.

Best laboratory practices for GBS screening have not changed substantially since 2010. In the 2020 guidelines, culture remains the backbone of GBS detection. Methods for the identification of candidate GBS isolates from culture have expanded over the last decade. While biochemical testing and latex agglutination are acceptable for GBS identification, protein-based identification is ideal for laboratories with a matrix-assisted laser desorption ionization–time of flight mass spectrometry (MALDI-TOF) instrument. Advantageously, MALDI-TOF provides rapid and specific identification of *Streptococcus agalactiae* and enables its differentiation from group B antigen-agglutinating strains of *Streptococcus halichoeri* and *Streptococcus pseudoporcinus*, whose clinical significance in EOD is poorly understood (4). However, these identification methods, including MALDI-TOF, require an isolated organism, increasing the time to results compared to identification directly from enrichment broth.

GBS is predictably susceptible to penicillin and ampicillin, the primary antibiotics for intrapartum prophylaxis. For women with severe penicillin allergy, clindamycin is the preferred agent. However, *Streptococcus agalactiae* rates of resistance to macrolides and lincosamides can exceed 40%, necessitating AST on GBS isolated from women who

Citation Filkins L, Hauser JR, Robinson-Dunn B, Tibbetts R, Boyanton BL, Jr, Revell P, on behalf of the American Society for Microbiology Clinical and Public Health Microbiology Committee, Subcommittee on Laboratory Practices. 2021. American Society for Microbiology provides 2020 guidelines for detection and identification of group B *Streptococcus*. *J Clin Microbiol* 59:e01230-20. <https://doi.org/10.1128/JCM.01230-20>.

Editor Sandra S. Richter, bioMérieux

Copyright © 2020 American Society for Microbiology. All Rights Reserved.

Address correspondence to Laura Filkins, Laura.Filkins@utsouthwestern.edu.

Accepted manuscript posted online 28 October 2020

Published 17 December 2020

are severely allergic to penicillin (5). Unless AST is automatically performed on all isolates, regardless of penicillin allergy status, a mechanism for clear communication to the laboratory must be established to identify which patients' isolates require AST.

Nucleic acid amplification tests (NAATs) performed after organism enrichment in broth medium are acceptable for the detection and identification of GBS. Indeed, recent studies report comparable or enhanced GBS sensitivity by NAATs after enrichment compared to culture (6–8). However, we recommend that all laboratories performing NAATs maintain culture procedures to permit AST for appropriate isolates.

Finally, the 2020 laboratory guidelines do not endorse direct-from-specimen GBS NAATs. Without enrichment, GBS NAATs suffer from poor sensitivity and a low negative predictive value (9, 10). Recent reports suggest that direct-from-specimen NAATs may inform additional prophylaxis in some situations but should not replace antenatal testing and risk-based prophylaxis (11). Evidence-based studies are required to determine when, how, and if direct-from-specimen NAATs are clinically beneficial. Currently, best-practice laboratory testing for GBS from screening specimens necessitates broth enrichment prior to NAAT or subculture.

REFERENCES

1. Filkins L, Hauser J, Robinson-Dunn B, Tibbetts R, Boyanton B, Revell P, on behalf of the American Society for Microbiology Clinical and Public Health Microbiology Committee, Subcommittee on Laboratory Practices. 2020. Guidelines for detection and identification of group B *Streptococcus*. American Society for Microbiology, Washington, DC. <https://asm.org/Guideline/Guidelines-for-the-Detection-and-Identification-of>.
2. Committee on Obstetric Practice. 2019. Prevention of early-onset group B streptococcal disease in newborns: ACOG committee opinion, number 782. *Obstet Gynecol* 134:e19–e40. <https://doi.org/10.1097/AOG.0000000000003334>.
3. Martin JA, Osterman MJK, Kirmeyer SE, Gregory ECW. 2015. Measuring gestational age in vital statistics data: transitioning to the obstetric estimate. *Natl Vital Stat Rep* 64:1–20.
4. Salimnia H, Robinson-Dunn B, Gundel A, Campbell A, Mitchell R, Taylor M, Fairfax MR. 2019. Suggested modifications to improve the sensitivity and specificity of the 2010 CDC-recommended routine *Streptococcus agalactiae* screening culture for pregnant women. *J Clin Microbiol* 57:e00446–19. <https://doi.org/10.1128/JCM.00446-19>.
5. Metcalf BJ, Chochua S, Gertz RE, Jr, Hawkins PA, Ricaldi J, Li Z, Walker H, Tran T, Rivers J, Mathis S, Jackson D, McGee L, Beall B, Active Bacterial Core Surveillance Team. 2017. Short-read whole genome sequencing for determination of antimicrobial resistance mechanisms and capsular serotypes of current invasive *Streptococcus agalactiae* recovered in the USA. *Clin Microbiol Infect* 23:574.e7–574.e14. <https://doi.org/10.1016/j.cmi.2017.02.021>.
6. Shin JH, Pride DT. 2019. Comparison of three nucleic acid amplification tests and culture for detection of group B *Streptococcus* from enrichment broth. *J Clin Microbiol* 57:e01958–18. <https://doi.org/10.1128/JCM.01958-18>.
7. Hernandez DR, Wolk DM, Walker KL, Young S, Dunn R, Dunbar SA, Rao A. 2018. Multicenter diagnostic accuracy evaluation of the Luminex Aries real-time PCR assay for group B *Streptococcus* detection in Lim broth-enriched samples. *J Clin Microbiol* 56:e01768–17. <https://doi.org/10.1128/JCM.01768-17>.
8. Buchan BW, Faron ML, Fuller D, Davis TE, Mayne D, Ledebroer NA. 2015. Multicenter clinical evaluation of the Xpert GBS LB assay for detection of group B *Streptococcus* in prenatal screening specimens. *J Clin Microbiol* 53:443–448. <https://doi.org/10.1128/JCM.02598-14>.
9. Alfa MJ, Sepehri S, De Gagne P, Helawa M, Sandhu G, Harding GKM. 2010. Real-time PCR assay provides reliable assessment of intrapartum carriage of group B *Streptococcus*. *J Clin Microbiol* 48:3095–3099. <https://doi.org/10.1128/JCM.00594-10>.
10. Plainvert C, El Alaoui F, Tazi A, Joubrel C, Anselem O, Ballon M, Frigo A, Branger C, Mandelbrot L, Goffinet F, Poyart C. 2018. Intrapartum group B *Streptococcus* screening in the labor ward by Xpert GBS real-time PCR. *Eur J Clin Microbiol Infect Dis* 37:265–270. <https://doi.org/10.1007/s10096-017-3125-2>.
11. Ramesh Babu S, McDermott R, Farooq I, Le Blanc D, Ferguson W, McCallion N, Drew R, Eogan M. 2018. Screening for group B *Streptococcus* (GBS) at labour onset using PCR: accuracy and potential impact—a pilot study. *J Obstet Gynaecol* 38:49–54. <https://doi.org/10.1080/01443615.2017.1328490>.