



DATA NOTE

The EU one-stop-shop collection of publicly available information on COVID-19 *in vitro* diagnostic medical devices [version 1; peer review: 2 approved]

Mauro Petrillo ¹, Maddalena Querci¹, Olga Tkachenko², Ioana-Raluca Siska², Enrico Ben¹, Alexandre Angers-Loustau ^{3,4}, Alessia Bogni¹, Antonino Brunetto ⁵, Marco Fabbri¹, Linda Garland⁶, Antoon Lievens⁶, Amalia Munoz⁶, Valentina Paracchini¹, Danilo Pietretti¹, Antonio Puertas-Gallardo ¹, Barbara Raffael¹, Eleonora Sarno⁷, Virginie Tregogat⁶, Fabrizio Zaro⁸, Guy Van den Eede⁶

¹European Commission, Joint Research Centre (JRC), Ispra, Italy

²European Commission, Directorate General for Health and Food Safety (SANTE), Brussels, Belgium

³Past affiliation: European Commission, Joint Research Centre (JRC), Ispra, Italy

⁴Current affiliation: European Commission, European Publication Office, Luxembourg, Luxembourg

⁵GFT Italia s.r.l., Milan, Italy

⁶European Commission, Joint Research Centre (JRC), Geel, Belgium

⁷Past affiliation (until 30-06-2020): European Commission, Joint Research Centre (JRC), Ispra, Italy

⁸Engineering Ingegneria Informatica S.p.A, Taino, Italy

V1 First published: 03 Nov 2020, 9:1296
<https://doi.org/10.12688/f1000research.27308.1>

Latest published: 03 Nov 2020, 9:1296
<https://doi.org/10.12688/f1000research.27308.1>

Abstract

The *JRC COVID-19 In Vitro Diagnostic Devices and Test Methods Database*, aimed to collect in a single place all publicly available information on performance of CE-marked *in vitro* diagnostic medical devices (IVDs) as well as *in house* laboratory-developed devices and related test methods for COVID-19, is here presented. The database, manually curated and regularly updated, has been developed as a follow-up to the Communication from the European Commission “Guidelines on *in vitro* diagnostic tests and their performance” of 15 April 2020 and is freely accessible at <https://covid-19-diagnostics.jrc.ec.europa.eu/>.

Keywords

In-vitro diagnostics, covid-19, sars-cov-2, detection method

Open Peer Review

Reviewer Status  


Invited Reviewers

	1	2
version 1		
03 Nov 2020	report	report

1. **Angelo Facchiano** , National Research Council, Avellino, Italy

2. **Christopher Viljoen** , University of the Free State, Bloemfontein, South Africa

Any reports and responses or comments on the article can be found at the end of the article.



This article is included in the [Disease Outbreaks](#) gateway.



This article is included in the [Science Policy Research](#) gateway.



This article is included in the [Coronavirus](#) collection.

Corresponding author: Mauro Petrillo (mauro.petrillo@ec.europa.eu)

Author roles: **Petrillo M:** Conceptualization, Writing – Original Draft Preparation; **Querci M:** Conceptualization, Data Curation, Project Administration, Supervision; **Tkachenko O:** Conceptualization, Writing – Review & Editing; **Siska IR:** Conceptualization, Writing – Review & Editing; **Ben E:** Data Curation, Software, Supervision, Visualization, Writing – Review & Editing; **Angers-Loustau A:** Data Curation, Software; **Bogni A:** Data Curation; **Brunetto A:** Software, Visualization; **Fabbri M:** Data Curation, Formal Analysis, Validation; **Garlant L:** Data Curation; **Lievens A:** Data Curation; **Munoz A:** Data Curation, Formal Analysis, Validation; **Paracchini V:** Data Curation; **Pietretti D:** Data Curation; **Puertas-Gallardo A:** Data Curation, Software; **Raffael B:** Data Curation, Validation; **Sarno E:** Data Curation; **Tregogat V:** Data Curation, Validation; **Zaro F:** Software, Visualization; **Van den Eede G:** Conceptualization, Data Curation, Project Administration, Supervision, Writing – Review & Editing

Competing interests: No competing interests were disclosed.

Grant information: This work was supported by the Joint Research Centre.

Copyright: © 2020 Petrillo M *et al.* This is an open access article distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

How to cite this article: Petrillo M, Querci M, Tkachenko O *et al.* **The EU one-stop-shop collection of publicly available information on COVID-19 *in vitro* diagnostic medical devices [version 1; peer review: 2 approved]** F1000Research 2020, 9:1296 <https://doi.org/10.12688/f1000research.27308.1>

First published: 03 Nov 2020, 9:1296 <https://doi.org/10.12688/f1000research.27308.1>

Introduction

The Communication from the Commission “Guidelines on *in vitro* diagnostic tests and their performance”¹, published on 15 April 2020, states the following under the *Further Actions Needed* section: “*The Commission, supported by the ECDC, health technology assessment experts and in vitro diagnostics competent authorities, will assist Member States with a centralised overview of available information on test performance and act as a single point of contact for management of this information. Taking stock of the state of the art on a regular basis will support Member States’ informed decisions on national testing strategies, as well as support the continuous development of devices by manufacturers.*”

As an initial step in collecting performance information of devices and *in house* methods, to address the above need, European Commission services (Directorate-General for Health and Food Safety [DG SANTE], Directorate-General Joint Research Centre [DG JRC], Directorate-General for Research and Innovation [DG RTD]), together with the European Centre for Disease Prevention and Control (ECDC), several experts from *in vitro* diagnostics competent authorities and from the European Network for Health Technology Assessment (EUnetHTA)¹, published the working document “Current performance of COVID-19 test methods and devices and proposed performance criteria”² on 16 April 2020.

The JRC capitalised on its expertise in knowledge management to conduct the literature review as part of this work and, as a follow-up action to the need identified in the Communication, committed to make the information broadly accessible and to update the compilation as new data become available.

The outcome of these actions is the **JRC COVID-19 In Vitro Diagnostic Devices and Test Methods Database** presented here, a single place collection of all publicly available information on performance of CE-marked *in vitro* diagnostic medical devices (IVDs) as well as *in house* laboratory-developed devices and related test methods for COVID-19. The database is freely accessible at <https://covid-19-diagnostics.jrc.ec.europa.eu/>.

Methods

Information retrieval

The initial information is gathered following the strategy for documentation on test methods and devices indicated in Section 3 “Methodology used” of the *Current performance of COVID-19 test methods and devices and proposed performance criteria - Working document of Commission services*².

¹ EUnetHTA (European network for Health Technology Assessment) Joint Action 3 (2016 – 2021) is the scientific and technical component of EU cooperation on HTA. The Joint Action is co-funded by the EU Health Programme and includes government appointed organisations and a large number of relevant regional agencies and not-for-profit organisations that produce or contribute to HTA in Europe.

Verification of information

The retrieved information is manually verified and curated. In particular:

- For CE-marked devices, only the information that the manufacturer has chosen to make publicly available is included in the database. Full information on the manufacturer’s performance evaluation of the device is contained in the technical documentation required by EU legislation³. As manufacturer technical documentation is usually not publicly available and is therefore not included, a form is available to stakeholders to update and integrate the information: manufacturers are invited to submit performance information on new devices, which are not yet listed in the database, or to provide data not available to the authors at the time of the last update. The submitted information, once verified against the source provider, is taken into consideration for updating the database.
- Performance details (as retrieved from manufacturers’ web pages) are provided only for devices commercially available with CE-IVD mark. Products labelled as for ‘research use only’ or ‘under development’ as well as products fulfilling regulatory frameworks other than the one in place in the EU are listed for information only. The correctness of information, such as performance data of the listed devices, has not been confirmed by checking raw experimental data or full technical documentation of the manufacturer nor by own laboratory verification or by any clinical validation studies.
- A team of JRC experts regularly takes care of updating the Scientific Literature section, which includes scientific articles reporting about the use and performances of devices and related test methods. This section regards also the so-called *in house*² or laboratory developed devices, which are used in healthcare institutions but are not commercially available. These are not included in the list of devices, as they are not easily identifiable by name, but scientific publications on such devices and related test methods are included in the section.
- Information on publicly available results of validation studies and on publications by national/regional health technology assessment bodies and Joint Action EUnetHTA⁴ is included. These reports are produced in response to requests from national authorities and support national policy and decision-making on testing for COVID-19.

In silico NAAT methods simulations

In silico simulations of NAAT methods, on available high-quality and full length sequenced SARS-CoV-2 genomes

² The *in house* devices are those referred to in Art.1(5) of Directive 98/79/EC: “*devices manufactured and used only within the same health institution and on the premises of their manufacture or used on premises in the immediate vicinity without having been transferred to another legal entity*”. They are exempt from requirements of Directive 98/79/EC but may be subject to national law.

(provided as collaboration with GISAID⁵) are performed by using ThernucleotideBLAST⁶. The tool is run on each assembled genome with default parameters, except for the following ones: $-e=30$, $-E=20$, $--max-gap=4$, $--max-mismatch=4$. Putative not redundant target amplicons are then extracted to build the final target dataset that is added to the database.

At the time of writing, high-quality (i.e. <1% Ns and <0.05% unique amino acid mutations) and full length (i.e. >29,000 bp) sequenced SARS-CoV-2 genomes have been manually downloaded from the GISAID website. To simplify the sharing of information and the related analysis, a formal agreement is under finalisation between GISAID and JRC.

Results

The database is structured in a way aimed to:

- facilitate sharing of information among researchers, in line with the principle of FAIR sharing;
- link the devices, and their features, to scientific articles that reported their use and performance, in order to have a clear, transparent and fully open data source where users can look at and make the right choice in the selection of the device or method to use;
- verify and monitor if natural occurring SARS-CoV-2 mutations give rise to “false negative”.

Data records are currently organised and available in three main sections:

1. **COVID-19 *in vitro* diagnostic medical devices** - This section gives access to publicly available *in vitro* diagnostic medical devices for COVID-19 and is being updated periodically. Additional performance (as retrieved from manufacturers’ web pages) is provided only for devices or kits commercially available with CE-IVD mark.

Direct link to this dataset page, where data can be downloaded in CSV format, is: <https://covid-19-diagnostics.jrc.ec.europa.eu/devices>.

2. **Scientific literature on COVID-19 test methods and devices** - This section allows browsing on performance of test methods and devices for COVID-19 diagnostics retrieved from selected scientific articles and is being updated periodically.

Direct link to this dataset page, where data can be downloaded in CSV format, is: <https://covid-19-diagnostics.jrc.ec.europa.eu/literature>.

3. **SARS-CoV-2 target regions** - This section is an inventory of PCR-based nucleic acid amplification tests (NAATs) used by laboratories and *in silico* simulations of these assays on available high-quality and full length sequenced SARS-CoV-2 genomes (in collaboration with GISAID⁵). The latter are pre-computed values corresponding to the extent of matching of the primers and probes from NAAT database methods

against high-quality, full length genomic sequences, which are made available by GISAID to enable this analysis.

In addition, the database has pointers to the EUnetHTA publications repository. EUnetHTA Joint Action 3 is the scientific and technical component of EU cooperation on HTA. It builds on years of long standing collaboration between HTA agencies, financially supported by the EU since early 2000. The current Joint Action (2016 – 2021) is co-funded by the EU Health Programme and includes government appointed organisations and a large number of relevant regional agencies and not-for-profit organisations that produce or contribute to HTA in Europe (i.e. 86 organisations from 26 Member States plus Norway, Serbia, Switzerland, Ukraine and United Kingdom). Responding to the needs of policy makers and public health authorities for researched, timely and reliable information, EUnetHTA has launched its COVID-19-related repository of publications and outputs. The repository gathers publications by EUnetHTA and by HTA organisations on testing methods and devices, but also on treatment options, and other public health measures relevant to COVID-19. Direct link to this page is: <https://covid-19-diagnostics.jrc.ec.europa.eu/eunetha>.

Finally, a **Submit your device** section provides forms for manufacturers and NAAT developers to submit information on devices not yet listed in the database or to provide performance data not available to the authors at the time of the last update. The submitted information, once verified against the source provider, is taken into consideration for updating the database.

At the time of writing, the database includes 869 devices and 382 selected articles, easily linked to each other, as shown in [Figure 1](#). Considering the rapidly evolving situation in relation to the development and commercialisation of diagnostic devices for COVID-19, the completeness of the information is limited to the time of the last update as indicated in the database for each individual item. For this reason, manufacturers are invited to submit and update the information through the corresponding “Submit your device” section’s form.

With respect to NAAT methods, currently eight NAAT methods are available in the database: seven are from the World Health Organization (WHO) support to COVID-19 and are *in-house* PCR protocols assays posted online on the WHO website⁷ while one was developed by the JRC in the context of the production of EURM-019, a universal positive control material to be used in the testing of SARS-CoV-2 by RT-PCR (described in reference 8). The possibility to evaluate whether variations occurring overtime in the viral genomes can affect PCR-based detection methods is fundamental to guarantee reliability of the detection. The aims of the **SARS-CoV-2 target regions** section in the database are therefore both to detect potential target regions of the NAAT methods and to highlight possible differences with the expected reference sequence that might affect the performance of NAATs. As of 01 May 2020, about 10% of the retrieved genomes were found to be not detectable *in silico* by at least one of the eight




Figure 1. Device-article relations. This graph visualises the relations (grey lines) among devices (coloured according to the corresponding detection principle) and the scientific articles (black spots) in which they are mentioned. The large “corona” near the centre and the group near the top right corner are clustered around two documents that refer to 172⁹ and 23¹⁰ devices, respectively. It is quite evident that publications are specialised with respect to devices detection principle (i.e. the used technology): few publications (only three) mention different detection principles’ devices. Please note that devices without references are not represented. The size of the symbol representing the device is proportional to the corresponding number of references (min 1, max 8).

NAAT methods, as shown in [Figure 2](#). As an example, a case of variation potentially affecting NAAT methods detectability is shown in [Figure 3](#).

As for devices, considering the rapidly evolving situation in relation to COVID-19, NAAT methods developers are invited to submit updated information or new developed ones through the corresponding “Submit your device” section’s form.

Discussion

According to our knowledge, it is the first time that a repository provides and links information on COVID-19 *in vitro* diagnostic medical devices and the scientific articles testing and reporting their performance.

To date, more than 100 requests, either for addition of new *in vitro* diagnostic medical devices or for updating the information, have been received through the forms complied by manufacturers, as demonstration of the strong need of

having a structured data sharing point for this kind of information. However, it is important to highlight that this resource does NOT represent a list of devices approved or authorised for use either by the European Commission or by Member States’ national authorities. There is no central approval system for *in vitro* diagnostic medical devices in the EU. The currently applicable legislation for placing COVID-19 diagnostic devices on the market in the EU is Directive 98/79/EC on *in vitro* diagnostic medical devices³. Under the current Directive³, for COVID-19 devices that are designed for professional use, the manufacturer may affix the CE-mark to

³ From 26 May 2022, the Directive will be replaced by Regulation (EU) 2017/746 on *in vitro* diagnostic medical devices¹¹. The Regulation lays down a transitional period starting on the date of its entry into force (May-2017) during which the conformity of *in vitro* diagnostic medical devices can be assessed either under the Regulation or under the Directive. More information on the regulatory framework, together with a number of guidance documents, may be found on the Commission’s dedicated medical devices website https://ec.europa.eu/health/md_sector/overview_en.




Figure 2. NAAT methods target detectability. As of 01 May 2020, more than 16,000 high quality (<1% Ns and < 0.05% unique amino acid mutations) and full length (> 29,000 bp) viral genomes were made available by GISAID⁵. By performing *in silico* PCR simulations with eight NAAT methods, the graph shows that about 10% of the 16,000 selected genomes (cumulative numbers) were found to be not detectable *in silico* by at least one of the eight NAAT methods. Blue area represents the (cumulative) number of genomes *in silico* detected by all eight NAAT methods; red area represents the (cumulative) number of those not detected by at least one NAAT method.

*

GACCCCAAAATCAGCGAAATGCAC^{*}TCCGCATTACGTTTGGTGGACCCTCAGATTCAACTGGCAGTAACCAGA

Region	#	Country	#	First seen	Last seen
Asia	45	South Korea	2	26-Feb-20	27-Feb-20
		Taiwan	4	4-Mar-20	30-Mar-20
		Saudi Arabia	1	7-Mar-20	7-Mar-20
		Philippines	10	8-Mar-20	28-Mar-20
		India	1	16-Mar-20	16-Mar-20
		Singapore	23	17-Mar-20	15-Apr-20
		China	2	18-Mar-20	19-Mar-20
		Malaysia	2	18-Mar-20	20-Mar-20
Oceania	80	Australia	80	5-Mar-20	7-Apr-20
America	10	Canada	3	10-Mar-20	12-Mar-20
		USA	6	11-Mar-20	8-Apr-20
		Brazil	1	10-Mar-20	10-Mar-20
Africa	2	Gambia	1	21-Mar-20	21-Mar-20
		South Africa	1	31-Mar-20	31-Mar-20
Europe	3	United Kingdom	2	26-Mar-20	1-Apr-20
		Slovenia	1	29-Mar-20	29-Mar-20

Figure 3. Genomic variations that can affect NAAT methods target detectability. As of 01 May 2020, more than 16,000 high quality (<1% Ns and < 0.05% unique amino acid mutations) and full length (> 29,000 bp) viral genomes were made available by GISAID⁵. An example of identified genomic variations that can affect NAAT methods target detectability is here reported and it regards the first method developed by the US Centers for Disease Control and Prevention and supported by WHO in 7. 1) On the top, the target sequence (in red the left primer, in green right primer, and in azure the probe annealing region) that has been found in several genomes. With respect to the expected reference target, a variation (C/T transition, marked by *) is present at the beginning of the probe annealing region that can potentially affect the detection method performance. 2) Below, the table with details about the genomes where this variation of the target sequence has been found: the two columns with # indicate the number of genomes with such a sequence, aggregated by region and country, respectively. In the table, the dates of first and last appearance of the variation are reported for each country, in order to provide information on the “spread” of the variation.

the product after having ensured the compliance of the device with the Directive and drawn up a declaration of conformity. For COVID-19 devices that are designed for use by lay persons (self-tests), the manufacturer must also apply to a third party body called a notified body that will do additional verification and issue a certificate. As a consequence, the JRC should not be deemed responsible for the validity of such data.

As recently reported by Guglielmi¹², “*even once a test is working beautifully in the lab, it still faces an arduous journey to mass usage. The first challenge is to verify performance, because quality can vary*”. We believe that the here presented database contributes to the required prompt handling of the on-going public health COVID-19 situation. It makes sharing information easier and helps medical professionals, scientists and laboratories understand and properly assess COVID-19 medical devices quickly. The overview of the market may also help manufacturers develop and improve their own devices, in line with the current EU legislation.

Data availability

The database is freely accessible at <https://covid-19-diagnostics.jrc.ec.europa.eu/>. CSVs of the ‘COVID-19 *in vitro* diagnostic medical devices’ and ‘Scientific literature on COVID-19 test methods and devices’ databases at the time of publication have been archived on Zenodo repository (see below). However, the database is updated very frequently, and the latest data can be downloaded in these sections (described in Results) as CSVs.

Zenodo: JRC COVID-19 In Vitro Diagnostic Devices and Test Methods Database. <http://doi.org/10.5281/zenodo.4117601>¹³.

This project contains the following underlying data:

- covid-19-methods.csv (CSV of COVID-19 *in vitro* diagnostic medical devices database at time of publication)

- covid-19-literature.csv (CSV of Scientific literature on COVID-19 test methods and devices database at time of publication)

Data are available under the terms of the [Creative Commons Attribution 4.0 International license](https://creativecommons.org/licenses/by/4.0/) (CC-BY 4.0).

The ‘SARS-Cov2 target regions’ database is enabled using third party data from GISAID, subject to GISAID’s [Terms and Conditions](https://www.gisaid.org/terms/). The database is freely accessible at <https://covid-19-diagnostics.jrc.ec.europa.eu/amplicons> and GISAID data are available directly from GISAID at <https://www.gisaid.org/>. Access to GISAID data requires registration and agreement to the conditions for use at <https://www.gisaid.org/registration/register/>.

Code availability

Code used to perform scientific literature mining is freely available at https://github.com/ec-jrc/JRC_COV19_IVD-DEVs-TEMs_DB.

Archived code at time of publication: <https://doi.org/10.5281/zenodo.4084600>¹⁴.

License: [BSD-3-Clause](https://opensource.org/licenses/BSD-3-Clause)

Acknowledgements

We gratefully acknowledge the Authors and the Originating laboratories where the clinical specimen or virus isolate was first obtained and the Submitting laboratories, where sequence data have been generated and submitted to GISAID, on which this research is based.

We acknowledge all the manufacturers and repository website owners of the devices that have been analysed; all the authors, publishers and literature repositories of the papers that have been analysed.

References

1. **Communication from the Commission Guidelines on COVID-19 *in vitro* diagnostic tests and their performance 2020/C 122 I/01C/2020/2391.**
[Reference Source](#)
2. **Current performance of COVID-19 test methods and devices and proposed performance criteria - Working document of Commission services.**
[Reference Source](#)
3. **Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on *in vitro* diagnostic medical devices and repealing directive 98/79/EC and Commission Decision 2010/227/EU, OJ L 117, 5.5. 2017; 176–332.**
[Reference Source](#)
4. See <https://eunetha.eu/covid-19-diagnostics/>
5. Shu Y, McCauley J: **GISAID: Global initiative on sharing all influenza data - from vision to reality.** *Euro Surveill.* 2017; **22**(13): 30494.
[PubMed Abstract](#) | [Publisher Full Text](#) | [Free Full Text](#)
6. Gans JD, Wolinsky M: **Improved assay-dependent searching of nucleic acid sequence databases.** *Nucleic Acids Res.* 2008; **36**(12): e74.
[PubMed Abstract](#) | [Publisher Full Text](#) | [Free Full Text](#)
7. At time of writing available at <https://www.who.int/publications/m/item/molecular-assays-to-diagnose-covid-19-summary-table-of-available-protocols>
8. The method is available at <https://crm.jrc.ec.europa.eu/p/EURM-019>
9. Gonzalez JM, Shelton JW, Diaz-Vallejo M, *et al.*: **Immunological assays for SARS-CoV-2: an analysis of available commercial tests to measure antigen and antibodies.** 2020.
[Publisher Full Text](#)
10. Cheng MP, Papenburg J, Desjardins M, *et al.*: **Diagnostic Testing for Severe Acute Respiratory Syndrome-Related Coronavirus 2: A Narrative Review.** *Ann Intern Med.* 2020; **172**(11): 726–734.
[PubMed Abstract](#) | [Publisher Full Text](#) | [Free Full Text](#)
11. **Directive 98/79/EC of the European Parliament and of the Council of 27**

October 1998 on *in vitro* diagnostic medical devices, OJ L 331, 7.12. 1998; 1–37.

[Reference Source](#)

12. Guglielmi G: **The explosion of new coronavirus tests that could help to end the pandemic.** *Nature*. 2020; **583**(7817): 506–509.
[PubMed Abstract](#) | [Publisher Full Text](#)

13. Petrillo M: **JRC COVID-19 *In Vitro* Diagnostic Devices and Test Methods Database (Version v1) [Data set].** *Zenodo*. 2020.
<http://www.doi.org/10.5281/zenodo.4117601>

14. Gallardo AP, gabbaan: **ec-jrc/JRC_COV19_IVD-DEVs-TEMs_DB: First release of JRC_COVID-19_Literature_search (Version v1.0.0).** *Zenodo*. 2020.
<http://www.doi.org/10.5281/zenodo.4084600>

Open Peer Review

Current Peer Review Status:  

Version 1

Reviewer Report 01 February 2021

<https://doi.org/10.5256/f1000research.30175.r75227>

© 2021 Viljoen C. This is an open access peer review report distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.



Christopher Viljoen 

Human Molecular Biology Unit, School of Biomedical Sciences, University of the Free State, Bloemfontein, South Africa

- The Abstract would be strengthened with an introductory sentence regarding context to COVID-19 and the need for such information.
- The Introduction needs a short contextualized introduction to the need for this type of information in a COVID-19 pandemic.
- The Introduction could also indicate how the Commission has responded to the need for the information as indicated above.
- The paragraph on COVID-19 mutations and test methods needs to include some context to the problem.
- The Discussion - 1st paragraph - should indicate "in vitro diagnostic medical devices and assays".
- The article would benefit from some concluding comments reiterating the importance of having access to the information presented in the database.

Is the rationale for creating the dataset(s) clearly described?

Partly

Are the protocols appropriate and is the work technically sound?

Yes

Are sufficient details of methods and materials provided to allow replication by others?

Yes

Are the datasets clearly presented in a useable and accessible format?

Yes

Competing Interests: No competing interests were disclosed.**Reviewer Expertise:** Human Molecular Biology**I confirm that I have read this submission and believe that I have an appropriate level of expertise to confirm that it is of an acceptable scientific standard.**

Reviewer Report 06 January 2021

<https://doi.org/10.5256/f1000research.30175.r75984>

© 2021 Facchiano A. This is an open access peer review report distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

**Angelo Facchiano** 

Institute of Food Science, National Research Council, Avellino, Italy

The manuscript presents a significant work aimed to collect publicly available information on test methods for COVID-19. Due to the amount of devices and methods available, the database described offers a useful tool to be addressed in the field. I find of special interest the example on genomic variations that can affect detectability.

The manuscript is well presented and written.

I find in the manuscript only minor points to be addressed:

- As a reader, I would find in the article an example of the output obtained when the database is searched.
- I would not find explanation of abbreviations in the abstract if not really used there.
- abbreviations should be explained at the first appearance (check for NAAT; IVD - if removed from abstract, FAIR - explain or add a reference; and so on).
- under keywords, check upper/lower case.

Is the rationale for creating the dataset(s) clearly described?

Yes

Are the protocols appropriate and is the work technically sound?

Yes

Are sufficient details of methods and materials provided to allow replication by others?

Yes

Are the datasets clearly presented in a useable and accessible format?

Yes

Competing Interests: No competing interests were disclosed.

Reviewer Expertise: Bioinformatics, biochemistry, molecular biology

I confirm that I have read this submission and believe that I have an appropriate level of expertise to confirm that it is of an acceptable scientific standard.

The benefits of publishing with F1000Research:

- Your article is published within days, with no editorial bias
- You can publish traditional articles, null/negative results, case reports, data notes and more
- The peer review process is transparent and collaborative
- Your article is indexed in PubMed after passing peer review
- Dedicated customer support at every stage

For pre-submission enquiries, contact research@f1000.com

F1000Research