

1 **Supplementary Information for:**

2

3 **A robotic platform for fluidically-linked human body-on-chips**
4 **experimentation**

5 Richard Novak^{1*}, Miles Ingram^{1*}, Susan Clauson^{1*}, Debarun Das^{2*}, Aaron Delahanty^{1*}, Anna Herland^{1,*},
6 Ben M. Maoz^{1,3*}, Sauveur S. F. Jeanty^{1**}, Mahadevabharath R. Somayaji^{2*}, Morgan Burt¹, Elizabeth
7 Calamari¹, Angeliki Chalkiadaki¹, Alexander Cho², Youngjae Choe¹, David Benson Chou^{1,4}, Michael
8 Cronic¹, Stephanie Dauth^{1,3}, Toni Divic¹, Jose Fernandez-Alcon^{1*}, Thomas Ferrante¹, John Ferrier^{1,3},
9 Edward A. FitzGerald¹, Rachel Fleming¹, Sasan Jalili-Firoozinezhad^{1,5}, Thomas Grevesse^{1,3}, Josue A.
10 Goss^{1,3}, Tiama Hamkins-Indik¹, Olivier Henry¹, Chris Hinojosa^{1*}, Tessa Huffstater¹, Kyung-Jin Jang^{1*}, Ville
11 Kujala^{1,3*}, Lian Leng^{1*}, Robert Mannix^{1,6}, Yuka Milton¹, Janna Nawroth^{1,3*}, Bret A. Nestor¹, Carlos F. Ng¹,
12 Blakely O'Connor^{1,3}, Tae-Eun Park¹, Henry Sanchez¹, Josiah Sliz^{1*}, Alexandra Sontheimer-Phelps^{1,7}, Ben
13 Swenor¹, Guy Thompson II^{1*}, George J. Touloumes^{1,3}, Zachary Tranchemontagne¹, Norman Wen^{1*},
14 Moran Yadid^{1,3}, Anthony Bahinski^{1,8}, Geraldine A. Hamilton^{1*}, Daniel Levner^{1*}, Oren Levy¹, Andrzej
15 Przekwas², Rachele Prantil-Baun¹, Kevin K. Parker^{1,3}, and Donald E. Ingber^{1,3,6}

16

17 ¹Wyss Institute for Biologically Inspired Engineering at Harvard University, Boston, MA 02115, USA; ²CFD
18 Research Corporation, Huntsville, AL, USA 35806; ³Harvard John A. Paulson School of Engineering and
19 Applied Sciences, Harvard University, Cambridge, MA 02138, USA; ⁴Department of Pathology,
20 Massachusetts General Hospital, Boston, MA 02115, USA; ⁵Department of Bioengineering and iBB -
21 Institute for Bioengineering and Biosciences, Instituto Superior Técnico, Universidade de Lisboa, Lisboa,
22 Portugal; ⁶Vascular Biology Program and Department of Surgery, Boston Children's Hospital and Harvard
23 Medical School, Boston, MA 02115, USA; ⁷Department of Biology, University of Freiburg, Freiburg,
24 Germany; ⁸GlaxoSmithKline, 1250 S. Collegeville Rd., Collegeville, PA 19426, USA;

25

26 * Authors contributed equally

27 + Current address: Emulate, Inc., 27 Drydock Avenue, Boston, MA 02210, USA

28

29

30

Stage Calibration

The deck layout exists in three spaces: the physical layout, the CAD assembly and the web-based GUI. The calibration protocol developed here relates the positional error between the physical and CAD layouts. The measured positional errors can then be used to calculate linear scaling and offset values that, when applied in the web-based GUI, significantly reduce alignment error between the liquid handler and the individual deck components. A given x, y, z coordinate as described by the CAD layout is calibrated to the physical layout by multiplying the spatial coordinates by a calibration matrix and adding an offset matrix. A schematic of this calibration is shown in **Fig. 23b**. In a perfectly aligned system, the calibration matrix would simply be an identity matrix (Equation 1) and all offset values would be zero (Equation 2). In order to use this calibration procedure, robot alignment error must be determined by installing four precision ground post subassemblies (Thorlabs P series 1.5" diameter post, MiSUMi M6 10mm diameter locating pin) used as calibration references and a ruby-probe tip (Renishaw 1mm ball end 10 mm length, M4 thread) attached to the end of the liquid handler with a custom tapped 20 mm long spacer (M4 female to ¼-28 female). Two matrices are generated; matrix 'P', the theoretical, CAD-based, distance between the four measurement references (equation 3), and matrix 'Q', the physically measured distance between the four calibration references (equation 4). The web-based GUI features an interactive aid to complete the measurement process; upon completion, the calibration posts are removed. The process of determining the calibration and offset matrix is shown in equations 5 and 6, respectively. Once the M and O matrices have been determined, they are permanently stored within the client code and are applied to all motion system coordinates as directed by the web-based GUI. The calibration procedure is run when the system is assembled and re-calibrated upon software re-install, if any major changes are made to the deck layout, or at regularly scheduled maintenance intervals.

$$M = \begin{bmatrix} 1 & 0 & 0 \\ 0 & 1 & 0 \\ 0 & 0 & 1 \end{bmatrix} \quad [\text{Eq. 1}]$$

$$O = \begin{bmatrix} 0 & 0 & 0 \end{bmatrix} \quad [\text{Eq. 2}]$$

$$P = \begin{bmatrix} (B_{x1} - A_{x1}) & (B_{y1} - A_{y1}) & (B_{z1} - A_{z1}) \\ (C_{x1} - A_{x1}) & (C_{y1} - A_{y1}) & (C_{z1} - A_{z1}) \\ (D_{x1} - A_{x1}) & (D_{y1} - A_{y1}) & (D_{z1} - A_{z1}) \end{bmatrix} \quad [\text{Eq. 3}]$$

$$Q = \begin{bmatrix} (B_{x2} - A_{x2}) & (B_{y2} - A_{y2}) & (B_{z2} - A_{z2}) \\ (C_{x2} - A_{x2}) & (C_{y2} - A_{y2}) & (C_{z2} - A_{z2}) \\ (D_{x2} - A_{x2}) & (D_{y2} - A_{y2}) & (D_{z2} - A_{z2}) \end{bmatrix} \quad [\text{Eq. 4}]$$

$$M = [Q][P]^{-1} \quad [\text{Eq. 5}]$$

$$\begin{bmatrix} O_x & O_y & O_z \end{bmatrix} = \begin{bmatrix} A_{x2} & A_{y2} & A_{z2} \end{bmatrix} - \begin{bmatrix} A_{x1} & A_{y1} & A_{z1} \end{bmatrix} [M] \quad [\text{Eq. 6}]$$

SI Table 1. Organ Chip types and associated properties of the microfluidic features.

Organ Chip	Membrane	Channel Width (mm)	Apical Channel Height (mm)	Basal Channel Height (mm)	Channel Length (mm)
Gut	PDMS 7 μ m pores	1	1	0.2	91
Liver	PDMS 7 μ m pores	1	1	0.2	24
Kidney	PET 0.4 μ m pores	1	0.1	0.1	24
Heart	PC 5 μ m pores	2.5	0.25	1.2	25
Lung	PDMS 7 μ m pores	1	1	0.2	24
Skin	PDMS 7 μ m pores	5.5 (apical) 1 (basal)	5	0.3	24
BBB	PET 0.4 μ m pores	1	1	0.2	24
Brain	PC 5 μ m pores	2.5	0.1	1.2	25

PDMS: poly (dimethyl siloxane); PET: polyester terephthalate; PC: polycarbonate

72 **SI Table 2 – Cell culture parameters**


73

Organ	Chip Type	ECM Coating	Apical		Basal	
			Cell Type	Seeding Density	Cell Type	Seeding Density
Gut	Long Tall Channel, Stretchable, PDMS Membrane	Matrigel & Collagen I	Caco2 BBE	100,000 cells/cm ²	HUVEC	100,000 cells/cm ²
Liver	Standard Tall Channel, Stretchable, PDMS Membrane	Collagen I	Primary Human Hepatocytes	250,000 cells/cm ²	LSEC	100,000 cells/cm ²
Kidney	Standard tall channel, PET membrane	Collagen IV & Laminin	Human Renal Proximal Tubule	120,000 cells/cm ²	Human Glomerular Microvascular Endothelial	100,000 cells/cm ²
Heart	Dual channel PC, PET Membrane	Fibronectin	Human iPSC Cardiomyocytes (Cor4U, Axiogenesis)	200,000 cell/cm ²	HUVEC	100,000 cells/cm ²
Lung	Standard Tall Channel, Stretchable, PDMS Membrane	Fibronectin & Collagen I	A549 Human Lung Carcinoma Cell Line	200,000 cells/cm ²	HUVEC	60,000 cells/cm ²
BBB	Standard tall channel, PET membrane	Collagen IV & Fibronectin	Primary Human Astrocytes	70,000 cells/cm ²	Primary Human Brain Microvascular Endothelial	90,000 cells/cm ²
			Primary Human Pericytes	30,000 cells/cm ²		
Brain	Dual channel PC, PET Membrane	PDL & Laminin	Human Hippocampal Neural Stem Cells (hNSCs)	100,000 cells/cm ²	N/A	
Skin	Oval Open-top, Stretchable, PDMS Membrane	Fibronectin & Collagen I	Human Primary Adult Dermal Fibroblasts	82,000 cells/cm ²	Human Primary Dermal Microvascular Endothelial	100,000 cells/cm ²
			Human Primary Neonatal Epidermal Keratinocytes	800,000 cells/cm ²		

74

75

76




77

78

79 **Supplementary Figure S1.** Cross sectional views of Organ Chip CAD for all 8 Organ
 80 Chips used in this study. In all chips, two parallel fluidic channels are separated by a
 81 porous semi-permeable membrane.

82




Supplementary Figure S2. Characterization of Interrogator operation with respect to temperature during continuous fluid handling (A) indicates minimal impact of electronics on incubator thermal stability. Thermal stability was also not affected over the course of an Organ Chip linking study (B). Drops in temperature indicate incubator door opening to exchange plates and reagents. The peristaltic pump flow rate (C) remained within 10% of the nominal setpoint for at least one week of operation.

Formatted: Font: Bold


Formatted: Font: Bold

Formatted: Font: Bold




92 **Supplementary Figure S2S3.** Hydra software screen shots highlighting key features of
93 device setup, component customization, and experimental design. All Interrogator
94 programming is performed using this graphical Hydra software interface.


Landing
Page




Connect
to Device




Create New
Experiment



Setup
Fluidics





Setup
Serial
Dilution



Run
Experiment




99 **Supplementary Figure S3S4.** Hydra software screen shots of the rapid workflow
100 involved in designing and executing a typical experiment
101
102
103
104
105
106
107



108
109
110 **Supplementary Figure S4S5.** Schematic of Hydra (A) system communication structure
111 and (B) network software architecture.
112
113
114
115

Formatted: Font: Bold
Formatted: Font: Bold



Supplementary Figure S5S6. Rendering of the phase plate (A) used to generate phase contrast-like illumination in Organ Chips while enabling scanning across their channel length. Each Organ Chip cartridge location is directly below each phase slot feature (B), which was designed based on standard condenser phase ring geometry (C). An epifluorescence module (D) was designed to enable fluorescence imaging in the Interrogator (E).

Formatted: Font: Bold

Formatted: Font: Bold



Formatted: Font: Bold

Formatted: Font: Bold

Formatted: Font: Bold


a.

GUT




b.

LIVER




c.

KIDNEY




d.

HEART




e.

LUNG




f.

BBB




g.


BOC



h.

SKIN





Supplementary Figure S6S7. Organ Chip-specific readouts prior to linking demonstrating organ maturity. **(A)** Gut Chip: Cascade blue® (596 Da) permeability, n=7, CYP3A4 activity, and bright field imaging of villi, scale 50 μ m. **(B)** Liver Chip: albumin production and CYP3A4 activity, n=4 and n=6, influence of liver endothelial cells on Liver Chip function, n=2 **(C)** Kidney Chip: albumin reabsorption, n=3 **(D)** Heart Chip: beat rate and orientation order parameter (OOP¹¹), n=4 **(E)** Lung Chip: Texas Red™ (3 kDa) permeability, n=10 **(F)** BBB Chip permeability effect of Cyclic adenosine monophosphate (cAMP) and comparison to rat *in vivo* data : Inulin-FITC (2-5 kDa, listed as 4 kDa for comparison with literature *in vivo* data) p=0.002, n=3, Dextran-Cascade blue® (10 kDa) p= 0.024, n=3, Dextran-Texas Red™ (70 kDa) permeability, n=4, *in vivo* data¹² **(G)** Brain Chip: ratio of neuronal cells and glutamine:glutamate ratio, n=9 and n=8 **(H)** Skin Chip: Cascade Blue (596Da) permeability, n=3 at left and H&E stained section through chip showing squamous epidermis overlying a thick epidermis (bar, 50 μ m).

Formatted: Font: Bold

Formatted: Font: Bold

Formatted: Font: Bold

Formatted: Font: Bold


Formatted: Font: Bold

Formatted: Font: Bold




Formatted: Font: Bold

Formatted: Font: Bold

Formatted: Font: Symbol



Supplementary Figure S87. Length of culture time required for Organ Chips to reach a mature morphological and functional state. After the indicated time the Organ Chips were used for linkage experiments.



156
157
158 **Supplementary Figure S98.** Computational Microphysiological Model of a generic
159 Organ Chip **(A)**. Each Organ Chip model is spatially resolved in both stream-wise and
160 cross-stream direction (top PDMS, top channel media, epithelial cell layer, membrane,
161 endothelial cell layer, bottom channel media, bottom PDMS). **(B)** A quantitative reduced
162 order 3D model (Q3D) enables tracking drug dynamics through linked Organ Chips (Gut
163 to Liver shown here). Two time instances show initial transport of the gut microchannel
164 at 1.5 hr and at 6.1hr. The media from “vascular” collecting reservoirs is periodically
165 transferred to downstream linked Organ Chips using the Interrogator. **(C)** Detailed view
166 of the 3D spatially resolved Organ Chip model showing the geometrical/mesh resolution
167 in the cross-stream and stream-wise (axial) directions.
168
169
170


Formatted: Font: Bold

Formatted: Font: Bold


Formatted: Font: Bold

171

a.



b.




172
173
174
175
176
177
178
179
180


Supplementary Figure S910. Computational Model of the Lung Chip (A) Reduced-Order Model of Lung Chip in liquid-liquid interface culture, 4 compartments in axial direction, (B) High-Fidelity Model of drug bolus movement and apical→basal mass transfer through Lung Chip membrane over time.

Formatted: Font: Bold

Formatted: Font: Bold



Supplementary Figure S110. Computational Model of the Skin Chip: Reduced-Order Model of Skin Chip, 4 compartments in axial direction



Supplementary Figure S124. Computational model of the Heart Chip **(A)** Reduced-Order Model of MTF Heart Chip (6 compartments in axial direction), **(B)** High-fidelity models of MTF Heart Chip predict increased apical→basolateral mass transfer of high MW tracer when increasing chip length from 8 mm (top) to 12 mm (bottom).

Formatted: Font: Bold

Formatted: Font: Bold

194
195

a.

Parameter	Value	Unit
Step Pulse	10	μs
Step Idle Delay	100	ms
Step Size (X & Y, Z)	20.051, 200.000	step/mm
Max Velocity (X & Y, Z)	200, 56	mm/s
Acceleration (X & Y, Z)	400, 160	mm/s ²
Motion Envelope (X & Y, Z)	310, 135	mm

b.

```
$0=10 (step pulse, usec)
$1=100 (step idle delay, msec)
$2=0 (step port invert mask:00000000)
$3=3 (dir port invert mask:00000011)
$4=0 (step enable invert, bool)
$5=1 (limit pins invert, bool)
$6=0 (probe pin invert, bool)
$10=1 (status report mask:00000001)
$11=0.020 (junction deviation, mm)
$12=0.010 (arc tolerance, mm)
$13=0 (report inches, bool)
$20=1 (soft limits, bool)
$21=1 (hard limits, bool)
$22=1 (homing cycle, bool)
$23=7 (homing dir invert mask:00000111)
$24=60.000 (homing feed, mm/min)
$25=2000.000 (homing seek, mm/min)
$26=100 (homing debounce, msec)
$27=1.000 (homing pull-off, mm)
$100=20.051 (x, step/mm)
$101=20.051 (y, step/mm)
$102=200.000 (z, step/mm)
$110=12000.000 (x max rate, mm/min)
$111=12000.000 (y max rate, mm/min)
$112=3400.000 (z max rate, mm/min)
$120=400.000 (x accel, mm/sec^2)
$121=400.000 (y accel, mm/sec^2)
$122=160.000 (z accel, mm/sec^2)
$130=310.000 (x max travel, mm)
$131=310.000 (y max travel, mm)
$132=135.000 (z max travel, mm)
```

196
197
198
199
200
201
202
203

Supplementary Figure S132. Motion system parameters **(A)** and **GRBL-stepper motor control shield** settings **(B)** that define the characteristic behavior of the system.

Formatted: Font: Bold
Formatted: Font: Bold

204 **Supplementary Movies**

205

206

207 **SI Movie 1**

208

209 **Overview of Interrogator system components and Organ Chip linking.**

210

211

212 **SI Movie 2**

213

214

215 **Microscope module movie of gut chip stretching.**

216

217

218 **SI Movie 3**

219

220 **Heart Chip beating after being linked for 3 weeks on the Interrogator platform.**

221

222

223 **SI Movie 4**

224

225 **High-fidelity diffusion model of the Heart Chip.**

226

227

228 **SI Movie 5**

229

230 **CoBi Q-3D Model of Gut-Liver Chip Linking**

231