



Since January 2020 Elsevier has created a COVID-19 resource centre with free information in English and Mandarin on the novel coronavirus COVID-19. The COVID-19 resource centre is hosted on Elsevier Connect, the company's public news and information website.

Elsevier hereby grants permission to make all its COVID-19-related research that is available on the COVID-19 resource centre - including this research content - immediately available in PubMed Central and other publicly funded repositories, such as the WHO COVID database with rights for unrestricted research re-use and analyses in any form or by any means with acknowledgement of the original source. These permissions are granted for free by Elsevier for as long as the COVID-19 resource centre remains active.

sis pulmonares consolidativas. Dado que la consolidación de mayor tamaño, localizada en el lóbulo inferior derecho, presentaba amplio contacto con la pleura periférica, se optó por realizar una biopsia con guía ecográfica. Se administraron 2,4 ml de contraste ecográfico (SonoVue, Rovi, Pozuelo de Alarcón, Madrid, España) y se objetivaron dos zonas diferenciadas según el patrón de captación: un área periférica con captación tardía (más allá de los 6 segundos postinyección) y lavado precoz (desaparición de la captación de contraste a los pocos segundos de iniciarse la captación) y el resto de la consolidación con captación precoz (antes de 6 segundos) y homogénea y lavado tardío (más allá de un minuto) (fig. 1B-F). Ante estos hallazgos se decidió biopsiar la primera de las zonas descritas con aguja fina de 22G con el resultado de metástasis alveolar de adenocarcinoma pancreático (fig. 1G). Debido a que el resto de la consolidación mostró características sugestivas de proceso neumónico se decidió realizar fibrobroncoscopia y toma de muestras que demostraron una infección fúngica (fig. 1H).

El caso expuesto se incluye dentro de este grupo, pues existía una consolidación neumónica preexistente con un área hipoeoica en su interior. Tras el estudio con contraste ecográfico la consolidación presentó una captación tardía y homogénea con lavado tardío a excepción del área central y periférica cuyo comportamiento mostró una captación tardía pero lavado precoz, sugestivo de malignidad. El patrón de captación objetivado permitió además guiar la biopsia percutánea hacia la zona deseada más sospechosa para mejorar la rentabilidad diagnóstica en la obtención de muestras.

La biopsia con guía ecográfica es una alternativa a la realizada mediante TC para las lesiones pulmonares periféricas o pleurales<sup>1,3</sup> y alcanza una efectividad y rentabilidad diagnóstica similar a la obtenida con TC<sup>5</sup>. Adicionalmente, los procedimientos percutáneos guiados por ecografía ofrecen ciertas ventajas tales como la monitorización del procedimiento en tiempo real, la ausencia de radiación, menores costes y duración del procedimiento, con ratios de complicaciones similares o menores a los de la biopsia guiada

por TC<sup>1,2</sup>. En muchos casos el estudio con ecografía torácica con contraste ayuda a comprender la naturaleza de la lesión a estudio y guiar, en caso que fuere necesario, la toma de muestras hacia áreas de interés evitando focos necróticos<sup>1,2,6</sup> o hacia las zonas de mayor sospecha de malignidad como en el caso presentado.

## Bibliografía

- Caremani M, Bencì A, Lapini L, Tacconi D, Caremani A, Ciccotosto C, et al. Contrast enhanced ultrasonography (CEUS) in peripheral lung lesions: A study of 60 cases. *J Ultrasound*. 2008;11:89–96.
- Dong Y, Mao F, Wang WP, Ji ZB, Fan PL. Value of contrast-enhanced ultrasound in guidance of percutaneous biopsy in peripheral pulmonary lesions. *Biomed Res Int*. 2015;2015:531507.
- Sartori S, Postorivo S, di Vece F, Ermili F, Tassinari D, Tombesi P. Contrast-enhanced ultrasonography in peripheral lung consolidations: What's its actual role? *World J Radiol*. 2013;28:372–80.
- García-Ortega A, Briones-Gómez A, Fabregat S, Martínez-Tomás R, Martínez-García MÁ, Cases E. Benefit of chest ultrasonography in the diagnosis of peripheral thoracic lesions in an Interventional Pulmonology Unit. *Arch Bronconeumol*. 2016;52:244–9.
- Lee MH, Lubner MG, Hinshaw JL, Pickhardt PJ. Ultrasound guidance versus CT guidance for peripheral lung biopsy: Performance according to lesions size and pleural contact. *AJR Am J Roentgenol*. 2018;210:W110–7.
- Wang S, Yang W, Zhang H, Xu Q, Yan K. The role of contrast-enhanced ultrasound in selection indication and improving diagnosis for transthoracic biopsy in peripheral pulmonary and mediastinal lesions. *Biomed Res Int*. 2015;2015:231782.

Gemma Isus Olivé, Alfredo Páez Carpio, Daniel Martínez e Ivan Vollmer\*

Servicio de Radiodiagnóstico (CDIC), Hospital Clínic, Barcelona, España

\* Autor para correspondencia.

Correo electrónico: vollmer@clinic.cat (I. Vollmer).

<https://doi.org/10.1016/j.arbres.2021.01.033>

0300-2896/ © 2021 SEPAR. Publicado por Elsevier España, S.L.U. Todos los derechos reservados.

## High O<sub>2</sub> Flow Rates Required to Achieve Acceptable FiO<sub>2</sub> in CPAP-Treated Patients With Severe Covid-19: A Clinically Based Bench Study



### Se requieren altos índices de flujo de O<sub>2</sub> para alcanzar una FiO<sub>2</sub> aceptable en los pacientes con covid-19 tratados con CPAP: un estudio experimental basado en la clínica

Dear Editor,

During the Covid-19 pandemic, a lack of ventilatory equipment in intensive care units (ICU), or patient comorbidities meant that some patients received non invasive, continuous positive airway pressure (CPAP), with the highest fraction of inspired oxygen (FiO<sub>2</sub>) as a ceiling treatment.<sup>1–3</sup> In Le Havre hospital (France), between September and December 2020, around 30 patients were treated with bilevel home devices in CPAP mode using vented oronasal masks. In contrast with bench studies that reported high FiO<sub>2</sub> with oxygen flow rates <30 L/min in optimal experimental conditions (i.e., low minute ventilation, good pulmonary compliance and no leakage),<sup>4,5</sup> many patients required O<sub>2</sub> flow rates >70 L/min to maintain oxygen saturation (SpO<sub>2</sub>) ≥90%. We hypothesised that due to the high respiratory demand of patients with severe Covid-19, a high O<sub>2</sub> flow rate would be required to reach adequate pressure levels with CPAP, thereby substantially increasing the oxygen flow rate needed to achieve a high FiO<sub>2</sub>. We carried out a bench

study to measure the oxygen flow rates needed to reach high FiO<sub>2</sub> levels, using a pulmonary model that reproduced the characteristics of patients with severe Covid-19.

First, we extracted clinical data from 10 consecutive patients who were included in an ongoing clinical trial (EURO-CPAP CT2220141 approved by our institutional review board) from the CPAP built-in software. This case series is presented for descriptive purposes to illustrate our hypothesis. Then, we performed several bench experiments. We used a CPAP device (AirSense 10 AutoSet, ResMed, San Diego, CA, USA) set at 10 cmH<sub>2</sub>O, connected to an artificial lung (ASL 5000, Ingmar Medical, Pittsburgh, PA, USA). A low-resistance antibacterial filter (Clear-Guard, Intersurgical, Wokingham, UK), an oxygen inlet connector, a 15 mm single limb circuit, a Whisper Swivel II Exhalation port (Philips Respironics, Murrysville, PA, USA) that simulated intentional leaks, and the artificial lung were placed in series and connected to the CPAP device. A standard 4-mm diameter leak port was placed between the exhalation port and the artificial lung to mimic unintentional leakage. Oxygen was delivered into the system using two O<sub>2</sub> flow meters (EASY MED-O<sub>2</sub>, Air Liquide Healthcare, Paris, France) that could both deliver up to 50 L/min. Three different mechanical lung conditions were simulated by modulating the resistance (R) and the compliance (C) of the artificial lung, corresponding to the following experimental models:

A. Normal: R = 5 cmH<sub>2</sub>O/l s and C = 60 ml/cmH<sub>2</sub>O.

**Table 1**

Descriptive data from 10 consecutive patients treated with CPAP and O<sub>2</sub>. Leakage, Tidal Volume (Vt) and Respiratory Rate (RR) were collected from built-in CPAP software during the first 48 h of CPAP treatment.

Patient	Age (years)	Sex	BMI (kg/m <sup>2</sup> )	CRF	CF	AH	DT2	PaO <sub>2</sub> /FiO <sub>2</sub> at CPAP initiation	CPAP level (cm H <sub>2</sub> O)	Median leakage (l/min)	Mean Vt (ml)	Mean RR (cycle per minute)	O <sub>2</sub> flow <sup>a</sup> (l/min)	SpO <sub>2</sub> (%)
1	87	M	20.8	–	–	–	–	122	10	43.5 <sup>a</sup>	490	29	70	96
2	72	M	33.7	–	–	–	Yes	120	7	0 <sup>b</sup>	374.5	34.5	15	93
3	61	M	41.3	–	–	Yes	Yes	151	10	5 <sup>b</sup>	409.5	23	15	92
4	83	M	29.6	–	Yes	Yes	–	82	8	6 <sup>b</sup>	490	38.5	50	95
5	85	M	27.3	–	Yes	Yes	–	143	10	35 <sup>a</sup>	927	30.5	30	94
6	70	M	31.9	–	Yes	Yes	Yes	101	10	35 <sup>a</sup>	1108	38	100	92
7	80	M	33.3	Yes	–	Yes	Yes	58	7	0 <sup>b</sup>	938	34	100	89
8	77	M	23.6	Yes	Yes	Yes	Yes	67	12	7 <sup>b</sup>	972.5	38.5	50	93
9	55	M	21	–	Yes	–	–	134	8	43.5 <sup>a</sup>	473	24	30	98
10	84	M	26.1	–	Yes	Yes	–	90	10	MD	625	39.5	100	90

Body mass index: BMI; chronic respiratory failure: CRF; cardiac failure: CF; arterial hypertension: AH; type 2 diabetes: DT2; continuous positive airway pressure: CPAP; partial O<sub>2</sub> pressure in arterial blood: PaO<sub>2</sub>; fraction of inspired oxygen: FiO<sub>2</sub>.

<sup>a</sup> Total leakage.

<sup>b</sup> Unintentional leakage.

B. Restrictive:  $R = 5 \text{ cmH}_2\text{O/l/s}$  and  $C = 30 \text{ ml/cmH}_2\text{O}$ .

C. Normal with by-pass leak using a T-connector as recommended to reduce the risk of droplet aerosolization.<sup>6</sup>

Each model was run at breathing frequencies of 22, 30 and 35 cycles per minute (cpm), associated with inspiratory airway pressure drops of 2.5–6 cmH<sub>2</sub>O at 100 ms (P0.1) to reach tidal volumes of around 600 ml. The values of the inspiratory effort settings chosen for the simulation in the ASL5000 were based on published clinical values.<sup>7–11</sup> Each condition was run with and without unintentional leakage. FiO<sub>2</sub> was measured with 5, 15, 30, 50, 60 and 80 L/min of O<sub>2</sub>, and the experiment was stopped when 100% FiO<sub>2</sub> was reached. The sensor used was a paramagnetic oxygen transducer located in the ASL 5000 piston chamber which allowed a fast measurement with a response time below 350 ms. Ten respiratory cycles were analysed after the FiO<sub>2</sub> had reached an asymptote for each model at each O<sub>2</sub> flow rate. Normality and homogeneity were assessed using Shapiro-Wilk and Levene tests, respectively. Welch's ANOVA was used to compare respiratory rates at each O<sub>2</sub> flow rate. A Games-Howell post-hoc procedure was applied for multiple comparisons, with 22 cpm as the reference condition. Significance was set at 0.05. Statistical analysis was performed using IBM SPSS Statistics V25.

Data from 10 consecutive patients with Covid-19 treated with CPAP and O<sub>2</sub> are presented in Table 1. Six patients needed O<sub>2</sub> flow rates  $\geq 50 \text{ L/min}$  to maintain SpO<sub>2</sub> > 90%. Minimum O<sub>2</sub> flow rate was 15 L/min. All patients had severe infection (i.e., PaO<sub>2</sub>/FiO<sub>2</sub> ratio  $\leq 150$  at CPAP initiation).

Fig. 1 shows the relationship between the O<sub>2</sub> flow rate delivered and the resulting FiO<sub>2</sub> in the simulation. Without unintentional leaks, a flow rate of at least 30 L/min was necessary to reach FiO<sub>2</sub> > 95%. FiO<sub>2</sub> differed significantly for each respiratory rate at each O<sub>2</sub> flow rate with and without unintentional leaks ( $p < 0.001$ ). With unintentional leaks, O<sub>2</sub> flow rates of up to 80 l/min were required to reach 100% FiO<sub>2</sub> (Fig. 1A–C).

The results of this experimental bench study confirmed our clinical observation that high O<sub>2</sub> flow rates are required to achieve a high FiO<sub>2</sub> in patients with severe Covid-19 using CPAP and vented masks. An O<sub>2</sub> flow rate of at least 30 L/min, and up to 80 L/min, was required to achieve 100% FiO<sub>2</sub> for all three experimental models. Reduced compliance and unintentional leakage dramatically reduced FiO<sub>2</sub>. The arbitrary compliance of 30 ml/cmH<sub>2</sub>O used in model B is similar to that of patients with severe COVID-19.<sup>12</sup> Such a low level of compliance would reduce the FiO<sub>2</sub> with home CPAP and a vented mask, further confirming the need for high O<sub>2</sub> flow rates. In our clinical cohort, three patients required an O<sub>2</sub> flow rate

above 80 L/min: this very high rate causes two concerns. Firstly, we did not know the lung compliance of these patients. Secondly, the O<sub>2</sub> flow rate required to achieve a given FiO<sub>2</sub> is likely to be underestimated with the experimental model compared to the clinical situation, since no gas exchange occurs with an artificial lung. For instance, in a patient, with a 21% FiO<sub>2</sub>, the partial pressure of O<sub>2</sub> in the inspired air will be 160 mm Hg, while in the expired air (at sea level) it will be 120 mm Hg due to oxygen absorption. Thus, in vivo conditions even higher O<sub>2</sub> flows than those reported in this study may be required.<sup>13</sup>

Another important point is that many teams use helmets to ventilate patients with Covid-19.<sup>2,14</sup> Helmets could be more advantageous than oronasal masks because of the lower rate of unintentional leakage; this would allow lower oxygen flow rates to be used to achieve sufficient FiO<sub>2</sub> levels for patients with severe hypoxemia.

Studies have shown that patients who cannot be intubated require a high FiO<sub>2</sub> to maintain their SpO<sub>2</sub> above 92%,<sup>3,15</sup> and that poor oxygenation is strongly associated with mortality.<sup>3,16</sup> Our study provides valuable data to guide clinicians in the use of O<sub>2</sub> when using home CPAP in patients with severe Covid-19 who cannot be intubated. The results showed that these patients can be successfully treated with home CPAP devices, or NIV devices in CPAP mode, with vented masks, if a high O<sub>2</sub> flow rate is provided. However, close monitoring of unintentional leaks is crucial to ensure a sufficient FiO<sub>2</sub> for the treatment of severe hypoxemia in these patients.

## Authors' contribution

All authors participated in the study design, research, and manuscript preparation. Dr. Marius Lebret is the guarantor of the content of the manuscript, including the data and analysis.


## Funding

Funding was provided by Asten Santé.

## Conflict of interest

ML reports receiving personal fees from Air Liquide Medical Systems and non-financial support from NOMICS over the past 3 years, for work unrelated to the work presented here. ML is a part-time employee of Air Liquide Medical System (Med2lab).

YC, CM and GP work as consultants for Air Liquide Medical System.



**Fig. 1.** Dose-response relationship between oxygen (O<sub>2</sub>) flow rates delivered into the circuit, and fraction of inspired oxygen (FiO<sub>2</sub>) as a function of respiratory rate, presence/absence of non-intentional leakage, lung mechanics and leak port locations. Each Plot represents the mean FiO<sub>2</sub> of the 10 cycles analysed. For each FiO<sub>2</sub>, the standard deviation was systematically below 0.3 and is therefore not represented on the graph. Panel A displays the results for the conventional intentional leak model, in which the unintentional leakage port was placed in series between the exhalation valve (intentional leak port) and the ASL 5000. Panel B displays the results for the same model as in panel A, but with restrictive lung mechanics. Panel C displays the results for by-pass leakage ports using an antibacterial filter between the ASL 5000 and both intentional and non-intentional leak ports. The blue lines represent the relationship between O<sub>2</sub> flow and FiO<sub>2</sub> without unintentional leaks, and the green lines represent the relationship between O<sub>2</sub> flow and FiO<sub>2</sub> with unintentional leaks.

## Acknowledgements

Dr. Marius Lebret is the guarantor of the content of the manuscript, including the data and analysis.

The authors would like to thank Asten Santé for providing funding the trial. The authors would also like to thank Johanna Robertson for English editing.

## Bibliografía

1. Alviset S, Riller Q, Aboab J, Dilworth K, Billy P-A, Lombardi Y, et al. Continuous Positive Airway Pressure (CPAP) face-mask ventilation is an easy and cheap option to manage a massive influx of patients presenting acute respiratory failure during the SARS-CoV-2 outbreak: a retrospective cohort study. *PLOS ONE*. 2020;15:e0240645, <http://dx.doi.org/10.1371/journal.pone.0240645>.
2. Aliberti S, Radovanovic D, Billi F, Sotgiu G, Costanzo M, Pilocane T, et al. Helmet CPAP in patients with COVID-19 pneumonia: a multicentre cohort study. *Eur Respir J*. 2020;56, <http://dx.doi.org/10.1183/13993003.01935-2020>.
3. Bellani G, Grasselli G, Cecconi M, Antolini L, Borelli M, Giacomi F, et al. Noninvasive ventilatory support of covid-19 patients outside the intensive care units (WARD-COVID). *Ann Am Thorac Soc*. 2021, <http://dx.doi.org/10.1513/AnnalsATS.202008-1080OC>. Published Online First: 4 January 2021.
4. Landry SA, Mann DL, Djumas L, Messineo L, Terrill PI, Thomson LD, et al. Laboratory performance of oronasal CPAP and adapted snorkel masks to entrain oxygen and CPAP. *Respirology*. 2020;25:1309–12, <http://dx.doi.org/10.1111/resp.13922>.
5. Mechanical Ventilators for Non-Invasive Ventilation: Principles of Technology and Science – Nova Science Publishers [Internet]. [acceso 10 Feb 2021]. Disponible sur: <https://nam03.safelinks.protection.outlook.com/?url=https%3A%2F%2Fnovapublishers.com%2Fshop%2Fmechanical-ventilators-for-non-invasive-ventilation-principles-of-technology-and-science%2F&data=04%7C01%7CM.Valero%40elsevier.com%7C830dea6bd8124dd4e81908d915246021%7C9274ee3f94254109a27f9fb15c10675d%7C0%7C0%7C637564068248117156%7CUnknown%7CTWFpbGZsb3d8eyJWjoiMC4wLjAwMDAiLCjQjoiV2luMzliLCJB Til6lk1haWwiLCJVCI6Mn0%3D%7C1000&data=L90Hd%2FCMQBXKOMQuOrPb1HBy8m5UdqwnbWRxbOap2go%3D&reserved=0>
6. Rabec C, Gonzalez-Bermejo J, Mercy M, Grassion L, Pontier S, Patout M, et al. Respiratory support in patients with COVID-19 (outside intensive care unit). A position paper of the Respiratory Support and Chronic Care Group of the French Society of Respiratory Diseases. *Respir Med Res*. 2020;78:100768, <http://dx.doi.org/10.1016/j.resmer.2020.100768>.
7. ATS/ERS Statement on Respiratory Muscle Testing. *Am J Respir Crit Care Med*. 2002;166:518624.
8. Fresnel E, Muir J-F, Letellier C. Realistic human muscle pressure for driving a mechanical lung. *EPJ Nonlinear Biomed Phys*. 2014;2:7, <http://dx.doi.org/10.1140/epjnbp/s40366-014-0007-8>.
9. Murciano D, Aubier M, Bussi S, Derenne JP, Pariente R, Milic-Emili JAT Comparison of esophageal, tracheal, and mouth occlusion pressure in patients with chronic obstructive pulmonary disease during acute respiratory failure. *Am Rev Respir Dis*. 1982;126:837–41, <http://dx.doi.org/10.1164/arrd.1982.126.5.837>.
10. de Oca MM, Celli BR. Mouth occlusion pressure, CO<sub>2</sub> response and hypercapnia in severe chronic obstructive pulmonary disease. *Eur Respir J*. 1998;12:666–71, <http://dx.doi.org/10.1183/09031996.98.12030666>.
11. Baydur A. Respiratory muscle strength and control of ventilation in patients with neuromuscular disease. *Chest*. 1991;99:330–8, <http://dx.doi.org/10.1378/chest.99.2.330>.
12. Pandya A, Kaur NA, Sacher D, O'Corragain O, Salerno D, Desai P, et al. Ventilatory mechanics in early vs late intubation in a cohort of coronavirus disease 2019 patients with ARDS: a single center's experience. *CHEST*. 2021;159:653–6, <http://dx.doi.org/10.1016/j.chest.2020.08.2084>.
13. Hughes JMB. Review series: lung function made easy: assessing gas exchange. *Chron Respir Dis*. 2007;4:205–14, <http://dx.doi.org/10.1177/1479972307084446>.
14. Radovanovic D, Rizzi M, Pini S, Saad M, Chiumello DA, Santus PAT Helmet CPAP to treat acute hypoxic respiratory failure in patients with COVID-19: a management strategy proposal. *J Clin Med*. 2020;9:1191, <http://dx.doi.org/10.3390/jcm9041191>.
15. Vaschetto R, Barone-Adesi F, Racca F, Pissaia C, Maestrone C, Colombo D, et al. Outcomes of COVID-19 patients treated with continuous positive airway pressure outside the intensive care unit. *ERJ Open Res*. 2021;7, <http://dx.doi.org/10.1183/23120541.00541-2020>.
16. COVID-ICU Group on behalf of the REVA Network and the COVID-ICU Investigators. Clinical characteristics and day-90 outcomes of 4244 critically ill adults with COVID-19: a prospective cohort study. *Intensive Care Med*. 2021;47:60–73, <http://dx.doi.org/10.1007/s00134-020-06294-x>.

Marius Lebret <sup>a,\*</sup>, Emeline Fresnel <sup>b</sup>, Guillaume Prieur <sup>c,d</sup>, Jean Quieffin <sup>c</sup>, Johan Dupuis <sup>e</sup>, Bouchra Lamia <sup>c,f,g</sup>, Yann Combret <sup>c,d</sup>, Clément Medrinal <sup>c,h,i</sup>

<sup>a</sup> ERFPS, Rouen University Hospital, Physiotherapy School Rouen University, France

<sup>b</sup> KerNel Biomedical, Bois-Guillaume, France

<sup>c</sup> Pulmonology Department, Le Havre Hospital, Avenue Pierre Mendes France, Montivilliers, France

<sup>d</sup> Research and Clinical Experimentation Institute (IREC), Pulmonology, ORL and Dermatology, Louvain Catholic University, Brussels, Brussels Capital Region, Belgium

<sup>e</sup> Asten Santé, Isneauville, France

<sup>f</sup> UPRES EA3830 – GRHV, Institute for Research and Innovation in Biomedicine (IRIB), UNIROUEN, Rouen University Hospital, Rouen, Normandie, France

<sup>g</sup> Intensive Care Unit, Respiratory Section, Rouen University Hospital, Rouen, Normandie, France

<sup>h</sup> Université Paris-Saclay, UVSQ, Erphan, Versailles, France

<sup>i</sup> Saint Michel School of Physiotherapy, Paris-Saclay University, Paris, Île-de-France, France

\* Corresponding author.

E-mail address: [mariuslebret@gmail.com](mailto:mariuslebret@gmail.com) (M. Lebret).

<https://doi.org/10.1016/j.arbres.2021.04.018>

0300-2896/ © 2021 SEPAR. Published by Elsevier España, S.L.U. All rights reserved.