

Article title: Pragmatic study of a thromboprophylaxis algorithm in critically ill patients with sars-cov-2 infection

Journal name: Journal of Thrombosis and Thrombolysis

Author names: Maurizio Fattorutto. Yves Bouckaert. Jonathan Brauner. Stéphane Franck. Fabrice Bouton. Danielle Heuse. Charlotte Bouckaert. Arnaud Bruyneel.

Corresponding author: Maurizio Fattorutto, Department of Anesthesiology, Centre Hospitalier Universitaire Tivoli, Avenue Max Buset 34, 7100 La Louvière, Belgium. E-mail address: maurizio.fattorutto@chu-tivoli.be

Laboratory tests

Blood samples were taken from an arterial line and, to ensure clinically accurate sampling, a discard volume of three times the dead space was recommended as well as returning the dead space blood to the patient. Hemostasis tests included platelet count, prothrombin time (PT, expressed as International Normalized Ratio), activated partial thromboplastin time (APTT), fibrinogen level, D-dimer level, antithrombin III (ATIII) activity, anti-factor Xa (anti-FXa) activity, and ROTEM-EXTEM™. PT and APTT assays were measured on the ACL TOP 500 (Werfen, Instrumentation Laboratory, Bedford, MA, USA) with the reagents RECOMBIPLASTIN 2G (Werfen, Instrumentation Laboratory, Bedford, MA, USA) and ACTIN FS (SIEMENS HEALTHCARE DIAGNOSTICS PRODUCTS GMBH, Erlangen, Germany), respectively. Fibrinogen level, D-dimer level, anti-FXa, and ATIII activity were measured on the ACL TOP 500 using the HEMOSIL FIBRINOGEN C test (Werfen, Instrumentation Laboratory, Bedford, MA, USA), D-DIMER HS 500 (Werfen, Instrumentation Laboratory, Bedford, MA, USA), HEMOSIL LIQUID ANTI-FXa (Werfen, Instrumentation Laboratory, Bedford, MA, USA), and HEMOSIL LIQUID ANTITHROMBIN (Werfen, Instrumentation Laboratory, Bedford, MA, USA), respectively. The determination of platelet count was performed using the SYSMEX XN 2000 analyzer (SYSMEX, Hoellaart, Belgium) using the PLT FLUOROCCELL and WNR FLUOROCCELL reagents (SYSMEX, Kobe, Japan) according to the cytometric method.

The normal limits, as well as the analysis thresholds, were validated by our laboratory as follows: platelets ($150-400 \times 10^9/L$), INR (0.8-1.3), APTT (20.0-33.0 s), fibrinogen (2.0-4.0 g/L), ATIII (80%-120%), anti-FXa (≤ 2 IU-FXa/mL), D-dimer (upper limit of normal 500 ng/mL FEU ; upper limit of analysis 5000 ng/mL FEU), MCF-EXTEM™ (50-72 mm), C-reactive protein (0.4-12 mg/L), white blood cell count ($4.0-11.0 \times 10^3$ cells/ μ L), neutrophil count ($1.9-8.0 \times 10^3$ cells/ μ L), lymphocyte count ($0.9-5.2 \times 10^3$ cells/ μ L). D-dimer levels were measured on the ACL TOP 500 using D-DIMER HS 500 (Werfen, Instrumentation Laboratory, Bedford, MA, USA).