




Manogepix, the Active Moiety of the Investigational Agent Fosmanogepix, Demonstrates *In Vitro* Activity against Members of the *Fusarium oxysporum* and *Fusarium solani* Species Complexes

Hamid Badali,^a Hoja P. Patterson,^a Carmita J. Sanders,^a Barbara Mermella,^a Connie F. C. Gibas,^a Ashraf S. Ibrahim,^{b,c} Karen J. Shaw,^d  Nathan P. Wiederhold^a

^aThe University of Texas Health Science Center at San Antonio, San Antonio, Texas, USA

^bThe Lundquist Institute at Harbor-UCLA Medical Center, Torrance, California, USA

^cDavid Geffen School of Medicine at UCLA, Los Angeles, California, USA

^dHearts Consulting, San Diego, California, USA

ABSTRACT We evaluated the *in vitro* activity of manogepix against *Fusarium oxysporum* and *Fusarium solani* species complex (FOSC and FSSC, respectively) isolates per CLSI document M38 broth microdilution methods. Manogepix demonstrated activity against both FOSC (MEC [minimum effective concentration] range, ≤ 0.015 to $0.03 \mu\text{g/ml}$; MIC₅₀ range, ≤ 0.015 to $0.125 \mu\text{g/ml}$) and FSSC (MEC, $\leq 0.015 \mu\text{g/ml}$; MIC₅₀, ≤ 0.015 to $0.25 \mu\text{g/ml}$). Amphotericin B was also active (MIC, 0.25 to $4 \mu\text{g/ml}$), whereas the triazoles (MIC, 1 to $>16 \mu\text{g/ml}$) and micafungin (MEC, $\geq 8 \mu\text{g/ml}$) had limited activity.

KEYWORDS *Fusarium oxysporum*, *Fusarium solani*, manogepix, *in vitro* activity, minimum effective concentration, fosmanogepix, *Fusarium*, susceptibility

Fusarium spp. can cause a wide range of infections in humans, including keratitis and onychomycosis in immunocompetent individuals. *Fusarium* can also cause invasive and disseminated disease in immunocompromised hosts, including patients with neutropenia and hematological malignancies, hematopoietic stem cell transplant recipients, and those with severe T-cell deficiencies, and is associated with marked morbidity and mortality (1, 2). Human infections can be caused by members of eight different *Fusarium* spp. complexes, and those that are commonly seen as causing disease include members of the *Fusarium oxysporum* and *Fusarium solani* species complexes (FOSC and FSSC, respectively) (3). Although response rates have improved over the last 25 years with the use of voriconazole and amphotericin B lipid formulations, clinical outcomes are still suboptimal (4, 5), and both antifungals have unfavorable side-effect profiles. Manogepix (APX001A, Amplyx Pharmaceuticals, San Diego, CA; formerly E1210), the active moiety of the prodrug fosmanogepix (APX001), is a novel antifungal that targets inositol acyltransferase Gwt1, an enzyme in the glycosylphosphatidylinositol (GPI) anchor biosynthesis pathway. Inhibition of this enzyme prevents the maturation of GPI-anchored proteins (6). Both *in vitro* and *in vivo* activity have been demonstrated against *Candida* spp. (with the exception of *Candida krusei*), *Cryptococcus*, and *Coccidioides* spp., as well as molds, such as *Aspergillus* and *Scedosporium* spp. and *Rhizopus arrhizus* (6–18). *In vitro* activity has also been reported against a limited number of *Fusarium* isolates (8, 19), and this has translated to *in vivo* efficacy in murine models of disseminated fusariosis (18, 20). We evaluated the *in vitro* activity of manogepix against a larger number of clinical isolates of FOSC and FSSC.

Clinical isolates of FOSC ($n=49$) and FSSC ($n=19$) in the collection of the Fungus Testing Laboratory at the University of Texas Health Science Center at San Antonio were

Citation Badali H, Patterson HP, Sanders CJ, Mermella B, Gibas CFC, Ibrahim AS, Shaw KJ, Wiederhold NP. 2021. Manogepix, the active moiety of the investigational agent fosmanogepix, demonstrates *in vitro* activity against members of the *Fusarium oxysporum* and *Fusarium solani* species complexes. Antimicrob Agents Chemother 65:e02343-20. <https://doi.org/10.1128/AAC.02343-20>.

Copyright © 2021 American Society for Microbiology. All Rights Reserved.

Address correspondence to Nathan P. Wiederhold, wiederholdn@uthscsa.edu.

Received 5 November 2020

Returned for modification 17 December 2020

Accepted 6 March 2021

Accepted manuscript posted online 15 March 2021

Published 18 May 2021

TABLE 1 MEC and MIC ranges, MEC/MIC₅₀ and MEC/MIC₉₀, and GM MEC/MIC values for manogepix, amphotericin B, posaconazole, isavuconazole, voriconazole, and micafungin against *F. oxysporum* species complex and *F. solani* species complex isolates

	MGX ^a (μg/ml)			MIC (μg/ml) of ^b :				
Antifungal and parameter	MEC	MIC ₅₀	MIC	AMB	PSC	ISC	VRC	MFG MEC (μg/ml)
FOSC (n = 49)								
Range	≤0.015–0.03	≤0.015–0.125	>8	1–4	1 to >16	>16	4–16	≥8
MEC/MIC ₅₀	≤0.015	≤0.015	>8	2	4	16	8	>8
MEC/MIC ₉₀	≤0.015	0.125	>8	2	>16	>16	8	>8
GM MEC/MIC	≤0.015	0.021	>8	1.59	6.11	>16	6.94	>8
FSSC (n = 19)								
Range	≤0.015	≤0.015–0.25	<0.015 to >8	0.25–2	4 to >16	>16	2 to >16	≥8
MEC/MIC ₅₀	≤0.015	≤0.015	>8	1	>16	>16	>16	>8
MEC/MIC ₉₀	≤0.015	≤0.015	>8	2	>16	>16	>16	>8
GM MEC/MIC	≤0.015	0.017	>8	1.16	>16	>16	>16	>8

^aMGX, manogepix.^bAMB, amphotericin B; PSC, posaconazole; ISC, isavuconazole; VRC, voriconazole; MFG, micafungin.

used. Each isolate had previously been confirmed to the species complex level by combined phenotypic characteristics and DNA sequence analysis of the translation elongation factor 1- α (*TEF1 α*) and the RNA polymerase II second largest subunit (*RPB2*), as previously described (21). Antifungal susceptibility testing was performed by broth microdilution methods as described in CLSI document M38 (22), with RPMI 1640 (0.165 M MOPS [morpholinepropanesulfonic acid], pH 7.0, without bicarbonate) as the growth medium. Stock solutions of manogepix (Amplyx); amphotericin B, posaconazole, and voriconazole (Sigma); and isavuconazole and micafungin (Astellas) were prepared in DMSO (dimethyl sulfoxide), with further dilutions prepared in RPMI. For manogepix, activity was measured as the minimal effective concentration (MEC) and MICs at two endpoints: (i) an ~50% reduction in visual growth compared to the growth control, as allowed by the CLSI M38 standard for certain antifungals against filamentous fungi (i.e., fluconazole, ketoconazole, and 5-flucytosine) (22), and (ii) complete inhibition of growth, both of which were measured after 48 h of incubation at 35°C. The MEC is now the standard endpoint used to measure the *in vitro* activity of manogepix against filamentous fungi (6, 8). Similarly, the MEC endpoint was used for micafungin (22, 23). For amphotericin B, posaconazole, isavuconazole, and voriconazole, the MIC after 48 h of incubation was the endpoint used per CLSI recommendations (22). MEC/MIC ranges, MEC/MIC₅₀, MEC/MIC₉₀, and geometric mean (GM) MEC/MIC values were determined.

Manogepix demonstrated *in vitro* activity against FOSC and FSSC isolates when the MEC and MIC₅₀ endpoints were used (Table 1 and Tables S1 and S2 in the supplemental material). Against the FOSC isolates, the manogepix MEC range was ≤0.015 to 0.03 μg/ml, which was similar to the range with the MIC₅₀ endpoint (≤0.015 to 0.125 μg/ml). The GM MEC and MIC values were ≤0.015 and 0.021 μg/ml, respectively, and only 9 of the 49 FOSC isolates tested had a manogepix MIC₅₀ value higher than the lowest concentration tested (0.015 μg/ml) (see Fig. S1 in the supplemental material). Similar results were observed against FSSC isolates. Here, the MEC and MIC₅₀ ranges for manogepix were ≤0.015 and ≤0.015 to 0.25 μg/ml, respectively, and the GM MEC and MIC values were ≤0.015 and 0.017 μg/ml, respectively. In contrast, when the MIC endpoint was used, manogepix appeared to have reduced or no *in vitro* activity against *Fusarium* isolates at the highest concentration tested (Table 1; Supplemental Fig. S1).

Of the clinically available antifungals tested, amphotericin B demonstrated activity, with MIC ranges of 1 to 4 and 0.25 to 4 μg/ml against FOSC and FSSC, respectively. In contrast, the azoles demonstrated limited activity overall, with the MIC ranges for posaconazole, isavuconazole, and voriconazole falling between 1 and >16 μg/ml. In addition, the GM MIC values were higher than bloodstream concentrations that can consistently and safely be achieved with these antifungals. No activity was observed with micafungin at the highest concentration tested (MEC, ≥8 μg/ml against all isolates).

The *in vitro* activity of manogepix against *Fusarium* isolates observed in this study is consistent with that previously published by others. Against a limited number of *Fusarium* isolates ($n = 10$) from multiple species complexes, Pfaller et al. (19) reported the manogepix MECs to range between ≤ 0.008 and $8 \mu\text{g/ml}$ when determined by CLSI methods. In an earlier study by the same group using both CLSI and EUCAST methods that included a larger number of isolates ($n = 67$), the MECs against various *Fusarium* spp. ranged between 0.008 and $0.5 \mu\text{g/ml}$, and essential agreement between the two methods was reported to be 94% to 96.7% (8). In contrast, a recent study reported differences in manogepix activity against *Fusarium* isolates between the EUCAST and CLSI methods, with more potent activity when measured using the CLSI method (24). In the current study and others, the MEC was the endpoint chosen for manogepix *in vitro* activity because, similar to the echinocandins, this agent inhibits hyphal extension but does not completely inhibit growth (6). Here, we also report the MIC_{50} and MIC. The MIC_{50} endpoint demonstrated good agreement with the MEC value, although correlation analysis was not possible because many values for both endpoints were equal to the lowest concentration tested. The MIC is an inappropriate endpoint for manogepix against filamentous fungi, because this would suggest little to no *in vitro* activity against *Fusarium* spp. and other molds, contrary to *in vivo* efficacy model data. Previous studies have reported improved outcomes with manogepix treatment in murine models of invasive fusariosis, and these *in vivo* results are in agreement with the manogepix *in vitro* susceptibility result measured using either the MIC_{50} or the MEC as the susceptibility endpoint (18, 20).

In conclusion, manogepix demonstrated *in vitro* activity against FOSC and FSSC isolates. Clinical studies are ongoing to determine the efficacy and safety of fosmanogepix in patients with invasive fungal infections, and it is important for clinical laboratories to use the correct susceptibility endpoint for determination of *in vitro* activity for manogepix.

SUPPLEMENTAL MATERIAL

Supplemental material is available online only.

SUPPLEMENTAL FILE 1, PDF file, 0.3 MB.

ACKNOWLEDGMENTS

Manogepix powder was provided by Amplyx Pharmaceuticals, Inc., San Diego, CA.

N.P.W. has received research support to the UT Health San Antonio from Astellas, bioMérieux, Cepheid, and F2G and has served on advisory boards for Mayne Pharma. A.S.I. has received research support to the Lundquist Institute from Amplyx, Astellas, and Cidara and served on advisory boards for Amplyx, Astellas, Cidara, and Navigen. K.J.S. was previously an employee of Amplyx and is now an independent consultant at Hearts Consulting Group, LLC.

REFERENCES

1. Nucci F, Nouer SA, Capone D, Anaissie E, Nucci M. 2015. Fusariosis. *Semin Respir Crit Care Med* 36:706–714. <https://doi.org/10.1055/s-0035-1562897>.
2. Nucci M, Marr KA, Queiroz-Telles F, Martins CA, Trabasso P, Costa S, Voltarelli JC, Colombo AL, Imhof A, Pasquini R, Maiolino A, Souza CA, Anaissie E. 2004. Fusarium infection in hematopoietic stem cell transplant recipients. *Clin Infect Dis* 38:1237–1242. <https://doi.org/10.1086/383319>.
3. O'Donnell K, Sarver BA, Brandt M, Chang DC, Noble-Wang J, Park BJ, Sutton DA, Benjamin L, Lindsley M, Padhye A, Geiser DM, Ward TJ. 2007. Phylogenetic diversity and microsphere array-based genotyping of human pathogenic *Fusaria*, including isolates from the multistate contact lens-associated U.S. keratitis outbreaks of 2005 and 2006. *J Clin Microbiol* 45:2235–2248. <https://doi.org/10.1128/JCM.00533-07>.
4. Nucci M, Marr KA, Vehreschild MJGT, de Souza CA, Velasco E, Cappellano P, Carlesse F, Queiroz-Telles F, Sheppard DC, Kindo A, Cesaro S, Hamerschlak N, Solza C, Heinz WJ, Schaller M, Atalla A, Arikian-Akdagli S, Bertz H, Galvão Castro C, Herbrecht R, Hoenigl M, Härter G, Hermansen NEU, Josting A, Pagano L, Salles MJC, Mossad SB, Ogunc D, Pasqualotto AC, Araujo V, Troke PF, Lortholary O, Comely OA, Anaissie E. 2014. Improvement in the outcome of invasive fusariosis in the last decade. *Clin Microbiol Infect* 20:580–585. <https://doi.org/10.1111/1469-0691.12409>.
5. McCarthy MW, Katragkou A, Iosifidis E, Roilides E, Walsh TJ. 2018. Recent advances in the treatment of scedosporiosis and fusariosis. *J Fungi (Basel)* 4:73. <https://doi.org/10.3390/jof4020073>.
6. Miyazaki M, Horii T, Hata K, Watanabe NA, Nakamoto K, Tanaka K, Shirotori S, Murai N, Inoue S, Matsukura M, Abe S, Yoshimatsu K, Asada M. 2011. *In vitro* activity of E1210, a novel antifungal, against clinically important yeasts and molds. *Antimicrob Agents Chemother* 55:4652–4658. <https://doi.org/10.1128/AAC.00291-11>.
7. Pfaller MA, Hata K, Jones RN, Messer SA, Moet GJ, Castanheira M. 2011. *In vitro* activity of a novel broad-spectrum antifungal, E1210, tested against *Candida* spp. as determined by CLSI broth microdilution method. *Diagn Microbiol Infect Dis* 71:167–170. <https://doi.org/10.1016/j.diagmicrobio.2011.05.001>.
8. Castanheira M, Duncanson FP, Diekema DJ, Guarro J, Jones RN, Pfaller MA. 2012. Activities of E1210 and comparator agents tested by CLSI and

- EUCAST broth microdilution methods against *Fusarium* and *Scedosporium* species identified using molecular methods. *Antimicrob Agents Chemother* 56:352–357. <https://doi.org/10.1128/AAC.05414-11>.
9. Pfaller MA, Duncanson F, Messer SA, Moet GJ, Jones RN, Castanheira M. 2011. *In vitro* activity of a novel broad-spectrum antifungal, E1210, tested against *Aspergillus* spp. determined by CLSI and EUCAST broth microdilution methods. *Antimicrob Agents Chemother* 55:5155–5158. <https://doi.org/10.1128/AAC.00570-11>.
 10. Pfaller MA, Watanabe N, Castanheira M, Messer SA, Jones RN. 2011. Pre-clinical development of antifungal susceptibility test methods for the testing of the novel antifungal agent E1210 versus *Candida*: comparison of CLSI and European Committee on Antimicrobial Susceptibility Testing methods. *J Antimicrob Chemother* 66:2581–2584. <https://doi.org/10.1093/jac/dkr342>.
 11. Gebremariam T, Alkhazraji S, Alqarihi A, Jeon HH, Gu Y, Kapoor M, Shaw KJ, Ibrahim AS. 2018. APX001 is effective in the treatment of murine invasive pulmonary aspergillosis. *Antimicrob Agents Chemother* 63:e01713-18. <https://doi.org/10.1128/AAC.01713-18>.
 12. Viriyakosol S, Kapoor M, Okamoto S, Covell J, Soltow QA, Trzoss M, Shaw KJ, Fierer J. 2018. APX001 and other Gwt1 inhibitor prodrugs are effective in experimental *Coccidioides immitis* pneumonia. *Antimicrob Agents Chemother* 63:e01715-18. <https://doi.org/10.1128/AAC.01715-18>.
 13. Zhao M, Lepak AJ, Marchillo K, Vanhecker J, Sanchez H, Ambrose PG, Andes DR. 2019. APX001 pharmacokinetic/pharmacodynamic target determination against *Aspergillus fumigatus* in an *in vivo* model of invasive pulmonary aspergillosis. *Antimicrob Agents Chemother* 63:10. <https://doi.org/10.1128/AAC.02372-18>.
 14. Wiederhold NP, Najvar LK, Fothergill AW, McCarthy DI, Bocanegra R, Olivo M, Kirkpatrick WR, Everson MP, Duncanson FP, Patterson TF. 2015. The investigational agent E1210 is effective in treatment of experimental invasive candidiasis caused by resistant *Candida albicans*. *Antimicrob Agents Chemother* 59:690–692. <https://doi.org/10.1128/AAC.03944-14>.
 15. Shaw KJ, Schell WA, Covell J, Duboc G, Giamberardino C, Kapoor M, Moloney M, Soltow QA, Tenor JL, Toffaletti DL, Trzoss M, Webb P, Perfect JR. 2018. *In vitro* and *in vivo* evaluation of APX001A/APX001 and other Gwt1 inhibitors against *Cryptococcus*. *Antimicrob Agents Chemother* 62:10. <https://doi.org/10.1128/AAC.00523-18>.
 16. Alkhazraji S, Gebremariam T, Alqarihi A, Lin L, Gu Y, Shaw K, Ibrahim A. 2019. APX001 protects immunosuppressed mice from scedosporiosis, abstr 00734. Abstr 29th European Congress on Clinical Microbiology and Infectious Diseases, Amsterdam, Netherlands, 13 to 16 April 2019.
 17. Zhao Y, Lee MH, Paderu P, Lee A, Jimenez-Ortigosa C, Park S, Mansbach RS, Shaw KJ, Perlin DS. 2018. Significantly improved pharmacokinetics enhances *in vivo* efficacy of APX001 against echinocandin- and multidrug-resistant *Candida* isolates in a mouse model of invasive candidiasis. *Antimicrob Agents Chemother* 62:e00425-18. <https://doi.org/10.1128/AAC.00425-18>.
 18. Alkhazraji S, Gebremariam T, Alqarihi A, Gu Y, Mamouei Z, Singh S, Wiederhold NP, Shaw KJ, Ibrahim AS. 2019. Fosmanogepix (APX001) is effective in the treatment of immunocompromised mice infected with invasive pulmonary scedosporiosis or disseminated fusariosis. *Antimicrob Agents Chemother* 64:e01735-19. <https://doi.org/10.1128/AAC.01735-19>.
 19. Pfaller MA, Huband MD, Flamm RK, Bien PA, Castanheira M. 2019. *In vitro* activity of APX001A (manogepix) and comparator agents against 1,706 fungal isolates collected during an international surveillance program in 2017. *Antimicrob Agents Chemother* 63:e00840-19. <https://doi.org/10.1128/AAC.00840-19>.
 20. Hata K, Horii T, Miyazaki M, Watanabe NA, Okubo M, Sonoda J, Nakamoto K, Tanaka K, Shirotori S, Murai N, Inoue S, Matsukura M, Abe S, Yoshimatsu K, Asada M. 2011. Efficacy of oral E1210, a new broad-spectrum antifungal with a novel mechanism of action, in murine models of candidiasis, aspergillosis, and fusariosis. *Antimicrob Agents Chemother* 55:4543–4551. <https://doi.org/10.1128/AAC.00366-11>.
 21. Al-Hatmi AM, Van Den Ende AH, Stielow JB, Van Diepeningen AD, Seifert KA, McCormick W, Assabgui R, Grafenhan T, De Hoog GS, Levesque CA. 2016. Evaluation of two novel barcodes for species recognition of opportunistic pathogens in *Fusarium*. *Fungal Biol* 120:231–245. <https://doi.org/10.1016/j.funbio.2015.08.006>.
 22. Clinical and Laboratory Standards Institute. 2017. Reference method for broth dilution antifungal susceptibility testing of filamentous fungi; approved standard—3rd ed. CLSI document M38. Clinical and Laboratory Standards Institute, Wayne, PA.
 23. Arikan S, Lozano-Chiu M, Paetznick V, Rex JH. 2001. *In vitro* susceptibility testing methods for caspofungin against *Aspergillus* and *Fusarium* isolates. *Antimicrob Agents Chemother* 45:327–330. <https://doi.org/10.1128/AAC.45.1.327-330.2001>.
 24. Rivero-Menendez O, Cuenca-Estrella M, Alastruay-Izquierdo A. 2019. *In vitro* activity of APX001A against rare moulds using EUCAST and CLSI methodologies. *J Antimicrob Chemother* 74:1295–1299. <https://doi.org/10.1093/jac/dkz022>.